

AN ASSESSMENT OF
LASER VELOCIMETRY IN
HYPERSONIC FLOW

145752

p- 227

FINAL REPORT
NAS 2-12853

(NASA-CR-190895) AN ASSESSMENT OF
LASER VELOCIMETRY IN HYPERSONIC
FLOW Final Report (Complere)
227 p

N93-21556

Unclas

G3/02 0145752

COMPLERE INC.
P.O. Box 1697
PALO ALTO, CA
SEPTEMBER 1992

TABLE OF CONTENTS

| Section | Title | Page |
|---------|-----------------------|------|
| 1.0 | Abstract. | 2 |
| 2.0 | Introduction. | 2 |
| 3.0 | Experimental Details. | 3 |
| 4.0 | Test Results. | 4 |
| 5.0 | Concluding Remarks. | 6 |
| 6.0 | References. | 7 |

LIST OF TABLES

| Table | Title | Page |
|-------|-------------------|------|
| 1 | Seeder Operation. | 8 |
| 2 | Optics Details. | 10 |

LIST OF FIGURES

| Figure | Title | Page |
|--------|---|------|
| 1 | Schematic of the 3.5 ft. HWT LDV Seeder System. | 12 |
| 2 | Laser Doppler Velocimeter Signals. | 13 |
| 3 | Comparison of Probe and Laser Velocimeter Data. | 14 |
| 4 | Velocity Fluctuations across the Zero Pressure Gradient Boundary Layer. | 14 |
| 5 | Comparison of Axial Velocity Fluctuations. | 14 |
| 6 | Turbulent Velocity Cross-Correlation Coefficient. | 15 |
| 7 | Laser Velocimeter Measurements Across an Oblique Shock Wave. | 15 |
| 8 | Particle Response in Hypersonic Flow. | 16 |

SBIR RIGHTS NOTICE

Information contained in this report is subject to the SBIR Rights Notice (June 1987).

1.0 ABSTRACT

Although extensive progress has been made in computational fluid mechanics, reliable flight vehicle designs and modifications still cannot be made without recourse to extensive wind tunnel testing. Future progress in the computation of hypersonic flow fields is restricted by the need for a reliable mean flow and turbulence modeling data base which could be used to aid in the development of improved empirical models for use in numerical codes. Currently, there are few compressible flow measurements which could be used for this purpose. In this report, the results of experiments designed to assess the potential for laser velocimeter measurements of mean flow and turbulent fluctuations in hypersonic flow fields are presented. Details of a new laser velocimeter system which was designed and built for this test program are described.

2.0 INTRODUCTION

Current hypersonic flow field instrumentation is insufficient to meet present and future ground test requirements. Measurements are required to establish the basic physical mechanisms and turbulence models required for reliable prediction of transitional and turbulent hypersonic flow fields.

In recent years, experimental methods in lower speed regimes have also made significant advances due primarily to the availability of high power lasers. Their introduction has enabled the field of laser velocimetry to expand from low speed, small scale, closely controlled laboratory applications to the measurement of compressible flows in large scale wind tunnels (Ref. 1). The advent of the laser velocimeter allows us to measure velocity fluctuations directly in a linear, non-intrusive manner. Of particular value is the capability it offers to measure some of the compressible turbulent shear stresses, since this is an impractical task with hot wires (Ref. 2).

However, before laser velocimetry can be extended to hypersonic flow,

some basic questions must be addressed. The primary question is that of particle size requirement for reliable response combined with adequate Mie scattering. Practical assessments must therefore be made of flow seeding capability and the potential for laser velocimetry in hypersonic flows.

3.0 EXPERIMENTAL DETAILS

The laser velocimeter investigation was conducted in the NASA Ames 3.5 Ft. Hypersonic Wind Tunnel. In this facility, high-pressure air flows through a pebble bed heater and then through an open jet test section to lower pressure spheres. The tests were conducted at a nominal freestream Mach number of 7 and a freestream Reynolds number of 3 million per foot. The test model used in this study was a 10° cone-ogive-cylinder which was 79 inches long and 8 inches in diameter (Ref. 3). Measurements were made in the local freestream above the model and in the zero pressure gradient boundary layer flow on the cylindrical portion of the model. Measurements were also made across an oblique shock wave generated by the introduction of a 20 deg. flare installed 55 inches from the nose. A seed particle generator and injectors were designed and installed in the facility. The particles were injected through a thermocouple port into the plenum just upstream of the throat. A schematic of the seeding system and seeder operational procedures are given in Figure 1 and Table 1 respectively. Seed particles and seed mixtures detailed in Ref. 4 were used during the tests. The two component, forward scatter, fringe mode laser velocimeter system, which was used for the flow field measurements, utilized the 4880 and 5145 Angstrom lines of an argon-ion laser. Details of the optical system are presented in Table 2. Details of the traverse control and data acquisition systems are described in Ref. 5.

4.0 TEST RESULTS

Initially, measurements were confined to the local freestream until seeder mass flow rates and procedures were optimized for data rate and signal to noise ratio. Figure 2 shows examples of signal quality in *wind off* and *wind on* situations. Clearly, signal quality, visibility and fringe crossings were adequate in the freestream hypersonic flow. On occasion, data rates of more than 100,000 per second were measured. Figure 3 shows the mean boundary layer flow results along with the mean profile measurements which were obtained from previous conventional probe measurements (Ref. 6). Although, as expected, the signal to noise ratio decreased close to the wall, the good agreement between the two measurement methods confirms the seed particle response for mean velocity measurements in the zero pressure gradient boundary layer.

The results of a more stringent test of the particle response and the laser velocimeter measurements are shown in Figure 4 where the zero pressure gradient axial and vertical turbulence measurements are presented. These data show similarities in levels and trends to previous incompressible test results. The streamwise turbulence component has a pronounced maximum close to the wall whereas the vertical component, which is approximately half the axial value, is relatively flat in the wall region. These similarities are not altogether surprising since previous hot wire turbulence convection velocity measurements (Ref. 3) showed that the relative velocity between the disturbances and the local mean flow was always subsonic which allows the turbulent bursts to propagate as they would in an incompressible flow.

The axial component measurements are also compared with Klebanoff's incompressible results and previous hot wire hypersonic measurements in Figure 5. There is reasonably good agreement between the hypersonic laser velocimeter and incompressible hot wire data when normalized by the wall friction velocity. This is in contrast to previous hot wire compressible flow results, reviewed in Ref. 6, which show a monotonic decrease with increasing Mach number. However, all these past hot wire results have been evaluated

assuming zero pressure fluctuations which we would expect to become more important with increasing Mach number (Ref. 2). The turbulent velocity cross correlations are presented in Figure 6, which shows the variation of the turbulent velocity correlation coefficient across the boundary layer. The maximum value of approximately -0.4 is in close agreement with incompressible shear layer observations.

The most stringent test of particle response was made by perturbing the flow and measuring the particle velocity variation across an oblique shock wave and shear layer generated by the introduction of a 20 deg. flare. Unfortunately, these attempts to determine particle response were complicated by the proximity of the shock to the shear layer on the flare and by shock boundary/layer interaction instabilities. The results of a scan taken 2 inches above the model surface are presented in Figure 7 which shows the measured mean streamwise velocity and flow angularity distributions through the shock and shear layer region compared with conical flow theory and shadowgraph measurements of the shock location. The location of the measured mean velocity gradient is in good agreement with the shadowgraph shock location and the velocity change across the shock is comparable to conical flow predictions until the shear layer is encountered. The flow angularity measurements are consistent with conical flow predictions and the experimental flare angle. These comparisons indicate adequate particle response since some of the velocity and flow angularity gradient discrepancies across the shock are probably caused by small scale, time dependent oscillations of the shock wave about its mean location. Indeed, attempts to measure particle response across the 30 deg. shock wave were unsuccessful as the increased tunnel blockage led to excessive flow field instabilities and extensive shock motions.

The velocity probability density distributions, shown in Figure 8, are narrow in the freestream ahead of the shock where the turbulence level is low and wider in the more turbulent region within the shock layer. They are clearly bimodal in the region of the time averaged shock location. These bimodal distributions are of most interest as they give a clear indication of particle

response in hypersonic flow. The bimodal distributions shown in Figure 8 are due to shock wave fluctuations around its mean location. Thus, if the particles follow the flow, the two, bimodal peaks should be a measure of the velocity change across the shock. Since, when the instantaneous shock location is upstream of the focal volume, particles will register the lower velocity behind the shock and, when the focal volume is upstream of the instantaneous shock location, the higher freestream velocity will be recorded.

The shift from the dominant freestream peak ahead of the shock as the probe volume is traversed towards the model, is a measure of the probability of shock passage through the focal volume. The location of the most symmetrical bimodal distribution is the most likely, time-average shock location. Thus, from these velocity probability density distributions we can determine the particle velocity change across the shock and estimate the mean shock location above the plate. These results compare well with theoretical velocity change predictions and optical observations of the mean shock location.

These measurements can also be used to assess seed particle response and dynamics in hypersonic flow. Using the measured velocity change and calculated transit time through the shock wave region, seed particle response characteristics can be calculated. These calculations show that the measured particle response is equivalent to that of a 0.3 micron, specific gravity 1.0 sphere which undergoes a deceleration of almost seven million times the acceleration due to gravity; ie. 7 Mg. This size and acceleration is consistent with hypersonic modifications to the Stoke's drag law. Since, in hypersonic flow the particle drag coefficient is inversely proportional to the particle Reynolds and Knudsen numbers.

5.0 CONCLUDING REMARKS

Diagnostic tools are available to attempt the measurement of turbulent hypersonic flows, an area where comprehensive studies are lacking. Comparisons of new laser velocimeter turbulence measurements with previous hot wire results indicates that past data reduction assumptions can result in

significant measurement errors in hypersonic flows. It is felt that these new test results are the most convincing evidence to date of particle response in hypersonic flow. They clearly show that attempts to assess seed particle response must involve detailed studies of the velocity probability distributions. Particle response assessments inferred from conventional time-averaged velocity measurements could well be flawed by their failure to account for the hidden, adverse effects of large-scale, time-dependent mean flow variations which, on closer examination, may well manifest themselves in the velocity probability density distributions. Clearly, extensive work is still needed to establish a reliable data base for turbulence modeling and to define the reliable ranges of laser anemometer application.

6.0 REFERENCES

1. Owen, F. K., *Application of Laser Velocimetry to Unsteady Flows in Large Scale, High Speed Wind Tunnels*, ICIASF '83 Record, IEEE Publication 83CH1954-7, 1983.
2. Owen, F. K. and Fiore, A. W., *Turbulent Boundary Layer Measurement Techniques*, AFWAL-TR-86-3031, 1986.
3. Owen, F. K. and Horstman, C. C., *On the Structure of Hypersonic Turbulent Boundary Layers*, J. Fluid Mech., Vol. 53, pt. 4, p. 611, 1972.
4. Seegmiller, H. L., *Seeding Subsonic, Transonic and Supersonic Flows with 0.5 Micron Polystyrene Spheres*, NASA CP 2393, 1985.
5. Complane Inc., *The NASA Ames 3.5-Ft. Hypersonic Wind Tunnel Laser Velocimeter System*, Contract Report 92-0401, April 1992.
6. Owen, F. K., Horstman, C. C., and M. I. Kussoy, *Mean and Fluctuating Flow Measurements of a Fully Developed, Non-adiabatic Hypersonic Boundary Layer*, J. Fluid Mech., Vol. 70, pt. 4, p. 393, 1975.

Table 1. Seeder Operation.

LDV SEEDER OPERATION : LIQUID SEED

All manual and solenoid valves closed except #10 (drain). Liquid seed nozzle installed in heater port. Flex hose lines attached to liquid seed tank. Control valves backed off.

FILL TANK

- 1) Open MV-5 (manual).
- 2) Open MV-4 (manual).
- 3) Fill tank with liquid seed until it flows out of drain valve.
- 4) Close MV-5, MV-4, and drain #10.

PRESSURIZE LIQUID SEED TANK AND SEED LINE TO HEATER

- 5) Open valve to 3000 PSIA nitrogen bank.
- 6) Open MV-3 (manual) and CV-3 (solenoid).
- 7) Adjust FCV-3 (control valve) until pressure is above tunnel total pressure. Use pressure gage PT-3 to adjust. N.B. FCV-3 valve will be incrementally adjusted to achieve optimum seeding during run.
- 8) Open MV-2 (manual).
- 9) Open CV-2 (solenoid).

PRESSURIZE 3000 PSIA AIR LINE TO HEATER

- 10) Open valve to 3000 PSIA air line. Open CV-1 (solenoid).
- 11) Adjust FRV-1 and FCV-1 (control) so pressure read by gage PT-1 is 100 PSI above tunnel total pressure.
- 12) Close CV-1 (solenoid) and open MV-1 (manual).

START BLOWDOWN

IN CONTROL ROOM

- 13) When tunnel conditions are met, open CV-1 (solenoid).
- 14) Adjust control valves FCV-3 and FCV-1 so that seed pressure and air pressure are higher than tunnel operating pressure. Use LDV Data Acquisition System plus oscilloscope to determine optimum seeding.

END BLOWDOWN

IN CONTROL ROOM

- 15) Back off FCV-3 (control). Close CV-3 (solenoid).
- 16) Back off FCV-1 and FRV-1 (control). Close CV-1 (solenoid).

BY HEATER

- 17) Close MV-1 and MV-3 (manual).
- 18) Close valves to 3000 nitrogen and air lines.
- 19) Open MV-5 (manual) to relieve pressure in seeder tank to heater flue.
- 20) When flow stops, close MV-5, open drain #10, remove and/or replace seeder nozzle and seed filter (#6).

Table 1. Seeder Operation Continued.

LDV SEEDER OPERATION : DRY SEED

All manual and solenoid valves closed except #10 (drain). Flex hose lines attached to dry seed tank. Control valves backed off. Liquid seed filter screen removed to prevent clogging.

FILL TANK

- 1) Unscrew dry seed filler cap from dry seed tank.
- 2) Fill with dry seed, and replace filler cap.
- 3) Hook up to flex hose lines.
- 4) Close drain valve (#10).
- 5) Follow steps 5 to 18 of liquid seed operation then complete the following steps 6 & 7.
- 6) Open MV-6 (manual) to relieve pressure in seeder tank to heater flue.
- 7) When flow stops, close MV-6, open drain #10, remove and/or replace seeder nozzle.

Table 2. Optical Details.

| <u>Parameter</u> | <u>Symbol or Equation</u> | <u>Value</u> | <u>Units</u> |
|--------------------------------------|--|---------------|--------------|
| Wavelength | Lambda | 5145 | Å |
| Focal Length (transmitting lens) | Ft | 0.7620 | meters |
| Focal Length (receiving lens) | Fr | 0.7620 | meters |
| Focal Length (lens to fiber) | Ff | 0.7620 | meters |
| Aperture Diameter at Fiber | Df | 0.0006 | meters |
| Receiving Side Lens Diameter | Ld | 0.1524 | meters |
| Beam to Receiving Lens Gap | Gap | 0.0889 | meters |
| Beam Separation at Transmitting Lens | Bt | 0.007 938 | meters |
| Beam Diameter at Lens | Dl | 0.002 200 | meters |
| Convergence Full Angle | $Tf = 2*ATAN(Bt/2/Ft)$ | 0.597 | degrees |
| Convergence Half Angle | $Th = 1*ATAN(Bt/2/Ft)$ | 0.298 | degrees |
| Fringe Spacing | $X = Lambda / (2*SIN(Th))$ | 0.000 049 | meters |
| Number of Fringes in Probe Volume | $Npv = Dpv/X$ | 5 | --- |
| Off Axis Collecting Angle | $Tc = Th + ATAN(Gap/Fr) + ATAN(Ld/2/Fr)$ | 12.7 | degrees |
| Beam Diameter at Waist | $Dw = 4*Lambda*Ft/Pi/Dl$ | 0.000 227 | meters |
| Beam Diameter at Probe Volume | $Dpv = Dw/COS(Th)$ | 0.000 227 | meters |
| Length of Probe Volume | $Lpv = Dw/SIN(Th)$ | 0.043 565 | meters |
| Probe Volume Effective Length | $Vl = Df*Fr/Ff/SIN(Tc)$ | 0.002 737 | meters |
| Macrodyne Frequency | $Fmac = Fringes*Clock / (Bin*2^{(Range-0)})$ | --- | Hz |
| Bragg Frequency | Fbrag | 40 000 000 | Hz |
| Mixing Frequency | Fmix | 0 | Hz |
| Sign of Macrodyne Frequency | Smac | -1 | --- |
| Sign of Bragg Frequency | Sbrag | 1 | --- |
| Sign of Mixing Frequency | Smix | 1 | --- |
| Counter Clock Rate | Clock | 1 000 000 000 | Hz |
| Fringes Counted | Fringes | 8 | --- |
| Total Frequency | $Ftotal = Smac*Fmac + Sbrag*Fbrag + Smix*Fmix$ | --- | Hz |
| Velocity | $Velocity = X*Ftotal$ | --- | m/s |
| Velocity Resolution | Resolution | --- | m/s |
| Time in Focal Volume | $T = ABS(Dpv/Velocity)$ | --- | s |
| Number of Fringes Seen | $Ns = Fmac*T$ | --- | --- |
| Power at fiber exit (nominal) | Power | 0.3 | Watts |

Table 2. Optical Details Continued.

| <u>Parameter</u> | <u>Symbol or Equation</u> | <u>Value</u> | <u>Units</u> |
|--------------------------------------|---|---------------|--------------|
| Wavelength | λ | 4880 | Å |
| Focal Length (transmitting lens) | F_t | 0.7620 | meters |
| Focal Length (receiving lens) | F_r | 0.7620 | meters |
| Focal Length (lens to fiber) | F_f | 0.7620 | meters |
| Aperture Diameter at Fiber | D_f | 0.0006 | meters |
| Receiving Side Lens Diameter | L_d | 0.1524 | meters |
| Beam to Receiving Lens Gap | Gap | 0.0889 | meters |
| Beam Separation at Transmitting Lens | B_t | 0.007 938 | meters |
| Beam Diameter at Lens | D_l | 0.002 200 | meters |
| Convergence Full Angle | $T_f = 2 \cdot \text{ATAN}(B_t/2/F_t)$ | 0.597 | degrees |
| Convergence Half Angle | $T_h = 1 \cdot \text{ATAN}(B_t/2/F_t)$ | 0.298 | degrees |
| Fringe Spacing | $X = \lambda / (2 \cdot \text{SIN}(T_h))$ | 0.000 047 | meters |
| Number of Fringes in Probe Volume | $N_{pv} = D_{pv}/X$ | 5 | --- |
| Off Axis Collecting Angle | $T_c = T_h + \text{ATAN}(\text{Gap}/F_r) + \text{ATAN}(L_d/2/F_r)$ | 12.7 | degrees |
| Beam Diameter at Waist | $D_w = 4 \cdot \lambda \cdot F_t \cdot P_1/D_l$ | 0.000 215 | meters |
| Beam Diameter at Probe Volume | $D_{pv} = D_w / \text{COS}(T_h)$ | 0.000 215 | meters |
| Length of Probe Volume | $L_{pv} = D_w / \text{SIN}(T_h)$ | 0.041 321 | meters |
| Probe Volume Effective Length | $V_l = D_f \cdot F_r / F_f / \text{SIN}(T_c)$ | 0.002 737 | meters |
| Macrodyne Frequency | $F_{mac} = \text{Fringes} \cdot \text{Clock} / (\text{Bin} \cdot 2^{(\text{Range}-0)})$ | --- | Hz |
| Bragg Frequency | F_{brag} | 40 000 000 | Hz |
| Mixing Frequency | F_{mix} | 0 | Hz |
| Sign of Macrodyne Frequency | S_{mac} | -1 | --- |
| Sign of Bragg Frequency | S_{brag} | 1 | --- |
| Sign of Mixing Frequency | S_{mix} | 1 | --- |
| Counter Clock Rate | Clock | 1 000 000 000 | Hz |
| Fringes Counted | Fringes | 8 | --- |
| Total Frequency | $F_{total} = S_{mac} \cdot F_{mac} + S_{brag} \cdot F_{brag} + S_{mix} \cdot F_{mix}$ | --- | Hz |
| Velocity | $\text{Velocity} = X \cdot F_{total}$ | --- | m/s |
| Velocity Resolution | Resolution | --- | m/s |
| Time in Focal Volume | $T = \text{ABS}(D_{pv}/\text{Velocity})$ | --- | s |
| Number of Fringes Seen | $N_s = F_{mac} \cdot T$ | --- | --- |
| Power at fiber exit (nominal) | Power | 0.3 | Watts |






Figure 1. Schematic of the 3.5 ft. HWT LDV Seeder System.



a.) Wind off.



b.) Wind on.

Figure 2. Laser Doppler Velocimeter Signals
(Vertical Velocity Component).




Figure 3. Comparison of Probe and Laser Velocimeter Data.




Figure 4. Velocity Fluctuations across the Zero Pressure Gradient Boundary Layer.




Figure 5. Comparison of Axial Velocity Fluctuations.




Figure 6. Turbulent Velocity Cross-Correlation Coefficient.




Figure 7. Laser Velocimeter Measurements Across an Oblique Shock Wave.




Figure 8. Particle Response in Hypersonic Flow.

LASER VELOCIMETER DATA
ACQUISITION SYSTEM
TO
SUN SPARC STATION
S11W 16 BIT
PARALLEL INTERFACE
DOCUMENTATION

COMPLERE INC.
December 1992

SBIR RIGHTS NOTICE

Final Report Contract Number: NAS 2-12853.

Information contained in this documentation is subject to the **SBIR Rights Notice** (June 1987).

LVDAS to SUN 16 Bit Parallel Interface.

TABLE OF CONTENTS

| Section | Title | Page |
|----------------|--|-------------|
| 1.0 | Introduction. | 2 |
| 2.0 | Interface Cable. | 2 |
| 2.1 | LVDAS to SUN Interface Cable. | 3 |
| 2.2 | SUN High Density Connector. | 4 |
| 2.3 | LVDAS Circular Connector. | 5 |
| 2.4 | Handshake Timing Diagram for Transfer of Data from SUN to LVDAS. | 6 |
| 2.5 | Handshake Timing Diagram for Transfer of Data from LVDAS to SUN. | 7 |
| 3.0 | Data Acquisition Commands. | 8 |
| 3.1 | "CS" Command: Sample All Channels with Coincidence. | 8 |
| 3.2 | "SC" Command: Sample One Channel. | 12 |
| 3.3 | "SA" Command: Sample All Channels. | 15 |

LIST OF FIGURES

| Figure | Title | Page |
|---------------|--|-------------|
| 1. | LVDAS to SUN Interface Cable Schematic Drawing. | 3 |
| 2. | SUN High Density Connector Pin Locations. | 4 |
| 3. | LVDAS Circular Connector Pin Locations. | 5 |
| 4. | Handshake Timing Diagram for Transfer of Data from SUN to LVDAS. | 6 |
| 5. | Handshake Timing Diagram for Transfer of Data from LVDAS to SUN. | 7 |

SBIR RIGHTS NOTICE

Final Report Contract Number: NAS 2-12853.

Information contained in this documentation is subject to the SBIR Rights Notice (June 1987).

1.0 **INTRODUCTION**

This documentation describes the LVDAS to SUN interface as well as the data acquisition commands that control the flow of data between the two devices. Section 2 of this documentation provides a detailed schematic drawing of the interface cable, a drawing showing the SUN high density connector pin locations, a drawing showing the LVDAS circular connector pin locations, and timing diagrams for the transfer of data between the two devices.

Section 3 of this documentation provides a detailed description of the data acquisition commands sent to the LVDAS to control the flow of data between the two devices. The types of data, quantity of data, the data acquisition time, and the data formats are also described in Section 3.

The LVDAS can acquire up to 10,000 coincident data sets. Each data sets is composed of 10 words where the word size is 16 bits or 2 bytes. Therefore, the total buffer size is $10,000 \times 10 \times 2$ which is equal to 200,000 bytes.

2.0 **INTERFACE CABLE**

The interface between the Laser Velocimeter Data Acquisition System (LVDAS) and the SUN Sparc Station Computer is a 16 bit parallel general purpose input / output interface. The interface cable shown in Figure 1 consists of a standard cable (SUN EDT Part Number: CAB-A-25) with the terminating connectors on one end removed and replaced with a 55 pin circular connector (Cannon Part Number: MS3470W22-55P). The 80 pin high density connector attaches to the single slot interface card (SUN Part Number: S11W / S16D) within the SUN computer. The pin locations for the high density connector are shown in Figure 2. The 55 pin circular Cannon connector attaches to the Parallel I/O port at the back of the LVDAS. The pin locations of the circular connector are shown in Figure 3.

The timing diagram in Figure 4 shows the handshake sequence for transferring commands or data from the SUN computer to the LVDAS. The timing diagram in Figure 5 shows the handshake sequence for transferring data from the LVDAS to the SUN computer.

2.1 LVDAS to SUN Interface Cable.




Figure 1. LVDAS to SUN Interface Cable Schematic Drawing.

2.2 SUN High Density Connector.




Figure 2. SUN High Density Connector Pin Locations.

2.3 LVDAS Circular Connector.





Figure 3. LVDAS Circular Connector Pin Locations.

SUN --> LVDAS




Figure 4. Handshake Timing Diagram for Transfer of Data from SUN Computer to LVDA's.

2.5 Handshake Timing Diagram for Transfer of Data from LVDAS to SUN.




Figure 5. Handshake Timing Diagram for Transfer of Data from LVDAS to SUN Computer.

3.0 DATA ACQUISITION COMMANDS

This section provides a detailed description of the data acquisition commands and parameters sent to the LVDAS to control the flow of data between the two devices.

Commands sent to the LVDAS tell the LVDAS to perform a specific task.

Parameters sent to the LVDAS specify the conditions under which the data acquisition is to take place. Parameters, depending on the command, might include the desired data acquisition time, the desired coincidence time, the inter-arrival and coincidence time exponents, the desired coincidence channel selection, and the desired number of coincident data set samples.

The types of data returned to the computer, depending on the command, might include the inter-arrival time, coincidence time and status, valid data indication, digital frequency data, and digitized analog voltage data.

3.1 “CS” Command: Sample All Channels with Coincidence.

The “CS” command will acquire a finite number of coincident data sets over a finite acquisition time. The following commands, parameters, and data are transferred between the LVDAS and the computer:

| <u>WORD</u> | <u>SYMBOL</u> | <u>DESCRIPTION</u> | <u>DIRECTION</u> | <u>LENGTH</u> | <u>TYPE</u> |
|-------------|---------------|-----------------------------|-------------------|---------------|-------------|
| 1 | Cmnd | “CS” command | Computer to LVDAS | 1 | Command |
| 2&3 | DAtime | Desired acquisition time | Computer to LVDAS | 2 | Parameter |
| 4&5 | DCtime | Desired coincidence time | Computer to LVDAS | 2 | Parameter |
| 6 | ATexp | Inter-arrival time exponent | Computer to LVDAS | 1 | Parameter |
| 7 | CTexp | Coincidence time exponent | Computer to LVDAS | 1 | Parameter |
| 8 | Cmask | Coincidence mask | Computer to LVDAS | 1 | Parameter |
| 9 | DNsam | Desired number of samples | Computer to LVDAS | 1 | Parameter |
| 10 | RNsam | Realized number of samples | LVDAS to Computer | 1 | Parameter |
| 11 | Data0 | Inter-arrival time | LVDAS to Computer | 1 | Data |
| 12 | Data1 | Coincidence time | LVDAS to Computer | 1 | Data |
| 13 | Data2 | Coincidence status | LVDAS to Computer | 1 | Data |
| 14 | Data3 | Not used | LVDAS to Computer | 1 | Data |
| 15 | Data4 | Data valid | LVDAS to Computer | 1 | Data |
| 16 | Data5 | Digital channel #1 raw data | LVDAS to Computer | 1 | Data |
| 17 | Data6 | Digital channel #2 raw data | LVDAS to Computer | 1 | Data |
| 18 | Data7 | Digital channel #3 raw data | LVDAS to Computer | 1 | Data |
| 19 | Data8 | Analog channel #1 raw data | LVDAS to Computer | 1 | Data |
| 20 | Data9 | Not used | LVDAS to Computer | 1 | Data |

The data words 11 through 20 above are repeated **RNsam** times.

The range (min & max), units, and format for the above commands, parameters, and data are shown below:

| <u>SYMBOL</u> | <u>MIN</u> | <u>MAX</u> | <u>UNITS</u> | <u>FORMAT</u> |
|---------------|------------|---------------|--------------|-------------------------|
| Cmnd | "CS" | - | none | 2 ASCII Bytes |
| DAtime | 0 | 4,294,967,295 | 100ns | Unsigned 32 bit integer |
| DCtime | 0 | 4,294,967,295 | 100ns | Unsigned 32 bit integer |
| ATexp | 0 | 16 | none | Unsigned 16 bit integer |
| CTexp | 0 | 16 | none | Unsigned 16 bit integer |
| Cmask | 1 | 7 | none | Unsigned 16 bit integer |
| DNsam | 0 | 10,000 | none | Unsigned 16 bit integer |
| RNsam | 0 | 10,000 | none | Unsigned 16 bit integer |
| Data0 | 0 | 65,535 | ns* | Unsigned 16 bit integer |
| Data1 | 0 | 65,535 | ns* | Unsigned 16 bit integer |
| Data2 | 0 | 15 | none | Unsigned 16 bit integer |
| Data3 | 0 | 0 | none | Unsigned 16 bit integer |
| Data4 | 1 | 1 | none | Unsigned 16 bit integer |
| Data5 | 0 | 65,535 | Hz* | Unsigned 16 bit integer |
| Data6 | 0 | 65,535 | Hz* | Unsigned 16 bit integer |
| Data7 | 0 | 65,535 | Hz* | Unsigned 16 bit integer |
| Data8 | -32,768 | 32,767 | volts* | Signed 16 bit integer |
| Data9 | 65,535 | 65,535 | none | Signed 16 bit integer |

The data words whose units are noted by a * are encoded. Their values in the specified units can be calculated using the raw encoded data.

The command word **Cmnd** (=CS) tells the LVDAS that the computer will want to acquire laser velocimeter data with coincidence. The maximum desired acquisition time **DAtime** and desired coincidence time **DCtime** are specified in 100 ns counts and each is sent to the LVDAS as two 16 bit words concatenated into one 32 bit unsigned integer. For example, counts of 50000, 10000000, and 600000000 would yield times of 5 milliseconds, 1 second, and 1 minute respectively.

The inter-arrival time exponent **ATexp** and coincidence time exponent **CTexp** are used to modify the inter-arrival and coincidence times. The LVDAS measures these times with a resolution of 100 ns and an unsigned integer data size of 32 bits. The 32 bit inter-arrival time is shifted right by the number of bits specified by the inter-arrival time exponent **ATexp**. The 32 bit coincidence time is shifted right by the number of bits specified by the coincidence time exponent **CTexp**. The resulting 16 bit words are later sent to the computer.

The coincidence mask **Cmask** determines the desired coincidence criterion. The least significant three bits individually select the digital channels #1, #2, and #3 for coincidence. Valid coincidence masks are as follows:

| Coincidence Mask (decimal) | Coincidence Mask (binary) | Channel #3 (selected) | Channel #2 (selected) | Channel #1 (selected) |
|----------------------------------|---------------------------------|-----------------------------|-----------------------------|-----------------------------|
| 0 | 0000 0000 0000 0000 | NO | NO | NO |
| 1 | 0000 0000 0000 0001 | NO | NO | YES |
| 2 | 0000 0000 0000 0010 | NO | YES | NO |
| 3 | 0000 0000 0000 0011 | NO | YES | YES |
| 4 | 0000 0000 0000 0100 | YES | NO | NO |
| 5 | 0000 0000 0000 0101 | YES | NO | YES |
| 6 | 0000 0000 0000 0110 | YES | YES | NO |
| 7 | 0000 0000 0000 0111 | YES | YES | YES |

The desired number of samples **DNsam** specifies the number of coincident data sets to be acquired within the previously specified desired data acquisition time **DAtime**. The data acquisition commences when **DNsam** is received by the LVDAS.

The data acquisition terminates when one of two conditions occur. The first terminating condition is that **DNsam** coincident data sets are realized before the allocated data acquisition time **DAtime** expires. In this case, the desired **DNsam** and realized **RNsam** number of samples are the same. The second terminating condition is that **DNsam** coincident data sets are not realized before the allocated data acquisition time **DAtime** expires. In this case, the realized number of samples **RNsam** may be less than the desired number of samples **DNsam**. In both terminating conditions, this value (**RNsam**) is then sent from the LVDAS to the computer to indicate data acquisition completion and to also indicate the size of the data array to be subsequently transferred to the computer.

Each coincident data set consists of ten 16 bit words. **RNsam** indicates the number of acquired coincident data sets. Therefore, there will be $10 \times \text{RNsam}$ words sent from the LVDAS to the computer. The computer's data array should be dimensioned accordingly. The 10 words will include the inter-arrival and coincidence times, the coincidence status and data valid words, as well as the digital and analog raw data words.

The inter-arrival time **Data0** and coincidence time **Data1** raw data words can be converted to the actual inter-arrival time **IAtime** and realized coincidence time **RCtime** in seconds using the following equations:

$$\text{IAtime} = \text{Data0} * (2^{\text{ATexp}}) / (10^7) \quad \text{seconds}$$

$$\text{RCtime} = \text{Data1} * (2^{\text{CTexp}}) / (10^7) \quad \text{seconds}$$

The coincidence status **Status** and data valid **Valid** words, **Data2** and **Data3** respectively, indicate the channels that have new data in the data set and the validity of the data. If **Valid**=0 then the data set does not contain valid data. If **Valid**=1 then the data set does contain valid data. The least significant four **Status** bits individually indicate whether

of not new data has been acquired on the digital and analog channels:

| Status Word (decimal) | Status Word (binary) | Analog Ch #1 (new) | Digital Ch #3 (new) | Digital Ch #2 (new) | Digital Ch #1 (new) |
|--------------------------|-------------------------|-----------------------|------------------------|------------------------|------------------------|
| 0 | 0000 0000 0000 0000 | NO | NO | NO | NO |
| 1 | 0000 0000 0000 0001 | NO | NO | NO | YES |
| 2 | 0000 0000 0000 0010 | NO | NO | YES | NO |
| 3 | 0000 0000 0000 0011 | NO | NO | YES | YES |
| 4 | 0000 0000 0000 0100 | NO | YES | NO | NO |
| 5 | 0000 0000 0000 0101 | NO | YES | NO | YES |
| 6 | 0000 0000 0000 0110 | NO | YES | YES | NO |
| 7 | 0000 0000 0000 0111 | NO | YES | YES | YES |
| 8 | 0000 0000 0000 1000 | YES | NO | NO | NO |
| 9 | 0000 0000 0000 1001 | YES | NO | NO | YES |
| 10 | 0000 0000 0000 1010 | YES | NO | YES | NO |
| 11 | 0000 0000 0000 1011 | YES | NO | YES | YES |
| 12 | 0000 0000 0000 1100 | YES | YES | NO | NO |
| 13 | 0000 0000 0000 1101 | YES | YES | NO | YES |
| 14 | 0000 0000 0000 1110 | YES | YES | YES | NO |
| 15 | 0000 0000 0000 1111 | YES | YES | YES | YES |

The raw data words **Data5**, **Data6**, and **Data7** contain the digital data from the Macrodyne laser velocimeter counter signal processors. These digital data can be converted into frequencies using the following equations:

Mantissa1 = Bits 0 to 9 of **Data5**
Mantissa2 = Bits 0 to 9 of **Data6**
Mantissa3 = Bits 0 to 9 of **Data7**

Exponent1 = Bits 10 to 13 of **Data5**
Exponent2 = Bits 10 to 13 of **Data6**
Exponent3 = Bits 10 to 13 of **Data7**

Fringes1 : If bit 14 of **Data5**=0 then **Fringes1**=16 else **Fringes1**=8
Fringes2 : If bit 14 of **Data6**=0 then **Fringes2**=16 else **Fringes2**=8
Fringes3 : If bit 14 of **Data7**=0 then **Fringes3**=16 else **Fringes3**=8

Period1 = **Mantissa1** * (2^{Exponent1}) / (10⁹) (seconds)
Period2 = **Mantissa2** * (2^{Exponent2}) / (10⁹) (seconds)
Period3 = **Mantissa3** * (2^{Exponent3}) / (10⁹) (seconds)

Frequency1 = **Fringes1** / **Period1** (Hz)
Frequency2 = **Fringes2** / **Period2** (Hz)
Frequency3 = **Fringes3** / **Period3** (Hz)

The following equation is used to convert the raw data word **Data8** into a voltage:

Analog = **Data8** * 5 / 32768 (volts)

3.2 “SC” Command: Sample One Channel.

The “SC” command will acquire 1000 data samples from one channel. The following commands, parameters, and data are transferred between the LVDAS and the computer:

| <u>WORD</u> | <u>SYMBOL</u> | <u>DESCRIPTION</u> | <u>DIRECTION</u> | <u>LENGTH</u> | <u>TYPE</u> |
|-------------|----------------|--------------------|-------------------|---------------|-------------|
| 1 | Cmnd1 | “DT” command | Computer to LVDAS | 1 | Command |
| 2 | Cmnd2 | “SC” command | Computer to LVDAS | 1 | Command |
| 3 | Channel | Channel Number | Computer to LVDAS | 1 | Parameter |
| 4 | Cmnd3 | “ET” command | Computer to LVDAS | 1 | Command |
| 5 | Cmnd4 | “RM” command | Computer to LVDAS | 1 | Command |
| 6&7 | First | Memory location | Computer to LVDAS | 2 | Parameter |
| 8&9 | Last | Memory location | Computer to LVDAS | 2 | Parameter |
| 10&11* | Data1 | Inter-arrival time | LVDAS to Computer | 2 | Data |
| 12* | Data2 | Channel number | LVDAS to Computer | 1 | Data |
| 13* | Data3 | Channel data | LVDAS to Computer | 1 | Data |

The data words 10 through 13 are repeated 1000 times.

The range (min & max), units, and format for the above commands, parameters, and data are shown below:

| <u>SYMBOL</u> | <u>MIN</u> | <u>MAX</u> | <u>UNITS</u> | <u>FORMAT</u> |
|----------------|--------------|---------------|--------------|-------------------------|
| Cmnd1 | “DT” | - | none | 2 ASCII Bytes |
| Cmnd2 | “SC” | - | none | 2 ASCII Bytes |
| Channel | 1 | 7 | none | Unsigned 16 bit integer |
| Cmnd3 | “ET” | - | none | 2 ASCII Bytes |
| Cmnd4 | “RM” | - | none | 2 ASCII Bytes |
| First | 08F00000 hex | - | none | Unsigned 32 bit integer |
| Last | 08F01F3F hex | - | none | Unsigned 32 bit integer |
| Data1 | 0 | 4,294,967,295 | 100ns | Unsigned 32 bit integer |
| Data2 | 0 | 6 | none | Unsigned 16 bit integer |
| Data3 | see text | see text | see text | see text |

The first command word **Cmnd1** (=DT) tells the LVDAS to disable internal timers which temporarily stops updating of the front panel displays. Sending out **Cmnd1** is optional. The second command word **Cmnd2** (=SC) tells the LVDAS that the computer will want to acquire laser velocimeter or analog data on one channel only. The data word **Channel** specifies the channel number for which data will be acquired. Valid channel

numbers are as follows:

| Channel Number | Channel Description | Generates Data Word | Generates Inter-Arrival Time Words |
|----------------|--------------------------|---------------------|------------------------------------|
| 1 | Digital channel #1 | YES | YES |
| 2 | Digital channel #2 | YES | YES |
| 3 | Digital channel #3 | YES | YES |
| 4 | Analog channel #1 | YES | YES |
| 6 | External trigger timer | NO | YES |
| 7 | Inter-arrival time timer | NO | YES |

The data acquisition commences when **Channel** is received by the LVDAS. The third command word **Cmnd3** (=ET) tells the LVDAS to enable internal timers which activates the updating of the front panel displays. After 1000 data samples have been acquired, then the third command word **Cmnd3** will be executed. The computer can now read back the data from the buffer's memory. Reading memory is initiated by sending the forth command word **Cmnd4** (=RM) and the two memory buffer parameters **First** and **Last**.

The LVDAS will respond by sending 4 words of data per sample to the computer. The first 2 words in **Data1** contain the inter-arrival time **IAtime**; the third word in **Data2** contains the channel number **Channel**; and the forth word in **Data3** contains the channel's data.

The two inter-arrival time raw data words in **Data1** can be converted to the actual inter-arrival time **IAtime** in seconds using the following equation:

$$\text{IAtime} = \text{Data1} / (10^7) \quad \text{seconds}$$

The type of data, its range (min & max), units, and format returned in **Data3** depend on which channel, specified by **Channel**, the data was acquired on. (Note: The LVDAS will return channel numbers minus one: 0..6; not 1..7).

| Channel Number | Channel Description | Generates Data Word | Generates Inter-Arrival Time Words |
|----------------|--------------------------|---------------------|------------------------------------|
| 0 | Digital Channel #1 | YES | YES |
| 1 | Digital Channel #2 | YES | YES |
| 2 | Digital Channel #3 | YES | YES |
| 3 | Analog Channel #1 | YES | YES |
| 5 | External Trigger Timer | NO | YES |
| 6 | Inter-Arrival Time Timer | NO | YES |

| <u>CHANNEL</u> | <u>MIN</u> | <u>MAX</u> | <u>UNITS</u> | <u>FORMAT</u> |
|----------------|------------|------------|--------------|-------------------------|
| 0 | 0 | 65,535 | Hz* | Unsigned 16 bit integer |
| 1 | 0 | 65,535 | Hz* | Unsigned 16 bit integer |
| 2 | 0 | 65,535 | Hz* | Unsigned 16 bit integer |
| 3 | -32,768 | 32,767 | volts* | Signed 16 bit integer |

The data words whose units are noted by a * are encoded. Their values in the specified units can be calculated using the raw encoded data.

If the data was acquired on channels 0 through 2, then the following equations should be used to convert the digital data from the Macrodyne laser velocimeter counter signal processors into frequencies:

$$\begin{aligned}
 \text{Mantissa} &= \text{Bits 0 to 9 of Data3} \\
 \text{Exponent} &= \text{Bits 10 to 13 of Data3} \\
 \text{Fringes} &: \text{ If bit 14 of Data3}=0 \text{ then Fringes}=16 \text{ else Fringes}=8 \\
 \text{Period} &= \text{Mantissa} * (2^{\text{Exponent}}) / (10^9) \quad (\text{seconds}) \\
 \text{Frequency} &= \text{Fringes} / \text{Period} \quad (\text{Hz})
 \end{aligned}$$

If the data was acquired on channel 3, then the following equation should be used to convert the raw data word into a voltage:

$$\text{Analog} = \text{Data3} * 5 / 32768 \quad (\text{volts})$$

Channels 4 through 6 produce an inter-arrival time but do not generate any meaningful data. Their data is ignored.

3.3 “SA” Command: Sample All Channel.

The “SA” command will acquire 1000 data samples from all channels. The 1000 samples will be spread out over all enabled channels. Channels with higher data rates will generate more samples than channels with lower data rates. The sum total of all samples will be 1000 samples. The following commands, parameters, and data are transferred between the LVDAS and the computer:

| <u>WORD</u> | <u>SYMBOL</u> | <u>DESCRIPTION</u> | <u>DIRECTION</u> | <u>LENGTH</u> | <u>TYPE</u> |
|-------------|---------------|--------------------|-------------------|---------------|-------------|
| 1 | Cmnd1 | “DT” command | Computer to LVDAS | 1 | Command |
| 2 | Cmnd2 | “SA” command | Computer to LVDAS | 1 | Command |
| 3 | Mask | Channel Mask | Computer to LVDAS | 1 | Parameter |
| 4 | Cmnd3 | “ET” command | Computer to LVDAS | 1 | Command |
| 5 | Cmnd4 | “RM” command | Computer to LVDAS | 1 | Command |
| 6&7 | First | Memory location | Computer to LVDAS | 2 | Parameter |
| 8&9 | Last | Memory location | Computer to LVDAS | 2 | Parameter |
| 10&11* | Data1 | Inter-arrival time | LVDAS to Computer | 2 | Data |
| 12* | Data2 | Channel number | LVDAS to Computer | 1 | Data |
| 13* | Data3 | Channel data | LVDAS to Computer | 1 | Data |

The data words 10 through 13 are repeated 1000 times.

The range (min & max), units, and format for the above commands, parameters, and data are shown below:

| <u>SYMBOL</u> | <u>MIN</u> | <u>MAX</u> | <u>UNITS</u> | <u>FORMAT</u> |
|---------------|--------------|---------------|--------------|-------------------------|
| Cmnd1 | “DT” | - | none | 2 ASCII Bytes |
| Cmnd2 | “SA” | - | none | 2 ASCII Bytes |
| Mask | 1 | 127 | none | Unsigned 16 bit integer |
| Cmnd3 | “ET” | - | none | 2 ASCII Bytes |
| Cmnd4 | “RM” | - | none | 2 ASCII Bytes |
| First | 08F00000 hex | - | none | Unsigned 32 bit integer |
| Last | 08F01F3F hex | - | none | Unsigned 32 bit integer |
| Data1 | 0 | 4,294,967,295 | 100ns | Unsigned 32 bit integer |
| Data2 | 0 | 6 | none | Unsigned 16 bit integer |
| Data3 | see text | see text | see text | see text |

The first command word **Cmnd1** (=DT) tells the LVDAS to disable internal timers which temporarily stops updating of the front panel displays. Sending out **Cmnd1** is

optional. The second command word **Cmnd2** (=SA) tells the LVDAS that the computer will want to acquire laser velocimeter and/or analog data on all channels. The data word **Mask** specifies the channel numbers for which data will be acquired. Each bit in the **Mask** enable data acquisition on the relevant channels. Channels whose **Mask** bit equals zero will be ignored. Channels whose **Mask** bit equals one will be serviced each time data becomes available.

| Channel Number | Mask Bit | Channel Description | Generates Data Word | Generates Inter-Arrival Time Words |
|----------------|----------|--------------------------|---------------------|------------------------------------|
| 1 | 0 | Digital channel #1 | YES | YES |
| 2 | 1 | Digital channel #2 | YES | YES |
| 3 | 2 | Digital channel #3 | YES | YES |
| 4 | 3 | Analog channel #1 | YES | YES |
| 6 | 5 | External trigger timer | NO | YES |
| 7 | 6 | Inter-arrival time timer | NO | YES |

The data acquisition commences when **Mask** is received by the LVDAS. The third command word **Cmnd3** (=ET) tells the LVDAS to enable internal timers which activates the updating of the front panel displays. After 1000 data samples have been acquired, then the third command word **Cmnd3** will be executed. The computer can now read back the data from the buffer's memory. Reading memory is initiated by sending the forth command word **Cmnd4** (=RM) and the two memory buffer parameters **First** and **Last**.

The LVDAS will respond by sending 4 words of data per sample to the computer. The first 2 words in **Data1** contain the channel inter-arrival time **IAtime**; the third word in **Data2** contains the channel number **Channel**; and the forth word in **Data3** contains the channel's data.

The channel inter-arrival times are the inter-arrival times of data samples acquired on the same channel. The average channel inter-arrival time for all samples acquired on a specific channel yield that channels data rate (rate=1/period). The two channel inter-arrival time raw data words in **Data1** can be converted to the actual inter-arrival time **IAtime** in seconds using the following equation:

$$\text{IAtime} = \text{Data1} / (10^7) \quad \text{seconds}$$

The type of data, its range (min & max), units, and format returned in **Data3** depend on which channel, specified by **Channel**, the data was acquired on. (Note: The LVDAS will return channel numbers minus one: 0..6; not 1..7).

| Channel Number | Channel Description | Generates Data Word | Generates Inter-Arrival Time Words |
|----------------|--------------------------|---------------------|------------------------------------|
| 0 | Digital Channel #1 | YES | YES |
| 1 | Digital Channel #2 | YES | YES |
| 2 | Digital Channel #3 | YES | YES |
| 3 | Analog Channel #1 | YES | YES |
| 5 | External Trigger Timer | NO | YES |
| 6 | Inter-Arrival Time Timer | NO | YES |

| <u>CHANNEL</u> | <u>MIN</u> | <u>MAX</u> | <u>UNITS</u> | <u>FORMAT</u> |
|----------------|------------|------------|--------------|-------------------------|
| 0 | 0 | 65,535 | Hz* | Unsigned 16 bit integer |
| 1 | 0 | 65,535 | Hz* | Unsigned 16 bit integer |
| 2 | 0 | 65,535 | Hz* | Unsigned 16 bit integer |
| 3 | -32,768 | 32,767 | volts* | Signed 16 bit integer |

The data words whose units are noted by a * are encoded. Their values in the specified units can be calculated using the raw encoded data.

If the data was acquired on channels 0 through 2, then the following equations should be used to convert the digital data from the Macrodyne laser velocimeter counter signal processors into frequencies:

$$\begin{aligned}
 \text{Mantissa} &= \text{Bits 0 to 9 of Data3} \\
 \text{Exponent} &= \text{Bits 10 to 13 of Data3} \\
 \text{Fringes} &: \text{ If bit 14 of Data3}=0 \text{ then Fringes}=16 \text{ else Fringes}=8 \\
 \text{Period} &= \text{Mantissa} * (2^{\text{Exponent}}) / (10^9) \quad (\text{seconds}) \\
 \text{Frequency} &= \text{Fringes} / \text{Period} \quad (\text{Hz})
 \end{aligned}$$

If the data was acquired on channel 3, then the following equation should be used to convert the raw data word into a voltage:

$$\text{Analog} = \text{Data3} * 5 / 32768 \quad (\text{volts})$$

Channels 4 through 6 produce an inter-arrival time but do not generate any meaningful data. Their data is ignored.

**THE NASA AMES
3.5-FT. HYPERSONIC WIND TUNNEL
LASER VELOCIMETER SYSTEM**

**CONTRACT REPORT
92-0401**

**COMPLERE INC.
P. O. BOX 1697
PALO ALTO, CA
APRIL 1992**

TABLE OF CONTENTS

| Title | Chapter |
|--|----------------|
| 3.5 Ft. HWT Optical System. | 1 |
| Laser Velocimeter Data Acquisition System (LVDAS). | 2 |
| Data Acquisition Computer Hardware and Software. | 3 |
| Traverse Control System (TCS8). | 4 |

LIST OF APPENDICES

| Title | Appendix |
|--|-----------------|
| Original Software Code Listing. | A |
| Revised Software Code Listing. | B |
| Data Reduction and Coordinate System Transformation Equations. | C |

SBIR RIGHTS NOTICE

Information contained in this report is subject to the SBIR Rights Notice (June 1987).

CHAPTER 1

3.5 FT HWT OPTICAL SYSTEM.

CHAPTER 1

3.5 FT HWT Optical System

TABLE OF CONTENTS

| Section | Title | Page |
|----------------|--|-------------|
| 1.0 | The Laser Doppler Velocimeter--Introduction. | 3 |
| 2.0 | Transmitting Color Separation System--Alignment. | 6 |
| 2.1 | Optics Enclosure. | 6 |
| 2.2 | Outside Steering Mirrors. | 7 |
| 2.3 | Polarization Rotator. | 7 |
| 2.4 | Acousto-Optic Modulator (Bragg Cell). | 7 |
| 2.5 | Dispersion Prisms. | 8 |
| 2.6 | Polarization Rotator. | 8 |
| 2.7 | Inside Steering Mirrors. | 8 |
| 2.8 | Iris Diaphragms. | 8 |
| 2.9 | Separation Prisms. | 9 |
| 2.10 | Final Steering Prisms. | 9 |
| 3.0 | Fiber Optic Link. | 9 |
| 3.1 | Laser to Fiber Couplers. | 9 |
| 3.2 | Laser to Fiber Coupler Alignment. | 9 |
| 3.3 | Preliminary Z-Axis Adjustment. | 12 |
| 3.4 | Polarization Axis Adjustment. | 12 |
| 3.5 | Plenum Feedthrough for Optical Fibers. | 14 |
| 4.0 | Fiber Collimators and Transmitting Optics. | 15 |
| 5.0 | Collecting Optics. | 16 |
| 6.0 | Collected Light Color Separation System. | 16 |

LIST OF FIGURES

| Figure | Title | Page |
|---------------|---|-------------|
| 1 | The 3.5 FT HWT LDV System. | 3 |
| 2 | Plenum Optics. | 4 |
| 3 | Transmitting Optics Color Separation System. | 5 |
| 4 | High Power Laser to Fiber Coupler. | 10 |
| 5 | Multi-mode Fiber Preadjustment. | 11 |
| 6 | Polarization Preserving Connector. | 13 |
| 7 | Plenum Feedthrough for Optical Fibers. | 14 |
| 8 | Collected Light Color Separation and Signal Detection System. | 16 |

1.0 THE LASER DOPPLER VELOCIMETER

The layout of the 3.5 FT HWT Laser Doppler Velocimeter is shown schematically in Fig. 1. Details of the plenum optics are shown in Fig. 2.




Figure 1 The 3.5 FT HWT 2-D LDV System.

Mean velocity and turbulence measurements are made with a dual-beam velocimeter utilizing a Bragg cell that enables moving interference fringes to be generated in the focal volume so that instantaneous velocity magnitude and direction measurements can be achieved from the frequency shift (f_D) around the incident and modulated laser beam interference frequency (f_0). i.e. $U = \lambda(f_D - f_0) / 2 \sin(\theta/2)$ where λ is the wavelength of the incident laser light.




Figure 2. Plenum Optics.

The transmitting optics color separation system (Fig. 3) is straightforward with a few unique features addressing the common problem of beam distortion or thermal blooming at higher laser powers. Frequency shifting is done before the color separation prisms, using a single acousto-optic modulator made of a selected flint glass, which can handle substantial laser power with minimal distortion. This is followed by color separation prisms, the first of which are made of fused silica for power handling capacity. A final prism of dense flint provides maximum angular displacement once the light has been dispersed into numerous beams. Final color selection is made using right angle prisms. The lines used for this application were 514.5 nm and 488 nm. Other laser lines could have been selected.




Figure 3 Transmitting Optics Separation System.

Pure fused-silica core single-mode polarization-preserving fibers are used for light transmission; two fibers per color. The use of optical fibers avoids the tedium of mirror-traverse

alignment. The pure fused silica core fibers are less susceptible to the progressive transmission losses which are found in other fibers. Polarization preserving fibers provide greater modal stability when the fibers are flexed or manipulated. For mechanical and thermal protection, the fibers are armored and contained within a conduit which is air cooled within the plenum. Upon exiting the fibers the light is collimated at 2.2 mm dia. with a separation of 60 mm. Adjustable rhomboid prisms reduce the beam separation to .3125 inches. The final focusing lens is 50.8 mm diameter and 750 mm focal length.

Forward-scattered light is collected with a 6 inch diameter, 30 inch focal length lens and focused into a 600 μm multi-mode optical fiber, which conducts it to the color separation and signal detection box through an air cooled conduit. For maximum throughput efficiency of the collected light color separator, a prism separation scheme is used rather than di-chroic filter and interference filters..

Experience has shown that accurate positioning is vital to a successful test program. Position is maintained by a custom designed eight axis capable traverse controller with micro-stepping drives, optical encoder feedback, and limit switch safety stops. Chapter 4 contains a detailed description of the traverse control system.

2.0 TRANSMITTING COLOR SEPARATION SYSTEM ALIGNMENT

2.1 Optics Enclosure

In order to steer the laser beam into the transmitting color separation box (Fig. 3), the best approach is often to use two steering mirrors between the laser and the box. This allows the beam to be fully manipulated without moving the laser or the box.

WARNING

**LASER SHOULD BE OPERATED AT
MINIMUM POWER DURING ALIGNMENT.**

**DO NOT STARE AT THE BEAM
OR DIFFUSE REFLECTIONS.**

WEAR APPROPRIATE PROTECTIVE EYEWEAR.

**PROJECT BEAMS ONTO A DULL, FLAT
BLACK, NON-FLAMMABLE SURFACE.**

Select a position for the color separation box on the optical table which supports the laser and bolt down the box baseplate using the clearance holes at the center of each side of the baseplate.

2.2 Outside Steering Mirrors

Move the polarization rotator and Bragg cell out of the way and open the iris diaphragm fully. Adjust the position and orientation of the input steering mirrors to direct the beam into the color separation box at a height of 4.25 inches from the top surface of the baseplate and 1.25 inches from the inside surface of the front plate. This can be checked by placing transparent rulers into the beam path at each end of the optics box.

2.3 Polarization Rotator

Put the polarization rotator into its holder and adjust it so that the beam travels through the center of the aperture. The polarization will be set after the Brewster angle dispersion prisms are put into position.

2.4 Acousto-Optic Modulator (Bragg Cell)

Set the precision micrometer adjustment of the the Bragg cell mount to the middle of its travel; about three full turns from the stops.

Slide the Bragg cell into the beam.

Connect the Bragg cell RF input to the inside front panel BNC feed through with a short length of RG-58 cable.

Connect the Bragg cell driver to the outside front panel BNC connector.

Loosen the gross movement set screw and position the Bragg cell so that the beam travels through the center of both the input and output apertures. Tighten the gross movement set screw.

Switch on the Bragg cell driver and turn up the drive power so that the diffracted beams can be seen.

Tilt the Bragg cell up and down to identify the first order diffracted beam. The unshifted beam is the one which remains when the Bragg cell driver is switched off. As the Bragg cell is tilted back and forth, the first order diffracted beam will appear above and then below the undiffracted beam. The Bragg cell should be set so that the undiffracted beam is on the bottom. Using the precision adjustment knob, set the Bragg cell tilt to put the maximum amount of light into the upper, first diffracted beam. A laser power meter can be used for this.

Adjust the Bragg cell drive power so that equal power is in both beams. Again, use a laser power meter for the greatest precision in this adjustment. This adjustment should be checked again at the measurement volume after the transmitting optics are completely set up. Coupling efficiency will vary from fiber to fiber. Also, the percentage of laser power diffracted into the shifted beam is wavelength dependent. Bragg cell drive power should be set to the best compromise, remembering that the ideal is equal power in each beam of each pair.

2.5 Dispersion Prisms

The two fused silica Brewster angle prisms are fixed to their mount. The incident beam should strike the first prism about in the middle. Increase the laser power just enough so that all beams are visible.

Rotate the prism pair while watching the refracted beams some distance past the prisms. As the prisms are rotated the refracted beams will be seen to move in one direction and then reverse. The position at which the beams reverse is the place to stop. Rotate the prisms just slightly to either side of that maximum deflection point. In one direction it will be seen that the beams are more circular than in the other. Fasten the mount at the point near the maximum deflection where the beams are circular.

Put the mounted flint prism into its holder and adjust it in the same manner as the Brewster angle prisms.

2.6 Polarization Rotator

The polarization rotator should be set for maximum transmission through the dispersion prisms. A laser power meter will provide the greatest precision in this adjustment.

2.7 Inside Steering Mirrors

The three steering mirrors should be placed so that the incident beams strike them in the middle. The three mirrors direct the beams along a path of sufficient length to allow adequate separation of the beams. The mean beam height above the baseplate should be 4.25 inches. Some adjustment of the outside steering mirrors may be required. If the outside steering mirrors are adjusted, the other optical components should be checked and readjusted as necessary. A transparent ruler placed in the beam path will show the undiffracted beams below 4.25 inches and the diffracted beams an equal distance above 4.25 inches. The third mirror should be adjusted to send the beams onto the middle of the separation prisms. Ensure that all beams enter and exit cleanly without striking any edges.

2.8 Iris Diaphragms

In normal operation, the iris diaphragm should be wide open and set close to the third steering mirror. The two iris diaphragms define the position of the incident laser beam. After the positions of all components ahead of the second iris are set, reduce laser power to a minimum and switch off the Bragg cell driver. The only remaining beam will normally be the undiffracted blue beam. Now close down the first iris, center it over the single beam, and fasten it down. Open the first diaphragm and close down the second one. Position the second iris so it is centered on the single beam and fasten it in position. Open the iris, switch on all beams and ensure that all beams travel through the open iris. The two iris diaphragms can now be used to pre-position the system if alignment is lost.

2.9 Separation Prisms

Two large right angle prisms are used to further separate the diffracted and undiffracted (Bragged and un-Bragged) beams and direct them to the final individual prisms. The junction of the two prisms should be set to 4.25 inches up from the baseplate. Ensure that all beams travel through the large prisms cleanly. This is a good place to clip any unwanted short wavelength beams.

2.10 Final Steering Prisms

Each beam is picked off by a small right angle prism and directed to the fiber launching optics. The lower beams are 3.5 inches and the upper 5.5 inches above the baseplate. It is important that the final prisms be positioned so that the beams are directed squarely into the launching optics. Each beam should be orthogonal to the front plate which holds the launching optics. Transparent plastic rulers may be used to make this setting.

3.0 FIBER OPTIC LINK

3.1 Laser to Fiber Coupler

The laser to fiber couplers are mounted on the outside of the front panel by 1"-32 mounting threads. The function of the couplers is to launch the laser beams into the fibers efficiently. Each coupler will focus its respective laser beam down to a small waist and maneuver the single mode fiber to the image plane of the lens system. The coupler (Fig. 4) is comprised of two baseplates, each with an axial bore, with an O-ring sandwiched between. Lateral (radial) movement across the beam is accomplished by adjustment of the three small socket-head screws which compress the outer baseplate against an O-ring. The clearance holes for the three socket head screws in the outer plate are slightly oversized to allow some lateral movement when the screws are loose. Final precision adjustment and stability is provided by another three screws which push against the inner baseplate in opposition to the compressing screws. Z axis adjustment along the beam waist is provided by a fine threaded adjustment, which is locked down with a set screw.

3.2 Laser to Fiber Coupler Alignment

The procedure described applies to each beam and coupler.

WARNING

**KEEP LASER POWER LOW UNTIL
ALL COUPLERS ARE ALIGNED TO
PREVENT BURNING FIBER CLADDING.**




Figure 4 High Power Laser to Fiber Coupler

Operating the laser at low power, first center the beam in the front plate hole using a transparent plastic ruler or template, then screw in the coupler. Project the expanding laser beam onto a flat black screen one or two feet away. Slightly loosen the screws which hold the inner and outer coupler plates together to allow lateral movement. Displace the coupler disc by hand while observing the projected beam on the opaque screen. Tighten the screws in the position where the beam intensity is centered and symmetric. This aligns the lens axis with the beam, increasing coupling efficiency.

Insert the 50 micron multi-mode fiber into the coupler and set the fiber output in a position to project light onto an opaque screen (Fig. 5). Identify the three screws which pull the coupler plates together. The tilt mechanism will displace the fiber core laterally relative to the image at the focal plane. As the focused image nears the center of the core, lower numerical aperture (N.A.) modes will be excited and more light will be concentrated into the center of the output spot (Fig. 5b). This will occur only under launching conditions where the N.A. of the focused rays is

smaller than the N.A. of the fiber. Using the small ball driver, adjust capscrew 1 while observing the output. The distribution of light in the output should change. Try to concentrate most of the light into the center of the spot by rotating the capscrew. Adjust capscrew 2 and continue to concentrate more and more light into the modes closest to the center of the spot. Repeat with capscrew 3. Sequentially adjust each screw for the maximum light coupling while steadily pulling the coupler tighter.

When the coupler is quite tight and adjusted for maximum coupling, replace the multi-mode fiber with a single-mode fiber and optimize for maximum coupling efficiency. An optical power meter should be used for the final adjustments. At the focal plane, where the fiber is located, there are multiple maxima due to diffraction (Airy discs). If the maximum light coupled is very low ($<10\%$ of input) it might be that a side order maxima is positioned on the fiber core. While watching the power meter, tilt each of the screws sufficiently to verify that the most powerful maximum is being coupled into the fiber.

If coupling efficiency is still low, the Z-axis may need some adjustment. The easiest way to verify that the Z-axis does need adjustment is to loosen the FC type fiber optic connector (FC connector) a turn and pull the fiber back from the focal plane. Then slowly tighten the FC




Figure 5 Multi-mode Fiber Pre-adjustment.

connector while watching the power meter to identify the position where the maximum coupling efficiency occurs. The aim is to peak coupling at the point when the FC connector is tight. If Z axis adjustment is necessary, alignment may be maintained by making exact 360 degree adjustments; that is, loosen or tighten the Z axis by exactly one turn. If alignment is lost or if more than a minor adjustment is necessary, the multi-mode fiber should be installed initially. If significant adjustment is necessary, the procedure described in the section on preliminary Z-axis adjustment should be followed.

3.3 Preliminary Z-Axis Adjustment

The Z-axis is pre-adjusted and this procedure should not normally be required. If adjustment becomes necessary, the procedure is to adjust the coupler as you would a collimator. Couple light into a fiber. Install the launching coupler on the output end of the fiber. Loosen the radial set screws which secure the Z-axis ferrule. Project the beam onto an opaque screen and adjust the Z-axis ferrule until you get a minimum diameter collimated beam at some distance. The lens is now positioned so that the fiber is at its focal plane. The coupler is now set up to launch a collimated input beam when used as a laser to fiber coupler.

3.4 Polarization Axis Adjustment

The single mode fibers provided with this system are highly birefringent, polarization-maintaining fibers, generally with two perpendicular principal axes. By maintaining polarization in the fibers we are able to match the polarization of the beam pairs in the measurement volume. Also, the output of properly aligned polarization-maintaining fiber will not fluctuate when the fibers are moved or manipulated. Polarization is maintained only when the polarization axis of the light is matched with that of the fiber. Improper alignment will cause the output polarization state of the fiber to oscillate between elliptical and linear polarization states. The polarization of the light can be rotated using a half wave plate placed ahead of the launching optics to match the orientation of the fiber or the fiber can be rotated to match the polarization orientation of the light.

Polarization alignment of fibers is measured by determining the extinction ratio of the output. First align the coupler for best coupling efficiency. Measure the fiber output through a polarizer with a light powermeter. Rotate the polarizer until maximum light transmission through the polarizer is achieved. Record this value. Rotate the polarizer until the minimum output is achieved and record the powermeter reading. Calculate the difference (extinction ratio) between the maximum and minimum readings in dB. Then rotate the knurled section of the fiber connector (Fig. 6) very slightly and repeat the procedure. Continue to rotate the fiber connector until the extinction ratio is maximum. Extinction ratios of 20 to 35 dB should be achieved. As mentioned above, a half wave rotator placed ahead of the launching optics in a rotary mount can easily set the light polarization to match the principal axis of the fiber. Placed in position temporarily, a half wave rotator can be used to help determine if the best polarization matching has been achieved or approximately how much adjustment is required. When the best extinction ratio has been

achieved, press or bend the fiber slightly, which may result in a small change in the power output after the polarizer. If the change is less than a few dBs, the polarization axes are aligned. Ideally, there should be no change. If the change is more than a few dBs, then rotate the fiber connector slightly until the required extinction ratio is achieved.

When the fiber is rotated, if an increase in insertion loss is noticed, it is due to fiber core / cladding concentricity problems. In lens style couplers, this could be compensated for by adjusting the angle between the incoming collimated beam and the receiver lens. Adjustment as described above in the previous section should be performed.




Figure 6 Polarization Preserving Connector.

Once the polarization axes have been properly set, the connectors may be glued to fix them in position. Loctite 290® for preassembled fasteners, Duco® cement or instant glue may be used,

depending on the degree of permanence desired. The glue should be put in the slot as shown at the bottom of Fig. 6.

3.5 Plenum Feedthrough for Optical Fibers

For additional protection, the armored fiber optic cable is contained within a conduit. Figure 7 shows how the fiber and conduit was plumbed through the plenum portholes. Inside the plenum, shop air was blown through the conduit to cool the fiber.




Figure 7 Plenum Feedthrough for Optical Fibers.

4.0 FIBER COLLIMATORS AND TRANSMITTING OPTICS

A fiber collimator is very similar in design to a laser to fiber coupler. Collimator length is proportional to the output beam diameter desired. The output beam is collimated by releasing the set screw(s) and adjusting the Z-axis, which is the distance from the fiber output to the collimating lens. This is best done before the collimator is attached to the mounting plate. Direct the beam some distance away and adjust the Z-axis until the best collimation is achieved.

With the output of each collimator set, they are attached to the mounting plate with three screws and an o-ring. To set the four beams mutually parallel, set the mounting plate in its mounting ring in a stand on the optics table. Set up a 60 mm template at the same height several feet away. Adjust each collimator to steer its beam onto the appropriate spot. The three mounting screws on each collimator should be quite snug. When adjustment is complete, tighten the three locking screws while checking that alignment is maintained. Repeat the procedure for the other mounting plate.

Fasten the collimator plate into the 100 mm I.D. mounting ring, then to the rail. The rhomboid beam separation reducer should be fastened to the rail ahead of the collimator plate. The beam separation used for this application was .3125 inches. Ensure clean passage of the laser beams through the rhomboids. Lastly, the focusing lens should be attached to the rail.

WARNING

WEAR LASER SAFETY EYEWEAR.

BEWARE OF SPECULAR REFLECTIONS.

Using a Polaroid filter, check that polarization of each beam pair is matched. Output polarization is set by rotating the knurled ring on the FC connector (Fig. 6). The fiber axis and beam polarization at the fiber input should have been matched already, as described above in the section on polarization axis adjustment.

With a clear plastic ruler, check that all beams cross at the focal point of the converging lens. Be sure to use appropriate eye safety precautions. Place a microscope objective or eyepiece in the beams to project the crossover point. Move the eyepiece axially along the beams. Each beam should be waisting at the focal volume. Use the Z-axis adjustment of each collimator to set the beam waist. Notice that the tightness of the FC connector at the collimator affects the collimation and thus the beam waist somewhat.

Each beam pair should be crossing fully without shearing and all pairs should cross together. For fine adjustment use the three screws on each collimator which press the plates apart.

5.0 COLLECTING OPTICS

The scattered light collection lens should be positioned for optimum forward scatter collection. A thin scattering center such as a piece of tape should be placed near the center of the proposed scan and the collecting optics driven and adjusted to focus the collected light into the multi-mode fiber. Direct the output of the multi-mode fiber onto a flat black surface. When the fiber is at the focus of the collecting optics, the center of the projected fiber output is illuminated. If the fiber is not located at the focus, there will be a bright ring around the outer edge of the illuminated area.




Figure 8 Collected Light Color Separation and Signal Detection System.

6.0 COLLECTED LIGHT COLOR SEPARATION SYSTEM

Figure 8 illustrates the layout for the collected light color separation and signal detection system. The multi-mode fiber collimator is adjusted in a similar fashion to the single mode unit. To set beam collimation, the set screw on the shaft is loosened and the distance from fiber end to collimating lens is varied while viewing the beam size at some distance. The collimator is then

screwed into the color separation box using the threaded 1"-32 hole. The collimated beam is directed through a pair of dense flint dispersion prisms to a right angle prism, which reverses the beam to travel through the dispersion prisms a second time at a slightly lower level. Another right angle prism directs the diverging colors through sets of light baffles to the final steering prisms which send each color to its photo-multiplier tube.

The dispersion prisms are fixed to their mount. During initial alignment, the dispersion prism pair, together with the direction reversing right angle prism, were positioned to provide the greatest dispersion with unclipped beams. If a component is knocked out of alignment, the best course is usually to leave the other components undisturbed and replace and adjust that component until proper orientation is again achieved.

CHAPTER 2

LASER VELOCIMETER DATA ACQUISITION SYSTEM.

CHAPTER 2

Laser Velocimeter Data Acquisition System.

TABLE OF CONTENTS

| Section | Title | Page |
|----------------|---|-------------|
| 1.0 | Introduction. | 3 |
| 1.1 | LDV Signal Conditioning Instrumentation Signals. | 3 |
| 2.0 | Laser Velocimeter Data Acquisition System Description. | 6 |
| 2.1 | Analog Data Description. | 8 |
| 2.2 | Digital Data Description. | 9 |
| 3.0 | Old Model Macrodynes. | 10 |
| 3.1 | Data Format. | 10 |
| 3.2 | Frequency Mode. | 10 |
| 3.3 | Velocity Mode. | 11 |
| 3.4 | Macrodyne Front Panel Digital Output Pinouts. | 11 |
| 3.5 | Interface Cable Schematic and Handshake Timing Diagram. | 11 |
| 3.6 | Data Reduction. | 11 |
| 3.7 | Period Calculation. | 11 |
| 3.8 | Frequency Calculation. | 12 |
| 3.9 | Velocity Mode Frequency Calculation. | 12 |
| 4.0 | New Model Macrodynes. | 13 |
| 4.1 | Data Format. | 13 |
| 4.2 | Frequency Mode. | 14 |
| 4.3 | Velocity Mode. | 14 |
| 4.4 | Macrodyne Front Panel Digital Output Pinouts. | 14 |
| 4.5 | Interface Cable Schematic and Handshake Timing Diagram. | 14 |
| 4.6 | Data Reduction. | 14 |
| 4.7 | Period Calculation. | 14 |
| 4.8 | Frequency Calculation. | 15 |
| 4.9 | Velocity Mode Frequency Calculation. | 15 |

LIST OF FIGURES

| Figure | Title | Page |
|---------------|---|-------------|
| 1. | Laser Velocimeter System Configuration. | 4 |
| 2. | Two-Component LDV Signal Conditioning. | 5 |
| 3. | Laser Velocimeter Data Acquisition System. | 7 |
| 4. | Macrodyne Digital Output Port Pinouts. | 16 |
| 5. | Old Macrodyne to LVDAS Interface Cable Schematic Drawing. | 17 |
| 6. | New Macrodyne to LVDAS Interface Cable Schematic Drawing. | 18 |
| 7. | Macrodyne to LVDAS Interface Handshake Timing Diagram. | 19 |

1.0 INTRODUCTION

The NASA Ames Research Center 3.5 Foot Hypersonic Wind Tunnel Laser Doppler Velocimeter System provides the capability to acquire and process simultaneous analog data and two-component Laser Doppler Velocimeter (LDV) data. The system consists of the following five sub-systems:

1. LDV Signal Conditioning Instrumentation.
2. LDV Counter Signal Processor Instrumentation.
3. Laser Velocimeter Data Acquisition System (LVDAS).
4. Data Acquisition, Data Reduction, and Data Presentation Computer System.
5. Traverse Control System (TCS8).

This document will discuss the theory of operation of the LVDAS and the LDV Counter Signal Processors as well as provide sources of documentation drawings for these instruments. The manner in which they are connected to and interact with each of the other optical and electronic sub-systems listed above will also be included. Figure 1 shows the configuration setup of the Laser Doppler Velocimeter system. This shows how the LVDAS fits into the complete system.

1.1 LDV Signal Conditioning Instrumentation.

Figure 2 shows the signal conditioning that is applied to the two-component Laser Doppler Velocimeter signals, the tunnel static temperature signal, and other optional analog voltage signals. The tunnel static temperature voltage output is fed directly to the first analog input channel of the LVDAS. Other optional analog voltage outputs can also be fed directly to one of the analog inputs of the LVDAS.

The LDV signal conditioning instrumentation is composed of the following elements:

1. RF Amplifier.
2. Frequency Filter.
3. Macrodyne LDV Counter Processors.

The voltage outputs of each RF amplifier are fed directly to the inputs of the Macrodyne LDV Counter Processors. The 16bit digital outputs of the two Macrodyne LDV Counter Processors are connected directly to the digital inputs of the LVDAS. The LVDAS is described in Section 2 while the Macrodyne LDV Counter Processors are described in Sections 3 and 4.




Figure 1. Laser Velocimeter System Configuration.




Figure 2. Two-Component LDV Signal Conditioning.

2.0 LASER VELOCIMETER DATA ACQUISITION SYSTEM.

New applications in laser velocimetry have brought about the need for a more advanced laser velocimeter data acquisition system. These new applications require high data rates that are not hindered by on-line time dependent data sorting and real time graphic data presentation. The new Laser Velocimeter Data Acquisition System (LVDAS) was designed specifically to meet these advanced requirements.

The Laser Velocimeter Data Acquisition System (LVDAS) provides the capability to acquire, process, and present real time digital data and analog data. The digital, for LDV systems, is typically the output of LDV signal processors. The analog data, for hypersonic wind tunnels, might include the raw signals containing tunnel temperature and/or pressure. The output of a filter amplifier whose input comes from flow sensors, such as a hot wire, might also be acquired with the LDV data. Additional analog data might originate from such sources as temperature probes, position sensors, etc. A functional schematic diagram of the LVDAS is shown in Figure 3. The LVDAS acquires simultaneous digital data, analog data, and time information data. The data are sampled, multiplexed, buffered, and then transferred to the facility's host computer for further data reduction, analysis, and presentation.

The digital data are sampled in a manner which ensures that the required coincidence time criterion is met. This is achieved by comparing 32bit 10MHz time of arrival counters for each of the digital channels. If data arrives on all of the selected digital channels within the coincidence time, then the analog channels are immediately sampled and converted. Otherwise, the digital data are rejected and the process is repeated. The 16 bit word parallel input ports are provided to accept the digital output of LDV counter processors and/or other instrumentation. High data acquisition rates are achieved by providing a separate latched input for each laser velocimeter digital input and a separate converter for each temperature, pressure, hot wire sensor or other analog input. The system will allow for data acquisition rates of approximately 100,000 samples per second simultaneously on each of the laser velocimeter and analog inputs.

A 32 bit time of day (TOD) 10MHz counter is used to tag arrival times to acquired digital LDV data as they become available on each of digital inputs. When a data valid "sync" pulse is sensed for a particular channel, the LVDAS latches the current TOD into a 32 bit time of arrival register (TOA). A separate TOA register is available for each digital input, so that particle arrival times of measured velocity information for U, V, and W can be monitored for coincidence. The latched times of arrivals have a resolution of 100 ns and maximum time of over 7 minutes.

The coincidence control logic allows for up to 3 channel coincidence. The coincidence time can be adjustable to any resolution or duration within the capability of the time of arrival registers. The coincidence time is adjustable from 100 ns to 1 s. In addition to the laser velocimeter inputs, three additional data words are generated internally. They are the inter-arrival time, the coincidence time, and status words. The inter-arrival and coincidence time is provided by a clock whose resolution is 100 ns and the maximum elapsed time is over 7 minutes. The status word contains information about coincidence which indicates whether or not valid data have been acquired.




Figure 3. Laser Velocimeter Data Acquisition System.

When coincident criteria are met, the analog inputs can be sampled and converted to provide concurrent data with the digital data. A single time of arrival is latched for all each of the analog to digital inputs, since they are all sampled and converted simultaneously. Additionally, a time of arrival is latched for external events, if they occur. These might be derived from such sources as oscillating models or model surfaces, rotating helicopter blades, rotating engine fans, or flow sensors.

All of the acquired digital velocity data with corresponding time of arrival data can be processed and stored even if coincidence is not required. However, if coincident data are required, then the arrival time of the various channels can be conditionally accepted if they all occur within a finite window of time. These coincident events can then be assigned inter-arrival times, which represent elapsed time since the previous event.

During data acquisition, it is important that the user obtain some visual feedback about the data being acquired. This is necessary so that the user can make informed decisions about both the quality and quantity of data received. The user is either reassured about the quality of the data or can make alterations and improvements in technique while "on line". To help achieve this, the instantaneous velocities are used to generate real time histograms from which probability density distributions are determined for all velocity components.

Additionally, the laser velocimeter data acquisition system has the capability of reducing the raw laser velocimeter data. Each laser velocimeter output contains the information required to calculate the instantaneous velocities. From the instantaneous velocity determinations, the average velocities, turbulence levels, and the turbulence cross correlations are all to be calculated.

All digital Macrodyne data, optional digital data, analog to digital data, and time of arrival data can be sent by the LVDAS to other computers via two serial and two parallel input/output ports. One parallel port will be used for the LDV system's data acquisition computer while the other can be used by the facility host computer. The serial ports can be used by PC type computers such as IBMs or MACs.

2.1 Analog Data Description.

The analog inputs are provided to accept differential voltages from such sources as hot wires, temperature probes, pressure probes, and other such sensors and/or instrumentation. The inputs are differential inputs and accept ± 5 volts.

The inputs are sampled and converted at a rate 200KHz with 16 bit resolution. The converted raw data (16 bit word) are encoded into signed two's complement binary format. Examples of the raw data to voltage conversion is shown on the next page.

| msb | BINARY WORD | | | | lsb | INTEGER | VOLTAGE |
|------|-------------|------|------|------|-----|---------|----------|
| 0111 | 1111 | 1111 | 1111 | 1111 | | 32767 | +4.99985 |
| 0111 | 1111 | 1111 | 1111 | 1110 | | 32766 | +4.99969 |
| ↓ | ↓ | ↓ | ↓ | ↓ | | ↓ | ↓ |
| 0000 | 0000 | 0000 | 0000 | 0001 | | 1 | +0.00015 |
| 0000 | 0000 | 0000 | 0000 | 0000 | | 0 | 0.00000 |
| 1111 | 1111 | 1111 | 1111 | 1111 | | -1 | -0.00015 |
| ↓ | ↓ | ↓ | ↓ | ↓ | | ↓ | ↓ |
| 1000 | 0000 | 0000 | 0000 | 0001 | | -32767 | -4.99985 |
| 1000 | 0000 | 0000 | 0000 | 0000 | | -32768 | -5.00000 |

One Bit Resolution: 0.00015

The signed binary two's complement integer word can be converted to a real precision floating point voltage using the following equation.

$$A_v = \frac{5A_R}{2^{15}} \quad \text{volts}$$

Where A_R is the analog raw voltage and A_v is the converted analog voltage.

2.2 Digital Data Description.

Three digital inputs are provided to accept 16bit digital data from such sources as LDV counter signal processors and/or instrumentation. The format for the Macrodyne Counter Signal Processors, which are being used with this system, will depend on whether or not the old or new model Macrodynes are used. The old models provide 16 bits of frequency information on a DB type 25 socket connector. The new models provide 18 bits of frequency information on a DB type 37 socket connector. The data formats and cable schematics for the old models are included in Section 3. The data formats and cable schematics for the new models are included in Section 4. (Note: The new model Macrodynes were delivered with the two-component LDV system. Therefore, the description of the "New Model Macrodynes" in section 4.0 and the cable schematic shown in Figure 6 would apply to this system.)

3.0 OLD MODEL MACRODYNES.

This chapter describes the data format of the old model Macrodyne's digital output port. Also described are the equations necessary to convert the raw data into frequency which represents the rate of fringe crossings of the Doppler burst. Additionally, a detailed description and schematic drawing is provided for the Macrodyne to LVDAS interface cable (see Figure 5). Figure 7 is a timing diagram showing the handshaking sequence of the control lines. Figure 7 also shows the timing sequence of the data lines. This indicates when the data become valid and then later latched.

3.1 Data Format.

The old model Macrodyne LDV counter signal processors provide the digital frequency output in the following 16 bit format:

| MSB | | | | | | | | | | | | | | | LSB | |
|------|-----|----|----|----|----|----|----|----|----|----|----|----|----|----|-----|--|
| 15 | 14 | 13 | 12 | 11 | 10 | 9 | 8 | 7 | 6 | 5 | 4 | 3 | 2 | 1 | 0 | |
| TIME | 5/8 | X3 | X2 | X1 | X0 | D9 | D8 | D7 | D6 | D5 | D4 | D3 | D2 | D1 | D0 | |

The mantissa (D9..D0) is contained within the lower 10 bits of the 16 bit word while the exponent (X3..X0) is contained within bits 10 through 13. The number of fringes measured is defined as 8 if 5/8=1 or 16 if 5/8=0. The time bit tells whether the mantissa and exponent represent a period (TIME=0) or a velocity (TIME=1). The period of the Doppler frequency is measured by the number of clock pulses (C=500MHz or 1000MHz) that are required to sense 8 or 16 fringe crossings.

The one bit 5/8 bit is set by the front panel 5/8 - 10/16 switch. It specifies whether 8 or 16 fringes are to be measured by the counter processor. The one bit TIME bit is set by the front panel Time/Velocity switch. It specifies one of two encoding schemes used by the Macrodynes to represent the frequency data.

The logic level for each of the old model Macrodynes varies with different units. Each unit may have a mixture of positive true or negative true logic for the TIME, 5/8, X3..X0, and D9..D0 data pins on the digital output port. Typically, when one orders multiple units, they all come configured with the same logic. But, units ordered at a later date may have a different mixture of positive true and negative true logic on the digital output port data pins.

3.2 Frequency Mode.

If TIME=0 then the data represent the time that had elapsed while a specified number of fringe crossings were observed by the counter processor. The elapsed time is encoded into a 10 bit mantissa (D9..D0) and 4 bit exponent (X3..X0). Both the mantissa and the exponent are in unsigned binary format.

3.3 Velocity Mode.

If TIME=1 then the data represents the frequency of the observed doppler burst. The frequency information is encoded into a 14 bit concatenation of the 4 bit exponent (X3..X0) and the 10 bit mantissa (D9..D0). The resulting 14 bit frequency word (X3..X0 D9..D0) is in unsigned binary format.

3.4 Macrodyne Front Panel Digital Output Pinouts.

The pinout assignment for the old model Macrodynes is shown in the first column of Figure 4.

3.5 Interface Cable Schematic and Handshake Timing Diagram.

Figure 5 shows a detailed schematic drawing for the Macrodyne to LVDAS 16bit parallel data interface cable. Figure 7 is a timing diagram of the handshake processes that happen each time data are transferred from the Macrodyne to the LVDAS.

3.6 Data Reduction.

The following sections describe how to convert the raw data into useful period or frequency data. The raw data are encoded into the 5/8 fringe count bit, 4 bit X3..X0 period exponent, and 10 bit D9..D0 period mantissa.

3.7 Period Calculation

The time T for the selected number of fringes can be calculated using the following equation. (Note: T is the time for the entire measured burst of 8 or 16 fringes.)

$$T = M 2^{(E-2)} \quad \text{ns}$$

Where M is the mantissa bits D9 through D0 and E is the exponent bits X3 through X0. To determine the doppler period for 8 fringes (T_8) and for 16 fringes (T_{16}) the following equations would apply. (Note: Both T_8 and T_{16} are the average time for only one fringe of the entire measured burst of 8 or 16 fringes.)

$$T_8 = \frac{M2^{(E-2)}}{2^3 10^9} \quad \text{s}$$

$$T_{16} = \frac{M2^{(E-2)}}{2^4 10^9} \quad \text{s}$$

3.8 Frequency Calculation

The doppler frequency can be calculated using one of the following equations depending on whether 8 or 16 fringes were measured.

$$F_8 = \frac{1}{T_8} = \frac{2^3 10^9}{M2^{(E-2)}} \quad \text{Hz}$$

$$F_{16} = \frac{1}{T_{16}} = \frac{2^4 10^9}{M2^{(E-2)}} \quad \text{Hz}$$

3.9 Velocity Mode Frequency Calculation

The Macrodyne manuals provide no information as to the conversion of the raw velocity mode data into useful frequencies of velocities. Therefore, no equations are provided here for raw data to frequency conversion. This data format is not used at the present time.

4.0 NEW MODEL MACRODYNES

This chapter describes the data format of the new model Macrodyne's digital output port. Also described are the equations necessary to convert the raw data into frequency which represents the rate of fringe crossings of the Doppler burst. Additionally, a detailed description and schematic drawing is provided for the Macrodyne to LVDAS interface cable (see Figure 6). Figure 7 is a timing diagram showing the handshaking sequence of the control lines. Figure 7 also shows the timing sequence of the data lines. This indicates when the data become valid and then later latched.

4.1 Data Format.

The new model Macrodyne LDV counter processors provide frequency information in a similar format to the old models. New model Macrodyne LDV counter signal processors provide the digital frequency output in the following 18 bit format:

| | | | | | | | | | | | | | | | | | |
|------|-----|----|----|----|----|-----|-----|----|----|----|----|----|----|----|----|----|-----|
| MSB | | | | | | | | | | | | | | | | | LSB |
| 17 | 16 | 15 | 14 | 13 | 12 | 11 | 10 | 9 | 8 | 7 | 6 | 5 | 4 | 3 | 2 | 1 | 0 |
| TIME | 5/8 | X3 | X2 | X1 | X0 | D11 | D10 | D9 | D8 | D7 | D6 | D5 | D4 | D3 | D2 | D1 | D0 |

The mantissa (D11..D0) is contained within the lower 12 bits of the 18 bit word while the exponent (X3..X0) is contained within bits 12 through 15. The number of fringes measured is defined as 8 if 5/8=1 or 16 if 5/8=0. The time bit tells whether the mantissa and exponent represent a period (TIME=0) or a velocity (TIME=1). The period of the Doppler frequency is measured by the number of clock pulses (C=500MHz or 1000MHz) that are required to sense 8 or 16 fringe crossings.

The major difference is that the addition of two more mantissa bits (D11 and D10) to provide increased dynamic range for a fixed exponent. The equations for determining T, T8, T16, F8, and F16 would be identical to those illustrated in the previous section. However, D0 and D1 are ignored, since the current LVDAS digital interface provides for 16 input lines. Therefore, the following 16 bit format is actually transmitted to the LVDAS.

| | | | | | | | | | | | | | | | |
|------|-----|----|----|----|----|-----|-----|----|----|----|----|----|----|----|-----|
| MSB | | | | | | | | | | | | | | | LSB |
| 15 | 14 | 13 | 12 | 11 | 10 | 9 | 8 | 7 | 6 | 5 | 4 | 3 | 2 | 1 | 0 |
| TIME | 5/8 | X3 | X2 | X1 | X0 | D11 | D10 | D9 | D8 | D7 | D6 | D5 | D4 | D3 | D2 |

The one bit 5/8 bit is set by the front panel 5/8 - 10/16 switch. It specifies whether 8 or 16 fringes are to be measured by the counter processor. The one bit TIME bit is set by the front panel Time/Velocity switch. It specifies one of two encoding schemes used by the Macrodynes to represent the frequency data.

4.2 Frequency Mode.

If TIME=0 then the data represent the time that had elapsed while a specified number of fringe crossings were observed by the counter processor. The elapsed time is encoded into a 12 bit mantissa (D11..D0) and 4 bit exponent (X3..X0). Both the mantissa and the exponent are in unsigned binary format.

4.3 Velocity Mode.

If TIME=1 then the data represents the frequency of the observed doppler burst. The frequency information is encoded into a 16 bit concatenation of the 4 bit exponent (X3..X0) and the 12 bit mantissa (D11..D0). The resulting 16 bit frequency word (X3..X0 D11..D0) is in unsigned binary format.

4.4 Front Panel Digital Output Pinouts.

The pinout assignment for the old model Macrodyne is shown in the second column of Figure 4.

4.5 Interface Cable Schematic and Handshake Timing Diagram.

Figure 6 shows a detailed schematic drawing for the Macrodyne to LVDAS 16bit parallel data interface cable. Figure 7 is a timing diagram of the handshake processes that happen each time data are transferred from the Macrodyne to the LVDAS.

4.6 Date Reduction.

The following sections describe how to convert the raw data into useful period or frequency data. The raw data are encoded into the 5/8 fringe count bit, 4 bit X3..X0 period exponent, and 12 bit D11..D0 period mantissa.

With the deletion of the D1 and D0 the mantissa M would be represented by D11..D2 instead of D11..D0 and equations for T, T₈, T₁₆, F₈, and F₁₆ would be modified as shown in Sections 4.7 and 4.8.

4.7 Period Calculation

$$T = M 2^{(E-0)} \quad \text{ns}$$

$$T_8 = \frac{M2^{(E-0)}}{2^3 10^9} \quad \text{s}$$

$$T_{16} = \frac{M2^{(E-0)}}{2^4 10^9} \quad \text{s}$$

4.8 Frequency Calculation.

$$F_8 = \frac{1}{T_8} = \frac{2^3 10^9}{M2^{(E-0)}} \quad \text{Hz}$$

$$F_{16} = \frac{1}{T_{16}} = \frac{2^4 10^9}{M2^{(E-0)}} \quad \text{Hz}$$

4.9 Velocity Mode Frequency Calculation.

The Macrodyne manuals provide no information as to the conversion of the raw velocity mode data into useful frequencies of velocities. Therefore, no equations are provided here for raw data to frequency conversion. This data format is not used at the present time.

| | 3002 (OLD) FRONT | 3002 (NEW) FRONT | 300X REAR | 3003 FRONT |
|----|---------------------|---------------------|--------------|---------------|
| 1 | -D0 | D00 | -D0 | -D0 |
| 2 | -D2 | D02 | -D2 | -D2 |
| 3 | -D4 | D04 | -D4 | -D4 |
| 4 | -D6 | D06 | -D6 | -D6 |
| 5 | -D8 | D08 | -D8 | -D8 |
| 6 | -X0 | D10 | -X0 | -X0 |
| 7 | -X2 | X0 | -X2 | -X2 |
| 8 | -5/8 | X2 | -5/8 | -5/8 |
| 9 | SYNC | -HOLD OFF | -SYNC | -SYNC |
| 10 | -INHIBIT | SYNC | -SprInh | -INHIBIT |
| 11 | +5V | TIME | -RESET | -CC1 |
| 12 | -TIME | NotUsed | 10MHz | -CompAcc2 |
| 13 | GROUND | NotUsed | -CC2 | -CompAcc8 |
| 14 | -D1 | GROUND | -D1 | -D1 |
| 15 | -D3 | GROUND | -D3 | -D3 |
| 16 | -D5 | NotUsed | -D5 | -D5 |
| 17 | -D7 | NotUsed | -D7 | -D7 |
| 18 | -D9 | NotUsed | -D9 | -D9 |
| 19 | -X1 | +5V | -X1 | -X1 |
| 20 | -X3 | D01 | -X3 | -X3 |
| 21 | GROUND | D03 | GROUND | GROUND |
| 22 | GROUND | D05 | GROUND | GROUND |
| 23 | GROUND | D07 | GROUND | GROUND |
| 24 | GROUND | D09 | GROUND | -CompAcc1 |
| 25 | GROUND | D11 | GROUND | -CompAcc4 |
| 26 | | X1 | | |
| 27 | | X3 | | |
| 28 | | -ExtRes | | |
| 29 | | AnalogOut | | |
| 30 | | -EIGHT | | |
| 31 | | NotUsed | | |
| 32 | | GROUND | | |
| 33 | | GROUND | | |
| 34 | | GROUND | | |
| 35 | | NotUsed | | |
| 36 | | NotUsed | | |
| 37 | | NotUsed | | |

Figure 4. Macrodyne Digital Output Port Pinouts.




Figure 5. Old Macrodyne to LVDAS Interface Cable Schematic Drawing.




Figure 6. New Macrodyne to LVDAS Interface Cable Schematic Drawing.




Figure 7. Macrodyne to LVDAS Interface Handshake Timing Diagram.

CHAPTER 3

DATA ACQUISITION COMPUTER HARDWARE AND SOFTWARE.

CHAPTER 3

Data Acquisition Computer Hardware and Software.

TABLE OF CONTENTS

| Section | Title | Page |
|----------------|--|-------------|
| 1.0 | Data Acquisition Computer Hardware Description. | 2 |
| 1.1 | Hewlett-Packard Series 9000 Model 375 Computer. | 2 |
| 1.2 | System Interface Board. | 2 |
| 1.3 | General Purpose Input/Output High Speed Interface. | 2 |
| 1.4 | Integral Hard Disk and Floppy Disk Drives. | 2 |
| 1.5 | Paint Jet Printer. | 2 |
| 2.0 | Data Acquisition Computer Software Description. | 4 |
| 2.1 | Instantaneous Velocities, Voltages, and Temperatures. | 4 |
| 2.2 | Velocity, Voltage, and Temperature Averages. | 5 |
| 2.3 | Velocity, Voltage, and Temperature Standard Deviations. | 6 |
| 2.4 | Velocity, Voltage, and Temperature Cross Correlations. | 8 |
| 3.0 | HP 98622A GPIO Card to LVDAS 16bit Parallel I/O Interface. | 10 |
| 4.0 | HP 9000 Model 375 to LVDAS RS-232 Serial I/O Interface. | 13 |
| 5.0 | HP 9000 Model 375 to TCS8 RS-232 Serial I/O Interface. | 14 |

LIST OF FIGURES

| Figure | Title | Page |
|---------------|---|-------------|
| 1. | Laser Velocimeter System Configuration. | 3 |
| 2. | HP Series 9000 Model 3XX to LVDAS Interface Cable. | 10 |
| 3. | LVDAS Interface Cable Round Connector Pinouts. | 11 |
| 4. | Interface Cable Labels. | 12 |
| 5. | HP Series 9000 Model 375 to LVDAS Serial Interface Cable. | 13 |
| 6. | HP Series 9000 Model 375 to TCS8 Serial Interface Cable. | 14 |

1.0 DATA ACQUISITION COMPUTER HARDWARE DESCRIPTION.

A simplified schematic drawing of the computer hardware is shown in Figure 1. The Figure also shows the how computer hardware interconnects with the Traverse Control System (TCS8) and the Laser Velocimeter Data Acquisition System (LVDAS). The data acquisition, data reduction, and data presentation computer system hardware was comprised of the following elements:

1. A Hewlett-Packard Series 9000 Model 375 Computer.
2. A System Interface Board.
3. A General Purpose Input/Output High Speed Interface.
4. Integral Hard Disk and Floppy Disk Drives.
5. Paint Jet Printer.

1.1 Hewlett-Packard Series 9000 Model 375 Computer.

The HP Series 9000 Model 375 computer was used to control the traverse system, acquire LDV data, perform data reduction and analysis, present the reduced data in graphical form, and to store the raw and reduced data on hard disk.

1.2 System Interface Board.

The system interface board possesses multiple serial and parallel interfaces. A normal IEEE-488 HPIB interface is used to send data to the Paint Jet printer. A high speed IEEE-488 HPIB interface is used to read data from and write data to the integral 40MByte Hard Disk and Floppy Disk Drives. The RS-232 serial interface is used to send commands as well as to send and receive position information from the Traverse Control System (TCS8).

1.3 General Purpose Input/Output High Speed Interface.

The General Purpose Input/Output (GPIO) High Speed Parallel Interface is used to send commands to the Laser Velocimeter Data Acquisition System (LVDAS). The LVDAS subsequently transmits back LDV data over this GPIO interface to the HP 9000-375 Computer.

1.4 Integral Hard Disk and Floppy Disk Drives.

The hard disk is partitioned into volumes. One volume contains system related files and the data acquisition program. The system files include the BASIC operating system and initialization programs that configure the computer, CRT display, and keyboard. The data acquisition program, which also resides on this volume, is automatically loaded and executed as part of the computers "power up" sequence. Another volume is used to store raw and reduced data for archival purposes and for future data reduction and analysis.

1.5 Paint Jet Printer.

The Paint Jet printer is used for listing programs and to print reduced data in tabular form. Additionally, graphs are "dumped" to provide a hard copy of histogram and profile plots.




Figure 1. Laser Velocimeter System Configuration.

2.0 DATA ACQUISITION COMPUTER SOFTWARE DESCRIPTION.

The software used to control the traverse system and acquire tunnel data is listed in Appendix A and Appendix B of this report. Appendix A contains a catalog of the back-up floppy disk containing the system files and the data acquisition program named "3.5'HWT91". The system file "SYSB60" contains the BASIC 6.0 operating system. The "AUTOST" program is loaded and executed as part of the computer's "power up" sequence. This "AUTOST" program sets default values for the CRT and keyboard and then automatically loads and executes the "3.5'HWT91" program. Appendix A also contains a hardcopy listing of this "3.5'HWT91" program. This program is the original program that was used to acquire data during the hypersonic wind tunnel testing.

Appendix B is essentially a revised version with documentation of the program in Appendix A. The documentation is integrated into the code and includes information on how to boot the system and software. How to operate the menu driven software is also described. The documentation for the main program, each sub-routine, and each sub-program includes a description of the software, its purpose, and a list of variables with definitions. Appendix B contains a catalog of the back-up floppy disk containing the system files and the data acquisition program named "3.5'HWT92". The system file "SYSB60" contains the BASIC 6.0 operating system. The "AUTOST" program is loaded and executed as part of the computer's "power up" sequence. This "AUTOST" program sets default values for the CRT and Keyboard and then automatically loads and executes the "3.5'HWT92" program. Appendix B also contains a hardcopy listing of this "3.5'HWT92" program.

The following parts of Section 2 of this chapter contain a brief description of the data reduction applied to the raw data acquired by the "3.5'HWT" programs. A more complete set of documentation on the data reduction as well as coordinate system transformations can be found in Appendix C. Other topics are documented within the software code listing itself (refer to Appendix B for this software code listing.)

2.1 **Instantaneous Velocities, Voltages, and Temperatures.**

The following are the instantaneous velocities (U_1 , V_1) which are derived from the digital data outputted from the Macrodyne counter signal processors. All velocities are measured in meters/second (m/s).

U_1 : Instantaneous Streamwise Velocity.

V_1 : Instantaneous Vertical Velocity.

The following are the instantaneous voltages (A_1 , B_1) for the first two analog channels.
All analog inputs are measured in volts (v).

A_1 : Instantaneous Voltage on Analog Channel #1.

B_1 : Instantaneous Voltage on Analog Channel #2.

The following is the stagnation temperature (T_1) which is inputted from the first analog channel (A_1). The temperatures are measured in degrees Rankine (°R).

T_1 : Instantaneous Stagnation Temperature.

2.2 Velocity, Voltage, and Temperature Averages.

The instantaneous velocities (U_1 , V_1), the instantaneous voltages (A_1 , B_1), and the instantaneous stagnation temperatures (T_1) are summed so that the average velocities (\bar{U} , \bar{V}), the average voltages (\bar{A} , \bar{B}), and the average stagnation temperature (\bar{T}) can be calculated. All velocities are measured in meters/second (m/s), all analog inputs are measured in volts (v), and the stagnation temperature is measured in degrees Rankine (°R).

\bar{U} : Average Velocity (Streamwise).

\bar{V} : Average Velocity (Vertical).

\bar{A} : Average Voltage (Analog Channel #1).

\bar{B} : Average Voltage (Analog Channel #2).

\bar{T} : Average Stagnation Temperature.

$$\bar{U} = \frac{\sum_{i=1}^n [U_i]}{n} \quad \text{m/s}$$

$$\bar{V} = \frac{\sum_{i=1}^n [V_i]}{n} \quad \text{m/s}$$

$$\bar{A} = \frac{\sum_{i=1}^n [A_i]}{n} \quad v$$

$$\bar{B} = \frac{\sum_{i=1}^n [B_i]}{n} \quad v$$

$$\bar{T} = \frac{\sum_{i=1}^n [T_i]}{n} \quad ^\circ R$$

2.3 Velocity, Voltage, and Temperature Standard Deviations.

The velocity (U' , V'), voltage (A' , B'), and stagnation temperature (T') standard deviations are defined as shown here:

U' : Velocity Standard Deviation (Streamwise).

V' : Velocity Standard Deviation (Vertical).

A' : Voltage Standard Deviation (Analog Channel #1).

B' : Voltage Standard Deviation (Analog Channel #2).

T' : Stagnation Temperature Standard Deviation.

The following equations can be used to calculate the velocity (U' , V'), voltage (A' , B'), and stagnation temperature (T') standard deviations:

$$U' = \sqrt{\frac{\sum_{i=1}^n [U_i - \bar{U}]^2}{n}} \quad m/s$$

$$V' = \sqrt{\frac{\sum_{i=1}^n [V_i - \bar{V}]^2}{n}} \quad m/s$$

$$A' = \sqrt{\frac{\sum_{i=1}^n [A_i - \bar{A}]^2}{n}} \quad \text{v}$$

$$B' = \sqrt{\frac{\sum_{i=1}^n [B_i - \bar{B}]^2}{n}} \quad \text{v}$$

$$T' = \sqrt{\frac{\sum_{i=1}^n [T_i - \bar{T}]^2}{n}} \quad ^\circ\text{R}$$

The above equations are simplified to produce the following equations. The instantaneous velocities (U_i, V_i), voltages (A_i, B_i), and temperatures (T_i) are summed so that velocity (U', V'), voltage (A', B'), and stagnation temperature (T') standard deviations can be calculated. All velocity standard deviations are measured in meters/second (m/s), all voltage standard deviations are measured in volts (v), and all temperatures are measured in degrees Rankine ($^\circ\text{R}$).

$$U' = \sqrt{\frac{\sum_{i=1}^n [U_i^2]}{n} - \bar{U}^2} \quad \text{m/s}$$

$$V' = \sqrt{\frac{\sum_{i=1}^n [V_i^2]}{n} - \bar{V}^2} \quad \text{m/s}$$

$$A' = \sqrt{\frac{\sum_{i=1}^n [A_i^2]}{n} - \bar{A}^2} \quad \text{v}$$

$$B' = \sqrt{\frac{\sum_{i=1}^n [B_i^2]}{n} - \bar{B}^2} \quad \text{v}$$

$$T' = \sqrt{\frac{\sum_{i=1}^n [T_i^2]}{n} - \bar{T}^2} \quad ^\circ R$$

The equations are simplified to these forms so that the software can compute summations of the instantaneous velocities, voltages, and temperature as well as the summations of their squares within the same software loop. This eliminates the need to calculate the difference values ($U_i - \bar{U}$, $V_i - \bar{V}$, $A_i - \bar{A}$, $B_i - \bar{B}$, $T_i - \bar{T}$). Also, the need to calculate the averages before the squared summations is removed.

2.4 Velocity, Voltage, and Temperature Cross Correlations.

The velocity:velocity shear stress ($\overline{U'V'}$), velocity:voltage cross correlations ($\overline{U'A'}$, $\overline{V'A'}$), and voltage:voltage cross correlations ($\overline{A'B'}$) are defined as shown here:

$\overline{U'V'}$: Velocity:Velocity Shear Stress.

$\overline{U'A'}$: Velocity:Voltage Cross Correlation.

$\overline{V'A'}$: Velocity:Voltage Cross Correlation.

$\overline{A'B'}$: Voltage:Voltage Cross Correlation.

The following equations can be used to calculate the shear stress and the cross correlations ($\overline{U'V'}$, $\overline{U'A'}$, $\overline{V'A'}$, $\overline{A'B'}$):

$$\overline{U'V'} = \frac{\sum_{i=1}^n [(U_i - \bar{U})(V_i - \bar{V})]}{n} \quad m^2/s^2$$

$$\overline{U'A'} = \frac{\sum_{i=1}^n [(U_i - \bar{U})(A_i - \bar{A})]}{n} \quad mv/s$$

$$\overline{V' A'} = \frac{\sum_{i=1}^n [(V_i - \bar{V})(A_i - \bar{A})]}{n} \quad \text{mv/s}$$

$$\overline{A' B'} = \frac{\sum_{i=1}^n [(A_i - \bar{A})(B_i - \bar{B})]}{n} \quad \text{v}^2$$

The above equations are simplified to produce the following equations. Summations of the instantaneous velocity and voltage (U_i , V_i , A_i , B_i) products are summed so that velocity:velocity shear stress ($\overline{U' V'}$), velocity:voltage cross correlations ($\overline{U' A'}$, $\overline{V' A'}$), and voltage:voltage cross correlation ($\overline{A' B'}$) can be calculated. All velocity:velocity shear stresses are measured in meters²/second² (m²/s²). All velocity:voltage cross correlations are measured in meters•volts/second (mv/s). All voltage:voltage cross correlations are measured in volts² (v²).

$$\overline{U' V'} = \frac{\sum_{i=1}^n [U_i V_i]}{n} - \bar{U} \bar{V} \quad \text{m}^2/\text{s}^2$$

$$\overline{U' A'} = \frac{\sum_{i=1}^n [U_i A_i]}{n} - \bar{U} \bar{A} \quad \text{mv/s}$$

$$\overline{V' A'} = \frac{\sum_{i=1}^n [V_i A_i]}{n} - \bar{V} \bar{A} \quad \text{mv/s}$$

$$\overline{A' B'} = \frac{\sum_{i=1}^n [A_i B_i]}{n} - \bar{A} \bar{B} \quad \text{v}^2$$

The equations are simplified to this form so that the software can compute summations of the instantaneous velocities and voltages (U_i , V_i , A_i , B_i) as well as the summations of their products within the same software loop. This eliminates the need to calculate the difference values ($U_i - \bar{U}$, $V_i - \bar{V}$, $A_i - \bar{A}$, $B_i - \bar{B}$). Also, the need to calculate the averages before the product summations is removed.

3.0 HP 98622A GPIO CARD TO LVDAS 16BIT PARALLEL I/O INTERFACE.




Figure 1. Laser Velocimeter System Configuration.





Figure 3. LVDAS Interface Circular Connector Pinout Positions.





Figure 4. HP Series 9000 Model 3xx to LVDAS Interface Cable Labels.

4.0 HP 9000 MODEL 375 TO LVDAS RS-232 SERIAL I/O INTERFACE.

HP Series 9000 Model 375
to
LVDAS
Serial Interface Cable



Cable Label





Figure 5. HP Series 9000 Model 375 to LVDAS Serial Interface Cable.

5.0 HP 9000 MODEL 3XX TO TCS8 RS-232 SERIAL I/O INTERFACE.

HP Series 9000 Model 375 to TCS8 Serial Interface Cable



Cable Label




Figure 6. HP Series 9000 Model 375 to TCS8 Serial Interface Cable.

CHAPTER 4

TRAVERSE CONTROL SYSTEM.

CHAPTER 4

Traverse Control System.

TABLE OF CONTENTS

| Section | Title | Page |
|----------------|--|-------------|
| | The Traverse Control System. | 3 |
| 1.00 | | 4 |
| 1.01 | The TCS8. | 4 |
| 1.02 | The Motor Drive System. | 4 |
| 1.03 | Positioning Resolution. | 6 |
| 1.04 | Tunnel Penetration of Traverse Cables. | 7 |
| | Front Panel Description of the TCS8. | 7 |
| 2.00 | | 7 |
| 2.01 | Position Display Windows. | 7 |
| 2.02 | Power Key. | 7 |
| 2.03 | Jog Control Keys. | 8 |
| 2.04 | Scroll Keys. | 8 |
| 2.05 | Command Windows. | 8 |
| 2.06 | Execute Key. | 8 |
| 2.07 | Data Window. | 8 |
| 2.08 | Stop Key. | 8 |
| 2.09 | Status Window. | 8 |
| 2.10 | Numeric Key Pad. | 9 |
| | Local Command Descriptions of the TCS8. | 10 |
| 3.00 | | 10 |
| 3.01 | Move to Zero. | 10 |
| 3.02 | Move Absolute. | 11 |
| 3.03 | Move Relative. | 11 |
| 3.04 | Jog Mode. | 12 |
| 3.05 | Set Counts Per Unit. | 12 |
| 3.06 | Set Counts Per Revolution. | 13 |
| 3.07 | Set Position. | 13 |
| 3.08 | Set Velocity. | 14 |
| 3.09 | Set Acceleration. | 14 |
| 3.10 | Set Currents On. | 14 |
| 3.11 | Set Currents Off. | 15 |
| 3.12 | Set Inits On. | 15 |
| 3.13 | View Counts Per Unit. | 15 |
| 3.14 | View Counts Per Revolution. | 16 |
| 3.15 | View Velocity. | 16 |
| 3.16 | View Acceleration. | 16 |
| 3.17 | View Init. | 17 |
| 3.18 | View Currents. | 17 |
| 3.19 | View Plus Limit Switches. | 17 |
| 3.20 | View Minus Limit Switches. | 18 |
| 3.21 | View Home Switches. | 18 |
| 3.22 | View Stall Indication. | 18 |
| 3.23 | Init Default. | 18 |
| 3.24 | Init Drive On. | 18 |

| Section | Title | Page |
|----------------|---|-------------|
| 3.25 | Init Drive Off. | 19 |
| 3.26 | COM1/COM2 Baud Rate. | 19 |
| 3.27 | COM1/COM2 Bits Per Character. | 20 |
| 3.28 | COM1/COM2 Parity. | 20 |
| 3.29 | COM1/COM2 Stop Bits. | 20 |
| 3.30 | COM1/COM2 Handshake. | 21 |
| 4.00 | Serial Interface Command Descriptions of the TCS8. | 22 |
| 4.01 | Change Serial Configuration. | 23 |
| 4.02 | Move to Absolute Position. | 24 |
| 4.03 | Move Relative to Current Position. | 24 |
| 4.04 | Set Acceleration. | 25 |
| 4.05 | View Acceleration. | 25 |
| 4.06 | Set Velocity. | 26 |
| 4.07 | View Velocity. | 26 |
| 4.08 | Set Encoder Counts Per Unit of Travel. | 27 |
| 4.09 | View Encoder Counts Per Unit of Travel. | 27 |
| 4.10 | Set Counts Per Motor Revolution. | 28 |
| 4.11 | View Counts Per Motor Revolution. | 28 |
| 4.12 | Set Position. | 29 |
| 4.13 | View Position. | 29 |
| 4.14 | Set Current to Motor Windings. | 30 |
| 4.15 | View Current to Motor Windings. | 30 |
| 4.16 | Set Initialization of Indexer/Drivers. | 31 |
| 4.17 | View Initialization of Indexer/Drivers. | 31 |
| 5.00 | Traverse Control System Cables. | 32 |

LIST OF FIGURES

| Figure | Title | Page |
|---------------|---|-------------|
| 1 | NASA Ames 3.5' HWT Traverse Control System. | 3 |
| 2 | Schematic of Motor Drive System Back Panel. | 5 |
| 3 | Schematic of TCS8 Back Panel. | 5 |
| 4 | Tunnel Penetration Plate. | 6 |
| 5 | The Front Panel. | 7 |
| 6 | Motor Drive System to Compumotor Stepper Motor Cable. | 32 |
| 7 | Motor Drive System to Dynamics Research Encoder Cable. | 33 |
| 8 | Motor Drive System to Linear Industries Limit Switches Cable. | 34 |
| 9 | Motor Drive System to TCS8 Encoder Signals Cable. | 35 |
| 10 | HP Series 9000 Model 375 to TCS8 Serial Cable. | 36 |
| 11 | TCS8 to Motor Drive System Serial Cable. | 37 |

1.00 THE TRAVERSE CONTROL SYSTEM

The traverse control system is made up of four sub-systems, see Fig. 1. The first sub-system is the main data taking computer (host computer). The second sub-system, the TCS8 (Traverse Control System 8 Axis), receives high level traverse commands from the host computer. The full duplex serial communications that links these two sub-systems allows the host computer to monitor the position and status of each axis in the system, see Section 4.00 Serial Interface Command Descriptions of the TCS8. The TCS8 can also function as a "stand alone" traverse controller. Through the use of the TCS8's front panel, an operator can execute all of the commands that the host computer can in addition the operator can control all axes in jog mode, see Section 2.00 Front Panel Descriptions of the TCS8 and Section 3.00 Local Command Descriptions of the TCS8. The third sub-system, the MDS (Motor Drive System), is controlled solely by the TCS8. The TCS8 translates the high level commands from the host computer and its front panel into low level indexer commands, see The Compumotor AX Drive User Manual. The TCS8 also receives encoder pulses from the traverses via the MDS. This allows the TCS8 to display real time position information on its front panel. The fourth and final sub-system of the traverse control system is the slide, motor, encoder, and limit switches that make up each axis. A drawing of each cable which is used to connect the traverse control system is included in Section 5.00 Traverse Control System Cables.




Figure 1 NASA Ames 3.5' HWT Traverse Control System.

1.01 The TCS8

The TCS8 is a microprocessor controlled system designed to interface an operator with a traverse system. The operator can utilize the TCS8 through the front panel, see Section 2.00 TCS8's Front Panel Descriptions and Section 3.00 TCS8's Local Command Descriptions, and/or with one or two host computers over serial interfaces, see Section 4.00 TCS8's Serial Interface Command Descriptions. The TCS8 stores all the critical parameters of motion, for each of the eight axes that it controls, in non-volatile memory. The critical parameters of motion being: position, encoder counts per unit travel, encoder counts per motor revolution, velocity, and acceleration. All of these parameters may be viewed, set, and saved. The TCS8 has three modes of motion; absolute, relative, and jog. With absolute movements, the operator specifies the final location. With relative movements, a distance is specified. With jogged movements, the operator presses a jog key on the front panel of the TCS8 until the desired location is obtained.

1.02 The Motor Drive System

There are four indexer/drivers used in this system. The TCS8 communicates with the indexers in the MDS's over a closed loop serial daisy chain. The 4/8 switch is located on the back panel of the MDS and must be set to 4, see Fig. 2. This figure also shows the location of all the motor, limit, and encoder connections. Channels X1, X2, Y1, and Y2 of the TCS8 control axis 1 through 4 on the first MDS. The TCS8 Encoders connector on the back of each MDS has a corresponding connector on the back of the TCS8, see Fig. A3 Schematic of TCS8 Back Panel. The interconnecting cable is detailed in Section 5.00 Traverse Control System Cables.

1.03 Positioning Resolution

The indexer/drivers that are used in the MDS can drive the motors at 12,800 steps/revolution. The encoders used on each axis are 100 pulses/revolution with quadrature encoding. Quadrature encoding adds a factor of 4 to the number of pulses/revolution to make this number 400 pulses/revolution. The final factor in the product of the resolution of an axis is the number of threads/inch of the lead screw. All of the axes of the traverse system have lead screws of 10 threads/inch. Thus, the positioning resolution of the axes with a 10 threads/inch lead screw is 0.00025 inches.




Figure 2 Schematic of Motor Drive System Back Panel.




Figure 3 Schematic of TCS8 Back Panel.

1.04 Tunnel Penetration of Traverse Cables

The traverse slides for the 3.5' LDV System are located inside the pressurized test chamber and the traverse electronics are located outside, in the control room (see Fig 1. of Chapter 1). The traverse cables are fed through an existing access port on the north-east side of the test chamber. A special plate was designed (see Fig. 4) to replace an existing one. Bulkhead cable clamps are used to seal around the cables. When tightened down, these cable clamps compress a rubber grommet to create a seal. The four encoder cables are fed through one cable clamp and the four limit switch cables through another. The four motor cables, which are a larger diameter, are each fed through their own cable clamp.




Figure 4 Tunnel Penetration Plate.

2.00 FRONT PANEL DESCRIPTION OF THE TCS8




Figure 5 The Front Panel.

2.01 Position Display Windows.

There are eight windows corresponding to the eight axes that the TCS8 is capable of controlling. The position of each axis is continuously updated by monitoring its encoder, and displayed in a fixed format of a sign, two digits, a decimal point, and four digits.

2.02 Power Key.

The power key is used to store the current configuration to non-volatile memory before turning off power to the TCS8. Pressing the power key turns the displays off and saves the current configuration. Pressing it again turns the displays back on. This key can be used to implement a screen saver function.

2.03 Jog Control Keys.

These keys are used to control up to eight axes in a jog mode. The mode (slaved, one's only, or two's only) can be set through the jog menu. When the operator presses a jog key, the respective axis will begin to move. The direction that the axis moves is determined by the operator pressing either a plus or minus jog key. A plus jog key will turn the lead screw in a clockwise direction (away from the motor), a minus jog key will turn it in the counter-clockwise direction (towards the motor). By releasing the jog key, the operator stops motion on that axis. Motion will also stop if the axis reaches the limit for the direction it is moving, or if the indexer determines that the axis has stalled.

2.04 Scroll Keys.

These keys are used to scroll items through the MENU, COMMAND, and CHANNEL windows. All of the menus, their commands, and channel variations will be detailed in Section 3.00.

2.05 Command Windows.

These three windows (MENU, COMMAND, and CHANNEL) are used, in tandem with their respective scroll keys, to formulate a command to be executed by the TCS8.

2.06 Execute Key.

This key is used to execute the command currently formulated in the MENU, COMMAND, and CHANNEL windows.

2.07 Data Window.

Many of the TCS8's commands require some added data, e.g. the distance to move. Data for these commands are entered from the numeric key pad on the lower right of the TCS8 into the DATA window. Only a valid real number can be entered into the DATA window. If the operator enters an invalid real number, the character that is invalid will flash until the operator presses backspace or a valid character.

2.08 Stop Key.

The stop key, when pressed, will stop motion on all axes. The TCS8 will not lose track of the position of any axis. A move command started by the host computer and stopped by the stop key will finish normally with the position being reported. The position reported is the instantaneous position when the stop key was pressed. The final position of the axis being moved could be different than what was reported, thus the host computer should read the position again after a panic stop.

2.09 Status Window.

The STATUS window reflects the result of all commands. For commands that are not instantaneous, this window displays a busy status and then when the command completes it displays a ready status. The results of all view commands are displayed in the STATUS window. The STATUS window also displays the activity over the COM interfaces. For example, when the command for viewing position is sent over the COM1 interface, the STATUS window will display "COM1 VP" and when the command is completed the window will display "COM1 vp".

2.10 Numeric Key Pad.

The numeric key pad is used to enter a number into the data window. The user may backspace in the window or clear (shift-backspace) the window.

3.00 LOCAL COMMAND DESCRIPTIONS OF THE TCS8

This section describes the command set that can be executed from the front panel of the TCS8. Using the up and down keys under the MENU, COMMAND, and CHANNEL windows, the operator can formulate a command and then execute it by pressing the EXECUTE key. Some commands require extra information to be entered into the DATA window through the use of the numeric key pad. Each description includes a list of related commands that should be referred to in order to enhance the operator's understanding of the command. Also, where applicable, the default setting is given.

| Section | Title | Page |
|---------|-------------------------------|------|
| 3.01 | Move to Zero. | 10 |
| 3.02 | Move Absolute. | 10 |
| 3.03 | Move Relative. | 11 |
| 3.04 | Jog Mode. | 11 |
| 3.05 | Set Counts Per Unit. | 12 |
| 3.06 | Set Counts Per Revolution. | 12 |
| 3.07 | Set Position. | 13 |
| 3.08 | Set Velocity. | 13 |
| 3.09 | Set Acceleration. | 14 |
| 3.10 | Set Currents On. | 14 |
| 3.11 | Set Currents Off. | 14 |
| 3.12 | Set Inits On. | 15 |
| 3.13 | View Counts Per Unit. | 15 |
| 3.14 | View Counts Per Revolution. | 15 |
| 3.15 | View Velocity. | 16 |
| 3.16 | View Acceleration. | 16 |
| 3.17 | View Init. | 16 |
| 3.18 | View Currents. | 17 |
| 3.19 | View Plus Limit Switches. | 17 |
| 3.20 | View Minus Limit Switches. | 17 |
| 3.21 | View Home Switches. | 18 |
| 3.22 | View Stall Indication. | 18 |
| 3.23 | Init Default. | 18 |
| 3.24 | Init Drive On. | 18 |
| 3.25 | Init Drive Off. | 19 |
| 3.26 | COM1/COM2 Baud Rate. | 19 |
| 3.27 | COM1/COM2 Bits Per Character. | 20 |
| 3.28 | COM1/COM2 Parity. | 20 |
| 3.29 | COM1/COM2 Stop Bits. | 20 |
| 3.30 | COM1/COM2 Handshake. | 21 |

3.01 Move to Zero.

MENU: MOVE

COMMAND: TO ZERO

CHANNELS: ALL, ONE'S, TWO'S, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The MOVE TO ZERO command is an easy way to move some or all of the axes to the zero position. This command can also be accomplished with the MOVE ABSOLUTE command and a zero in the DATA window. Before using this command, selected axes must be initialized with the INIT DRIVE ON command. This command can be canceled by pressing the STOP key. When the STOP key is pressed, all axes will stop immediately. If an axis encounters a limit before reaching zero, the rest of its movement is aborted.

RELATED COMMANDS: MOVE ABSOLUTE, MOVE RELATIVE, INIT Drive ON.

3.02 Move Absolute.

MENU: MOVE

COMMAND: ABSOLUTE

CHANNELS: X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The MOVE ABSOLUTE command requires a position to be entered in the DATA window. This position and the current position of the axis is used to calculate the relative distance the axis must move. Before using this command, selected axes must be initialized with the INIT DRIVE ON command. This command can be canceled by pressing the STOP key. When the STOP key is pressed, all axes will stop immediately. If an axis encounters a limit before reaching the position entered in the DATA window, the rest of its movement is aborted.

RELATED COMMANDS: MOVE TO ZERO, MOVE RELATIVE, INIT Drive ON

3.03 Move Relative.

MENU: MOVE

COMMAND: RELATIVE

CHANNELS: X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The MOVE RELATIVE command requires a distance to be entered in the DATA window. This position is used to calculate the relative distance the axis must move. Before using this command, selected axes must be initialized with the INIT DRIVE ON command. This command can be canceled by pressing the STOP key. When the STOP key is pressed, all axes will stop immediately. If an axis encounters a limit before moving the distance entered in the DATA window, the rest of its movement is aborted.

RELATED COMMANDS: MOVE TO ZERO, MOVE ABSOLUTE, INIT Drive ON

3.04 Jog Mode.

MENU: JOG

COMMAND: MODE

CHANNELS: SLAVED, ONE'S, TWO'S

DESCRIPTION: The JOG MODE command sets the way the JOG keys operate. When SLAVED is the setting, both the one and two axis of the X, Y, Z, or A coordinate will move the same amount. When ONE'S is the setting, only the one axes of the X, Y, Z, or A coordinate will move. And finally, when TWO'S is the setting, only the two axes of the X, Y, Z, or A coordinate will move. The current mode is marked with an asterisk. After setting the jog mode, jogged movements can be made using the jog control keys. As with the other movement commands, the axis or axes that are to be jogged must be initialized with the INIT Drive ON command.

RELATED COMMANDS: INIT Drive ON

DEFAULT: SLAVED

3.05 Set Counts Per Unit.

MENU: SET

COMMAND: CPU

CHANNELS: ALL, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET CPU command allows the user to change the counts per unit travel. The CPU for an axis is determined by multiplying the encoder resolution (counts/revolution) by the lead screw resolution (revolutions/unit of travel). A units conversion can be added here to change, for example, from inches to centimeters. When the CPU for an axis is changed, the position is automatically converted. This command requires a value to be entered in the DATA window.

RELATED COMMANDS: SET CPR, SET POSITION

| | | | | |
|-----------------|----|------|----|------|
| DEFAULT: | X1 | 4000 | X2 | 4000 |
| | Y1 | 4000 | Y2 | 4000 |
| | Z1 | 4000 | Z2 | 4000 |
| | A1 | 4000 | A2 | 4000 |

3.06 Set Counts Per Revolution.

MENU: SET

COMMAND: CPR

CHANNELS: ALL, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET CPR command allows the user to change the encoder counts per motor revolution. The CPR for an axis is determined by dividing the encoder resolution (counts/revolution) by the lead screws resolution (revolutions/unit of travel). The encoder counts per motor revolution entered in the DATA window, must be a positive integer.

RELATED COMMANDS: SET CPU

| | | | | |
|-----------------|----|-----|----|-----|
| DEFAULT: | X1 | 400 | X2 | 400 |
| | Y1 | 400 | Y2 | 400 |
| | Z1 | 400 | Z2 | 400 |
| | A1 | 400 | A2 | 400 |

3.07 Set Position.

MENU: SET

COMMAND: POSITION

CHANNELS: ALL, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET POSITION command allows the user to change the current position of an axis without moving the axis. For example, the present position may be identified as zero or four inches or some other value. The new position must be entered in the DATA window before executing the command.

RELATED COMMANDS: SET CPU

3.08 Set Velocity.

MENU: SET

COMMAND: VELOCITY

CHANNELS: ALL, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET VELOCITY command allows the user to change the maximum speed at which an axis will travel. The range of valid velocities is 0.002 to 50.000 revolutions per second. The default is 5 revs/sec. Some stepper motor configurations will stall above a certain speed. To verify that a stall occurred, use the VIEW STALL command. When a stall happens, reduce the current velocity setting and continue normal operations. The new velocity must be entered in the DATA window before executing the command.

RELATED COMMANDS: SET ACCEL.

| | | | | |
|-----------------|----|-------|----|-------|
| DEFAULT: | X1 | 5.000 | X2 | 5.000 |
| | Y1 | 5.000 | Y2 | 5.000 |
| | Z1 | 5.000 | Z2 | 5.000 |
| | A1 | 5.000 | A2 | 5.000 |

3.09 Set Acceleration.

MENU: SET

COMMAND: ACCEL.

CHANNELS: ALL, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET ACCEL. command allows the user to change the maximum acceleration for an axis. The range of valid accelerations is 0.01 to 999.99 revolutions per second per second. The default is 5 revs/sec/sec. The new acceleration must be entered in the DATA window before executing the command.

RELATED COMMANDS: SET VELOCITY

| | | | | |
|-----------------|----|------|----|------|
| DEFAULT: | X1 | 5.00 | X2 | 5.00 |
| | Y1 | 5.00 | Y2 | 5.00 |
| | Z1 | 5.00 | Z2 | 5.00 |
| | A1 | 5.00 | A2 | 5.00 |

3.10 Set Currents On.

MENU: SET

COMMAND: CrntsOn

CHANNELS: ALL, ONE'S, TWO'S, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET CrntsOn command allows the user to turn the motor currents on. The motor current must be on for an axis to be moved. The information in the DATA window is ignored.

RELATED COMMANDS: SET CrntsOff

3.11 Set Currents Off.

MENU: SET

COMMAND: CrntsOff

CHANNELS: ALL, ONE'S, TWO'S, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET CrntsOff command allows the user to power down motors when they will not be used for long periods of time. The information in the DATA window is ignored.

RELATED COMMANDS: SET CrntsOn

3.12 Set Inits On.

MENU: SET

COMMAND: INITS ON

CHANNELS: ALL, ONE'S, TWO'S, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The SET INITS ON command allows the user to initialize the indexers without turning on the power to the motors. This command gives the indexers their velocity, acceleration, and counts per motor revolution information. The indexers must have this information before any movement can occur. This information needs only to be given once after powering up the system. The information in the DATA window is ignored.

RELATED COMMANDS: INIT Drive ON

3.13 View Counts Per Unit.

MENU: VIEW

COMMAND: Cnt/Unit

CHANNELS: X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The VIEW Cnt/Unit command displays the current setting of the encoder counts per unit travel parameter for the selected axis in the STATUS window. The information in the DATA window is ignored.

RELATED COMMANDS: SET CPU

3.14 View Counts Per Revolution.

MENU: VIEW

COMMAND: Cnt/MRev

CHANNELS: X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The VIEW Cnt/MRev command displays the current setting of the encoder counts per motor revolution parameter for the selected axis in the STATUS window. The information in the DATA window is ignored.

RELATED COMMANDS: SET CPR

3.15 View Velocity.

MENU: VIEW

COMMAND: VELOCITY

CHANNELS: X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The VIEW VELOCITY command displays the current setting of the velocity parameter for the selected axis in the STATUS window. The information in the DATA window is ignored.

RELATED COMMANDS: SET VELOCITY

3.16 View Acceleration.

MENU: VIEW

COMMAND: ACCEL.

CHANNELS: X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The VIEW ACCEL. command displays the current setting of the acceleration parameter for the selected axis in the STATUS window. The information in the DATA window is ignored.

RELATED COMMANDS: SET ACCEL.

3.17 View Init.

MENU: VIEW

COMMAND: INIT

CHANNELS: none

DESCRIPTION: The VIEW INIT command uses the STATUS window to display a one(initialized) or a zero(uninitialized) for each axis. The STATUS window has eight characters; left to right respectively reflecting the status of X1, X2 ..., A1, A2. The information in the DATA window is ignored.

RELATED COMMANDS: SET INITs, INIT Drive ON

3.18 View Currents.

MENU: VIEW

COMMAND: CURRENTS

CHANNELS: none

DESCRIPTION: The VIEW CURRENTS command uses the STATUS window to display a one(current on) or a zero(current off) for each axis. The STATUS window has eight characters; left to right respectively reflecting the status of X1, X2 ..., A1, A2. The information in the DATA window is ignored.

RELATED COMMANDS: SET CrntsOn, SET CrntsOff, INIT Drive ON, INIT Drive OFF

3.19 View Plus Limit Switches.

MENU: VIEW

COMMAND: Plus LMT

CHANNELS: none

DESCRIPTION: The VIEW Plus LMT command uses the STATUS window to display a one(on limit) or a zero(not on limit) for each axis. The STATUS window has eight characters; left to right respectively reflecting the status of X1, X2 ..., A1, A2. The plus limit switches are located at the positive movement end of travel. The information in the DATA window is ignored.

RELATED COMMANDS: VIEW Minus LMT, VIEW HOME

3.20 View Minus Limit Switches.

MENU: VIEW

COMMAND: Minus LMT

CHANNELS: none

DESCRIPTION: The VIEW Minus LMT command uses the STATUS window to display a one(on limit) or a zero(not on limit) for each axis. The STATUS window has eight characters; left to right respectively reflecting the status of X1, X2 ..., A1, A2. The minus limit switches are located at the negative movement end of travel. The information in the DATA window is ignored.

RELATED COMMANDS: VIEW Plus LMT, VIEW HOME

3.21 View Home Switches.

MENU: VIEW

COMMAND: HOME

CHANNELS: none

DESCRIPTION: The VIEW HOME command uses the STATUS window to display a one(on limit) or a zero(not on limit) for each axis. The STATUS window has eight characters; left to right respectively reflecting the status of X1, X2 ..., A1, A2. The home limit switch can be adjusted by the user for application specific tasks. The information in the DATA window is ignored.

RELATED COMMANDS: VIEW Plus LMT, VIEW Minus LMT

3.22 View Stall Indication.

MENU: VIEW

COMMAND: STALL

CHANNELS: none

DESCRIPTION: The VIEW STALL command uses the STATUS window to display a one(stalled) or a zero(not stalled) for each axis. The STATUS window has eight characters; left to right respectively reflecting the status of X1, X2 ..., A1, A2. A stall is indicated when the indexer is making a move and the amount of pulses send to the motor does not match the corresponding number of pulses received from the encoder. A stall can occur if the velocity or acceleration is set to high, the encoder counts per motor revolution are set incorrectly, or the axis is physically jammed. The information in the DATA window is ignored.

RELATED COMMANDS: none

3.23 Init Default.

MENU: INIT

COMMAND: DEFAULT

CHANNELS: none

DESCRIPTION: The INIT DEFAULT command restores the initial factory defaults (CPU, CPR, VELOCITY, ACCELERATION,BAUD RATE,BITS/CHAR,PARITY,STOP BITS,HANDSHAKE) of the TCS8. After executing this command, execute the command INIT Drive ON to initialize the indexers. The information in the DATA window is ignored.

RELATED COMMANDS: SET CPU, SET CPR, SET VELOCITY, SET ACCEL.

3.24 Init Drive On.

MENU: INIT

COMMAND: Drive ON

CHANNELS: ALL, ONE'S, TWO'S, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The INIT Drive ON command initializes the selected axes for movement. This command does the same thing as SET INITS ON except that it also turns on the current to the motors. The information in the DATA window is ignored.

RELATED COMMANDS: SET CPU, SET CPR, SET VELOCITY, SET ACCEL., SET CntrsOn, SET CntrsOff, INIT DEFAULT

3.25 Init Drive Off.

MENU: INIT

COMMAND: Drive OFF

CHANNELS: ALL, ONE'S, TWO'S, X1&X2, Y1&Y2, Z1&Z2, A1&A2, X1, X2, Y1, Y2, Z1, Z2, A1, A2

DESCRIPTION: The INIT Drive OFF command is an alias for SET CntrsOff.

RELATED COMMANDS: SET CntrsOff

3.26 COM1/COM2 Baud Rate.

MENU: COM1/COM2

COMMAND: BaudRate

CHANNELS: 19.2K, 9600, 4800, 2400, 1200, 300, 110

DESCRIPTION: The COM1/COM2 BaudRate commands set the baud rate for the selected communication channel. The information in the DATA window is ignored. The current baud rate is marked with an asterisk.

RELATED COMMANDS: none

DEFAULT: 9600

3.27 COM1/COM2 Bits Per Character.

MENU: COM1/COM2

COMMAND: Bit/Char

CHANNELS: SEVEN, EIGHT

DESCRIPTION: The COM1/COM2 Bit/Char command set the bits per character for the selected communication channel. The information in the DATA window is ignored. The current number of bits per character is marked with an asterisk.

RELATED COMMANDS: none

DEFAULT: EIGHT

3.28 COM1/COM2 Parity.

MENU: COM1/COM2

COMMAND: Parity

CHANNELS: NONE, EVEN, ODD

DESCRIPTION: The COM1/COM2 Parity command set the parity for the selected communication channel. The information in the DATA window is ignored. The current parity is marked with an asterisk.

RELATED COMMANDS: none

DEFAULT: EVEN

3.29 COM1/COM2 Stop Bits.

MENU: COM1/COM2

COMMAND: StopBits

CHANNELS: 1, 1.5, 2

DESCRIPTION: The COM1/COM2 StopBits command set the stop bits for the selected communication channel. The information in the DATA window is ignored. The current number of stop bits is marked with an asterisk.

RELATED COMMANDS: none

DEFAULT: 1

3.30 COM1/COM2 Handshake.

MENU: COM1/COM2

COMMAND: HandShak

CHANNELS: NO, YES

DESCRIPTION: The COM1/COM2 HandShak command set the handshake for the selected communication channel. The information in the DATA window is ignored. An asterisk marks whether there is handshaking or not.

RELATED COMMANDS: none

DEFAULT: YES

4.00 SERIAL INTERFACE COMMAND DESCRIPTIONS OF THE TCS8

This section describes the command set that can be executed through the serial interfaces of the TCS8. Each description includes a code section that outlines the characters that must be sent to execute the command. The vertical bar in this section is used as a separator and is not sent as part of the command code. The symbol "CRLF" stands for the two characters carriage return and line feed. Also where applicable, the default setting is given.

| Section | Title | Page |
|---------|---|------|
| 4.01 | Change Serial Configuration. | 23 |
| 4.02 | Move to Absolute Position. | 24 |
| 4.03 | Move Relative to Current Position. | 24 |
| 4.04 | Set Acceleration. | 25 |
| 4.05 | View Acceleration. | 25 |
| 4.06 | Set Velocity. | 26 |
| 4.07 | View Velocity. | 26 |
| 4.08 | Set Encoder Counts Per Unit of Travel. | 27 |
| 4.09 | View Encoder Counts Per Unit of Travel. | 27 |
| 4.10 | Set Counts Per Motor Revolution. | 28 |
| 4.11 | View Counts Per Motor Revolution. | 28 |
| 4.12 | Set Position. | 29 |
| 4.13 | View Position. | 29 |
| 4.14 | Set Current to Motor Windings. | 30 |
| 4.15 | View Current to Motor Windings. | 30 |
| 4.16 | Set Initialization of Indexer/Drivers. | 31 |
| 4.17 | View Initialization of Indexer/Drivers. | 31 |

4.01 Change Serial Configuration.

COMMAND: CHANGE SERIAL CONFIGURATION

CODE: CS COM;CATEGORY;ATTRIBUTE;

| | | |
|--------------------|-------------------|---|
| PARAMETERS: | COM: | 1/COM1 2/COM2 |
| | CATEGORY: | 0/BAUDRATE |
| | ATTRIBUTE: | 0/19.2K 1/9600 2/4800 3/2400 4/1200 5/300 6/110 |
| | CATEGORY: | 1(BITS PER CHARACTER) |
| | ATTRIBUTE: | 0/SEVEN 1/EIGHT |
| | CATEGORY: | 2(PARITY) |
| | ATTRIBUTE: | 0/NONE 1/EVEN 2/ODD |
| | CATEGORY: | 3(STOP BITS) |
| | ATTRIBUTE: | 0/ONE 1/ONE AND A HALF 2/TWO |
| | CATEGORY: | 4(HANDSHAKE) |
| | ATTRIBUTE: | 0/NO 1/YES |

DESCRIPTION: This command must be executed with extreme caution and forethought. If the user changes an attribute of the same COM port from which he is sending the command, he must change to that attribute on the host computer before sending the next command. The best way to change the serial configuration of a COM port is to utilize the front panel commands.

DEFAULT: 9600 baud, EIGHT bits/char, EVEN parity, ONE stop bit, handshaking YES

EXAMPLE: To change the baudrate of COM1 to 2400 the user must send CS1;0;3;

4.02 Move to Absolute Position.

COMMAND: MOVE TO ABSOLUTE POSITION AND REPORT FINAL POSITION

CODE: MA CHANNEL:POSITION,CHANNEL:POSITION,...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
1/X1
2/X2
3/Y1
4/Y2
5/Z1
6/Z2
7/A1
8/A2
POSITION: Real number free format

DESCRIPTION: This command moves selected channels to absolute positions.

EXAMPLES: To move all channels to zero the user may send MA0:0,CRLF or MA12345678:0,CRLF. To move channel X1 to zero the user must send MA1:0,CRLF. To move channels X1 and X2 to zero the user may send MA12:0,CRLF or MA1:0,2:0,CRLF or MA1:0,CRLF and MA2:0,CRLF.

4.03 Move Relative to Current Position.

COMMAND: MOVE TO RELATIVE DISTANCE AND REPORT FINAL POSITION

CODE: MR CHANNEL:DISTANCE,CHANNEL:DISTANCE,...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
1/X1
2/X2
3/Y1
4/Y2
5/Z1
6/Z2
7/A1
8/A2
POSITION: Real number free format

DESCRIPTION: This command moves selected channels relative distances.

EXAMPLES: To move all channels one unit the user may send MR0:1,CRLF or MR12345678:1,CRLF. To move channel X1 one unit the user must send MR1:1,CRLF. To move channels X1 and X2 one unit the user may send MR12:1,CRLF or MR1:1,2:1,CRLF or MR1:1,CRLF and MR2:1,CRLF.

4.04 Set Acceleration.

COMMAND: SET ACCELERATION

CODE: SA CHANNEL:ACCELERATION,CHANNEL:ACCELERATION,...\CRLF

| | | |
|--------------------|----------------------|--|
| PARAMETERS: | CHANNEL: | 0/ALL CHANNELS |
| | | 1/X1 |
| | | 2/X2 |
| | | 3/Y1 |
| | | 4/Y2 |
| | | 5/Z1 |
| | | 6/Z2 |
| | | 7/A1 |
| | | 8/A2 |
| | ACCELERATION: | Real number free format between 0.01 and 99.99 inclusive. |

DESCRIPTION: This command sets the acceleration for selected channels.

DEFAULT: All channels 5.00 revolutions/second/second

EXAMPLES: To set the acceleration for all channels to 4.00 revolutions/second/second the user may send SA0:4.00,CRLF or SA12345678:4.00,CRLF. To set the acceleration for channel X1 to 4.00 revolutions/second/second the user must send SA1:4.00,CRLF. To set the acceleration for channels X1 and X2 to 4.00 revolutions/second/second the user may send SA12:4.00,CRLF or SA1:4.00 ,2:4.00,CRLF or SA1:4.00,CRLF and SA2:4.00,CRLF.

4.05 View Acceleration.

COMMAND: VIEW ACCELERATION

CODE: VA CHANNEL\CHANNEL...\CRLF

| | | |
|--------------------|-----------------|----------------|
| PARAMETERS: | CHANNEL: | 0/ALL CHANNELS |
| | | 1/X1 |
| | | 2/X2 |
| | | 3/Y1 |
| | | 4/Y2 |
| | | 5/Z1 |
| | | 6/Z2 |
| | | 7/A1 |
| | | 8/A2 |

DESCRIPTION: This command views the acceleration for selected channels. The TCS8 transmits each of the accelerations requested back to the host computer separated by carriage return line feeds.

EXAMPLES: To view the acceleration for all channels the user may send VA0CRLF or VA12345678CRLF. To view the acceleration for channel X1 the user must send VA1CRLF. To view the acceleration for channels X1 and X2 the user may send VA12CRLF or VA1CRLF and VA2CRLF.

4.06 Set Velocity.

COMMAND: SET VELOCITY

CODE: SV CHANNEL:VELOCITY,CHANNEL:VELOCITY,...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
1/X1
2/X2
3/Y1
4/Y2
5/Z1
6/Z2
7/A1
8/A2
VELOCITY: Real number free format between
0.001 and 50.000 inclusive.

DESCRIPTION: This command sets the velocity for selected channels.

DEFAULT: All channels 5.000 revolutions/second

EXAMPLES: To set the velocity for all channels to 4.00 revolutions/second the user may send SV0:4.00,CRLF or SV12345678:4.00,CRLF. To set the velocity for channel X1 to 4.00 revolutions/second the user must send SV1:4.00,CRLF. To set the velocity for channels X1 and X2 to 4.00 revolutions/second the user may send SV12:4.00,CRLF or SV1:4.00 ,2:4.00,CRLF or SV1:4.00,CRLF and SV2:4.00,CRLF.

4.07 View Velocity.

COMMAND: VIEW VELOCITY

CODE: VV CHANNEL|CHANNEL...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
1/X1
2/X2
3/Y1
4/Y2
5/Z1
6/Z2
7/A1
8/A2

DESCRIPTION: This command views the velocity for selected channels. The TCS8 transmits each of the velocities requested back to the host computer separated by carriage return line feeds.

EXAMPLES: To view the velocity for all channels the user may send VV0CRLF or VV12345678CRLF. To view the velocity for channel X1 the user must send VV1CRLF. To view the velocity for channels X1 and X2 the user may send VV12CRLF or VV1CRLF and VV2CRLF.

4.08 Set Encoder Counts Per Unit of Travel.

COMMAND: SET ENCODER COUNTS PER UNIT TRAVEL

CODE: SU CHANNEL:CPU,CHANNEL:CPU,...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
1/X1
2/X2
3/Y1
4/Y2
5/Z1
6/Z2
7/A1
8/A2
CPU: Non-zero real number free format.

DESCRIPTION: This command sets the encoder counts per unit travel for selected channels.

DEFAULT: X1,X2,Y1,Y2,Z1,Z2,A1, and A2 4000 counts/inch

EXAMPLES: To set the encoder counts per unit travel for all channels to 5000 the user may send SU0:5000,CRLF or SU12345678:5000,CRLF. To set the encoder counts per unit travel for channel X1 to 5000 the user must send SU1:5000,CRLF. To set the encoder counts per unit travel for channels X1 and X2 to 5000 the user may send SU12:5000,CRLF or SU1:5000 ,2:5000,CRLF or SU1:5000,CRLF and SU2:5000,CRLF.

4.09 View Encoder Counts Per Unit of Travel.

COMMAND: VIEW ENCODER COUNTS PER UNIT TRAVEL

CODE: VU CHANNEL|CHANNEL...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
1/X1
2/X2
3/Y1
4/Y2
5/Z1
6/Z2
7/A1
8/A2

DESCRIPTION: This command views the encoder counts per unit travel for selected channels. The TCS8 transmits each of the encoder counts per unit travel requested back to the host computer separated by carriage return line feeds.

EXAMPLES: To view the encoder counts per unit travel for all channels the user may send VU0CRLF or VU12345678CRLF. To view the encoder counts per unit travel for channel X1 the user must send VU1CRLF. To view the encoder counts per unit travel for channels X1 and X2 the user may send VU12CRLF or VU1CRLF and VU2CRLF.

4.10 Set Counts Per Motor Revolution.

COMMAND: SET ENCODER COUNTS PER MOTOR REVOLUTION

CODE: SR CHANNEL:CPR,CHANNEL:CPR,...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
1/X1
2/X2
3/Y1
4/Y2
5/Z1
6/Z2
7/A1
8/A2
CPU: Non-zero integer free format.

DESCRIPTION: This command sets the encoder counts per motor revolution for selected channels.

DEFAULT: X1,X2,Y1,Y2,Z1,Z2 and A1,A2 400 counts/inch

EXAMPLES: To set the encoder counts per motor revolution for all channels to 500 the user may send SR0:500,CRLF or SR12345678:500,CRLF. To set the encoder counts per motor revolution for channel X1 to 500 the user must send SR1:500,CRLF. To set the encoder counts per motor revolution for channels X1 and X2 to 500 the user may send SR12:500,CRLF or SR1:500,2:500,CRLF or SR1:500,CRLF and SR2:500,CRLF.

4.11 View Counts Per Motor Revolution.

COMMAND: VIEW ENCODER COUNTS PER MOTOR REVOLUTION

CODE: VR CHANNEL|CHANNEL...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
1/X1
2/X2
3/Y1
4/Y2
5/Z1
6/Z2
7/A1
8/A2

DESCRIPTION: This command views the encoder counts per motor revolution for selected channels. The TCS8 transmits each of the encoder counts per motor revolution requested back to the host computer separated by carriage return line feeds.

EXAMPLES: To view the encoder counts per motor revolution for all channels the user may send VR0CRLF or VR12345678CRLF. To view the encoder counts per motor revolution for channel X1 the user must send VR1CRLF. To view the encoder counts per motor revolution for channels X1 and X2 the user may send VR12CRLF or VR1CRLF and VR2CRLF.

4.12 Set Position.

COMMAND: SET POSITION

CODE: SP CHANNEL:POSITION,CHANNEL:POSITION,...|CRLF

| | | |
|--------------------|------------------|----------------|
| PARAMETERS: | CHANNEL: | 0/ALL CHANNELS |
| | | 1/X1 |
| | | 2/X2 |
| | | 3/Y1 |
| | | 4/Y2 |
| | | 5/Z1 |
| | | 6/Z2 |
| | | 7/A1 |
| | | 8/A2 |
| | POSITION: | real number. |

DESCRIPTION: This command sets the position for selected channels.

EXAMPLES: To set the position for all channels to 1.5 the user may send SP0:1.5,CRLF or SP12345678:1.5,CRLF. To set the position for channel X1 to 1.5 the user must send SP1:1.5,CRLF. To set the position for channels X1 and X2 to 1.5 the user may send SP12:1.5,CRLF or SP1:1.5,2:1.5,CRLF or SP1:1.5,CRLF and SP2:1.5,CRLF.

4.13 View Position.

COMMAND: VIEW POSITION

CODE: VP CHANNEL|CHANNEL...|CRLF

| | | |
|--------------------|-----------------|----------------|
| PARAMETERS: | CHANNEL: | 0/ALL CHANNELS |
| | | 1/X1 |
| | | 2/X2 |
| | | 3/Y1 |
| | | 4/Y2 |
| | | 5/Z1 |
| | | 6/Z2 |
| | | 7/A1 |
| | | 8/A2 |

DESCRIPTION: This command views the position for selected channels. The TCS8 transmits each of the positions requested back to the host computer separated by carriage return line feeds.

EXAMPLES: To view the position for all channels the user may send VP0CRLF or VP12345678CRLF. To view the position for channel X1 the user must send VP1CRLF. To view the position for channels X1 and X2 the user may send VP12CRLF or VP1CRLF and VP2CRLF.

4.14 Set Current to Motor Windings.

COMMAND: SET CURRENT TO MOTOR WINDINGS

CODE: SC CHANNEL:ON/OFF,CHANNEL:ON/OFF,...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
1/X1
2/X2
3/Y1
4/Y2
5/Z1
6/Z2
7/A1
8/A2
ON/OFF: 1/ON
0/OFF

DESCRIPTION: This command sets the current to the motor windings for selected channels on or off.

EXAMPLES:

To set the current to the motor windings for all channels on the user may send SC0:1,CRLF or SC12345678:1,CRLF to set them off the user may send SC0:0,CRLF or SC12345678:0,CRLF.

4.15 View Current to Motor Windings.

COMMAND: VIEW CURRENT TO MOTOR WINDINGS

CODE: VC CHANNEL|CHANNEL...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
1/X1
2/X2
3/Y1
4/Y2
5/Z1
6/Z2
7/A1
8/A2

DESCRIPTION: This command views the current to the motor windings for selected channels. The TCS8 transmits each response of on/off (1/0) back to the host computer separated by carriage return line feeds.

EXAMPLES:

To view the current to the motor windings for all channels the user may send VC0CRLF or VC12345678CRLF

To view the current to the motor windings for channel X1 the user must send VC1CRLF

To view the current to the motor windings for channels X1 and X2 the user may send VC12CRLF or VC1CRLF and VC2CRLF.

4.16 Set Initialization of Indexer/Drivers.

COMMAND: SET INITIALIZATION OF INDEXER/DRIVERS

CODE: SI CHANNEL|CHANNEL...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
1/X1
2/X2
3/Y1
4/Y2
5/Z1
6/Z2
7/A1
8/A2

DESCRIPTION: This command sends the current value of the acceleration, velocity, and the encoder counts per motor revolution to the indexer/driver for the selected channels. This command must be sent before any move commands may be sent.

EXAMPLES:

To initialize all channels the user may send SI0CRLF or SI12345678CRLF

To initialize channel X1 the user must send SI1CRLF

To initialize channels X1 and X2 the user may send SI12CRLF or SI1CRLF and SI2CRLF

4.17 View Initialization of Indexer/Drivers.

COMMAND: VIEW INITIALIZATION OF INDEXER/DRIVERS

CODE: VI CHANNEL|CHANNEL...|CRLF

PARAMETERS: CHANNEL: 0/ALL CHANNELS
1/X1
2/X2
3/Y1
4/Y2
5/Z1
6/Z2
7/A1
8/A2

DESCRIPTION: This command returns "1" if the indexer/driver has been initialized since the TCS8 was turned on and "0" if it has not. The TCS8 transmits each of the responses back to the host computer separated by carriage return line feeds.

EXAMPLES:

To check the initialization of all channels the user may send VI0CRLF or VI12345678CRLF

To check the initialization of channel X1 the user must send VI1CRLF

To check the initialization of channels X1 and X2 the user may send VI12CRLF or VI1CRLF and VI2CRLF

5.00 Traverse Control System Cables.




Figure 6 Motor Drive System to Compumotor Stepper Motor Cable.




Figure 7 Motor Drive System to Dynamics Research Encoder Cable.




Figure 8 Motor Drive System to Linear Industries Limit Switches Cable.




Figure 9 Motor Drive System to TCS8 Encoder Signals Cable.




Figure 10 HP Series 9000 Model 375 to TCS8 Serial Cable.




Figure 11 TCS8 to Motor Drive System Serial Cable.

APPENDIX A

ORIGINAL SOFTWARE CODE LISTING.

APPENDIX A

Original Software Code Listing.

TABLE OF CONTENTS

| Item | Title | Page |
|-------------|--|-------------|
| 1 | Original Hard Disk Directory Catalog Listing.. | 2 |
| 2 | Original Program "3.5'HWT91" Hard Copy Code Listing. | 3 |

Hard Disk Directory Catalog Listing.

:CS80, 1400, 0, 0

VOLUME LABEL: B9826

| FILE NAME | PRO | TYPE | REC/FILE | BYTE/REC | ADDRESS | DATE | TIME |
|-----------|-----|-------|----------|----------|---------|-----------|-------|
| SYSB60 | | SYSTM | 3388 | 256 | 32 | 17-Jul-91 | 13:10 |
| CDUMP6 | | PROG | 44 | 256 | 3420 | 17-Jul-91 | 13:10 |
| BPLOT6 | | PROG | 40 | 256 | 3464 | 17-Jul-91 | 13:10 |
| AUTOST | | PROG | 10 | 256 | 3504 | 17-Jul-91 | 13:10 |
| ARRAY | | BDAT | 50 | 256 | 3515 | 17-Jul-91 | 13:11 |
| KEYS | | BDAT | 4 | 256 | 3566 | 17-Jul-91 | 13:11 |
| COPY | | PROG | 25 | 256 | 3570 | 17-Jul-91 | 13:11 |
| 3.5'HWT91 | | PROG | 372 | 256 | 3595 | 17-Jan-92 | 16:03 |

```

100 ! *****
110 ! *
120 ! *      Property of COMPLERE INC.
130 ! *      Proprietary software
140 ! *      Copyright June 18, 1991
150 ! *      Developed by: T. Kevin McDevitt
160 ! *
170 ! *      LASER DOPPLER VELOCIMETER TEST
180 ! *
190 ! *      3.5 FOOT HYPER SONIC WIND TUNNEL
200 ! *      NASA AMES RESEARCH CENTER
210 ! *
220 ! *****
230 !
240 DEG
250 !
260 OPTION BASE 1
270 COM /Data/ INTEGER Raw(1000,10),Valid(1000),REAL Table(0:32766),Ui(1000),Vi(1000),Wi(1000),Ai(1000),
    Bi(1000),Ii(1000),Ci(1000)
280 COM /Array/ Names(100,4)[10],Image$(100,4)[10],Units$(100,4)[10],REAL Array(100,4)
290 COM /Pos/ Pnames(25,1)[10],Pimage$(25,1)[10],Punits$(25,1)[10],REAL Pos(25,1),Npos
300 COM /Graph/ Wndw(9,4),Vwprt(9,4),Xdiv(9),Ydiv(9),Xlabel$(9)[80],Ylabel$(9)[80],Titles(9)[80],
    Ximage$(9)[80],Yimage$(9)[80],Legend$(9,5)[80]
310 COM Run,File,Paxis
320 !
330 DIM Menu$(5,8)[80],Systems$(20),Data$(20),Files$(50),L$(160)
340 INTEGER Gsave(1280,1024),At_exp,Ct_exp,Cmask,Nsam,N(10,3)
350 REAL Atime,Ctime,Sum(10,3),Symbols(5,0:20,3)
360 !
370 DIM Tcs2tun1(4,4),Tun2tcs1(4,4),Tun2mod(3,3),Tun2ldv(3,3),Tun(4),Tcs1(4)
380 DIM Tcs2tun2(4,4),Tun2tcs2(4,4),Mod2tun(3,3),Ldv2tun(3,3),Mod(4),Tcs2(4)
390 !
400 DIM Beam_spc(3),Foc1_len(3),Mea_sgn(3),Mix_frq(3),Mix_sgn(3),Frng_spc(3),Thata(3,3)
410 DIM Beam_sep(3),Wave_len(3),Brq_frq(3),Brq_sgn(3),Index(3),Coin(3)
420 !
430 PRINTER IS CRT
440 CLEAR SCREEN
450 GCLEAR
460 !
470 GOSUB Lvds_set_up
480 GOSUB File_set_up
490 GOSUB Tcs8_set_up
500 GOSUB Menu_set_up
510 GOSUB Grph_set_up
520 GOTO 580
530 CLEAR SCREEN
540 CALL Tcs8move(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Mod2tun(*),Tun2tcs1(*),Tun2tcs2(*),"Tx & Rx","LASER",
    "ABSOLUTE",2,.250)
550 CALL Tcs8move(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Mod2tun(*),Tun2tcs1(*),Tun2tcs2(*),"Tx & Rx","LASER",
    "ABSOLUTE",2,.251)
560 GOTO 540
570 !
580 Date=TIMEDATE
590 Time=Date
600 CALL Purge(System$,Data$)
610 Here: DISP TIMES(TIMEDATE),DATES(TIMEDATE)
620 GOTO Here
630 STOP
640 File_set_up: System$=":",1400,0,0"
650 Data$=":",1400,0,1"
660 LOAD KEY "KEYS"&System$
670 IF NOT INMEM("Gdump_colored") THEN LOADSUB ALL FROM "CDUMP6"&System$
680 IF NOT INMEM("Bload") THEN LOADSUB ALL FROM "BPLOT6"&System$
690 IF NOT INMEM("Bstore") THEN LOADSUB ALL FROM "BPLOT6"&System$
700 GOSUB Read_array
710 GOSUB Read_calc_fill
720 GOSUB Save_array
730 CLEAR SCREEN
740 RETURN
750 Grph_set_up: CALL Read_symbols(Symbols(*))
760 CALL Crt_init
770 CALL Setup_graph(Array(*),Image$(*),Paxis,Symbols(*))
780 RETURN
790 Menu_set_up: CALL Menu_read(Menu$(*))
800 CALL Menu_disp(Menu,Menu$(*))
810 GOSUB On_key
820 Busy=0
830 Ready=1
840 RETURN
850 Tcs8_set_up: CALL Tcs8init(@Tcs8)

```

```

860 CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
870 GOSUB Calc
880 GOSUB Fill
890 RETURN
900 Lvds_set_up: CALL Lvds_init(@Lvds)
910 CALL Table(Table(*))
920 RETURN
930 On_key: ON KEY 1 GOSUB Key1
940 ON KEY 2 GOSUB Key2
950 ON KEY 3 GOSUB Key3
960 ON KEY 4 GOSUB Key4
970 ON KEY 5 GOSUB Key5
980 ON KEY 6 GOSUB Key6
990 ON KEY 7 GOSUB Key7
1000 ON KEY 8 GOSUB Key8
1010 RETURN
1020 Key1: Key=1
1030 CALL Menu_status(Menu,Key,Busy,Menu$(*))
1040 ON Menu GOSUB M1k1,M2k1,M3k1,M4k1,M5k1,M6k1,M7k1
1050 CALL Menu_status(Menu,Key,Ready,Menu$(*))
1060 CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
1070 RETURN
1080 Key2: Key=2
1090 CALL Menu_status(Menu,Key,Busy,Menu$(*))
1100 ON Menu GOSUB M1k2,M2k2,M3k2,M4k2,M5k2,M6k2,M7k2
1110 CALL Menu_status(Menu,Key,Ready,Menu$(*))
1120 CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
1130 RETURN
1140 Key3: Key=3
1150 CALL Menu_status(Menu,Key,Busy,Menu$(*))
1160 ON Menu GOSUB M1k3,M2k3,M3k3,M4k3,M5k3,M6k3,M7k3
1170 CALL Menu_status(Menu,Key,Ready,Menu$(*))
1180 CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
1190 RETURN
1200 Key4: Key=4
1210 CALL Menu_status(Menu,Key,Busy,Menu$(*))
1220 ON Menu GOSUB M1k4,M2k4,M3k4,M4k4,M5k4,M6k4,M7k4
1230 CALL Menu_status(Menu,Key,Ready,Menu$(*))
1240 CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
1250 RETURN
1260 Key5: Key=5
1270 CALL Menu_status(Menu,Key,Busy,Menu$(*))
1280 ON Menu GOSUB M1k5,M2k5,M3k5,M4k5,M5k5,M6k5,M7k5
1290 CALL Menu_status(Menu,Key,Ready,Menu$(*))
1300 CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
1310 RETURN
1320 Key6: Key=6
1330 CALL Menu_status(Menu,Key,Busy,Menu$(*))
1340 ON Menu GOSUB M1k6,M2k6,M3k6,M4k6,M5k6,M6k6,M7k6
1350 CALL Menu_status(Menu,Key,Ready,Menu$(*))
1360 CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
1370 RETURN
1380 Key7: Key=7
1390 CALL Menu_status(Menu,Key,Busy,Menu$(*))
1400 ON Menu GOSUB M1k7,M2k7,M3k7,M4k7,M5k7,M6k7,M7k7
1410 CALL Menu_status(Menu,Key,Ready,Menu$(*))
1420 CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
1430 RETURN
1440 Key8: Key=8
1450 CALL Menu_status(Menu,Key,Busy,Menu$(*))
1460 ON Menu GOSUB M1k8,M2k8,M3k8,M4k8,M5k8,M6k8,M7k8
1470 CALL Menu_status(Menu,Key,Ready,Menu$(*))
1480 CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
1490 RETURN
1500 M1k1: Menu=2
1510 CALL Menu_disp(Menu,Menu$(*))
1520 RETURN
1530 M1k2: Menu=3
1540 CALL Menu_disp(Menu,Menu$(*))
1550 RETURN
1560 M1k3: KEY LABELS OFF
1570 PRINTER IS CRT;WIDTH 132
1580 DISP ""
1590 FOR L=1 TO 9
1600 PRINT TABXY(L,L);RPTS(" ",120)
1610 NEXT L
1620 PRINTER IS PRT
1630 PRINT USING "#,@"
1640 DUMP GRAPHICS
1650 PRINT USING "#,@"

```

```

1660      PRINTER IS CRT
1670      CALL Menu_disp(Menu,Menu$(*))
1680      RETURN
1690 M1k4:  CALL Enter_value("number of traverse positions",Npos,"K")
1700      REDIM Pos(Npos,1),Pname$(Npos,1),Pimage$(Npos,1),Punits$(Npos,1)
1710      MAT Pimage$= ("M4D.4D")
1720      MAT Punits$= ("in")
1730      FOR K=1 TO Npos
1740          Pname$(K,1)="Pos#"+VAL$(K)
1750      NEXT K
1760      GSTORE Gsave(*)
1770      CALL Change("VALUES",Pos(*),Pname$(*),Pimage$(*),Punits$(*))
1780      GLOAD Gsave(*)
1790      CALL Menu_disp(Menu,Menu$(*))
1800      RETURN
1810 M1k5:  GOSUB Read_calc_fill
1820      CALL Enter_value(CHRS$(NUM("X")+Paxis-1),Mod(Paxis),"K")
1830 M1k5a: ON KBD CALL Do_nothing
1840      DISP "Moving"
1850      Movement=Mod(Paxis)
1860      CALL Tcs8move(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Mod2tun(*),Tun2tcs1(*),Tun2tcs2(*),"Tx &
          Rx","MODEL","ABSOLUTE",Paxis,Movement)
1870      CALL Tcs8print(Mod(*),Tun(*),Tcs1(*),Tcs2(*))
1880      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
1890      GOSUB Calc
1900      GOSUB Fill
1910      DISP ""
1920      OFF KBD
1930      RETURN
1940 M1k6:  CALL Fix(Array(*),Name$(*),Image$(*),Units$(*))
1950      DISP "Press any key to TAKE DATA"
1960      CALL Rt_histo(@Lvdas,Symbols(*),1)
1970      Cmask=Coin(1)*1+Coin(2)*2+Coin(3)*4
1980      Nsam=MIN(Nreads,1000)
1990      Date=TIMEDATE
2000      Time=Date
2010      Atime=10
2020      CALL Lvdas_take(@Lvdas,Atime,Ctime,At_exp,Ct_exp,Cmask,Nsam)
2030      IF Nsam>1 THEN
2040          OUTPUT PRT;RPTS("=",140)
2050          CALL Data_reduce(At_exp,Ct_exp,Nsam)
2060          !CALL Data_histo(Array(*),Nsam)
2070          CALL Pt_histo(Symbols(*),Run,File,Mod(Paxis),Nsam)
2080          CALL Data_clip(Nsam,Umin,Umax,Vmin,Vmax)
2090          CALL Data_sum(Sum(*),N(*),Nsam)
2100          CALL Data_calc(N(*),Sum(*),U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,A1b1,U1a1,V1a1,W1a1)
2110          CALL
          Data_print(Paxis,Mod(Paxis),Nsam,"MHz",U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,A1b1,U1a1,V1
          a1,W1a1)
2120          CALL Data_fconvert(Array(*))
2130          CALL Data_aconvert(Gain)
2140          CALL Data_sum(Sum(*),N(*),Nsam)
2150          CALL Data_calc(N(*),Sum(*),U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,A1b1,U1a1,V1a1,W1a1)
2160          CALL
          Data_print(Paxis,Mod(Paxis),Nsam,"LDV",U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,A1b1,U1a1,V1
          a1,W1a1)
2170          CALL Data_trnsfrm(Ldv2tun(*),U,V,W,U1,V1,W1,U1v1,V1w1,W1u1)
2180          CALL
          Data_print(Paxis,Mod(Paxis),Nsam,"TUN",U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,A1b1,U1a1,V1
          a1,W1a1)
2190          CALL Data_trnsfrm(Tun2mod(*),U,V,W,U1,V1,W1,U1v1,V1w1,W1u1)
2200          CALL
          Data_print(Paxis,Mod(Paxis),Nsam,"MOD",U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,A1b1,U1a1,V1
          a1,W1a1)
2210          CALL Data_plot(Array(*),Symbols(*),6,Mod(Paxis),U,V,W,1/Uinf,N(1,1),N(2,1),N(3,1))
2220          CALL Data_plot(Array(*),Symbols(*),7,Mod(Paxis),U1,V1,W1,1/Uinf,N(1,2),N(2,2),N(3,2))
2230          CALL Data_plot(Array(*),Symbols(*),8,Mod(Paxis),U1v1,V1w1,W1u1,1/Uinf^2,N(1,3),N(2,3),N(3,3))
2240          CALL Data_plot(Array(*),Symbols(*),9,Mod(Paxis),Ttemp,Uinf,Uedge,1,N(4,1),1,1)
2250          OUTPUT PRT;RPTS("=",140)
2260          GOSUB Store_file
2270          File=File+1
2280      END IF
2290      RETURN
2300 M1k7:  Quit=0
2310      ON KBD GOSUB Quit
2320      FOR J=1 TO Npos
2330          Mod(Paxis)=Pos(J,1)
2340          GOSUB M1k5a
2350          GOSUB M1k6
2360          IF Quit THEN 2380

```

```

2370      NEXT J
2380      OFF KBD
2390      GOSUB On_key
2400      CALL Menu_disp(Menu,Menu$(*))
2410      RETURN
2420 M1k8:  DISP "Press any key to return to main menu"
2430      CALL Rt_histo(@Lvdas,Symbols{*},1)
2440      RETURN
2450 M2k1:  Menu=1
2460      CALL Menu_disp(Menu,Menu$(*))
2470      RETURN
2480 M2k2:  SELECT TRIM$(Menu$(Menu,Key){20})
2490      CASE "Tx & Rx"
2500          Menu$(Menu,Key){20}="Tx"
2510      CASE "Tx"
2520          Menu$(Menu,Key){20}="Rx"
2530      CASE "Rx"
2540          Menu$(Menu,Key){20}="Tx & Rx"
2550      END SELECT
2560      CALL Menu_disp(Menu,Menu$(*))
2570      RETURN
2580 M2k3:  SELECT TRIM$(Menu$(Menu,Key){20})
2590      CASE "MODEL"
2600          Menu$(Menu,Key){20}="TUNNEL"
2610      CASE "TUNNEL"
2620          Menu$(Menu,Key){20}="LASER"
2630      CASE "LASER"
2640          Menu$(Menu,Key){20}="MODEL"
2650      END SELECT
2660      CALL Menu_disp(Menu,Menu$(*))
2670      RETURN
2680 M2k4:  SELECT TRIM$(Menu$(Menu,Key){20})
2690      CASE "ABSOLUTE"
2700          Menu$(Menu,Key){20}="RELATIVE"
2710      CASE "RELATIVE"
2720          Menu$(Menu,Key){20}="ABSOLUTE"
2730      END SELECT
2740      CALL Menu_disp(Menu,Menu$(*))
2750      RETURN
2760 M2k5:  !
2770 M2k6:  !
2780 M2k7:  !
2790 M2k8:  Side$=TRIM$(Menu$(Menu,2){20})
2800      Coord$=TRIM$(Menu$(Menu,3){20})
2810      Mode$=TRIM$(Menu$(Menu,4){20})
2820      CALL Enter_value(Mode$&" Movement",Movement,"4D.5D")
2830      ON KBD CALL Do_nothing
2840      DISP "Moving"
2850      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
2860      CALL Tcs8move(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Mod2tun(*),Tun2tcs1(*),Tun2tcs2(*),Side$,Coord$,Mode$,Key-
        4,Movement)
2870      CALL Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
2880      DISP ""
2890      OFF KBD
2900      RETURN
2910 M3k1:  Menu=1
2920      CALL Menu_disp(Menu,Menu$(*))
2930      RETURN
2940 M3k2:  CALL Enter_value("Run",Run,"3D.2D")
2950      CALL Enter_value("File",File,"3D")
2960      RETURN
2970 M3k3:  CALL Enter_value("Number of Samples ",Nreads,"K")
2980      RETURN
2990 M3k4:  CALL Enter_string("Traverse Axis for Profile ",Paxis$,"K")
3000      SELECT Paxis$
3010      CASE "X"
3020          Paxis=1
3030      CASE "Y"
3040          Paxis=2
3050      CASE "Z"
3060          Paxis=3
3070      CASE "A"
3080          Paxis=4
3090      CASE ELSE
3100          GOTO M3k4
3110      END SELECT
3120      RETURN
3130 M3k5:  GOSUB Read_calc_fill
3140 M3k5a: !OUTPUT PRT USING "#,0,2/"
3150      OUTPUT PRT USING "#,2/"

```

```

3160 OUTPUT PRT USING "20X,K,/,";"TRAVERSE COORDINATE TRANSFORMATION MATRICIES"
3170 OUTPUT PRT USING "20X,K,/,4(13X,4(8D.5D),/),";"Transmitting side TCS8 to TUNNEL",Tcs2tun1(*)
3180 OUTPUT PRT USING "20X,K,/,4(13X,4(8D.5D),/),";"Receiving side TCS8 to TUNNEL",Tcs2tun2(*)
3190 OUTPUT PRT USING "20X,K,/,4(13X,4(8D.5D),/),";"Transmitting side TUNNEL to TCS8",Tun2tcs1(*)
3200 OUTPUT PRT USING "20X,K,/,4(13X,4(8D.5D),/),";"Receiving side TUNNEL to TCS8",Tun2tcs2(*)
3210 OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/),";"TUNNEL to MODEL",Tun2mod(*)
3220 OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/),";"MODEL to TUNNEL",Mod2tun(*)
3230 OUTPUT PRT USING "20X,K,/,";"VELOCITY COORDINATE TRANSFORMATION MATRICIES"
3240 OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/),";"LASER to TUNNEL",Ldv2tun(*)
3250 OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/),";"TUNNEL to LASER",Tun2ldv(*)
3260 OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/),";"TUNNEL to MODEL",Tun2mod(*)
3270 OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/),";"MODEL to TUNNEL",Mod2tun(*)
3280 OUTPUT PRT USING "#,@"
3290 RETURN
3300 M3k6: CALL Setup_graph(Array(*),Image$(*),Paxis,Symbols(*))
3310 RETURN
3320 M3k7: Menu=4
3330 CALL Menu_disp(Menu,Menu$(*))
3340 RETURN
3350 M3k8: Menu=5
3360 CALL Menu_disp(Menu,Menu$(*))
3370 RETURN
3380 M4k1: Menu=3
3390 CALL Menu_disp(Menu,Menu$(*))
3400 RETURN
3410 M4k2: GOSUB Read_array
3420 GOSUB Read_calc_fill
3430 RETURN
3440 M4k3: GOSUB Read_calc_fill
3450 GOSUB Save_array
3460 RETURN
3470 M4k4: GOSUB Read_calc_fill
3480 GOSUB Print_header
3490 RETURN
3500 M4k5: GSTORE Gsave(*)
3510 GOSUB Read_calc_fill
3520 CALL Change("VALUES",Array(*),Name$(*),Image$(*),Units$(*))
3530 GOSUB Read_calc_fill
3540 GLOAD Gsave(*)
3550 RETURN
3560 M4k6: GSTORE Gsave(*)
3570 GOSUB Read_calc_fill
3580 CALL Change("NAMES",Array(*),Name$(*),Image$(*),Units$(*))
3590 GOSUB Read_calc_fill
3600 GLOAD Gsave(*)
3610 RETURN
3620 M4k7: GSTORE Gsave(*)
3630 GOSUB Read_calc_fill
3640 CALL Change("UNITS",Array(*),Name$(*),Image$(*),Units$(*))
3650 GOSUB Read_calc_fill
3660 GLOAD Gsave(*)
3670 RETURN
3680 M4k8: GSTORE Gsave(*)
3690 GOSUB Read_calc_fill
3700 CALL Change("IMAGES",Array(*),Name$(*),Image$(*),Units$(*))
3710 GOSUB Read_calc_fill
3720 GLOAD Gsave(*)
3730 RETURN
3740 M5k1: Menu=3
3750 CALL Menu_disp(Menu,Menu$(*))
3760 RETURN
3770 M5k2: CALL Tcs8set("P",@Tcs8) ! View and set TCS8 Positions
3780 GRAPHICS ON
3790 CALL Menu_disp(Menu,Menu$(*))
3800 RETURN
3810 M5k3: CALL Tcs8set("U",@Tcs8) ! View and set TCS8 counts per Unit length
3820 GRAPHICS ON
3830 CALL Menu_disp(Menu,Menu$(*))
3840 RETURN
3850 M5k4: CALL Tcs8set("R",@Tcs8) ! View and set TCS8 counts per Revolution
3860 GRAPHICS ON
3870 CALL Menu_disp(Menu,Menu$(*))
3880 RETURN
3890 M5k5: CALL Tcs8set("V",@Tcs8) ! View and set TCS8 Velocities
3900 GRAPHICS ON
3910 CALL Menu_disp(Menu,Menu$(*))
3920 RETURN
3930 M5k6: CALL Tcs8set("A",@Tcs8) ! View and set TCS8 Accelerations
3940 GRAPHICS ON
3950 CALL Menu_disp(Menu,Menu$(*))

```

```

3960 RETURN
3970 M5k7: CALL Enter_value("Run",Run,"3D.2D")
3980 CALL Enter_value("File",File,"3D")
3990 FOR Run=Run TO Run
4000 CLEAR SCREEN
4010 FOR File=1 TO 100
4020 GOSUB Read_file
4030 IF File$="" THEN 4170
4040 CALL Data_reduce(At_exp,Ct_exp,Nsam)
4050 Vwprt(1,1)=50
4060 Vwprt(1,2)=225
4070 Vwprt(2,1)=275
4080 Vwprt(2,2)=450
4090 Vwprt(4,1)=500
4100 Vwprt(4,2)=675
4110 FOR G=1 TO 5
4120 Vwprt(G,3)=1025-65*File
4130 Vwprt(G,4)=1065-65*File
4140 NEXT G
4150 CALL Pt_histo(Symbols(*),Run,File,Mod(Paxis),Nsam)
4160 NEXT File
4170 DISP ""
4180 PRINTER IS PRT
4190 PRINT USING "#,@"
4200 DUMP GRAPHICS
4210 PRINT USING "#,@"
4220 PRINTER IS CRT
4230 NEXT Run
4240 CLEAR SCREEN
4250 CALL Setup_graph(Array(*),Image$(*),Paxis,Symbols(*))
4260 RETURN
4270 M5k8: CALL Enter_value("Run",Run,"3D.2D")
4280 CALL Enter_value("File",File,"3D")
4290 GOSUB Read_file
4300 IF File$="" THEN RETURN
4310 GOSUB Print_header
4320 CALL Setup_graph(Array(*),Image$(*),Paxis,Symbols(*))
4330 !BEEP
4340 !DISP "SWITCH SWITCH TO B AND THEN PRESS <CONTINUE>"
4350 !PAUSE
4360 FOR File=1 TO 100
4370 GOSUB Read_file
4380 IF File$="" THEN 4630
4390 Cmask=Coin(1)*1+Coin(2)*2+Coin(3)*4
4400 OUTPUT PRT;RPTS("=",140)
4410 CALL Data_reduce(At_exp,Ct_exp,Nsam)
4420 !CALL Data_xfer(@Tcs8,Run,File,Ui(*),Vi(*),Ai(*),Valid(*),Nsam)
4430 CALL Pt_histo(Symbols(*),Run,File,Mod(Paxis),Nsam)
4440 CALL Data_clip(Nsam,Umin,Umax,Vmin,Vmax)
4450 CALL Data_sum(Sum(*),N(*),Nsam)
4460 CALL Data_calc(N(*),Sum(*),U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,Alb1,U1a1,V1a1,W1a1)
4470 CALL
Data_print(Paxis,Mod(Paxis),Nsam,"MHZ",U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,Alb1,U1a1,V1a1,W1a1)
4480 CALL Data_fconvert(Array(*))
4490 CALL Data_aconvert(Gain)
4500 CALL Data_sum(Sum(*),N(*),Nsam)
4510 CALL Data_calc(N(*),Sum(*),U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,Alb1,U1a1,V1a1,W1a1)
4520 CALL
Data_print(Paxis,Mod(Paxis),Nsam,"LDV",U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,Alb1,U1a1,V1a1,W1a1)
4530 CALL Data_trnsfrm(Ldv2tun(*),U,V,W,U1,V1,W1,U1v1,V1w1,W1u1)
4540 CALL
Data_print(Paxis,Mod(Paxis),Nsam,"TUN",U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,Alb1,U1a1,V1a1,W1a1)
4550 CALL Data_trnsfrm(Tun2mod(*),U,V,W,U1,V1,W1,U1v1,V1w1,W1u1)
4560 CALL
Data_print(Paxis,Mod(Paxis),Nsam,"MOD",U,V,W,A,B,IO,CO,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,Alb1,U1a1,V1a1,W1a1)
4570 CALL Data_plot(Array(*),Symbols(*),6,Mod(Paxis),U,V,W,1/Uinf,N(1,1),N(2,1),N(3,1))
4580 CALL Data_plot(Array(*),Symbols(*),7,Mod(Paxis),U1,V1,W1,1/Uinf,N(1,2),N(2,2),N(3,2))
4590 CALL Data_plot(Array(*),Symbols(*),8,Mod(Paxis),U1v1,V1w1,W1u1,1/Uinf^2,N(1,3),N(2,3),N(3,3))
4600 CALL Data_plot(Array(*),Symbols(*),9,Mod(Paxis),Ttemp,Uinf,Uedge,1,N(4,1),1,1)
4610 OUTPUT PRT;RPTS("=",140)
4620 NEXT File
4630 GOSUB M1k3
4640 File=File-1
4650 GOSUB Read_file
4660 GOSUB Print_header
4670 GOSUB M3k5a

```

```

4680          PRINTER IS CRT
4690          !BEEP
4700          !DISP "SWITCH SWITCH TO A AND THEN PRESS <CONTINUE>"
4710          !PAUSE
4720          RETURN
4730 Quit:    Quit=1
4740          RETURN
4750 Print_header: PRINTER IS PRT;WIDTH 144
4760          PRINT USING "#,0,5(K)";CHR$(27) & "&k2S" & CHR$(27) & "&l9D"
4770          CALL Array_print(Array(*),Name$(*),Image$(*),Units$(*))
4780          PRINT USING "#,0,5(K)";CHR$(27) & "E"
4790          PRINTER IS CRT
4800          RETURN
4810 Read_calc_fill: GOSUB Read
4820          GOSUB Calc
4830          GOSUB Fill
4840          FOR X=1 TO SIZE(Array,2)
4850              FOR Y=1 TO SIZE(Array,1)
4860                  Array(Y,X)=PROUND(Array(Y,X),-15)
4870              NEXT Y
4880          NEXT X
4890          RETURN
4900 Copy_file:  FOR Run=1.01 TO 6.01
4910              IF Run=1.01 THEN F2=9
4920              IF Run=2.01 THEN F2=6
4930              IF Run=3.01 THEN F2=5
4940              IF Run=4.01 THEN F2=11
4950              IF Run=5.01 THEN F2=0
4960              IF Run=6.01 THEN F2=5
4970              FOR File=1 TO F2
4980                  Data$=":,1400,0,1"
4990                  GOSUB Read_file
5000                  Data$=":,1400,0,0"
5010                  GOSUB Store_file
5020                  Data$=":,1400,0,1"
5030              NEXT File
5040          NEXT Run
5050          RETURN
5060 Store_header: DISP "Storing Header"
5070          File$="R"&VAL$(Run)&Data$
5080          ON ERROR GOTO 5280
5090          ASSIGN @Data TO File$
5100          OFF ERROR
5110          FOR K=1 TO 10
5120              WAIT .2
5130              BEEP
5140          NEXT K
5150          CALL Enter_string("Over Write old file ",LS,"K")
5160          SELECT LS[1,1]
5170          CASE "Y","y"
5180              ASSIGN @Data TO *
5190              PURGE File$
5200              GOTO 5280
5210          CASE "N","n"
5220              CALL Enter_value("Run",Run,"3D.2D")
5230              CALL Enter_value("File",File,"3D")
5240              GOTO 5060
5250          CASE ELSE
5260              GOTO 5060
5270          END SELECT
5280          OFF ERROR
5290          Fsize=INT((3200+4000*3+128*4+72*4)/256*1.05+1)
5300          CREATE BDAT File$,Fsize
5310          ASSIGN @Data TO File$
5320          OUTPUT @Data;Array(*),Name$(*),Image$(*),Units$(*)
5330          OUTPUT @Data;Tun2tcs1(*),Tun2tcs2(*),Mod2tun(*),Tun2ldv(*)
5340          OUTPUT @Data;Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*),Ldv2tun(*)
5350          ASSIGN @Data TO *
5360          PROTECT "R"&VAL$(Run)&Data$,"TKM"
5370          RETURN
5380 Store_file:  GOSUB Calc
5390          GOSUB Fill
5400          IF File=1 THEN GOSUB Store_header
5410          DISP "Storing Data"
5420          File$="R"&VAL$(Run)&"F"&VAL$(File)&Data$
5430          ON ERROR GOTO 5630
5440          ASSIGN @Data TO File$
5450          OFF ERROR
5460          FOR K=1 TO 10
5470              WAIT .2

```



```

5480         BEEP
5490     NEXT K
5500     CALL Enter_string("Over Write old file ",L$, "K")
5510     SELECT L$(1,1)
5520     CASE "Y","y"
5530         ASSIGN @Data TO *
5540         PURGE File$
5550         GOTO 5630
5560     CASE "N","n"
5570         CALL Enter_value("Run",Run,"3D.2D")
5580         CALL Enter_value("File",File,"3D")
5590         GOTO 5380
5600     CASE ELSE
5610         GOTO 5380
5620     END SELECT
5630     OFF ERROR
5640     Fsize=INT((3200+Nsam*10*2+60+240)/256*1.05+1)
5650     CREATE BDAT File$,Fsize
5660     ASSIGN @Data TO File$
5670     OUTPUT @Data;Array(*),Raw(*),N(*),Sum(*)
5680     ASSIGN @Data TO *
5690     PROTECT "R"&VAL$(Run)&"F"&VAL$(File)&Data$, "TKM"
5700     RETURN
5710 Read_header:
5720     DISP "Reading Header"
5730     File$="R"&VAL$(Run)&"<TKM>"&Data$
5740     ON ERROR GOTO 5820
5750     ASSIGN @Data TO File$
5760     ENTER @Data;Array(*),Name$(*),Image$(*),Units$(*)
5770     CALL Fix(Array(*),Name$(*),Image$(*),Units$(*))
5780     ENTER @Data;Tun2tcs1(*),Tun2tcs2(*),Mod2tun(*),Tun2ldv(*)
5790     ENTER @Data;Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*),Ldv2tun(*)
5800     ASSIGN @Data TO *
5810     OFF ERROR
5820     RETURN
5830     File$=""
5840     RETURN
5850 Read_file:
5860     IF File=1 THEN GOSUB Read_header
5870     DISP "Reading Data"
5880     File$="R"&VAL$(Run)&"F"&VAL$(File)&"<TKM>"&Data$
5890     ON ERROR GOTO 6020
5900     ASSIGN @Data TO File$
5910     ENTER @Data;Array(*)
5920     CALL Fix(Array(*),Name$(*),Image$(*),Units$(*))
5930     GOSUB Read
5940     REDIM Raw(1:Nsam,1:10)
5950     ENTER @Data;Raw(*),N(*),Sum(*)
5960     ASSIGN @Data TO *
5970     OFF ERROR
5980     Date=Array(1,1)
5990     Time=Array(2,1)
6000     Run=Array(3,1)
6010     File=Array(4,1)
6020     RETURN
6030     OFF ERROR
6040     File$=""
6050     RETURN
6060 Fill:
6070     Array(1,1)=Date
6080     Array(1,2)=Mach
6090     Array(1,3)=Stemp
6100     Array(1,4)=Alpha1
6110     Array(2,1)=Time
6120     Array(2,2)=Temp
6130     Array(2,3)=Ttemp
6140     Array(2,4)=Alpha2
6150     Array(3,1)=Run
6160     Array(3,2)=Uedge
6170     Array(3,3)=Tt_mv
6180     Array(3,4)=Alpha3
6190     Array(4,1)=File
6200     Array(4,2)=Uinf
6210     Array(4,3)=Tt_raw
6220     Array(4,4)=Theta
6230     MAT Array(11:14,1)= Mod
6240     MAT Array(11:14,2)= Tun
6250     MAT Array(11:14,3)= Tcs1
6260     MAT Array(11:14,4)= Tcs2
6270     MAT Array(21,1:3)= Beam_spc
6280     MAT Array(22,1:3)= Foc1_len
6290     MAT Array(23,1:3)= Beam_sep

```

```

6280 MAT Array(24,1:3)= Wave_len      ! Wave length
6290 MAT Array(25,1:3)= Frng_spc      ! Fringe spacing
6300 MAT Array(26,1:3)= Brg_frq       ! Bragg frequency
6310 MAT Array(27,1:3)= Mix_frq       ! Mixing frequency
6320 MAT Array(28,1:3)= Mea_sgn       ! Sign of measured frequency in velocity equation
6330 MAT Array(29,1:3)= Brg_sgn       ! Sign of bragg frequency in velocity equation
6340 MAT Array(30,1:3)= Mix_sgn       ! Sign of mixing frequency in velocity equation
6350 MAT Array(31,1:3)= Coin          ! Coincidence criteria
6360 MAT Array(32:34,1:3)= Thata      ! Angles between measured (ABC) & tunnel (UVW) coordinate systems
6370 MAT Array(21:23,4)= Index        ! Index of refraction of for laser light (eg: Nair,Nglass,Nwater)
6380 Array(24,4)=Nreads              ! Number of desired samples
6390 Array(25,4)=Nsam                 ! Number of acquired samples
6400 Array(26,4)=Atime                ! Acquisition time
6410 Array(27,4)=Ctime                ! Coincidence time
6420 Array(28,4)=At_exp                ! Acquisition time exponent
6430 Array(29,4)=Ct_exp                ! coincidence time exponent
6440 Array(30,4)=Gain                 ! Tunnel Total Temperature Voltage Gain
6450 Array(31,4)=Paxis                ! Axis for plots
6460 Array(35,1)=Umin                 ! Frequency minimum for U calculation
6470 Array(35,2)=Vmin                 ! Frequency minimum for V calculation
6480 Array(35,3)=Wmin                 ! Frequency minimum for W calculation
6490 Array(36,1)=Umax                 ! Frequency maximum for U calculation
6500 Array(36,2)=Vmax                 ! Frequency maximum for V calculation
6510 Array(36,3)=Wmax                 ! Frequency maximum for W calculation
6520 Array(36,4)=Clip                 ! Clip
6530 RETURN
6540 Read: !Date=Array(1,1)           ! Date
6550 Date=TIMEDATE                    ! Date
6560 Mach=Array(1,2)                  ! Mach Number
6570 Stemp=Array(1,3)                 ! Stagnation Temperature (°R)
6580 Alpha=Array(1,4)                 ! Angle of Attack
6590 !Time=Array(2,1)                  ! Time
6600 Time=Date                        ! Time
6610 Temp=Array(2,2)                  ! Room Temperature (°F)
6620 Ttemp=Array(2,3)                 ! Total Temperautue (°R)
6630 Alpha2=Array(2,4)                ! Cone angle
6640 !Run=Array(3,1)                   ! Run Number
6650 Uedge=Array(3,2)                 ! Uedge
6660 Tt_mv=Array(3,3)                 ! Total Temperautue (mv)
6670 Alpha3=Array(3,4)                ! Roll angle
6680 !File=Array(4,1)                  ! File Number
6690 Uinf=Array(4,2)                  ! Freestrem Velocity
6700 Tt_raw=Array(4,3)                ! Total Temperautue (raw voltage w/gain)
6710 Theta=Array(4,4)                 ! Tx Side OffAxis Angle
6720 MAT Mod= Array(11:14,1)           ! Probe volume position in Model coordinates
6730 MAT Tun= Array(11:14,2)           ! Probe volume position in Tunnel coordinates
6740 MAT Tcs1= Array(11:14,3)          ! Tx side traverse position in Tcs8 coordinates
6750 MAT Tcs2= Array(11:14,4)          ! Rx side traverse position in Tcs8 coordinates
6760 MAT Beam_spc= Array(21,1:3)       ! Beam spacing at lens
6770 MAT Focl_len= Array(22,1:3)       ! Focal length
6780 MAT Beam_sep= Array(23,1:3)       ! Beam separation agnle in degrees (full angle)
6790 MAT Wave_len= Array(24,1:3)       ! Wave length
6800 MAT Frng_spc= Array(25,1:3)       ! Fringe spacing
6810 MAT Brg_frq= Array(26,1:3)        ! Bragg frequency
6820 MAT Mix_frq= Array(27,1:3)        ! Mixing frequency
6830 MAT Mea_sgn= Array(28,1:3)        ! Sign of measured frequency in velocity equation
6840 MAT Brg_sgn= Array(29,1:3)        ! Sign of bragg frequency in velocity equation
6850 MAT Mix_sgn= Array(30,1:3)        ! Sign of mixing frequency in velocity equation
6860 MAT Coin= Array(31,1:3)           ! Coincidence criteria
6870 MAT Thata= Array(32:34,1:3)        ! Angles between measured (ABC) & tunnel (UVW) coordinate systems
6880 MAT Index= Array(21:23,4)         ! Index of refraction of for laser light (eg: Nair,Nglass,Nwater)
6890 Nreads=Array(24,4)                ! Number of desired samples
6900 Nsam=Array(25,4)                  ! Number of acquired samples
6910 Atime=Array(26,4)                 ! Acquisition time
6920 Ctime=Array(27,4)                 ! Coincidence time
6930 At_exp=Array(28,4)                ! Acquisition time exponent
6940 Ct_exp=Array(29,4)                ! coincidence time exponent
6950 Gain=Array(30,4)                  ! Tunnel Total Temperature Voltage Gain
6960 Paxis=Array(31,4)                 ! Axis for plots
6970 Umin=Array(35,1)                  ! Frequency minimum for U calculation
6980 Vmin=Array(35,2)                  ! Frequency minimum for V calculation
6990 Wmin=Array(35,3)                  ! Frequency minimum for W calculation
7000 Umax=Array(36,1)                  ! Frequency maximum for U calculation
7010 Vmax=Array(36,2)                  ! Frequency maximum for V calculation
7020 Wmax=Array(36,3)                  ! Frequency maximum for W calculation
7030 Clip=Array(36,4)                 ! Clip
7040 RETURN
7050 Calc: FOR K=1 TO 3
7060 IF I=2 THEN
7070 Beam1=Theta+ATN(Beam_spc(K)/2/Focl_len(K))

```

```

7080      Beam2=Theta-ATN(Beam_spc(K)/2/Focl_len(K))
7090      ELSE
7100          Beam1=0+ATN(Beam_spc(K)/2/Focl_len(K))
7110          Beam2=0-ATN(Beam_spc(K)/2/Focl_len(K))
7120      END IF
7130      Beam1=ASN(Index(1)/Index(3)*SIN(Beam1))
7140      Beam2=ASN(Index(1)/Index(3)*SIN(Beam2))
7150      Beam_sep(K)=Beam1-Beam2
7160      Frng_spc(K)=Wave_len(K)/(2*SIN(Beam_sep(K)/2))/1000
7170  NEXT K
7180  MAT Array(23,1:3)= Beam_sep      ! Beam separation agnie in degrees (full angle)
7190  MAT Array(25,1:3)= Frng_spc      ! Fringe spacing
7200  CALL Ctm_tcs1(Tcs2tun1(*),Tun2tcs1(*))
7210  CALL Ctm_tcs2(Tcs2tun2(*),Tun2tcs2(*))
7220  CALL Ctm_ldv(Index(*),Thata(*),Tun2ldv(*),Ldv2tun(*))
7230  CALL Ctm_mod(Alpha1,Alpha2,Alpha3,Mod2tun(*),Tun2mod(*))
7240  CALL Lvdsas_sample_c(@Lvdsas,4,Table(*),Vave,Vsdv,Tave,Tsdv)
7250  Tt_raw=Vave
7260  Tt_mv=Tt_raw/Gain*1000
7270  CALL Temp(Mach,Tt_mv,Stemp,Ttemp)
7280  Uinf=Mach*49.0*SQR(Stemp)*.3048
7290  !Uinf=20.043*Mach*SQR((273+5/9*(Temp-32))/(1+.2*Mach^2))
7300  Uedge=Uinf
7310  Cmask=Coin(1)*1+Coin(2)*2+Coin(3)*4
7320  SELECT Paxis
7330  CASE 1
7340      Paxis$="X"
7350  CASE 2
7360      Paxis$="Y"
7370  CASE 3
7380      Paxis$="Z"
7390  CASE 4
7400      Paxis$="A"
7410  CASE ELSE
7420      Paxis=2
7430      Paxis$="Y"
7440      GOTO M3k4
7450  END SELECT
7460  IF Run=0 OR File=0 THEN
7470      CALL Enter_value("Run Number ",Run,"3D.2D")
7480      CALL Enter_value("File Number ",File,"3D")
7490      GOTO 7460
7500  END IF
7510  RETURN
7520 Read_array:
7530  ASSIGN @File TO "ARRAY"&System$
7540  ENTER @File;Array(*),Name$(*),Image$(*),Units$(*)
7550  ENTER @File;Tun2tcs1(*),Tun2tcs2(*),Mod2tun(*),Tun2ldv(*)
7560  ENTER @File;Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*),Ldv2tun(*)
7570  ASSIGN @File TO *
7580  OFF ERROR
7590  RETURN
7600  OFF ERROR
7610  ASSIGN @File TO *
7620  ON ERROR GOTO 7640
7630  PURGE "ARRAY"&System$
7640  OFF ERROR
7650  CALL Array_init(Name$(*),Array(*),Image$(*),Units$(*))
7660  CREATE BDAT "ARRAY"&System$,50
7670  GOSUB Save_array
7680  RETURN
7690 Save_array:
7700  ASSIGN @File TO "ARRAY"&System$
7710  OUTPUT @File;Array(*),Name$(*),Image$(*),Units$(*)
7720  OUTPUT @File;Tun2tcs1(*),Tun2tcs2(*),Mod2tun(*),Tun2ldv(*)
7730  OUTPUT @File;Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*),Ldv2tun(*)
7740  ASSIGN @File TO *
7750  RETURN
7760 Do_nothing:
7770  KS=KBDS
7780  SUBEND
7790 Menu:
7800 Menu_read:
7810  OPTION BASE 1
7820  DIM L$(80)
7830  FOR Menu=1 TO SIZE(Menu$,1)
7840      FOR Key=1 TO 8
7850          Menu$(Menu,Key)="M"&VAL$(Menu)&"K"&VAL$(Key)&" ":"
7860      NEXT Key
7870  NEXT Menu

```

```

7880      ON ERROR GOTO 7950
7890      WHILE 1=1
7900          READ LS
7910          Menu=VAL(LS(2,2))
7920          Key=VAL(LS(4,4))
7930          Menu$(Menu,Key)=LS
7940      END WHILE
7950      SUBEXIT
7960      DATA "M1K1: Menu2: Laser Alignment"
7970      DATA      "M2K1: Return to main menu"
7980      DATA      "M2K2: Sides      : Tx & Rx"
7990      DATA      "M2K3: Coordinates: MODEL"
8000      DATA      "M2K4: Mode      : ABSOLUTE"
8010      DATA      "M2K5: Move X"
8020      DATA      "M2K6: Move Y"
8030      DATA      "M2K7: Move Z"
8040      DATA      "M2K8: Move A"
8050      DATA "M1K2: Menu3: Pre Run"
8060      DATA      "M3K1: Return to MAIN menu"
8070      DATA      "M3K2: Enter Run & File Numbers"
8080      DATA      "M3K3: Enter Number of Samples"
8090      DATA      "M3K4: Select Traverse Axis for Profile"
8100      DATA      "M3K5: Print Coordinate Transformation Matrices"
8110      DATA      "M3K6: Setup Graphics"
8120      DATA      "M3K7: Menu4: Tunnel Conditions"
8130      DATA          "M4K1: Return to PRE RUN menu"
8140      DATA          "M4K2: Load Tunnel Conditions"
8150      DATA          "M4K3: Save Tunnel Conditions"
8160      DATA          "M4K4: Print Tunnel Conditions"
8170      DATA          "M4K5: Enter Tunnel Condition Data"
8180      DATA          "M4K6: Enter Tunnel Condition Names"
8190      DATA          "M4K7: Enter Tunnel Condition Units"
8200      DATA          "M4K8: Enter Tunnel Condition Images"
8210      DATA      "M3K8: Menu5: Traverse"
8220      DATA          "M5K1: Return to TRAVERSE menu"
8230      DATA          "M5K2: View & Set TCS8 Positions"
8240      DATA          "M5K3: View & Set TCS8 Units"
8250      DATA          "M5K4: View & Set TCS8 Revolution"
8260      DATA          "M5K5: View & Set TCS8 Velocity"
8270      DATA          "M5K6: View & Set TCS8 Acceleration"
8280      DATA "M1K3: Post Run (Dump Graphics)"
8290      DATA "M1K4: Set Auto Move Positions"
8300      DATA "M1K5: Move traverse"
8310      DATA "M1K6: Take data"
8320      DATA "M1K7: Auto move and take"
8330      DATA "M1K8: Display Histograms"
8340      SUBEND
8350 Menu_disp:  SUB Menu_disp(Menu,Menu$(*))
8360                  PRINTER IS CRT
8370                  PRINT CHR$(128);
8380                  IF Menu=0 THEN Menu=1
8390                  FOR Key=1 TO 8
8400                      Menu$(Menu,Key)=Menu$(Menu,Key)&RPTS(" ",50-LEN(Menu$(Menu,Key)))
8410                      PRINT TABXY(1,Key);Menu$(Menu,Key) [3]
8420                  NEXT Key
8430                  PRINT CHR$(128);
8440      SUBEND
8450 Menu_status: SUB Menu_status(Menu,Key,Pen,Menu$(*))
8460                  PRINTER IS CRT
8470                  PRINT TABXY(1,Key);CHR$(129-Pen);Menu$(Menu,Key) [3];CHR$(128)
8480                  WAIT .1
8490      SUBEND
8500 Enter:      !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
8510 Enter_value: SUB Enter_value(Name$,Value,Image$)
8520                  IF Name$="Date" OR Name$="Time" THEN SUBEXIT
8530                  DISP CHR$(129);
8540                  DISP USING 8550;Name$
8550                  IMAGE #,"Old ",K,"="
8560                  IF Image$<>" " THEN DISP USING "#, "&Image$;Value
8570                  IF Image$="" THEN DISP USING "#, K";Value
8580                  DISP USING 8590;Name$
8590                  IMAGE #,"      Enter new ",K
8600                  INPUT " ? ",Value
8610                  DISP CHR$(128);
8620      SUBEND
8630 Enter_string: SUB Enter_string(Name$,Value$,Image$)
8640                  DISP CHR$(129);
8650                  DISP USING 8660;Name$
8660                  IMAGE #,"Old ",K,"="
8670                  DISP USING "#, "&Image$;Value$

```

```

8680      DISP USING 8690;Name$
8690      IMAGE #,"      Enter new ",K
8700      INPUT " ? ",Value$
8710      DISP CHR$(128);
8720      SUBEND
8730      Array:
8740      Array_init:
8750      ON ERROR GOTO 8930
8760      READ Y
8770      FOR X=1 TO SIZE(Name$,2)
8780          READ Name$(Y,X),Array(Y,X),Image$(Y,X),Units$(Y,X)
8790          SELECT Image$(Y,X)
8800              CASE "0"
8810                  Image$(Y,X)="9D"
8820              CASE "1" TO "7"
8830                  After=VAL(Image$(Y,X))
8840                  Before=8-After
8850                  Image$(Y,X)=VAL$(Before)&"D."&VAL$(After)&"D"
8860              CASE "K"
8870              CASE "N"
8880              CASE ELSE
8890                  Image$(Y,X)="9D"
8900          END SELECT
8910      NEXT X
8920      GOTO 8760
8930      SUBEXIT
8940      ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
8950      DATA 1, Date , 0,0,"", Mach , 7.0,4,"", STemp , 0,0,"R", Alpha1 , 0,4,"
8960      DATA 2, Time , 0,0,"", Temp , 68.5,4,"F", TTemp , 0,0,"R", Alpha2 , 0,4,"
8970      DATA 3, Run , 5,2,"", Uedge , 1,4,m/s, Tt , 0,3,mv, Alpha3 , 0,4,"
8980      DATA 4, File , 0,0,"", Uinf , 1,4,m/s, Tt (raw), 0,3,v, Theta , 0,4,"
8990      ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
9000      DATA 11, Xmod , 0,4,in, Xtun , 0,4,in, Xltcs , 0,4,in, X2tcs , 0,4,in
9010      DATA 12, Ymod , 0,4,in, Ytun , 0,4,in, Yltcs , 0,4,in, Y2tcs , 0,4,in
9020      DATA 13, Zmod , 0,4,in, Ztun , 0,4,in, Zltcs , 0,4,in, Z2tcs , 0,4,in
9030      DATA 14, Amod , 0,4,in, Atun , 0,4,in, Altcs , 0,4,in, A2tcs , 0,4,in
9040      ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
9050      DATA 21, UBeamSpc,.3125,3,in, VBeamSpc,.3438,3,in, WBeamSpc,.3125,3,in, Index1 ,1.000,3,"
9060      DATA 22, UFocLen,30.00,3,in, VFocLen,30.00,3,in, WFocLen,30.00,3,in, Index2 ,1.000,3,"
9070      DATA 23, UBeamSep,0.000,3,"", VBeamSep,0.000,3,"", WBeamSep,0.000,3,"", Index3 ,1.000,3,"
9080      DATA 24, UWaveLen,514.5,3,nm, VWaveLen,488.0,3,nm, WWaveLen,476.5,3,nm, Nreads ,1000,0,"
9090      DATA 25, UFrngSpc,00.00,3,um, VFrngSpc,00.00,3,um, WFrngSpc,00.00,3,um, Nsam ,1000,0,"
9100      DATA 26, Ubrag ,40.00,4,MHz, Vbrag ,40.00,4,MHz, Wbrag ,40.00,4,MHz, Atime ,5,6,s
9110      DATA 27, Umix , 0,0,4,MHz, Vmix , 0,0,4,MHz, WMix , 0,0,4,MHz, Ctime ,1E-2,6,s
9120      DATA 28, UmeaSgn , -1,0,"", VmeaSgn , +1,0,"", WmeaSgn , +1,0,"", ATexp ,12,0,"
9130      DATA 29, UbrgSgn , +1,0,"", VbrgSgn , -1,0,"", WbrgSgn , -1,0,"", CTexp ,7,0,"
9140      DATA 30, UmixSgn , -1,0,"", VmixSgn , +1,0,"", WmixSgn , +1,0,"", Tt Gain ,100,0,"
9150      DATA 31, U coin , 1,0,"", V coin , 1,0,"", W coin , 0,0,"", Paxis ,2,0,"
9160      DATA 32, ThetaAU , 0,4,"", ThetaAV , 90,4,"", ThetaAW , 90,4,"", " , 0,0,"
9170      DATA 33, ThetaBU , 90,4,"", ThetaBV , 0,4,"", ThetaBW , 90,4,"", " , 0,0,"
9180      DATA 34, ThetaCU , 90,4,"", ThetaCV , 90,4,"", ThetaCW , 0,4,"", Nose ,139,1,cm
9190      DATA 35, UFreqMin, 8,4,MHz, VFreqMin, 25,4,MHz, WFreqMin, 10,4,MHz, " , 0,0,"
9200      DATA 36, UFreqMax, 32,4,MHz, VFreqMax, 55,4,MHz, WFreqMax, 70,4,MHz, Clip ,1,0,"
9210      ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
9220      DATA 41, Xmin1 , 0.00,0,"", Xmax1 , 100,0,"", Ymin1 , 0,0,"", Ymax1 , 100,0,"
9230      DATA 42, Xmin2 , 0.00,0,"", Xmax2 , 100,0,"", Ymin2 , 0,0,"", Ymax2 , 100,0,"
9240      DATA 43, Xmin3 , 0.00,0,"", Xmax3 , 100,0,"", Ymin3 , 0,0,"", Ymax3 , 100,0,"
9250      DATA 44, Xmin4 , -1,2,"", Xmax4 , 1,2,"", Ymin4 , 0,0,"", Ymax4 , 100,0,"
9260      DATA 45, Xmin5 , -1,2,"", Xmax5 , 1,2,"", Ymin5 , 0,0,"", Ymax5 , 100,0,"
9270      DATA 46, Xmin6 , -0.5,1,"", Xmax6 , 1.5,1,"", Ymin6 , 0,2,"", Ymax6 , 4,2,"
9280      DATA 47, Xmin7 , 0,1,"", Xmax7 , .5,1,"", Ymin7 , 0,2,"", Ymax7 , 4,2,"
9290      DATA 48, Xmin8 , -1,1,"", Xmax8 , 1,1,"", Ymin8 , 0,2,"", Ymax8 , 4,2,"
9300      DATA 49, Xmin9 , 0,0,"", Xmax9 , 2000,0,"", Ymin9 , 0,2,"", Ymax9 , 4,2,"
9310      DATA 51, Xmin1 , 935,0,pxl, Xmax1 , 1235,0,pxl, Ymin1 , 725,0,pxl, Ymax1 , 825,0,pxl
9320      DATA 52, Xmin2 , 935,0,pxl, Xmax2 , 1235,0,pxl, Ymin2 , 585,0,pxl, Ymax2 , 685,0,pxl
9330      DATA 53, Xmin3 , 935,0,pxl, Xmax3 , 1235,0,pxl, Ymin3 , 445,0,pxl, Ymax3 , 545,0,pxl
9340      DATA 54, Xmin4 , 935,0,pxl, Xmax4 , 1235,0,pxl, Ymin4 , 305,0,pxl, Ymax4 , 405,0,pxl
9350      DATA 55, Xmin5 , 935,0,pxl, Xmax5 , 1235,0,pxl, Ymin5 , 165,0,pxl, Ymax5 , 265,0,pxl
9360      DATA 56, Xmin6 , 75,0,pxl, Xmax6 , 325,0,pxl, Ymin6 , 525,0,pxl, Ymax6 , 825,0,pxl
9370      DATA 57, Xmin7 , 425,0,pxl, Xmax7 , 675,0,pxl, Ymin7 , 525,0,pxl, Ymax7 , 825,0,pxl
9380      DATA 58, Xmin8 , 75,0,pxl, Xmax8 , 325,0,pxl, Ymin8 , 165,0,pxl, Ymax8 , 465,0,pxl
9390      DATA 59, Xmin9 , 425,0,pxl, Xmax9 , 675,0,pxl, Ymin9 , 165,0,pxl, Ymax9 , 465,0,pxl
9400      DATA 61, Xdiv1 , 10,0,"", Ydiv1 , 4,0,"", Xdiv6 , 4,0,"", Ydiv6 , 8,0,"
9410      DATA 62, Xdiv2 , 10,0,"", Ydiv2 , 4,0,"", Xdiv7 , 5,0,"", Ydiv7 , 8,0,"
9420      DATA 63, Xdiv3 , 10,0,"", Ydiv3 , 4,0,"", Xdiv8 , 4,0,"", Ydiv8 , 8,0,"
9430      DATA 64, Xdiv4 , 4,0,"", Ydiv4 , 4,0,"", Xdiv9 , 4,0,"", Ydiv9 , 8,0,"
9440      DATA 65, Xdiv5 , 4,0,"", Ydiv5 , 4,0,"", " , 0,0,"", " , 0,0,"
9450      ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
9460      ! Delta , 0,4,"", Beta , 0,4,"", Cfreq , 0,0,HZ, Ofreq , 0,4,HZ
9470      ! " , 0,0,"", Ujet/Ue , 1,4,m/s, " , 0,0,"", " , 0,0,"

```

```

!      Y      *****X=1*****      *****X=2*****      *****X=3*****      *****X=4*****
9480      SUBEND
9490      SUB Array_print(Array(*),Name$(*),Image$(*),Units$(*))
9500 Array_print:      PRINT USING "#,5/"
9510      FOR Y=1 TO SIZE(Array,1)
9520      MAT SEARCH Array(Y,*),#LOC(<>0);L1
9530      MAT SEARCH Name$(Y,*),#LOC(<>"");L2
9540      IF L1+L2=0 AND L3=0 THEN 9790
9550      L3=L1+L2
9560      PRINT USING "#,28X"
9570      FOR X=1 TO SIZE(Array,2)
9580      SELECT Name$(Y,X)
9590      CASE ""
9600      PRINT USING "#,28X"
9610      CASE "Date"
9620      LS=DATE$(Array(Y,X))
9630      LS=LS[1,2]&LS[4,6]&LS[8,11]
9640      PRINT USING "#,10A,A,9A,X,3A,4X";TRIM$(Name$(Y,X)),"=",LS,Units$(Y,X)
9650      CASE "Time"
9660      LS=" " & TIME$(Array(Y,X))
9670      PRINT USING "#,10A,A,9A,X,3A,4X";TRIM$(Name$(Y,X)),"=",LS,Units$(Y,X)
9680      CASE ELSE
9690      IF Image$(Y,X)="" THEN Image$(Y,X)="9D"
9700      ON ERROR GOTO 9740
9710      PRINT USING "#,10A,A,"&Image$(Y,X)&"X,3A,4X";TRIM$(Name$(Y,X)),"=",Array(Y,X),Units$(Y,X)
9720      GOTO 9760
9730      OFF ERROR
9740      PRINT USING "#,10A,A,K,X,3A,4X";TRIM$(Name$(Y,X)),"=",Array(Y,X),Units$(Y,X)
9750      END SELECT
9760      NEXT X
9770      PRINT
9780      NEXT Y
9790      SUBEND
9800      !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
9810 Change:      SUB Change(Type$,Array(*),Name$(*),Image$(*),Units$(*))
9820 Change:      PRINTER IS CRT
9830      FOR Y=1 TO SIZE(Array,1)
9840      FOR Y1=Y TO SIZE(Array,1)
9850      FOR X=1 TO SIZE(Array,2)
9860      IF Name$(Y1,X)<>"" THEN 9920
9870      NEXT X
9880      NEXT Y1
9890      CLEAR SCREEN
9900      SUBEXIT
9910      FOR Y2=Y1 TO SIZE(Array,1)
9920      FOR X=1 TO SIZE(Array,2)
9930      IF Name$(Y2,X)<>"" THEN 9970
9940      NEXT X
9950      GOTO 9980
9960      NEXT Y2
9970      FOR Y2=Y2 TO SIZE(Array,1)
9980      FOR X=1 TO SIZE(Array,2)
9990      IF Name$(Y2,X)<>"" THEN 10030
10000      NEXT X
10010      NEXT Y2
10020      Y2=Y2-1
10030      CLEAR SCREEN
10040      CALL Display(Type$,Y1,Y2,Array(*),Name$(*),Image$(*),Units$(*))
10050      Done=0
10060      X=1
10070      Y=Y1
10080      ON KBD ALL,15 GOSUB Kbd
10090      IF NOT Done THEN Wait
10100 Wait:      OFF KBD
10110      CLEAR SCREEN
10120      Y=Y2
10130      NEXT Y
10140      SUBEXIT
10150      CALL Update(Type$,X,Y,Y1,Y2,Done,Array(*),Name$(*),Image$(*),Units$(*))
10160 Kbd:      RETURN
10170      SUBEND
10180      SUB Display(Type$,Y1,Y2,Array(*),Name$(*),Image$(*),Units$(*))
10190 Display:      FOR Y=Y1 TO Y2
10200      FOR X=1 TO SIZE(Array,2)
10210      CALL Select(Type$,X,Y,Y1,Y2,0,Array(*),Name$(*),Image$(*),Units$(*))
10220      NEXT X
10230      NEXT Y
10240      CALL Select(Type$,1,Y1,Y1,Y2,1,Array(*),Name$(*),Image$(*),Units$(*))
10250      SUBEND
10260      SUB Select(Type$,X,Y,Y1,Y2,C,Array(*),Name$(*),Image$(*),Units$(*))
10270 Select:

```

```

10280      PRINT CHR$(128+C);TABXY(26*X-24,15+Y-Y1+1);
10290      PRINT RPT$( " ",23);TABXY(26*X-24,15+Y-Y1+1);
10300      IF Name$(Y,X)="" AND Array(Y,X)=0 THEN 10500
10310      Img$=Image$(Y,X)
10320      Unt$=Units$(Y,X)
10330      IF Image$(Y,X)="" THEN Img$="K"
10340      IF Units$(Y,X)="" THEN Unt$=" "
10350      SELECT Type$
10360      CASE "VALUES"
10370          SELECT Name$(Y,X)
10380          CASE "Date"
10390          CASE "Time"
10400          CASE ELSE
10410              PRINT USING "#,10A,A,"&Img$&","X,3A";Name$(Y,X),":",Array(Y,X),Unt$
10420      END SELECT
10430      CASE "NAMES"
10440          PRINT USING "#,10A,A,8A";Name$(Y,X),":",Name$(Y,X)
10450      CASE "UNITS"
10460          PRINT USING "#,10A,A,8A";Name$(Y,X),":",Units$(Y,X)
10470      CASE "IMAGES"
10480          PRINT USING "#,10A,A,8A";Name$(Y,X),":",Image$(Y,X)
10490      END SELECT
10500      PRINT CHR$(128);
10510  SUBEND
10520 Update: SUB Update(Type$,X,Y,Y1,Y2,Done,Array(*),Name$(*),Image$(*),Units$(*))
10530      DISABLE
10540      KS=KBD$
10550      IF KS="" THEN 11010
10560      SELECT NUM(K$(1,1))
10570      CASE 27                                     ! ESC
10580          Done=1
10590      CASE 255
10600          CALL Select(Type$,X,Y,Y1,Y2,0,Array(*),Name$(*),Image$(*),Units$(*))
10610          SELECT NUM(K$(2,2))
10620          CASE 73,80                               ! Break,Stop
10630              PAUSE
10640          CASE 124                                   ! Menu
10650              Done=1
10660          CASE 38                                   ! Select
10670              CALL Select(Type$,X,Y,Y1,Y2,1,Array(*),Name$(*),Image$(*),Units$(*))
10680              SELECT Type$
10690              CASE "VALUES"
10700                  IF Name$(Y,X)="" THEN CALL Enter_string("Name for "&Name$(Y,X),Name$(Y,X),"K")
10710                  IF Image$(Y,X)="" THEN CALL Enter_string("Image for "&Name$(Y,X),Image$(Y,X),"K")
10720                  CALL Enter_value(Name$(Y,X),Array(Y,X),Image$(Y,X))
10730              CASE "NAMES"
10740                  CALL Enter_string("Name for "&Name$(Y,X),Name$(Y,X),"K")
10750              CASE "UNITS"
10760                  CALL Enter_string("Units for "&Name$(Y,X),Units$(Y,X),"K")
10770              CASE "IMAGES"
10780                  CALL Enter_string("Image for "&Name$(Y,X),Image$(Y,X),"K")
10790              END SELECT
10800              CALL Select(Type$,X,Y,Y1,Y2,0,Array(*),Name$(*),Image$(*),Units$(*))
10810              IF X=SIZE(Array,2) THEN Y=Y+1
10820              X=X+1
10830          CASE 60                                     ! Left
10840              X=X-1
10850          CASE 62                                     ! Right
10860              X=X+1
10870          CASE 94                                     ! Up
10880              Y=Y-1
10890          CASE 86                                     ! Down
10900              Y=Y+1
10910          CASE 92                                     ! First
10920              X=1
10930              Y=1
10940          END SELECT
10950          X=(X-1) MOD SIZE(Array,2)+1
10960          Y=(Y-Y1+1-1) MOD (Y2-Y1+1)+Y1
10970          IF X<1 THEN X=SIZE(Array,2)
10980          IF Y<Y1 THEN Y=Y2
10990          CALL Select(Type$,X,Y,Y1,Y2,1,Array(*),Name$(*),Image$(*),Units$(*))
11000      END SELECT
11010      ENABLE
11020      SUBEXIT
11030  SUBEND
11040 Table: !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
11050 Table: SUB Table(Table(*))
11060      OPTION BASE 1
11070      REAL Mantisa(0:1023),Time(0:1023),Freq(0:1023)

```

```

11080 IF Table(32766) THEN SUBEXIT
11090 FOR Bin=0 TO 1023
11100 Mantisa(Bin)=Bin
11110 NEXT Bin
11120 Mantisa(0)=1
11130 Min=0
11140 FOR Fringes=0 TO 1
11150 FOR Exponent=0 TO 15
11160 Max=Min+1023
11170 IF Max=32767 THEN
11180 Max=32766
11190 REDIM Mantisa(0:1022),Time(0:1022),Freq(0:1022)
11200 END IF
11210 DISP Fringes,Exponent
11220 MAT Time= Mantisa*(2^(Exponent-1)/500000000)
11230 MAT Freq= (2^(4-Fringes))/Time
11240 MAT Freq= Freq/(1000000)
11250 MAT Table(Min:Max)= Freq
11260 Min=Min+1024
11270 NEXT Exponent
11280 NEXT Fringes
11290 SUBEND
11300 Ctm:
11310 Ctm_ldv: SUB Ctm_ldv(Index(*),Theta1(*),Tun2ldv(*),Ldv2tun(*))
11320 OPTION BASE 1
11330 REAL Theta2(3,3)
11340 ! Correct Theta for angles in water
11350 MAT Theta2= Theta1
11360 ! Theta2(2,1)=ASN(Index(1)/Index(3)*SIN(Theta2(2,1)))
11370 ! Theta2(2,2)=ASN(Index(1)/Index(3)*SIN(Theta2(2,2)))+90
11380 ! Tun2Lvd converts tunnel coordinates to laser coordinates.
11390 Tun2ldv(1,1)=COS(Theta2(1,1))
11400 Tun2ldv(1,2)=COS(Theta2(1,2))
11410 Tun2ldv(1,3)=COS(Theta2(1,3))
11420 Tun2ldv(2,1)=COS(Theta2(2,1))
11430 Tun2ldv(2,2)=COS(Theta2(2,2))
11440 Tun2ldv(2,3)=COS(Theta2(2,3))
11450 Tun2ldv(3,1)=COS(Theta2(3,1))
11460 Tun2ldv(3,2)=COS(Theta2(3,2))
11470 Tun2ldv(3,3)=COS(Theta2(3,3))
11480 ! Ldv2tun converts laser coordinates to tunnel coordinates.
11490 MAT Ldv2tun= INV(Tun2ldv)
11500 SUBEND
11510 Ctm_mod: SUB Ctm_mod(Alpha,Alpha2,Alpha3,Mod2tun(*),Tun2mod(*))
11520 OPTION BASE 1
11530 REAL T1(3,3),T2(3,3),T3(3,3),Abc(3),Abc1(3),Abc2(3),Temp(3,3)
11540 ! Define 1st coordinate transformation matrix for Mod2tun.
11550 T1(1,1)=COS(Alpha)
11560 T1(1,2)=SIN(Alpha)
11570 T1(1,3)=0
11580 T1(2,1)=-SIN(Alpha)
11590 T1(2,2)=COS(Alpha)
11600 T1(2,3)=0
11610 T1(3,1)=0
11620 T1(3,2)=0
11630 T1(3,3)=1
11640 ! Define 2nd coordinate transformation matrix for Mod2tun.
11650 T2(1,1)=1
11660 T2(1,2)=0
11670 T2(1,3)=0
11680 T2(2,1)=0
11690 T2(2,2)=COS(-Alpha2)
11700 T2(2,3)=SIN(-Alpha2)
11710 T2(3,1)=0
11720 T2(3,2)=-SIN(-Alpha2)
11730 T2(3,3)=COS(-Alpha2)
11740 ! Define 3rd coordinate transformation matrix for Mod2tun.
11750 Abc1(1)=1
11760 Abc1(2)=0
11770 Abc1(3)=0
11780 MAT Abc2= T1*Abc1
11790 MAT Abc= T2*Abc2
11800 T3(1,1)=Abc(1)*Abc(1)*(1-COS(-Alpha3))+COS(-Alpha3)
11810 T3(1,2)=Abc(2)*Abc(1)*(1-COS(-Alpha3))+Abc(3)*SIN(-Alpha3)
11820 T3(1,3)=Abc(3)*Abc(1)*(1-COS(-Alpha3))-Abc(2)*SIN(-Alpha3)
11830 T3(2,1)=Abc(1)*Abc(2)*(1-COS(-Alpha3))-Abc(3)*SIN(-Alpha3)
11840 T3(2,2)=Abc(2)*Abc(2)*(1-COS(-Alpha3))+COS(-Alpha3)
11850 T3(2,3)=Abc(3)*Abc(2)*(1-COS(-Alpha3))+Abc(1)*SIN(-Alpha3)
11860 T3(3,1)=Abc(1)*Abc(3)*(1-COS(-Alpha3))+Abc(2)*SIN(-Alpha3)
11870 T3(3,2)=Abc(2)*Abc(3)*(1-COS(-Alpha3))-Abc(1)*SIN(-Alpha3)

```



```

11880      T3(3,3)=Abc(3)*Abc(3)*(1-COS(-Alpha3))+COS(-Alpha3)
11890      ! Mod2tun converts model coordinates to tunnel coordinates.
11900      MAT Temp= T2*T1
11910      MAT Mod2tun= T3*Temp
11920      ! Tun2mod converts tunnel coordinates to model coordinates.
11930      MAT Tun2mod= INV(Mod2tun)
11940
11950 Ctm_tcs1: SUBEND
11960 SUB Ctm_tcs1(Tcs2tun(*),Tun2tcs(*))
11970   OPTION BASE 1
11980   REAL Nair,Nglass,Nwater
11990   REAL Flonaxis,Floffaxis,Bsonaxis,Bsoffaxis
12000   REAL Theta(4),Onaxis,Offaxis
12010   REAL Xon,Yon,Xoff,Yoff,X1,Y1,Y2
12020   REAL Ba,Bb,Xc,Yc
12030   REAL X(4),Yposition,Thickness
12040   INTEGER Offa,Offb,Ona,Onb,Beam,I,J
12050   Offa=1
12060   Offb=2
12070   Ona=3
12080   Onb=4
12090   Flonaxis=19.25
12100   Floffaxis=19.25
12110   Bsonaxis=60/25.4
12120   Bsoffaxis=60/25.4
12130   Thickness=1.25
12140   Onaxis=0.
12150   Offaxis=45.0
12160   Nair=1.00
12170   Nglass=1.43
12180   Nwater=1.33
12190   Yposition=0
12200   GOSUB Findstart
12210   Y1=Yon
12220   X1=Xoff
12230   Y2=Yoff
12240   Yposition=1
12250   GOSUB Findstart
12260   Y2=Yon-Y1+Y2
12270   MAT Tun2tcs= IDN
12280   Tun2tcs(2,2)=- (Yon-Y1)
12290   Tun2tcs(4,2)=-SQRT((Xoff-X1)^2+(Yoff-Y2)^2)
12300   Tun2tcs(4,4)=0
12310   MAT Tcs2tun= INV(Tun2tcs)
12320   Tcs2tun(4,2)=0
12330   MAT Tun2tcs= IDN
12340   MAT Tcs2tun= IDN
12350   SUBEXIT
12360   Theta(Offa)=Offaxis+ATN(Bsoffaxis/(2*Floffaxis))
12370   Theta(Offb)=Offaxis-ATN(Bsoffaxis/(2*Floffaxis))
12380   Theta(Ona)=Onaxis+ATN(Bsonaxis/(2*Flonaxis))
12390   Theta(Onb)=Onaxis-ATN(Bsonaxis/(2*Flonaxis))
12400   FOR Beam=Offa TO Onb
12410     X(Beam)=-Yposition*TAN(ASN(Nair/Nwater*SIN(Theta(Beam))))-
12420     Thickness*TAN(ASN(Nair/Nglass*SIN(Theta(Beam))))
12430     NEXT Beam
12440     Ba=-Thickness-X(Offa)/TAN(Theta(Offa))
12450     Bb=-Thickness-X(Offb)/TAN(Theta(Offb))
12460     Xc=(Bb-Ba)/(1/TAN(Theta(Offa))-1/TAN(Theta(Offb)))
12470     Yc=Xc/TAN(Theta(Offb))+Bb
12480     Xoff=Xc-Floffaxis*SIN(Offaxis)
12490     Yoff=Yc-Floffaxis*COS(Offaxis)
12500     Ba=-Thickness-X(Ona)/TAN(Theta(Ona))
12510     Bb=-Thickness-X(Onb)/TAN(Theta(Onb))
12520     Xc=(Bb-Ba)/(1/TAN(Theta(Ona))-1/TAN(Theta(Onb)))
12530     Yc=Xc/TAN(Theta(Onb))+Bb
12540     Xon=Xc-Flonaxis*SIN(Onaxis)
12550     Yon=Yc-Flonaxis*COS(Onaxis)
12560     RETURN
12570 SUBEND
12580 SUB Ctm_tcs2(Tcs2tun(*),Tun2tcs(*))
12590   OPTION BASE 1
12600   REAL Nair,Nglass,Nwater
12610   REAL Floffaxis,Bsoffaxis
12620   REAL Theta(2),Offaxis
12630   REAL Xoff,Yoff,X1,Y1
12640   REAL Ba,Bb,Xc,Yc
12650   REAL X(2),Yposition,Thickness
12660   INTEGER Offa,Offb,Beam,I,J
12670   Offa=1
12680   Offb=2

```

```

12670      Floffaxis=19.5
12680      Bsoffaxis=60/25.4
12690      Thickness=1.25
12700      Offaxis=-13.2
12710      Nair=1.00
12720      Nglass=1.43
12730      Nwater=1.33
12740      Yposition=0
12750      GOSUB Findstart
12760      X1=Xoff
12770      Y1=Yoff
12780      Yposition=1
12790      GOSUB Findstart
12800      X2=Xoff
12810      Y2=Yoff
12820      Kx=(X2-X1)
12830      Ky=(Y2-Y1)
12840      MAT Tun2tcs= IDN
12850      Tun2tcs(1,2)=Kx
12860      Tun2tcs(2,2)=-Ky
12870      Tun2tcs(4,4)=0
12880      MAT Tcs2tun= INV(Tun2tcs)
12890      Tcs2tun(4,2)=0
12900      MAT Tun2tcs= IDN
12910      MAT Tcs2tun= IDN
12920      SUBEXIT
12930 Findstart:  Theta(Offa)=Offaxis+ATN(Bsoffaxis/(2*Floffaxis))
12940      Theta(Offb)=Offaxis-ATN(Bsoffaxis/(2*Floffaxis))
12950      FOR Beam=Offa TO Offb
12960          X(Beam)=-Yposition*TAN(ASN(Nair/Nwater*SIN(Theta(Beam))))-
          Thickness*TAN(ASN(Nair/Nglass*SIN(Theta(Beam))))
12970      NEXT Beam
12980      Ba=-Thickness-X(Offa)/TAN(Theta(Offa))
12990      Bb=-Thickness-X(Offb)/TAN(Theta(Offb))
13000      Xc=(Bb-Ba)/(1/TAN(Theta(Offa))-1/TAN(Theta(Offb)))
13010      Yc=Xc/TAN(Theta(Offb))+Bb
13020      Xoff=Xc-Floffaxis*SIN(Offaxis)
13030      Yoff=Yc-Floffaxis*COS(Offaxis)
13040      RETURN
13050      SUBEND
13060 Tcs8:      !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
13070 Tcs8init:  SUB Tcs8init(@Tcs8)
13080      REAL I(1:8),C(1:8)
13090      ASSIGN @Tcs8 TO 9;BYTE,FORMAT OFF,EOL ""
13100      CONTROL 9,0;1
13110      CONTROL 9,3;9600
13120      CONTROL 9,4;31
13130      CONTROL 9,12;IVAL("EF",16)
13140      CONTROL 9,13;9600
13150      CONTROL 9,14;31
13160      OUTPUT @Tcs8 USING "K, /";"VIO"
13170      ENTER @Tcs8 USING "8(K)";I(*)
13180      IF SUM(I)<>8 THEN OUTPUT @Tcs8 USING "K, /";"SIO"
13190      OUTPUT @Tcs8 USING "K, /";"VCO"
13200      ENTER @Tcs8 USING "8(K)";C(*)
13210      IF SUM(C)<>8 THEN OUTPUT @Tcs8 USING "K, /";"SCO:1,"
13220      !OUTPUT @Tcs8 USING "K, /";"SCO:0,"
13230      SUBEND
13240 Tcs8set:  SUB Tcs8set(C$,@Tcs8)
13250      OPTION BASE 1
13260      DIM View(8,1),Set(8,2),Name$(8,1)[10],Image$(8,1)[10],Units$(8,1)[10]
13270      OUTPUT @Tcs8 USING "K, /";"V"&C$&"0"
13280      ENTER @Tcs8 USING "8(K)";View(*)
13290      READ Name$(*)
13300      MAT Image$= ("6D.3D")
13310      DATA X1,X2,Y1,Y2,Z1,Z2,A1,A2
13320      FOR Channel=1 TO 8
13330          Set(Channel,1)=Channel
13340          SELECT C$
13350              CASE "P"
13360                  Name$(Channel,1)=Name$(Channel,1)&" (pos)"
13370                  Units$(Channel,1)="in"
13380              CASE "U"
13390                  Name$(Channel,1)=Name$(Channel,1)&" (cpi)"
13400                  Units$(Channel,1)="cnt"
13410              CASE "R"
13420                  Name$(Channel,1)=Name$(Channel,1)&" (cpr)"
13430                  Units$(Channel,1)="cnt"
13440              CASE "V"
13450                  Name$(Channel,1)=Name$(Channel,1)&" (vel)"

```

```

13460         Units$(Channel,1)="rev"
13470     CASE "A"
13480         Name$(Channel,1)=Name$(Channel,1)&" (acc)"
13490         Units$(Channel,1)="rev"
13500     CASE "+"
13510         Name$(Channel,1)=Name$(Channel,1)&" (+LS)"
13520         Units$(Channel,1)=" "
13530     CASE "-"
13540         Name$(Channel,1)=Name$(Channel,1)&" (-LS)"
13550         Units$(Channel,1)=" "
13560     CASE "S"
13570         Name$(Channel,1)=Name$(Channel,1)&" (STALL)"
13580         Units$(Channel,1)=" "
13590     CASE "H"
13600         Name$(Channel,1)=Name$(Channel,1)&" (HS)"
13610         Units$(Channel,1)=" "
13620     END SELECT
13630 NEXT Channel
13640 CALL Change("VALUES",View(*),Name$(*),Image$(*),Units$(*))
13650 SELECT CS
13660 CASE "P","U","R","V","A"
13670     MAT Set(*,2)= View(*,1)
13680     OUTPUT @Tcs8 USING 13690;"S"&CS,Set(*)
13690     IMAGE K,8(D,"",M6D.4D,""),/
13700 END SELECT
13710 SUBEND
13720 Tcs8read: SUB Tcs8read(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Tcs2tun1(*),Tcs2tun2(*),Tun2mod(*))
13730     OUTPUT @Tcs8 USING "K,/","VP0"
13740     ENTER @Tcs8 USING "8(K)";Tcs1(1),Tcs2(1),Tcs1(2),Tcs2(2),Tcs1(3),Tcs2(3),Tcs1(4),Tcs2(4)
13750     MAT Tun= Tcs2tun1*Tcs1
13760     REDIM Tun(1:3),Mod(1:3)
13770     MAT Mod= Tun2mod*Tun
13780     REDIM Tun(1:4),Mod(1:4)
13790     Mod(4)=0
13800     Tun(4)=0
13810     CALL Tcs8print(Mod(*),Tun(*),Tcs1(*),Tcs2(*))
13820 SUBEND
13830 Tcs8print: SUB Tcs8print(Mod(*),Tun(*),Tcs1(*),Tcs2(*))
13840     PRINT CHR$(128);
13850     PRINT TABXY(50,1);"
13860     PRINT TABXY(50,2);"          MOD          TUN          TCS1          TCS2          "
13870     PRINT TABXY(50,3);"
13880     PRINT TABXY(50,4);
13890     PRINT USING "#,K,4 (M3D.4D),X";" X:",Mod(1),Tun(1),Tcs1(1),Tcs2(1)
13900     PRINT TABXY(50,5);
13910     PRINT USING "#,K,4 (M3D.4D),X";" Y:",Mod(2),Tun(2),Tcs1(2),Tcs2(2)
13920     PRINT TABXY(50,6);
13930     PRINT USING "#,K,4 (M3D.4D),X";" Z:",Mod(3),Tun(3),Tcs1(3),Tcs2(3)
13940     PRINT TABXY(50,7);
13950     PRINT USING "#,K,4 (M3D.4D),X";" A:",Mod(4),Tun(4),Tcs1(4),Tcs2(4)
13960     PRINT TABXY(50,8);"
13970 SUBEND
13980 Tcs8move: SUB
13990     Tcs8move(@Tcs8,Mod(*),Tun(*),Tcs1(*),Tcs2(*),Mod2tun(*),Tun2tcs1(*),Tun2tcs2(*),Side$,Coor$,Mode$,K,Mov
14000     ement)
14010     OPTION BASE 1
14020     DIM LS[100]
14030     REAL Move(8,2),I(8),C(8)
14040     IF Mode$="RELATIVE" THEN
14050         MAT Mod= (0)
14060         MAT Tun= (0)
14070         MAT Tcs1= (0)
14080         MAT Tcs2= (0)
14090     END IF
14100     SELECT Coor$
14110     CASE "MODEL"
14120         Mod(K)=Movement
14130         REDIM Tun(1:3),Mod(1:3)
14140         MAT Tun= Mod2tun*Mod
14150         REDIM Tun(1:4),Mod(1:4)
14160         IF POS(Side$,"Tx") THEN MAT Tcs1= Tun2tcs1*Tun
14170         IF POS(Side$,"Rx") THEN MAT Tcs2= Tun2tcs2*Tun
14180     CASE "TUNNEL"
14190         Tun(K)=Movement
14200         IF POS(Side$,"Tx") THEN MAT Tcs1= Tun2tcs1*Tun
14210         IF POS(Side$,"Rx") THEN MAT Tcs2= Tun2tcs2*Tun
14220     CASE "LASER"
14230         IF POS(Side$,"Tx") THEN Tcs1(K)=Movement
14240         IF POS(Side$,"Rx") THEN Tcs2(K)=Movement
14250     END SELECT

```

```

14240      FOR Channel=1 TO 8
14250          Move(Channel,1)=Channel
14260      NEXT Channel
14270      Move(1,2)=Tcs1(1)
14280      Move(2,2)=Tcs2(1)
14290      Move(3,2)=Tcs1(2)
14300      Move(4,2)=Tcs2(2)
14310      Move(5,2)=Tcs1(3)
14320      Move(6,2)=Tcs2(3)
14330      Move(7,2)=Tcs1(4)
14340      Move(8,2)=Tcs2(4)
14350      SELECT Mode$
14360      CASE "ABSOLUTE"
14370          OUTPUT @Tcs8 USING 14380;"MA",3,4,Move(3,2)
14380          IMAGE K,1(D,D,":",K,""),/
14390          ENTER @Tcs8 USING "K";L$ ! Tcs1(2)
14400          Tcs1(2)=VAL(L$)
14410          ENTER @Tcs8 USING "K";L$ ! Tcs2(2)
14420          Tcs2(2)=VAL(L$)
14430      CASE "RELATIVE"
14440          OUTPUT @Tcs8 USING 14450;"MR",Move(*)
14450          IMAGE K,8(D,":",S2D.5D,""),/
14460          ENTER @Tcs8 USING "8(K)";Tcs1(1),Tcs2(1),Tcs1(2),Tcs2(2),Tcs1(3),Tcs2(3),Tcs1(4),Tcs2(4)
14470      END SELECT
14480  SUBEND
14490 Tcs8view: SUB Tcs8view(@Tcs8)
14500      OPTION BASE 1
14510      REAL View(8)
14520      CS="--+HS"
14530      CLEAR SCREEN
14540      PRINT TABXY(1,1);"      X1 X2 Y1 Y2 Z1 Z2 A1 A2"
14550      FOR I=1 TO 4
14560          OUTPUT @Tcs8 USING "K,/";"V"&CS[I,I]&"0"
14570          ENTER @Tcs8 USING "8(K)";View(*)
14580          PRINT USING "AA,5X,8(3D)";"V"&CS[I,I],View(*)
14590      NEXT I
14600      BEEP
14610      GOTO 14540
14620  SUBEND
14630 Graph: !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
14640 Dump: SUB Dump(G1,G2,Prt,Array(*),INTEGER Gs(*))
14650      OPTION BASE 1
14660      ALLOCATE INTEGER Ws(400,400)
14670      GSTORE Gs(*)
14680      KEY LABELS OFF
14690      OUTPUT Prt USING "#,@"
14700      FOR G=G1 TO G2
14710          Xmin=Array(G+40,1)
14720          Xmax=Array(G+40,2)
14730          Ymin=Array(G+40,3)
14740          Ymax=Array(G+40,4)
14750          Xpix1=Array(G+50,1)-75
14760          Xpix2=Array(G+50,2)+25
14770          Ypix1=Array(G+50,3)-50
14780          Ypix2=Array(G+50,4)+25
14790          VIEWPORT Xpix1/10.23,Xpix2/10.23,Ypix1/10.23,Ypix2/10.23
14800          WINDOW 0,1,0,1
14810          CALL Bstore(Ws(*),(Xpix2-Xpix1)+1,(Ypix2-Ypix1)+1,3,0,1)
14820          GCLEAR
14830          CLEAR SCREEN
14840          Xnew=100-Xpix1
14850          Ynew=400-Ypix1
14860          Xpix1=Xpix1+Xnew
14870          Xpix2=Xpix2+Xnew
14880          Ypix1=Ypix1+Ynew
14890          Ypix2=Ypix2+Ynew
14900          WINDOW 0,1,0,1
14910          VIEWPORT Xpix1/10.23,Xpix2/10.23,Ypix1/10.23,Ypix2/10.23
14920          CALL Bload(Ws(*),(Xpix2-Xpix1)+1,(Ypix2-Ypix1)+1,3,0,1)
14930          CALL Gdump_colored(CRT,Prt,"NORMAL",180,"OFF","DITHER")
14940          GLOAD Gs(*)
14950      NEXT G
14960      DEALLOCATE Ws(*)
14970  SUBEND
14980 Crt_init: SUB Crt_init
14990      PLOTTER IS CRT,"INTERNAL"
15000      AREA PEN 0
15010      PEN 1
15020      PRINTER IS CRT
15030      PRINTALL IS CRT

```

```

15040      KEY LABELS OFF
15050      SUBEND
15060 Read_symbols: SUB Read_symbols(Symbols(*))
15070          OPTION BASE 1
15080          REAL Symbol(20,3),Dot(2,3)
15090          READ Dot(*)
15100          FOR S=1 TO 5
15110              READ Noc
15120              REDIM Symbol(Noc,3)
15130              READ Symbol(*)
15140              MAT Symbols(S,1:Noc,*)=Symbol
15150              MAT Symbols(S,Noc+1:Noc+2,*)=Dot
15160              Symbols(S,0,1)=Noc+2
15170          NEXT S
15180 Dot:      DATA 4.5, 7.5,-2, 4.5, 7.5,-1
15190 Square:    DATA 5, 0.5, 3.5,-2, 8.5, 3.5,-1, 8.5,11.5,-1, 0.5,11.5,-1, 0.5,3.5,-1
15200 Octagon:   DATA 9, 0.5, 5.5,-2, 2.5, 3.5,-1, 6.5, 3.5,-1, 8.5, 5.5,-1, 8.5,9.5,-1, 6.5,11.5,-1, 2.5,11.5,-1
15210 Diamond:   DATA 5, -0.5, 7.5,-2, 4.5, 2.5,-1, 9.5, 7.5,-1, 4.5,12.5,-1, -0.5,7.5,-1
15220 Utriangle: DATA 4, 0.5, 4.5,-2, 8.5, 4.5,-1, 4.5,13.5,-1, 0.5, 4.5,-1
15230 Dtriangle: DATA 4, 0.5,10.5,-2, 8.5,10.5,-1, 4.5, 1.5,-1, 0.5,10.5,-1
15240      SUBEND
15250 Setup_graph: SUB Setup_graph(Array(*),Image$(*),Paxis,Symbols(*))
15260          OPTION BASE 1
15270          COM /Graph/
15280          Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),Legend$(*)
15290          MAT Wndw= Array(41:49,*)
15300          MAT Vwprt= Array(51:59,*)
15310          MAT Xdiv(1:5)= Array(61:65,1)
15320          MAT Xdiv(6:9)= Array(61:64,3)
15330          MAT Ydiv(1:5)= Array(61:65,2)
15340          MAT Ydiv(6:9)= Array(61:64,4)
15350          MAT Ximage$= Image$(41:49,1)
15360          MAT Yimage$= Image$(41:49,3)
15370          FOR G=1 TO 9
15380              READ G,Xlabel$(G)
15390              FOR I=1 TO SIZE(Legend$,2)
15400                  READ Legend$(G,I)
15410              NEXT I
15420              SELECT G
15430                  CASE 1 TO 5
15440                      Ylabel$(G)=" "
15450                  CASE 6 TO 9
15460                      Ylabel$(G)=CHR$(NUM("X")+Paxis-1)
15470              END SELECT
15480              CALL Set_up(G,Symbols(*))
15490          NEXT G
15500          SUBEXIT
15510          DATA 1, " " , " " , " " , " " , " " , " " , " " , " "
15520          DATA 2, " " , " " , " " , " " , " " , " " , " " , " "
15530          DATA 3, " " , " " , " " , " " , " " , " " , " " , " "
15540          DATA 4, " " , " " , " " , " " , " " , " " , " " , " "
15550          DATA 5, " " , " " , " " , " " , " " , " " , " " , " "
15560          DATA 6, "Velocities / Uinf" , "U" , "V" , "W" , " " , " " , " " , " "
15570          DATA 7, "RMS / Uinf" , "U1" , "V1" , "W1" , " " , " " , " " , " "
15580          DATA 8, "Shear Stress / Uinf^2" , "U1V1" , "V1W1" , "W1U1" , " " , " " , " " , " "
15590          DATA 9, "Tt:3R Uinf:m/s Uedge:m/s" , "Tt" , "Uinf" , "Uedge" , " " , " " , " " , " "
15600      SUBEND
15610 Set_up:      SUB Set_up(G,Symbols(*))
15620          OPTION BASE 1
15630          COM /Graph/
15640          Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),Legend$(*)
15650          DIM LS(80)
15660          ON ERROR CALL Error
15670          PLOTTER IS CRT,"INTERNAL"
15680          Black=-1
15690          White=1
15700          CSIZE 100*15/1023
15710          Xmin=Wndw(G,1)
15720          Xmax=Wndw(G,2)
15730          Ymin=Wndw(G,3)
15740          Ymax=Wndw(G,4)
15750          Xpix1=Vwprt(G,1)
15760          Xpix2=Vwprt(G,2)
15770          Ypix1=Vwprt(G,3)
15780          Ypix2=Vwprt(G,4)
15790          Xstep=(Xmax-Xmin)/Xdiv(G)
15800          Ystep=(Ymax-Ymin)/Ydiv(G)
15810          Xpixel=(Xmax-Xmin)/(Xpix2-Xpix1)
15820          Ypixel=(Ymax-Ymin)/(Ypix2-Ypix1)

```

```

15810 AREA PEN Black
15820 PEN White
15830 GOSUB Back_ground
15840 GOSUB Axes
15850 !GOSUB Grid
15860 GOSUB Plot_area
15870 CLIP OFF
15880 GOSUB Ylabel
15890 GOSUB Xlabel
15900 CALL Legend(G, Symbols(*))
15910 OFF ERROR
15920 SUBEXIT
15930 Back_ground: VIEWPORT (Xpix1-75)/10.23, (Xpix2+25)/10.23, (Ypix1-33)/10.23, (Ypix2+6)/10.23
15940 WINDOW -1.E+9, 1.E+9, -1.E+9, 1.E+9
15950 MOVE 0,0
15960 WINDOW 0,1,0,1
15970 MOVE 0,0
15980 RECTANGLE 1,1, FILL
15990 RETURN
16000 Axes: VIEWPORT (Xpix1-1)/10.23, (Xpix2+1)/10.23, (Ypix1-6)/10.23, (Ypix2+1)/10.23
16010 WINDOW Xmin, Xmax, 1,0
16020 AXES Xstep, 2, Xmin, 0, 1, 1, 1
16030 VIEWPORT (Xpix1-1)/10.23, (Xpix2+1)/10.23, (Ypix2+1)/10.23, (Ypix2+6)/10.23
16040 WINDOW Xmin, Xmax, 0, 1
16050 AXES Xstep, 2, Xmin, 0, 1, 1, 1
16060 VIEWPORT (Xpix1-6)/10.23, (Xpix1-1)/10.23, (Ypix1-1)/10.23, (Ypix2+1)/10.23
16070 WINDOW 1,0, Ymin, Ymax
16080 AXES 2, Ystep, 0, Ymin, 1, 1, 1
16090 VIEWPORT (Xpix2+1)/10.23, (Xpix2+6)/10.23, (Ypix1-1)/10.23, (Ypix2+1)/10.23
16100 WINDOW 0,1, Ymin, Ymax
16110 AXES 2, Ystep, 0, Ymin, 1, 1, 1
16120 RETURN
16130 Grid: VIEWPORT (Xpix1-1)/10.23, (Xpix2+1)/10.23, (Ypix1-1)/10.23, (Ypix2+1)/10.23
16140 WINDOW Xmin, Xmax, Ymin, Ymax
16150 LINE TYPE 4
16160 GRID Xstep, Ystep, Xmin, Ymin
16170 LINE TYPE 1
16180 RETURN
16190 Plot_area: VIEWPORT Xpix1/10.23, Xpix2/10.23, Ypix1/10.23, Ypix2/10.23
16200 WINDOW Xmin, Xmax, Ymin, Ymax
16210 RETURN
16220 Xlabel: LONG 5
16230 FOR X=Xmin TO Xmax+Xstep/100 STEP Xstep
16240 MOVE X, Ymin-12*Ypixel
16250 OUTPUT L$ USING Ximage$(G); X
16260 LABEL TRIMS(L$)
16270 NEXT X
16280 MOVE (Xmin+Xmax)/2, Ymin-25*Ypixel
16290 LABEL Xlabel$(G)
16300 RETURN
16310 Ylabel: LONG 8
16320 Len=0
16330 FOR Y=Ymin TO Ymax+Ystep/100 STEP Ystep
16340 MOVE Xmin-5*Xpixel, Y
16350 OUTPUT L$ USING Yimage$(G); Y
16360 LABEL TRIMS(L$)
16370 Len=MAX(Len, LEN(TRIMS(L$)))
16380 NEXT Y
16390 MOVE Xmin-(5+7*Len)*Xpixel, (Ymin+Ymax)/2
16400 LABEL Ylabel$(G)
16410 LONG 5
16420 RETURN
16430 SUBEND
16440 Legend: SUB Legend(G, Symbols(*))
16450 OPTION BASE 1
16460 COM /Graph/
16470 Wndw(*), Vwprt(*), Xdiv(*), Ydiv(*), Xlabel$(*), Ylabel$(*), Title$(*), Ximage$(*), Yimage$(*), Legend$(*)
16480 DIM Symbol(20, 3)
16490 VIEWPORT Vwprt(G, 1)/10.23, Vwprt(G, 2)/10.23, Vwprt(G, 3)/10.23, Vwprt(G, 4)/10.23
16500 WINDOW Vwprt(G, 1), Vwprt(G, 2), Vwprt(G, 3), Vwprt(G, 4)
16510 Black=-1
16520 White=1
16530 CSIZE 100*15/1023
16540 AREA PEN -1 ! Black
16550 PEN 1 ! White
16560 LONG 2
16570 Len=0
16580 FOR S=1 TO SIZE(Legend$, 2)
16590 Len=MAX(Len, LEN(Legend$(G, S)), Len)
NEXT S

```

```

16600         FOR S=1 TO SIZE(Legend$,2)
16610             IF LEN(Legend$(G,S))=0 THEN 16690
16620                 Noc=Symbols(S,0,1)
16630                 REDIM Symbol(Noc,3)
16640                 MAT Symbol= Symbols(S,1:Noc,*)
16650                 MOVE Vwprt(G,2)-7*Len-23,Vwprt(G,4)-15*S+5
16660                 SYMBOL Symbol(*),FILL,EDGE
16670                 MOVE Vwprt(G,2)-7*Len-10,Vwprt(G,4)-15*S+4
16680                 LABEL Legend$(G,S)
16690             NEXT S
16700         SUBEND
16710 Lvdas:      !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
16720 Lvdas_init: SUB Lvdas_init(@Gpio)
16730             ASSIGN @Gpio TO 12;WORD,FORMAT OFF,EOL ""
16740             OUTPUT @Gpio USING "#,AA";"HP"
16750         SUBEND
16760 Lvdas_sample_a: SUB Lvdas_sample_a(@Lvdas,Channel,Symbol(*))
16770             OPTION BASE 1
16780             INTEGER Gx,Gy,Data(1000,4),G(128,102),Iv(1000)
16790             DIM Ls[80],V(1000),Vv(1000),T(1000),Wndw(4),Vwprt(4)
16800             Black=-1
16810             White=1
16820             READ Wndw(*),Xdiv,Ydiv,Vwprt(*),Ximage$,Yimage$,Xlabel$,Ylabel$
16830             ! Xmin,Xmax,Ymin,Ymax,Xdiv,Ydiv,Xpix1,Xpix2,Ypix1,Ypix2,Ximage$,Yimage$,Xlabel$,Ylabel$
16840             DATA 0,.001,-5,5,10,10,75,1235,165,825,6D.4D,6D.3D,t(sec),V
16850             CALL Set_up(Wndw(*),Vwprt(*),Xdiv,Ydiv,Xlabel$,Ylabel$,Ximage$,Yimage$)
16860             GSTORE G(*)
16870             PEN White
16880             OUTPUT @Lvdas USING "#,AA";"DT"
16890             OUTPUT @Lvdas USING "#,AA,W";"SC",Channel
16900             OUTPUT @Lvdas USING "AA";"RM"
16910             OUTPUT @Lvdas USING "W,W";IVAL("08F2",16),IVAL("0000",16)
16920             OUTPUT @Lvdas USING "W,W";IVAL("08F2",16),IVAL("1F3F",16)
16930             ENTER @Lvdas USING "#,W";Data(*)
16940             OUTPUT @Lvdas USING "#,AA";"ET"
16950             MAT T= Data(*,2)
16960             MAT V= Data(*,4)
16970             MAT Vv= V*(5./2.^15)
16980             MAT Vv= V . V
16990             Ave=SUM(V)/1000
17000             Sdv=SQR(SUM(Vv)/1000-Ave*Ave)
17010             MAT SEARCH V(*),MIN;Min
17020             MAT SEARCH V(*),MAX;Max
17030             Dif=Max-Min
17040             GLOAD G(*)
17050             MOVE Xmin+10*Xpixel,Ymax-20*Ypixel
17060             LOG 2
17070             LABEL USING "5(M5D.4D)";Ave,Sdv,Min,Max,Dif
17080             Ave=Ave/5*2^15
17090             Sdv=Sdv/5*2^15
17100             Min=Min/5*2^15
17110             Max=Max/5*2^15
17120             Dif=Max-Min
17130             LABEL
17140             LABEL USING "2(M8D.1D),3(M10D)";Ave,Sdv,Min,Max,Dif
17150             Time=0
17160             LOG 5
17170             CLIP ON
17180             FOR I=1 TO 1000
17190                 PLOT Time,V(I)
17200                 SYMBOL Symbol(*),EDGE
17210                 MOVE Time,V(I)
17220                 PLOT Time,V(I)
17230                 Time=Time+T(I)*.0000001
17240             NEXT I
17250             GOTO 16880
17260             SUBEXIT
17270         SUBEND
17280 Lvdas_average: SUB Lvdas_average(Table(*),INTEGER Data(*),REAL Vave,Vsdv,Tave,Tsdv)
17290             OPTION BASE 1
17300             REAL V(1000),Vv(1000),T(1000),Tt(1000)
17310             N=SIZE(Data,1)
17320             REDIM V(N),Vv(N),T(N),Tt(N)
17330             Channel=Data(1,3)+1
17340             SELECT Channel
17350             CASE 1,2,3
17360                 FOR I=1 TO N
17370                     V(I)=Table(BINAND(32767,BINCOMP(Data(I,4))))
17380                 NEXT I
17390             CASE 4,5

```



```

18200     FOR K=1 TO N
18210         OUTPUT @Mac USING 18220;K,Valid(K),DROUND(U(K),5),DROUND(V(K),5),DROUND(A(K),5)
18220         IMAGE K," ",K," ",2(K," "),K,/
18230     NEXT K
18240     OUTPUT @Mac USING "@,/"
18250 SUBEND
18260 Data_clip: SUB Data_clip(INTEGER Nsam,REAL Umin,Umax,Vmin,Vmax)
18270     OPTION BASE 1
18280     COM /Data/ INTEGER Raw(*),Valid(*),REAL Table(*),U(*),V(*),W(*),A(*),B(*),I(*),C(*)
18290     DISP "Clipping Histograms"
18300     FOR K=1 TO Nsam
18310         MAT SEARCH U(*),LOC(<Umin);L,K
18320         IF L<Nsam THEN Valid(L)=0
18330         K=L
18340     NEXT K
18350     FOR K=1 TO Nsam
18360         MAT SEARCH U(*),LOC(>Umax);L,K
18370         IF L<Nsam THEN Valid(L)=0
18380         K=L
18390     NEXT K
18400     FOR K=1 TO Nsam
18410         MAT SEARCH V(*),LOC(<Vmin);L,K
18420         IF L<Nsam THEN Valid(L)=0
18430         K=L
18440     NEXT K
18450     FOR K=1 TO Nsam
18460         MAT SEARCH V(*),LOC(>Vmax);L,K
18470         IF L<Nsam THEN Valid(L)=0
18480         K=L
18490     NEXT K
18500     MAT U= U . Valid
18510     MAT V= V . Valid
18520     MAT W= W . Valid
18530     MAT A= A . Valid
18540     MAT B= B . Valid
18550     MAT I= I . Valid
18560     MAT C= C . Valid
18570 SUBEND
18580 Data_aconvert: SUB Data_aconvert(Gain)
18590     DISP "Converting Data"
18600     OPTION BASE 1
18610     COM /Data/ INTEGER Raw(*),Valid(*),REAL Table(*),U(*),V(*),W(*),A(*),B(*),I(*),C(*)
18620     N=SIZE(Raw,1)
18630     DIM Mv(1000),Mvn(1000),Amvn(1000),Sum(1000)
18640     REDIM Mv(N),Mvn(N),Amvn(1000),Sum(N)
18650     ! A0, A1, A2, A3, A4, A5, A6, A7
18660     DATA 150,257.10163,-28.16138,6.064559,-.792687,.05708673,-.002103462,.00003110036
18670     MAT Mv= A*(1000/Gain) ! Tt_mv=Tt_raw/Gain*1000
18680     MAT Sum= (0)
18690     MAT Mvn= (1)
18700     FOR K=0 TO 7
18710         READ An
18720         MAT Amvn= (An)*Mvn
18730         MAT Sum= Sum+Amvn
18740         MAT Mvn= Mvn . Mv
18750     NEXT K
18760     MAT A= Sum+(460)
18770 SUBEND
18780 Data_fconvert: SUB Data_fconvert(Array(*))
18790     OPTION BASE 1
18800     COM /Data/ INTEGER Raw(*),Valid(*),REAL Table(*),U(*),V(*),W(*),A(*),B(*),I(*),C(*)
18810     DIM Frng_spc(3),Brq_frq(3),Mix_frq(3),Mea_sgn(3),Brq_sgn(3),Mix_sgn(3)
18820     DISP "Converting Data"
18830     MAT Frng_spc= Array(25,1:3)
18840     MAT Brq_frq= Array(26,1:3)
18850     MAT Mix_frq= Array(27,1:3)
18860     MAT Mea_sgn= Array(28,1:3)
18870     MAT Brq_sgn= Array(29,1:3)
18880     MAT Mix_sgn= Array(30,1:3)
18890     MAT U= U*(Mea_sgn(1))
18900     MAT V= V*(Mea_sgn(2))
18910     MAT W= W*(Mea_sgn(3))
18920     MAT U= U+(Brq_sgn(1)*Brq_frq(1)+Mix_sgn(1)*Mix_frq(1))
18930     MAT V= V+(Brq_sgn(2)*Brq_frq(2)+Mix_sgn(2)*Mix_frq(2))
18940     MAT W= W+(Brq_sgn(3)*Brq_frq(3)+Mix_sgn(3)*Mix_frq(3))
18950     MAT U= U*(Frng_spc(1))
18960     MAT V= V*(Frng_spc(2))
18970     MAT W= W*(Frng_spc(3))
18980     MAT W= (0)
18990 SUBEND

```

```

19000 Data_sum: SUB Data_sum(Sum(*),INTEGER N(*),Nsam)
19010 OPTION BASE 1
19020 COM /Data/ INTEGER Raw(*),Valid(*),REAL Table(*),U(*),V(*),W(*),A(*),B(*),I(*),C(*)
19030 REAL Uu(1000),Vv(1000),Ww(1000),Aa(1000),Bb(1000),Ii(1000),Cc(1000)
19040 REAL Uv(1000),Vw(1000),Wu(1000),Ab(1000),Ua(1000),Va(1000),Wa(1000)
19050 REDIM Uu(Nsam),Vv(Nsam),Ww(Nsam),Aa(Nsam),Bb(Nsam),Ii(Nsam),Cc(Nsam)
19060 REDIM Uv(Nsam),Vw(Nsam),Wu(Nsam),Ab(Nsam),Ua(1000),Va(1000),Wa(1000)
19070 DISP "Summing Data"
19080 !
19090 MAT Uu= U . U
19100 MAT Vv= V . V
19110 MAT Ww= W . W
19120 MAT Aa= A . A
19130 MAT Bb= B . B
19140 MAT Uv= U . V
19150 MAT Vw= V . W
19160 MAT Wu= W . U
19170 MAT Ab= A . B
19180 MAT Ua= U . A
19190 MAT Va= V . A
19200 MAT Wa= W . A
19210 MAT Ii= I . I
19220 MAT Cc= C . C
19230 !
19240 Sum(1,1)=SUM(U)
19250 Sum(2,1)=SUM(V)
19260 Sum(3,1)=SUM(W)
19270 Sum(4,1)=SUM(A)
19280 Sum(5,1)=SUM(B)
19290 Sum(6,1)=SUM(I)
19300 Sum(7,1)=SUM(C)
19310 Sum(1,2)=SUM(Uu)
19320 Sum(2,2)=SUM(Vv)
19330 Sum(3,2)=SUM(Ww)
19340 Sum(4,2)=SUM(Aa)
19350 Sum(5,2)=SUM(Bb)
19360 Sum(6,2)=SUM(Ii)
19370 Sum(7,2)=SUM(Cc)
19380 Sum(1,3)=SUM(Uv)
19390 Sum(2,3)=SUM(Vw)
19400 Sum(3,3)=SUM(Wu)
19410 Sum(4,3)=SUM(Ab)
19420 Sum(5,3)=SUM(Ua)
19430 Sum(6,3)=SUM(Va)
19440 Sum(7,3)=SUM(Wa)
19450 MAT N= (SUM(Valid))
19460 N(3,1)=0
19470 N(5,1)=0
19480 N(3,2)=0
19490 N(5,2)=0
19500 N(2,3)=0
19510 N(3,3)=0
19520 N(4,3)=0
19530 N(6,3)=0
19540 N(7,3)=0
19550 SUBEND
19560 Data_calc: SUB Data_calc(INTEGER N(*),REAL
Sum(*),U,V,W,A,B,I,C,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,A1b1,U1a1,V1a1,W1a1)
DISP "Calculating Results"
!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
U=0
V=0
W=0
A=0
B=0
I=0
C=0
IF N(1,1) THEN U=Sum(1,1)/N(1,1)
IF N(2,1) THEN V=Sum(2,1)/N(2,1)
IF N(3,1) THEN W=Sum(3,1)/N(3,1)
IF N(4,1) THEN A=Sum(4,1)/N(4,1)
IF N(5,1) THEN B=Sum(5,1)/N(5,1)
IF N(6,1) THEN I=Sum(6,1)/N(6,1)
IF N(7,1) THEN C=Sum(7,1)/N(7,1)
!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
U1=0
V1=0
W1=0
A1=0
B1=0

```

```

19790      I1=0
19800      C1=0
19810      IF N(1,2) THEN U1=SQR(ABS(Sum(1,2)/N(1,2)-U*U))
19820      IF N(2,2) THEN V1=SQR(ABS(Sum(2,2)/N(2,2)-V*V))
19830      IF N(3,2) THEN W1=SQR(ABS(Sum(3,2)/N(3,2)-W*W))
19840      IF N(4,2) THEN A1=SQR(ABS(Sum(4,2)/N(4,2)-A*A))
19850      IF N(5,2) THEN B1=SQR(ABS(Sum(5,2)/N(5,2)-B*B))
19860      IF N(6,2) THEN I1=SQR(ABS(Sum(6,2)/N(6,2)-I*I))
19870      IF N(7,2) THEN C1=SQR(ABS(Sum(7,2)/N(7,2)-C*C))
19880      !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
19890      U1v1=0
19900      V1w1=0
19910      W1u1=0
19920      A1b1=0
19930      U1a1=0
19940      V1a1=0
19950      W1a1=0
19960      IF N(1,3) THEN U1v1=Sum(1,3)/N(1,3)-U*V
19970      IF N(2,3) THEN V1w1=Sum(2,3)/N(2,3)-V*W
19980      IF N(3,3) THEN W1u1=Sum(3,3)/N(3,3)-W*U
19990      IF N(4,3) THEN A1b1=Sum(4,3)/N(4,3)-A*B
20000      IF N(5,3) THEN U1a1=Sum(5,3)/N(5,3)-U*A
20010      IF N(6,3) THEN V1a1=Sum(6,3)/N(6,3)-V*A
20020      IF N(7,3) THEN W1a1=Sum(7,3)/N(7,3)-W*A
20030      !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
20040      SUBEND
20050 Data_trnsfrm: SUB Data_trnsfrm(REAL K(*),U,V,W,U1,V1,W1,U1v1,V1w1,W1u1)
20060      OPTION BASE 1
20070      REAL Vabc(3),Vuvw(3),Kabc(3,3),First
20080      REAL Ku(3,3),Kv(3,3),Kw(3,3),Ktu(3,3),Ktv(3,3),Ktw(3,3)
20090      REAL Kuu(3,3),Kvv(3,3),Kww(3,3),Kuv(3,3),Kvw(3,3),Kwu(3,3)
20100      REAL Kulul(3,3),Kvlv1(3,3),Kwlw1(3,3),Kulv1(3,3),Kvlw1(3,3),Kwlul(3,3)
20110      DISP "Transforming Results"
20120      Vabc(1)=U
20130      Vabc(2)=V
20140      Vabc(3)=W
20150      Kabc(1,1)=U1*U1
20160      Kabc(1,2)=U1v1
20170      Kabc(1,3)=W1u1
20180      Kabc(2,1)=U1v1
20190      Kabc(2,2)=V1*V1
20200      Kabc(2,3)=V1w1
20210      Kabc(3,1)=W1u1
20220      Kabc(3,2)=V1w1
20230      Kabc(3,3)=W1*W1
20240      MAT Vuvw= K*Vabc
20250      !OUTPUT PRT USING "6A,/,3(3(5DZ.5D),/),/";"K      ",K(*)
20260      U=Vuvw(1)
20270      V=Vuvw(2)
20280      W=Vuvw(3)
20290      FOR I=1 TO 3
20300          FOR J=1 TO 3
20310              Ku(I,J)=K(1,I)
20320              Kv(I,J)=K(2,I)
20330              Kw(I,J)=K(3,I)
20340          NEXT J
20350      NEXT I
20360      !OUTPUT PRT USING "6A,/,3(3(5DZ.5D),/),/";"Ku      ",Ku(*)
20370      !OUTPUT PRT USING "6A,/,3(3(5DZ.5D),/),/";"Kv      ",Kv(*)
20380      !OUTPUT PRT USING "6A,/,3(3(5DZ.5D),/),/";"Kw      ",Kw(*)
20390      MAT Ktu= TRN(Ku)
20400      MAT Ktv= TRN(Kv)
20410      MAT Ktw= TRN(Kw)
20420      !OUTPUT PRT USING "6A,/,3(3(5DZ.5D),/),/";"Ktu      ",Ktu(*)
20430      !OUTPUT PRT USING "6A,/,3(3(5DZ.5D),/),/";"Ktv      ",Ktv(*)
20440      !OUTPUT PRT USING "6A,/,3(3(5DZ.5D),/),/";"Ktw      ",Ktw(*)
20450      MAT Kuu= Ku . Ktu
20460      MAT Kvv= Kv . Ktv
20470      MAT Kww= Kw . Ktw
20480      MAT Kuv= Ku . Ktv
20490      MAT Kvw= Kv . Ktw
20500      MAT Kwu= Kw . Ktu
20510      !OUTPUT PRT USING "6A,/,3(3(5DZ.5D),/),/";"Kuu      ",Kuu(*)
20520      !OUTPUT PRT USING "6A,/,3(3(5DZ.5D),/),/";"Kvv      ",Kvv(*)
20530      !OUTPUT PRT USING "6A,/,3(3(5DZ.5D),/),/";"Kww      ",Kww(*)
20540      !OUTPUT PRT USING "6A,/,3(3(5DZ.5D),/),/";"Kuv      ",Kuv(*)
20550      !OUTPUT PRT USING "6A,/,3(3(5DZ.5D),/),/";"Kvw      ",Kvw(*)
20560      !OUTPUT PRT USING "6A,/,3(3(5DZ.5D),/),/";"Kwu      ",Kwu(*)
20570      MAT Kulul= Kuu . Kabc
20580      MAT Kvlv1= Kvv . Kabc

```

```

20590      MAT Kw1w1= Kww . Kabc
20600      MAT Kulv1= Kuv . Kabc
20610      MAT Kv1w1= Kvw . Kabc
20620      MAT Kw1u1= Kwu . Kabc
20630      U1u1=SUM(Kulul)
20640      V1v1=SUM(Kvlv1)
20650      W1w1=SUM(Kwlw1)
20660      U1v1=SUM(Kulv1)
20670      V1w1=SUM(Kvlw1)
20680      W1u1=SUM(Kwlul)
20690      U1=SQR(ABS(U1u1))
20700      V1=SQR(ABS(V1v1))
20710      W1=SQR(ABS(W1w1))
20720      SUBEND
20730 Data_print: SUB Data_print(Axis,Pos,INTEGER Nsam,CS,REAL
                U,V,W,A,B,I,C,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,A1b1,U1a1,V1a1,W1a1)
                IF CS="LDV" OR CS="TUN" THEN SUBEXIT
                DISP "Printing Results"
                ON ERROR CALL Error
                PRINTER IS PRT;WIDTH 144
                PRINT CHR$(27)&"&k2S"&CHR$(27)&"&l9D";
                LS=CHR$(NUM("X")+Axis-1)
                SELECT CS
                CASE "MHz","MOD"
                    PRINT USING 20920;LS,Pos,U,U1,U1v1
                    PRINT USING 20950;A,A1,A1b1,U1a1
                    PRINT USING 20930;"N",Nsam,V,V1,V1w1
                    PRINT USING 20960;B,B1,I1,V1a1
                    PRINT USING 20940;CS[1,3],W,W1,W1u1
                    PRINT USING 20970;C,I,C1,W1a1
                    IF CS<>"MOD" THEN PRINT
                END SELECT
                PRINTER IS CRT
                OFF ERROR
                IMAGE $,8X, A,"=",3D.4D,"   U=",5D.3D,"   U1=",5D.3D,"   U1V1=",8D.2D
                IMAGE $,8X, A,"=", 8D,"   V=",5D.3D,"   V1=",5D.3D,"   V1W1=",8D.2D
                IMAGE $,8X,3A,    7X,"   W=",5D.3D,"   W1=",5D.3D,"   W1U1=",8D.2D
                IMAGE     " A =" ,5D.3D,"   A1 =" ,5D.3D,"   A1B1=" ,6D.2D,"   U1A1=" ,7D.2D
                IMAGE     " B =" ,5D.3D,"   B1 =" ,5D.3D,"   IAT1=" ,6D.2D,"   V1A1=" ,7D.2D
                IMAGE     " CT=" ,5D.3D,"   IAT=" ,5D.3D,"   CT1 =" ,6D.2D,"   W1A1=" ,7D.2D
20980      SUBEND
20990 Data_plot: SUB Data_plot(Array(*),Symbols(*),G,Y,P1,P2,P3,Scale,INTEGER N1,N2,N3)
                OPTION BASE 1
                DIM Wndw(4),Vwprt(4),Symbol(20,3)
                DISP "Plotting Results"
                AREA PEN -1
                PEN 1
                MAT Wndw= Array(40+G,*)
                MAT Vwprt= Array(50+G,*)
                VIEWPORT Vwprt(1)/10.23,Vwprt(2)/10.23,Vwprt(3)/10.23,Vwprt(4)/10.23
                WINDOW Wndw(1),Wndw(2),Wndw(3),Wndw(4)
                CLIP ON
                FOR I=0 TO 2
                    IF I=0 AND N1=0 THEN 21300
                    IF I=1 AND N2=0 THEN 21300
                    IF I=2 AND N3=0 THEN 21300
                    Sy=I+1
                    Noc=Symbols(Sy,0,1)
                    REDIM Symbol(Noc,3)
                    MAT Symbol= Symbols(Sy,1:Noc,*)
                    SELECT I
                    CASE 0
                        X=P1*Scale
                    CASE 1
                        X=P2*Scale
                    CASE 2
                        X=P3*Scale
                    END SELECT
                    Xm=MIN(MAX(X,Wndw(1)),Wndw(2))
                    Ym=MIN(MAX(Y,Wndw(3)),Wndw(4))
                    MOVE Xm,Ym
                    SYMBOL Symbol(*),FILL,EDGE
                NEXT I
                SUBEND
21320 Histo: !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
21330 Rt_histo: SUB Rt_histo(@Lvdas,Symbols(*),Repeat)
                OPTION BASE 1
                COM /Graph/
                Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),Legend$(*)
                INTEGER Histo(1000,3),Nplots,Nbins,F1,F2,A1,A2

```

```

21370 REAL Nnew,Nold,N(5)
21380 OUTPUT @Lvdas USING "AA";"CA"
21390 FOR Channel=1 TO 5
21400     CALL Set_up(Channel,Symbols(*))
21410 NEXT Channel
21420 CALL Convert2words(.1*1000000,A1,A2) ! Atime=.1 seconds
21430 ON KBD GOSUB Hdone
21440 REPEAT
21450     FOR Channel=1 TO 5
21460         G=Channel
21470         SELECT Channel
21480             CASE 1,2
21490                 Min=0
21500                 Bin=20
21510                 Ww=2^Bin
21520                 Kw=1000000
21530                 CALL Convert2words(Min,F1,F2)
21540             CASE 4
21550                 Min=-5
21560                 Bin=10
21570                 Ww=2^Bin
21580                 F1=-1
21590                 F2=-32768
21600                 Kw=32768/5
21610             CASE ELSE
21620                 GOTO 21880
21630             END SELECT
21640 Hsend:   OUTPUT @Lvdas USING "AA,6(W)";"TH",F1,F2,Bin,A1,A2,Channel
21650 Henter:  ENTER @Lvdas USING "#,W";Nbins
21660         IF Nbins>0 THEN
21670             REDIM Histo(Nbins,3)
21680             ENTER @Lvdas USING "#,W";Histo(*)
21690         END IF
21700         ENTER @Lvdas USING "#,W";Nnew,Nold
21710 Hplot:    VIEWPORT Vwprr(G,1)/10.23,Vwprr(G,2)/10.23,Vwprr(G,3)/10.23,Vwprr(G,4)/10.23
21720         WINDOW Kw*Wndw(G,1),Kw*Wndw(G,2),Wndw(G,3),Wndw(G,4)
21730         Xpixel=Kw*(Wndw(Channel,2)-Wndw(Channel,1))/(Vwprr(Channel,2)-Vwprr(Channel,1))
21740         N1=N(Channel)
21750         N2=N(Channel)-Nold+Nnew
21760         N(Channel)=N(Channel)-Nold+Nnew
21770         FOR I=1 TO Nbins
21780             Old=MIN(Histo(I,3),Wndw(Channel,4))
21790             New=MIN(Histo(I,2),Wndw(Channel,4))
21800             AREA PEN SGN(New-Old)
21810             X1=Histo(I,1)*Ww+Min*Kw
21820             X2=Ww
21830             Y1=Old
21840             Y2=New-Old
21850             MOVE X1,Y1
21860             RECTANGLE X2-Xpixel,Y2,FILL
21870         NEXT I
21880     NEXT Channel
21890 UNTIL KBDS<>" OR NOT Repeat
21900 SUBEXIT
21910 Hdone:   Done=1
21920     RETURN
21930 SUBEND
21940 Histo:   !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
21950 Pt_histo: SUB Pt_histo(Symbols(*),Run,File,Pos,INTEGER Nsam)
21960     OPTION BASE 1
21970     COM /Data/ INTEGER Raw(*),Valid(*),REAL Table(*),U(*),V(*),W(*),A(*),B(*),I(*),C(*)
21980     COM /Graph/
21990     Wndw(*),Vwprr(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),Legend$(*)
22000     INTEGER Histo(0:100)
22010     REAL Data(1000)
22020     REDIM Data(Nsam)
22030     FOR Channel=5 TO 1 STEP -1
22040         G=Channel
22050         IF Channel=1 THEN MAT Data= U
22060         IF Channel=2 THEN MAT Data= V
22070         IF Channel=4 THEN MAT Data= A
22080         SELECT Channel
22090             CASE 1,2,4
22100 Hsort:    CALL Set_up(Channel,Symbols(*))
22110             Xmin=Wndw(Channel,1)
22120             Xmax=Wndw(Channel,2)
22130             Xwin=(Xmax-Xmin)/100
22140             MAT Data= Data-(Xmin)
22150             MAT Data= Data/((Xmax-Xmin)/100)
22160             MAT Histo= (0)

```

```

22160         FOR K=1 TO Nsam
22170             L=MAX(MIN(Data(K),100),0)
22180             Histo(L)=Histo(L)+1
22190         NEXT K
22200 Hplot:    VIEWPORT Vwprt(G,1)/10.23,Vwprt(G,2)/10.23,Vwprt(G,3)/10.23,Vwprt(G,4)/10.23
22210         WINDOW 0,100,Wndw(G,3),Wndw(G,4)
22220         Xpixel=(100-0)/(Vwprt(Channel,2)-Vwprt(Channel,1))
22230         MOVE 55,70
22240         IF G=2 THEN LABEL USING "2A,2D.2D";"R=",Run
22250         IF G=1 THEN LABEL USING "2A,2D.3D";"Y=",Pos
22260         IF G=4 THEN LABEL USING "2A,2D "; "F=",File
22270         FOR K=0 TO 100
22280             IF Histo(K) THEN
22290                 MOVE K-.5,0
22300                 AREA PEN SGN(1)
22310                 RECTANGLE 1-Xpixel,Histo(K),FILL
22320             END IF
22330         NEXT K
22340         END SELECT
22350     NEXT Channel
22360     SUBEXIT
22370 SUBEND
22380 Misc:    !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
22390 Convert2words: SUB Convert2words(Real,INTEGER High,Low)
22400             Hex$=DVAL$(Real,16)
22410             High=IVAL(Hex$[1,4],16)
22420             Low=IVAL(Hex$[5,8],16)
22430 SUBEND
22440 Temp:    SUB Temp(Mach,Mv,Ts,Tt)
22450             !      A0,      A1,      A2,      A3,      A4,      A5,      A6,      A7
22460             DATA 150,257.10163,-28.16138,6.064559,-.792687,.05708673,-.002103462,.00003110036
22470             Tt=0
22480             FOR I=0 TO 7
22490                 READ K
22500                 Tt=Tt+K*Mv^I
22510             NEXT I
22520             Tt=Tt+460
22530             Ts=.09259*Tt
22540             IF Mach<>7 THEN BEEP
22550             IF Mach<>7 THEN PAUSE
22560 SUBEND
22570 Error:  SUB Error
22580             BEEP
22590             DISP ERRMS
22600             OUTPUT PRT;ERRMS
22610             Prt=VAL(SYSTEM$( "PRINTER IS" ))
22620             PRINTER IS CRT
22630             PRINT TABXY(95,1);ERRMS
22640             PRINTER IS Prt
22650             ERROR SUBEXIT
22660 SUBEND
22670 Fix:    SUB Fix(Array(*),Name$(*),Image$(*),Units$(*))
22680             OPTION BASE 1
22690             Run=Array(3,1)
22700             SELECT INT(Run)
22710             CASE ELSE
22720                 Image$(3,1)="6D.2D"      ! Run
22730                 Array(3,1)=INT(Run)+.01  ! Run
22740                 Array(21,1)=.3125        ! UBeamSpc
22750                 Array(21,2)=.34375       ! VBeamSpc
22760                 Array(21,3)=.3125       ! WBeamSpc
22770                 Array(22,1)=30          ! UFoclLen
22780                 Array(22,2)=30          ! VFoclLen
22790                 Array(22,3)=30          ! WFoclLen
22800                 Array(23,1)=2*ATN(Array(21,1)/2/Array(22,1)) ! UFrngSpc
22810                 Array(23,2)=2*ATN(Array(21,2)/2/Array(22,2)) ! VFrngSpc
22820                 Array(23,3)=2*ATN(Array(21,3)/2/Array(22,3)) ! WFrngSpc
22830                 Array(28,1)=-1          ! UMeaSgn
22840                 Array(29,1)=1           ! UBrSgn
22850                 Array(30,1)=1           ! UMixSgn
22860                 Array(41,1)=0           ! Xmin1
22870                 Array(41,2)=100         ! Xmax1
22880                 Array(42,1)=0           ! Xmin2
22890                 Array(42,2)=100         ! Xmax2
22900                 Array(43,1)=0           ! Xmin3
22910                 Array(43,2)=100         ! Xmax3
22920                 Array(44,1)=-5          ! Xmin4
22930                 Array(44,2)=5          ! Xmax4
22940                 Array(45,1)=-5          ! Xmin5
22950                 Array(45,2)=5          ! Xmax5

```

```

22960      Array(61,1)=5          ! Xdiv1
22970      Array(62,1)=5          ! Xdiv2
22980      Array(63,1)=5          ! Xdiv3
22990      Array(64,1)=4          ! Xdiv4
23000      Array(65,1)=4          ! Xdiv5
23010      !
23020      Array(46,4)=4          ! Ymax6
23030      Array(47,4)=4          ! Ymax7
23040      Array(48,4)=4          ! Ymax8
23050      Array(49,4)=4          ! Ymax9
23060      Array(49,1)=0          ! Xmin9
23070      Array(49,2)=2000       ! Xmax9
23080      Array(61,4)=8          ! Ydiv6
23090      Array(62,4)=8          ! Ydiv7
23100      Array(63,4)=8          ! Ydiv8
23110      Array(64,4)=8          ! Ydiv9
23120      Array(64,3)=4          ! Xdiv9
23130      !
23140      Array(35,1)=8          ! UFreqMin
23150      Array(36,1)=40         ! UFreqMax
23160      Array(35,2)=20         ! VFreqMin
23170      Array(36,2)=55         ! VFreqMax
23180      Array(35,3)=10         ! WFreqMin
23190      Array(36,3)=70         ! WFreqMax
23200      Array(36,4)=1          ! Clip
23210      !
23220      Name$(35,1)="UFreqMin" ! UFreqMin
23230      Name$(36,1)="UFreqMax" ! UFreqMax
23240      Name$(35,2)="VFreqMin" ! VFreqMin
23250      Name$(36,2)="VFreqMax" ! VFreqMax
23260      Name$(35,3)="WFreqMin" ! WFreqMin
23270      Name$(36,3)="WFreqMax" ! WFreqMax
23280      Name$(36,4)="Clip"     ! Clip
23290      !
23300      Units$(35,1)="MHz"     ! UFreqMin
23310      Units$(36,1)="MHz"     ! UFreqMax
23320      Units$(35,2)="MHz"     ! VFreqMin
23330      Units$(36,2)="MHz"     ! VFreqMax
23340      Units$(35,3)="MHz"     ! WFreqMin
23350      Units$(36,3)="MHz"     ! WFreqMax
23360      Units$(36,4)="        ! Clip
23370      !
23380      Image$(35,1)="4D.4D"   ! UFreqMin
23390      Image$(36,1)="4D.4D"   ! UFreqMax
23400      Image$(35,2)="4D.4D"   ! VFreqMin
23410      Image$(36,2)="4D.4D"   ! VFreqMax
23420      Image$(35,3)="4D.4D"   ! WFreqMin
23430      Image$(36,3)="4D.4D"   ! WFreqMax
23440      Image$(36,4)="9D"       ! Clip
23450      !
23460      END SELECT
23470      SUBEND
23480      Scale:
23490      SUB Scale(G)
23500      OPTION BASE 1
23510      COM /Graph/
23520      Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),Legend$(*)
23530      VIEWPORT Vwprt(G,1)/10.23,Vwprt(G,2)/10.23,Vwprt(G,3)/10.23,Vwprt(G,4)/10.23
23540      WINDOW Wndw(G,1),Wndw(G,2),Wndw(G,3),Wndw(G,4)
23550      SUBEND
23560      Purge:
23570      SUB Purge(System$,Data$)
23580      OPTION BASE 1
23590      DIM FS(400)(80)
23600      MASS STORAGE IS Data$
23610      CAT TO FS(*) ; NAMES
23620      MAT SEARCH FS(*),LOC(="");N
23630      N=N-1
23640      IF N>0 THEN
23650          REDIM FS(N)
23660          FOR I=1 TO N
23670              IF FS(I)(1,4)="R.01" THEN
23680                  PURGE FS(I)&"<TKM>"
23690                  DISP FS(I)&"<TKM>"
23700              END IF
23710          NEXT I
23720      END IF
23730      END IF
23740      MASS STORAGE IS System$
23750      SUBEND

```

APPENDIX B

REVISED SOFTWARE CODE LISTING.

APPENDIX B

Revised Software Code Listing.

TABLE OF CONTENTS

| Item | Title | Page |
|-------------|---|-------------|
| 1 | Revised Hard Disk Directory Catalog Listing. | 2 |
| 2 | Revised Program "3.5'HWT92" Hard Copy Code Listing. | 3 |

Hard Disk Directory Catalog Listing.

:CS80, 1400, 0, 0

VOLUME LABEL: B9826

| FILE NAME | PRO | TYPE | REC/FILE | BYTE/REC | ADDRESS | DATE | TIME |
|-----------|-----|-------|----------|----------|---------|-----------|-------|
| SYSB60 | | SYSTM | 3388 | 256 | 32 | 17-Jul-91 | 13:10 |
| CDUMP6 | | PROG | 44 | 256 | 3420 | 17-Jul-91 | 13:10 |
| BPLOT6 | | PROG | 40 | 256 | 3464 | 17-Jul-91 | 13:10 |
| AUTOST | | PROG | 10 | 256 | 3504 | 17-Jul-91 | 13:10 |
| ARRAY | | BDAT | 50 | 256 | 3515 | 17-Jul-91 | 13:11 |
| KEYS | | BDAT | 4 | 256 | 3566 | 17-Jul-91 | 13:11 |
| COPY | | PROG | 25 | 256 | 3570 | 17-Jul-91 | 13:11 |
| 3.5'HWT92 | | PROG | 816 | 256 | 3967 | 29-Mar-92 | 13:45 |

```

100 Main: ----- 3.5'HWT92
105
110 NASA AMES RESEARCH CENTER Property of COMPLERE INC.
115 3.5 FOOT HYPERSONIC WIND TUNNEL Proprietary software
120 Laser Doppler Velocimeter Test Copyright March 29, 1992
125 Developed by: T. Kevin McDevitt
130 -----
135
140
145 PROGRAM DESCRIPTION:
150
155 This program provides the capability to acquire simultaneous Laser Doppler Velocimeter (LDV), Stagnation
160 Temperature, and analog voltage data at user selectable traverse controlled probe volume positions within
165 the hypersonic wind tunnel flow.
170 The Laser Velocimeter Data Acquisition System (LVDAS) is used to sample the LDV and analog voltage
175 data simultaneously with a coincidence criterion being applied to LDV incoming data. The LVDAS also generates
180 inter-arrival times and coincidence times.
185 The measured LDV data provide the necessary frequency information from which two components of flow
190 velocities can be determined. These velocities are measured directly in "TUNNEL" coordinates. A coordinate
195 system transformation is applied to these measured velocities to obtain velocities in and "MODEL" coordinates
200 if the model is at angles of attack, yaw, and/or roll.
205 The Traverse Control System (TCS8) is used to precisely move the LDV probe volume within the tunnel and
210 about the model. The TCS8 provides three axes, plus one auxiliary axis, of traverse capability for both the
215 transmitting (Tx) and receiving (Rx) side optical packages. The Tx and Rx side traverses can be moved
220 independently to achieve laser alignment or they can be moved together to maintain laser alignment.
225 The TCS8 will give the traverse positions in "TUNNEL" coordinates where one inch of commanded movement will
230 yield one inch of movement on the traverse slides. This will also yield one inch of movement of the probe
235 volume crossover point within the tunnel. However, the traverse positions in "TUNNEL" coordinates will differ
240 from positions in "MODEL" coordinates if the model is at angles of attack, yaw, and/or roll. Therefore, a
245 coordinate system transformation is applied to TCS8 positions to obtain positions in "MODEL" coordinates.
250 During data acquisition, real time histograms will be displayed of the LDV and analog data. After the
255 data have been acquired, the averages, standard deviations, and shear stresses will be calculated and displayed
260 in profile plots where the data are plotted versus traverse position. The reduced data are also sent to the
265 printer in tabular form. The reduced data as well as the raw data are stored along with the tunnel conditions
270 on the hard disc for archival purposes and also to allow for further data reduction, data plotting, and/or data
275 transfer to other computers.
280
285 PROGRAM OPERATION:
290
295 The following power up sequences should be completed before this program is run:
300 1. Turn on the "MDS" Motor Drive System boxes.
305 2. Turn on the "TCS8" Traverse Control System.
310 3. Turn on the "LVDAS" Laser Velocimeter Data Acquisition System.
315 4. Turn on the HP series 9000 model 375 computer.
320 This program will automatically be loaded and executed when the computer is turned on. If it is not loaded,
325 then you can type in the following commands to load and then execute it.
330 LOAD "3.5'HWT:,1400,0,0"
335 RUN
340 When the program is ready for user operation, it will display three things on the CRT. These are the main
345 menu, TCS8 traverse positions, and new sets of histogram and profile graphs. If they do not appear on the CRT
350 then the following actions should be performed to reinitialize the systems.
355 1. Press shift reset on the HP series 9000 model 375 computer's keyboard.
360 2. Press reset on the back of the TCS8.
365 3. Press reset on the front (or back) of the LVDAS.
370 4. LOAD "3.5'HWT:,1400,0,0"
375 5. RUN
380
385 PROGRAM VARIABLES:
390
395 Mass Storage Variables:
400
405 System$ Tells the program where to read/store system data related files.
410 Data$ Tells the program where to read/store raw and reduced data related files.
415 File$ File name for tunnel conditions data or raw and reduced data.
420
425 Menu Variables:
430
435 Menu$(*) String array where each element describes its corresponding menu subroutine's function.
440 Menu Used as an index to the string array Menu$(*). Indicates which of the menus has been
445 selected as the current menu.
450 Key Used as an index to the string array Menu$(*). Indicates which one of eight menu
455 subroutines in the menu is to be executed.
460 Busy Tells the Menu Status subprogram to display the current menu selection in red text.
465 Ready Tells the Menu Status subprogram to display the current menu selection in blue text.
470
475 Traverse Position Variables:
480
485 Tun1(*) TCS8 transmitting side traverse positions (X1,Y1,Z1,A1) in "TUNNEL" coordinates.
490 Tun2(*) TCS8 receiving side traverse positions (X2,Y2,Z2,A2) in "TUNNEL" coordinates.
495 Mod1(*) TCS8 transmitting side traverse positions (X1,Y1,Z1,A1) in "MODEL" coordinates.

```

```

500      !      Mod2(*)      TCS8 transmitting side traverse positions (X2,Y2,Z2,A2) in "MODEL" coordinates.
505      !      Side$       Indicates which sides are to be moved:
510      !                   Tx       : Transmitting side only.
515      !                   Rx       : Receiving side only.
520      !                   Tx & Rx  : Both sides together.
525      !      Coord$      Indicates which coordinate system the movement is to be made in:
530      !                   TUNNEL   : TUNNEL coordinates.
535      !                   MODEL    : MODEL coordinates.
540      !      Mode$       Indicates which movement mode is to be completed:
545      !                   RELATIVE: Movements are relative to current positions.
550      !                   ABSOLUTE: Movements are to absolute positions.
555      !      Movement    Indicates the desired movement for the selected axis.
560      !      Paxis       Specifies which axis is to be traversed for the profile. Also defines axis for plots.
565      !
570      ! Auto Move Traverse Position Variables:
575      !
580      !      Pos(*)       Array of pre-programmed auto move positions.
585      !      Pname$(*)    Names for the variables in Pos(*).
590      !      Pimage$(*)   Image formats for the variables in Pos(*).
595      !      Punits$(*)   Units for the variables in Pos(*).
600      !      Npos        Number of pre-programmed auto move positions in Pos(*).
605      !      Paxis       Specifies which axis is to be traversed for the profile. Also defines axis for plots.
610      !
615      ! Traverse Positions and Velocity Coordinate System Transformation Variables:
620      !
625      !      Alpha(*)     Angles of attack, yaw, and roll.
630      !                   Alpha(1) : Angle of Attack.
635      !                   Alpha(2) : Angle of Yaw.
640      !                   Alpha(3) : Angle of Roll.
645      !      Mod2tun(*)   Coordinate system transformation matrix for converting positions & velocities from MODEL to TUNNEL.
650      !      Tun2mod(*)   Coordinate system transformation matrix for converting positions & velocities from TUNNEL to MODEL.
655      !
660      ! Tunnel Condition Variables:
665      !
670      !      Array(*)     Array of tunnel conditions, laser parameters, graph scales, etc.
675      !      Name$(*)    Names for the variables in Array(*).
680      !      Image$(*)   Image formats for the variables in Array(*).
685      !      Units$(*)   Units for the variables in Array(*).
690      !
695      ! Misc. Tunnel Condition Variables:
700      !
705      !      Date         Date.
710      !      Time        Time.
715      !      Run         Run Number.
720      !      File        File Number.
725      !      Mach        Mach Number.
730      !      Re_ft       Re/Ft (Reynolds Number per Foot).
735      !      Uedge       Freestream Velocity (m/s).
740      !      Uinf        Freestream Velocity (m/s).
745      !      Stemp        Stagnation Temperature (deg R).
750      !      Ttemp       Total Temperature (deg R).
755      !      Tt_mv        Total Temperature data in gained millivolts.
760      !      Tt_raw       Total Temperature raw data in ungained volts.
765      !      Gain        Gain for total temperature raw analog data in ungained volts to gained millivolts conversion.
770      !
775      ! LVDAS Variables:
780      !
785      !      Table(*)     Lookup table of frequencies.
790      !      Atime        The maximum desired acquisition time (seconds).
795      !      Ctime        The maximum desired coincidence time (seconds).
800      !      At_exp       Exponent for inter-arrival times.
805      !      Ct_exp       Exponent for coincidence times.
810      !      Nreads       Number of desired samples.
815      !      Nsam         Number of acquired samples.
820      !      Coin(*)      Coincidence criteria.
825      !      Cmask       Coincidence mask for U,V,W selection.
830      !      Raw(*)      Array of raw data acquired from the LVDAS.
835      !
840      ! Instantaneous Velocity and Voltage Variables:
845      !
850      !      Ui(*)        Read from LVDAS as the instantaneous U frequency data, then converted into U velocities.
855      !      Vi(*)        Read from LVDAS as the instantaneous V frequency data, then converted into V velocities.
860      !      Wi(*)        Read from LVDAS as the instantaneous W frequency data, then converted into W velocities.
865      !      Ai(*)        Read from LVDAS as the instantaneous A voltage data.
870      !      Bi(*)        Read from LVDAS as the instantaneous B voltage data.
875      !      Ii(*)        Read from LVDAS as the raw inter-arrival time data, then converted into inter-arrival times.
880      !      Ci(*)        Read from LVDAS as the raw coincidence time data, then converted into coincidence times.
885      !      Valid(*)     Validation words. Initially all ones, then some set to zero during histogram clipping.
890      !
895      ! Histogram Clipping Variables:

```

```

900      !
905      !   Umin      The minimum acceptable U frequency (MHz).
910      !   Umax      The maximum acceptable U frequency (MHz).
915      !   Vmin      The minimum acceptable V frequency (MHz).
920      !   Vmax      The maximum acceptable V frequency (MHz).
925      !   Wmin      The minimum acceptable W frequency (MHz).
930      !   Wmax      The maximum acceptable W frequency (MHz).
935      !   Clip      Clip: 1 turn histogram clipping on; 0 turns it off.
940      !
945      ! Frequency to Velocity Conversion Variables:
950      !
955      !   Beam_spc(*) Beam spacing at lens.
960      !   Focl_len(*) Focal length.
965      !   Beam_sep(*) Beam separation angle in degrees (full angle).
970      !   Wave_len(*) Wave length.
975      !   Frng_spc(*) Fringe Spacings.
980      !   Brq_frq(*) Bragg Frequencies.
985      !   Mix_frq(*) Mixing Frequencies.
990      !   Mea_sgn(*) Measured Frequencies' Signs.
995      !   Brq_sgn(*) Bragg Frequencies' Signs.
1000     !   Mix_sgn(*) Mixing Frequencies' Signs.
1005     !
1010     ! Summation Variables:
1015     !
1020     !   Sumu      Summation of all of the valid U1.
1025     !   Sumv      Summation of all of the valid V1.
1030     !   Sumw      Summation of all of the valid W1.
1035     !   Suma      Summation of all of the valid A1.
1040     !   Sumb      Summation of all of the valid B1.
1045     !   Sumi      Summation of all of the valid I1.
1050     !   Sumc      Summation of all of the valid C1.
1055     !   Sumuu     Summation of all of the valid U1*U1.
1060     !   Sumvv     Summation of all of the valid V1*V1.
1065     !   Sumww     Summation of all of the valid W1*W1.
1070     !   Sumaa     Summation of all of the valid A1*A1.
1075     !   Sumbb     Summation of all of the valid B1*B1.
1080     !   Sumii     Summation of all of the valid I1*I1.
1085     !   Sumcc     Summation of all of the valid C1*C1.
1090     !   Sumuv     Summation of all of the valid U1*V1.
1095     !   Sumvw     Summation of all of the valid V1*W1.
1100     !   Sumwu     Summation of all of the valid W1*U1.
1105     !   Sumab     Summation of all of the valid A1*B1.
1110     !   Sumua     Summation of all of the valid U1*A1.
1115     !   Sumva     Summation of all of the valid V1*A1.
1120     !   Sumwa     Summation of all of the valid W1*A1.
1125     !   Suml      Number of valid samples for the above summations.
1130     !
1135     ! Reduced Data Variables:
1140     !
1145     !   N          Number of valid samples acquired.
1150     !   U          Average U frequency or velocity.
1155     !   V          Average V frequency or velocity.
1160     !   W          Average W frequency or velocity.
1165     !   A          Average A voltage.
1170     !   B          Average B voltage.
1175     !   I          Average inter-arrival time.
1180     !   C          Average coincidence time.
1185     !   U1         Standard deviation for U frequency or velocity.
1190     !   V1         Standard deviation for V frequency or velocity.
1195     !   W1         Standard deviation for W frequency or velocity.
1200     !   A1         Standard deviation for A voltage.
1205     !   B1         Standard deviation for B voltage.
1210     !   I1         Standard deviation for inter-arrival time.
1215     !   C1         Standard deviation for coincidence time.
1220     !   Ulv1       Velocity:Velocity Shear Stress.
1225     !   Vlw1       Velocity:Velocity Shear Stress.
1230     !   Wlul       Velocity:Velocity Shear Stress.
1235     !   Alb1       Voltage :Voltage Cross Correlation.
1240     !   Ulal       Velocity:Voltage Cross Correlation.
1245     !   Vlal       Velocity:Voltage Cross Correlation.
1250     !   Wlal       Velocity:Voltage Cross Correlation.
1255     !
1260     ! Data Plotting Symbol Variables:
1265     !
1270     !   Symbols(*)   Array of Symbol arrays. Each symbol array contains a distinct geometric symbol.
1275     !   Symbol(*)    Array of coordinates which when connected produce a distinct geometric symbol.
1280     !   Dot(*)       Array of coordinates which produce a dot. The dot symbol is added to all symbols.
1285     !   Noc          The number of coordinates in a symbol.
1290     !   Sy           Used to index the Symbols array.
1295     !

```

```

1300 ! Histogram and Profile Graph Variables:
1305 !
1310 ! Wndw(*)      Array containing the plots' scales.
1315 ! Vwprt(*)     Array containing the plots' CRT positions.
1320 ! Xdiv(*)      Array containing the number of X divisions for the plot's X axis.
1325 ! Ydiv(*)      Array containing the number of Y divisions for the plot's Y axis.
1330 ! Xlabels(*)   String array containing labels for the X axis.
1335 ! Ylabels(*)   String array containing labels for the Y axis.
1340 ! Title$(*)    String array containing labels for the plots.
1345 ! Ximage$(*)   String array containing image formats for the X axis labeling.
1350 ! Yimage$(*)   String array containing image formats for the Y axis labeling.
1355 ! Legend$(*)   String array containing labels for each symbol in a profile plot.
1360 ! G           Used as an index to the above arrays. Specifies one of nine plots.
1365 ! Gsave(*)     Used to save the entire graphics contents of the CRT.
1370 !
1375 ! Dimension the variables and arrays defined above.
1380 OPTION BASE 1
1385 COM /Pos/ Pname$(25,1)[10],Pimage$(25,1)[10],Punits$(25,1)[10],REAL Pos(25,1),Npos
1390 COM /Array/ Name$(100,4)[10],Image$(100,4)[10],Units$(100,4)[10],REAL Array(100,4)
1395 COM /Data1/ REAL Table(0:32766),INTEGER Raw(1000,10),Valid(1000)
1400 COM /Data2/ REAL U1(1000),V1(1000),W1(1000),A1(1000),B1(1000),I1(1000),C1(1000)
1405 COM /Data3/ REAL Puu(1000),Pvv(1000),Pww(1000),Paa(1000),Pbb(1000),Pii(1000),Pcc(1000)
1410 COM /Data4/ REAL Puv(1000),Pvw(1000),Pwu(1000),Pab(1000),Pua(1000),Pwa(1000)
1415 COM /Graph1/ Wndw(9,4),Vwprt(9,4),Xdiv(9),Ydiv(9),Xlabel$(9)[80],Ylabel$(9)[80]
1420 COM /Graph2/ Title$(9)[80],Ximage$(9)[80],Yimage$(9)[80],Legend$(9,5)[80]
1425 COM /Color1/ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta
1430 COM /Color2/ White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
1435 COM /Sum1/ REAL Sumu,Sumv,Sumw,Suma,Sumb,Sumi,Sumc,Suml
1440 COM /Sum2/ REAL Sumuu,Sumvv,Sumww,Sumaa,Sumbb,Sumii,Sumcc
1445 COM /Sum3/ REAL Sumuv,Sumvw,Sumwu,Sumab,Sumua,Sumva,Sumwu
1450 COM /Reduced/ N,U,V,W,A,B,I,C,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,A1a1,U1a1,V1a1,W1a1
1455 COM Run,File,Paxis
1460 DIM Menu$(6,8)[80],Systems[20],Data[20],File[50],L$[160],Kbd$[160]
1465 INTEGER Gsave(1280,1024),At_exp,Ct_exp,Cmask,Nsam
1470 REAL Atime,Ctime,Symbols(5,0:20,3)
1475 DIM Tun2mod(3,3),Mod2tun(3,3),Tun1(4),Tun2(4),Mod1(4),Mod2(4),Alpha(3)
1480 DIM Beam_spc(3),Focl_len(3),Mea_sgn(3),Mix_frq(3),Mix_sgn(3),Frng_spc(3)
1485 DIM Beam_sep(3),Wave_len(3),Brg_frq(3),Brg_sgn(3),Coin(3)
1490 DEG ! Perform trigonometric operations in degrees.
1495 !
1500 ! Perform any necessary setup and initialization routines.
1505 CALL Crt_init ! Clear the CRT and direct printed output to it.
1510 GOSUB Lvds_set_up ! Initialize the HP to LVDAS interface.
1515 GOSUB File_set_up ! Select mass storage devices for system and data files.
1520 GOSUB Tcs8_set_up ! Initialize the HP to TCS8 interface.
1525 GOSUB Menu_set_up ! Initialize the user driven menus and display the main menu.
1530 GOSUB Grph_set_up ! Initialize the CRT and plot the nine empty plots for profiles and histograms.
1535 !
1540 Here: ! The main program, while continually displaying the time of day, will wait here for menu key selection.
1545 Date=TIMEDATE
1550 Time=Date
1555 PRINT PEN Blue
1560 DISP CHR$(129);" ";TIM$(TIMEDATE);" ";DATE$(TIMEDATE);" ";CHR$(128)
1565 GOTO Here
1570 STOP
1575 On_key: ON KEY 1 GOSUB Key1 ! If the user function key #1 is ever pressed then execute the "Key1" subroutine.
1580 ON KEY 2 GOSUB Key2 ! If the user function key #2 is ever pressed then execute the "Key2" subroutine.
1585 ON KEY 3 GOSUB Key3 ! If the user function key #3 is ever pressed then execute the "Key3" subroutine.
1590 ON KEY 4 GOSUB Key4 ! If the user function key #4 is ever pressed then execute the "Key4" subroutine.
1595 ON KEY 5 GOSUB Key5 ! If the user function key #5 is ever pressed then execute the "Key5" subroutine.
1600 ON KEY 6 GOSUB Key6 ! If the user function key #6 is ever pressed then execute the "Key6" subroutine.
1605 ON KEY 7 GOSUB Key7 ! If the user function key #7 is ever pressed then execute the "Key7" subroutine.
1610 ON KEY 8 GOSUB Key8 ! If the user function key #8 is ever pressed then execute the "Key8" subroutine.
1615 RETURN
1620 Keys: ! Subroutine Key1,Key2,Key3,Key4,Key5,Key6,Key7,Key8 descriptions:
1625 ! When one of the special user function keys is pressed, the main program will execute one the
1630 ! following eight subroutines. Each of these subroutines performs essentially the same basic
1635 ! function in that it subsequently executes one of the menu subroutines. The particular menu
1640 ! subroutine to be executed will depend on the current menu selected and the current key pressed.
1645 ! Before the selected menu subroutine is executed, the corresponding menu entry at the top of
1650 ! the CRT is redisplayed in red text. This indicates that the menu selection has been
1655 ! acknowledged and that any resultant actions are still in progress. When the highlighted menu
1660 ! subroutine has completed the current TCS8 traverse positions will be read and updated on the CRT
1665 ! display. The corresponding menu entry displayed at the top of the CRT is redisplayed in blue
1670 ! text to indicate the completion of the menu subroutine. The user can then select another special
1675 ! function key.
1680 ! Variables:
1685 ! Menu Indicates which of the menus has been selected as the current menu.
1690 ! Key Indicates which one of eight menu subroutines in the menu is to be executed.
1695 ! Menu$(*) String array where each element describes its corresponding menu subroutine's function.

```

```

1700      !      Busy      Tells the Menu Status subroutine to display the current menu selection in red text.
1705      !      Ready     Tells the Menu Status subroutine to display the current menu selection in blue text.
1710 Key1:      Key=1
1715      CALL Menu_status(Menu,Key,Busy,Menu$(*))
1720      ON Menu GOSUB M1k1,M2k1,M3k1,M4k1,M5k1,M6k1,M7k1
1725      CALL Menu_status(Menu,Key,Ready,Menu$(*))
1730      CALL Tcs8read(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*))
1735      RETURN
1740 Key2:      Key=2
1745      CALL Menu_status(Menu,Key,Busy,Menu$(*))
1750      ON Menu GOSUB M1k2,M2k2,M3k2,M4k2,M5k2,M6k2,M7k2
1755      CALL Menu_status(Menu,Key,Ready,Menu$(*))
1760      CALL Tcs8read(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*))
1765      RETURN
1770 Key3:      Key=3
1775      CALL Menu_status(Menu,Key,Busy,Menu$(*))
1780      ON Menu GOSUB M1k3,M2k3,M3k3,M4k3,M5k3,M6k3,M7k3
1785      CALL Menu_status(Menu,Key,Ready,Menu$(*))
1790      CALL Tcs8read(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*))
1795      RETURN
1800 Key4:      Key=4
1805      CALL Menu_status(Menu,Key,Busy,Menu$(*))
1810      ON Menu GOSUB M1k4,M2k4,M3k4,M4k4,M5k4,M6k4,M7k4
1815      CALL Menu_status(Menu,Key,Ready,Menu$(*))
1820      CALL Tcs8read(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*))
1825      RETURN
1830 Key5:      Key=5
1835      CALL Menu_status(Menu,Key,Busy,Menu$(*))
1840      ON Menu GOSUB M1k5,M2k5,M3k5,M4k5,M5k5,M6k5,M7k5
1845      CALL Menu_status(Menu,Key,Ready,Menu$(*))
1850      CALL Tcs8read(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*))
1855      RETURN
1860 Key6:      Key=6
1865      CALL Menu_status(Menu,Key,Busy,Menu$(*))
1870      ON Menu GOSUB M1k6,M2k6,M3k6,M4k6,M5k6,M6k6,M7k6
1875      CALL Menu_status(Menu,Key,Ready,Menu$(*))
1880      CALL Tcs8read(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*))
1885      RETURN
1890 Key7:      Key=7
1895      CALL Menu_status(Menu,Key,Busy,Menu$(*))
1900      ON Menu GOSUB M1k7,M2k7,M3k7,M4k7,M5k7,M6k7,M7k7
1905      CALL Menu_status(Menu,Key,Ready,Menu$(*))
1910      CALL Tcs8read(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*))
1915      RETURN
1920 Key8:      Key=8
1925      CALL Menu_status(Menu,Key,Busy,Menu$(*))
1930      ON Menu GOSUB M1k8,M2k8,M3k8,M4k8,M5k8,M6k8,M7k8
1935      CALL Menu_status(Menu,Key,Ready,Menu$(*))
1940      CALL Tcs8read(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*))
1945      RETURN
1950 Menu1:      ! Descriptions of the "Main Menu" subroutines M1K1,...,M1K8:
1955      !      The eight subroutines M1K1,...,M1K8 together implement the "Main Menu". The following will be
1960      !      displayed at the top left of the CRT display when the "Main Menu" is selected:
1965      !
1970      !      M1K1: Laser Alignment
1975      !      M1K2: Pre Run
1980      !      M1K3: Post Run (Dump Graphics)
1985      !      M1K4: Set Auto Move Positions
1990      !      M1K5: Move traverse
1995      !      M1K6: Take data
2000      !      M1K7: Auto move and take
2005      !      M1K8: Display Histograms
2010      !
2015      !      M1K1 will change the current active menu from the "Main Menu" to the "Laser Alignment Menu".
2020      !      M1K2 will change the current active menu from the "Main Menu" to the "Pre Run Menu". M1K3 will
2025      !      transfer the graphics contents of the CRT to the printer. This provides a hard copy of the profile
2030      !      plots. M1K4 has the user enter predefined traverse positions for a profile plot. M1K5 moves the
2035      !      traverse to a user selectable position. M1K6 acquires LVDAS data at the current TCS8 traverse
2040      !      position. M1K7 acquires LVDAS data at each of the pre programed TCS8 traverse positions set up by
2045      !      M1K4. M1K8 repeatedly displays five channels of real time histograms until the user presses any
2050      !      key on the keyboard.
2055      !
2060 M1k1:      ! Change the current active menu from the "Main Menu" to the "Laser Alignment Menu".
2065      Menu=2
2070      CALL Menu_disp(Menu,Menu$(*))
2075      RETURN
2080 M1k2:      ! Change the current active menu from the "Main Menu" to the "Pre Run Menu".
2085      Menu=3
2090      CALL Menu_disp(Menu,Menu$(*))
2095      RETURN

```

```

2100 M1k3:      ! Transfer the graphics contents of the CRT to the printer. This provides a hard copy of the plots.
                KEY LABELS OFF      ! Turn off the key labels so that they won't be printed.
2105      PRINTER IS CRT;WIDTH 132
2110      DISP ""      ! Clear the CRT's display line so that they won't be printed.
2115      FOR L=1 TO 9      ! Clear the CRT's menu lines so that it won't be printed.
2120      PRINT TABXY(1,L);RPTS(" ",120)
2125      NEXT L
2130      PRINTER IS PRT      ! Direct printed output to the printer.
2135      GOSUB Print_header      ! Print the "header" tunnel conditions.
2140      CALL Dump      ! Dump the entire CRT's contents to the printer.
2145      PRINT USING "#,@"      ! Move to the top of the next page on the printer.
2150      PRINTER IS CRT      ! Direct printed output to the CRT.
2155      CALL Menu_disp(Menu,Menu$(*))      ! Redisplay the menus.
2160      RETURN
2165      ! Have the user enter predefined traverse positions for a profile plot.
2170 M1k4:      CALL Enter_value("number of traverse positions",Npos,"K")
2175      REDIM Pos(Npos,1),Pname$(Npos,1),Pimage$(Npos,1),Punits$(Npos,1)
2180      MAT Pimage$= ("M4D.4D")
2185      MAT Punits$= ("in")
2190      FOR K=1 TO Npos
2195          Pname$(K,1)="Pos"&VAL$(K)
2200      NEXT K
2205      GSTORE Gsave(*)
2210      CALL Change("VALUES",Pos(*),Pname$(*),Pimage$(*),Punits$(*))
2215      GLOAD Gsave(*)
2220      CALL Menu_disp(Menu,Menu$(*))
2225      RETURN
2230      ! Moves the traverse to a user selectable position.
2235 M1k5:      GOSUB Read_calc_fill
2240      CALL Tcs8read(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*))
2245      CALL Enter_value(CHR$(NUM("X")+Paxis-1),Mod1(Paxis),"K")
2250      ON KBD CALL Do_nothing
2255 M1k5a:      DISP "Moving"
2260      Movement=Mod1(Paxis)
2265      CALL Tcs8move(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*),"Tx & Rx","MODEL","ABSOLUTE",
2270          Paxis,Movement)
2275      CALL Tcs8read(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*))
2280      GOSUB Calc
2285      GOSUB Fill
2290      DISP ""
2295      OFF KBD
2300      RETURN
2305 M1k6:      ! Acquire LVDAS data at the current TCS8 traverse position.
2310      DISP "Press any key to TAKE DATA"
2315      CALL Rt_histo(@Lvdas,Symbols(*),1,Kbd$)
2320      IF POS(Kbd$,"Q") THEN RETURN
2325      Cmask=Coin(1)*1+Coin(2)*2+Coin(3)*4
2330      Nsam=MIN(Nreads,1000)
2335      Date=TIMEDATE
2340      Time=Date
2345      CALL Lvdas_take(@Lvdas,Atime,Ctime,At_exp,Ct_exp,Cmask,Nsam)
2350      IF Nsam>1 THEN
2355          GOSUB Process_data
2360          OUTPUT PRT USING "K,K";CHR$(27)&"&k2S"&CHR$(27)&"&19D",RPTS("=",140)
2365          GOSUB Store_file
2370          File=File+1
2375      END IF
2380      RETURN
2385 M1k7:      ! Acquire LVDAS data at each of the pre programmed TCS8 traverse positions set up by M1K4.
2390      Quit=0
2395      FOR J=1 TO Npos
2400          CALL Tcs8read(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*))
2405          Mod1(Paxis)=Pos(J,1)
2410          GOSUB M1k5a
2415          GOSUB M1k6
2420          IF POS(Kbd$,"Q") THEN 2430
2425      NEXT J
2430      GOSUB On_key
2435      CALL Menu_disp(Menu,Menu$(*))
2440      RETURN
2445 M1k8:      ! Repeatedly displays five channels of real time histograms until the user presses any key on the keyboard.
2450      DISP "Press any key to return to main menu"
2455      CALL Rt_histo(@Lvdas,Symbols(*),1,Kbd$)
2460      RETURN
2465 Menu2:      ! Descriptions of the "Laser Alignment Menu" subroutines M2K1,...,M2K8:
2470      ! The eight subroutines M2K1,...,M2K8 together implement the "Laser Alignment Menu". The
2475      ! following will be displayed at the top left of the CRT display when the "Laser Alignment Menu" is
2480      ! selected:
2485      !
2490      ! M2K1: Return to main menu

```



```

2495      !           M2K2: Sides      : Tx & Rx
2500      !           M2K3: Coordinates: MODEL
2505      !           M2K4: Mode      : ABSOLUTE
2510      !           M2K5: Move X
2515      !           M2K6: Move Y
2520      !           M2K7: Move Z
2525      !           M2K8: Move A
2530      !
2535      !           M2K1 will change the current active menu from the "Laser Alignment Menu" to the "Main Menu".
2540      !           M2K2 selects whether the transmitting, receiving, or both sides of the traverse are to be moved.
2545      !           M2K3 selects the TUNNEL or MODEL coordinate systems for traverse movements. M2K4 specifies
2550      !           movements to be relative to the current's position or to absolute positions. M2K5 has the user
2555      !           enter a movement for the X axis and then the movement is performed. M2K6 has the user enter
2560      !           a movement for the Y axis and then the movement is performed. M2K7 has the user enter a movement
2565      !           for the Z axis and then the movement is performed. M2K8 has the user enter a movement for the A
2570      !           axis and then the movement is performed.
2575      !
2580      ! Change the current active menu from the "Laser Alignment Menu" to the "Main Menu".
2585      Menu=1
2590      CALL Menu_disp(Menu,Menu$(*))
2595      RETURN
2600      ! Select whether the transmitting, receiving, or both sides of the traverse are to be moved.
2605      SELECT TRIMS(Menu$(Menu,Key)[20])
2610      CASE "Tx & Rx"
2615      Menu$(Menu,Key)[20]="Tx"
2620      CASE "Tx"
2625      Menu$(Menu,Key)[20]="Rx"
2630      CASE "Rx"
2635      Menu$(Menu,Key)[20]="Tx & Rx"
2640      END SELECT
2645      CALL Menu_disp(Menu,Menu$(*))
2650      RETURN
2655      ! Selects the TUNNEL or MODEL coordinate systems for traverse movements.
2660      SELECT TRIMS(Menu$(Menu,Key)[20])
2665      CASE "MODEL"
2670      Menu$(Menu,Key)[20]="TUNNEL"
2675      CASE "TUNNEL"
2680      Menu$(Menu,Key)[20]="MODEL"
2685      END SELECT
2690      CALL Menu_disp(Menu,Menu$(*))
2695      RETURN
2700      ! Specifies movements to be relative to the current's position or to absolute positions.
2705      SELECT TRIMS(Menu$(Menu,Key)[20])
2710      CASE "ABSOLUTE"
2715      Menu$(Menu,Key)[20]="RELATIVE"
2720      CASE "RELATIVE"
2725      Menu$(Menu,Key)[20]="ABSOLUTE"
2730      END SELECT
2735      CALL Menu_disp(Menu,Menu$(*))
2740      RETURN
2745      ! The subroutines M2K5 thru M2K8 all execute the same code. The code will have the user enter a
2750      ! movement for the X,Y,Z, or A depending on what the value of "Key" is. The user specified movement
2755      ! for the selected axis will then be performed.
2760      !
2765      Side$=TRIMS(Menu$(Menu,2)[20])
2770      Coord$=TRIMS(Menu$(Menu,3)[20])
2775      Mode$=TRIMS(Menu$(Menu,4)[20])
2780      CALL Enter_value(Mode$&" Movement",Movement,"4D.5D")
2785      ON KBD CALL Do_nothing
2790      DISP "Moving"
2795      CALL Tcs8read(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*))
2800      CALL Tcs8move(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*),Side$,Coord$,Mode$,Key-4,Movement)
2805      CALL Tcs8read(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*))
2810      DISP ""
2815      OFF KBD
2820      RETURN
2825      ! Descriptions of the "Pre Run Menu" subroutines M3K1,...,M3K8:
2830      ! The eight subroutines M3K1,...,M3K8 together implement the "Pre Run Menu". The following will
2835      ! be displayed at the top left of the CRT display when the "Pre Run Menu" is selected:
2840      !
2845      !           M3K1: Return to MAIN menu
2850      !           M3K2: Enter Run & File Numbers
2855      !           M3K3: Enter Number of Samples
2860      !           M3K4: Select Traverse Axis for Profile
2865      !           M3K5: Print Coordinate Transformation Matrices
2870      !           M3K6: Setup Graphics
2875      !           M3K7: Tunnel Conditions
2880      !           M3K8: Traverse
2885      !
2890      ! M3K1 will change the current active menu from the "Pre Run Menu" to the "Main Menu". M3K2 has

```

```

2895      ! the user enter a the Run and File numbers. A new run number should be assigned to each profile
2900      ! while a new file number is assigned to each set of data. M3K3 has the user enter the desired
2905      ! number of samples. M3K4 has the user select which axis to traverse in for the profiles. M3K5
2910      ! prints the coordinate system transformation matrices for both traverse positions and velocities.
2915      ! M3K6 creates a new set of empty plots for new profiles. M3K7 will change the current active menu
2920      ! from the "Pre Run Menu" to the "Tunnel Conditions Menu". M3K8 will change the current active menu
2925      ! from the "Pre Run Menu" to the "Traverse Menu".
2930      !
2935 M3k1:      ! Change the current active menu from the "Pre Run Menu" to the "Main Menu".
2940      Menu=1
2945      CALL Menu_disp(Menu,Menu$(*))
2950      RETURN
2955 M3k2:      ! Have the user enter a the Run and File numbers.
2960      CALL Enter_value("Run",Run,"3D.2D")
2965      CALL Enter_value("File",File,"3D")
2970      RETURN
2975 M3k3:      ! Have the user enter the desired number of samples.
2980      CALL Enter_value("Number of Samples ",Nreads,"K")
2985      RETURN
2990 M3k4:      ! Have the user select which axis to traverse in for the profiles.
2995      CALL Enter_string("Traverse Axis for Profile ",Paxis$,"K")
3000      SELECT Paxis$
3005      CASE "X"
3010          Paxis=1
3015      CASE "Y"
3020          Paxis=2
3025      CASE "Z"
3030          Paxis=3
3035      CASE "A"
3040          Paxis=4
3045      CASE ELSE
3050          GOTO M3k4
3055      END SELECT
3060      GOSUB Fill
3065      RETURN
3070 M3k5:      ! Prints the coordinate system transformation matrices for both traverse positions and velocities.
3075      GOSUB Read_calc_fill
3080      OUTPUT PRT USING "#,2/"
3085      OUTPUT PRT USING "20X,K,/";"TRAVERSE COORDINATE TRANSFORMATION MATRICES"
3090      OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/);"TUNNEL to MODEL",Tun2mod(*)
3095      OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/);"MODEL to TUNNEL",Mod2tun(*)
3100      OUTPUT PRT USING "20X,K,/";"VELOCITY COORDINATE TRANSFORMATION MATRICES"
3105      OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/);"TUNNEL to MODEL",Tun2mod(*)
3110      OUTPUT PRT USING "20X,K,/,3(13X,3(8D.5D),/);"MODEL to TUNNEL",Mod2tun(*)
3115      OUTPUT PRT USING "#,@"
3120      RETURN
3125 M3k6:      ! Display a new set of plots for new profiles.
3130      CALL Setup_graph(Array(*),Image$(*),Paxis,Symbols(*))
3135      RETURN
3140 M3k7:      ! Change the current active menu from the "Pre Run Menu" to the "Tunnel Conditions Menu".
3145      Menu=4
3150      CALL Menu_disp(Menu,Menu$(*))
3155      RETURN
3160 M3k8:      ! Change the current active menu from the "Pre Run Menu" to the "Traverse Menu".
3165      Menu=5
3170      CALL Menu_disp(Menu,Menu$(*))
3175      RETURN
3180 Menu4:      ! Descriptions of the "Tunnel Conditions Menu" subroutines M4K1,...,M4K8:
3185      ! The eight subroutines M4K1,...,M4K8 together implement the "Tunnel Conditions Menu". The
3190      ! following will be displayed at the top left of the CRT display when the "Tunnel Conditions Menu" is
3195      ! selected:
3200      !
3205      ! M4K1: Return to PRE RUN menu
3210      ! M4K2: Load Tunnel Conditions
3215      ! M4K3: Save Tunnel Conditions
3220      ! M4K4: Print Tunnel Conditions
3225      ! M4K5: Enter Tunnel Condition Data
3230      ! M4K6: Enter Tunnel Condition Names
3235      ! M4K7: Enter Tunnel Condition Units
3240      ! M4K8: Enter Tunnel Condition Images
3245      !
3250      ! M4K1 will change the current active menu from the "Tunnel Conditions Menu" to the "Pre Run
3255      ! Menu". M4K2 loads the old tunnel conditions from a file on the disk. M4K3 saves the current
3260      ! tunnel conditions to a file on the disk. M4K2 & M4K3 load and save default tunnel conditions from
3265      ! the file "ARRAY" on the hard disk. The default values are not related to any particular run number.
3270      ! M4K4 sends the current tunnel conditions to the printer. M4K5 has the user enter values for the
3275      ! tunnel condition variables. M4K6 has the user enter names for the tunnel condition variables.
3280      ! M4K7 has the user enter units for the tunnel condition variables. M4K8 has the user enter image
3285      ! formats for the tunnel condition variables.
3290      !

```

```

6495 Menu:
6500 Menu_read: SUB Menu_read(Menu$(*))
6505 ! Description:
6510 ! This subprogram reads in the menu descriptors for each entry of the five menus.
6515 ! Variables:
6520 ! Menu Used as an index to the string array Menu$(*).
6525 ! Key Used as an index to the string array Menu$(*).
6530 ! Menu$(*) String array where each element describes its corresponding menu subroutine's function.
6535 ! LS String use to read in the menu descriptor from the data statements.
6540 OPTION BASE 1
6545 DIM LS(80)
6550 ! Fill all of the menu entry's descriptions with "MxKx".
6555 FOR Menu=1 TO SIZE(Menu$,1)
6560 FOR Key=1 TO 8
6565 Menu$(Menu,Key)="M"&VAL$(Menu)&"K"&VAL$(Key)&": "
6570 NEXT Key
6575 NEXT Menu
6580 ON ERROR GOTO 6620 ! The following while loop will get error#36 when the data statements run out.
6585 ! For each menu and key, enter the menu entry's description.
6590 WHILE 1=1
6595 READ LS
6600 Menu=VAL$(LS(2,2))
6605 Key=VAL$(LS(4,4))
6610 Menu$(Menu,Key)=LS
6615 END WHILE
6620 SUBEXIT
6625 DATA "M1K1: Laser Alignment"
6630 DATA "M2K1: Return to main menu"
6635 DATA "M2K2: Sides : Tx & Rx"
6640 DATA "M2K3: Coordinates: MODEL"
6645 DATA "M2K4: Mode : ABSOLUTE"
6650 DATA "M2K5: Move X"
6655 DATA "M2K6: Move Y"
6660 DATA "M2K7: Move Z"
6665 DATA "M2K8: Move A"
6670 DATA "M1K2: Pre Run"
6675 DATA "M3K1: Return to MAIN menu"
6680 DATA "M3K2: Enter Run & File Numbers"
6685 DATA "M3K3: Enter Number of Samples"
6690 DATA "M3K4: Select Traverse Axis for Profile"
6695 DATA "M3K5: Print Coordinate Transformation Matrices"
6700 DATA "M3K6: Setup Graphics"
6705 DATA "M3K7: Tunnel Conditions"
6710 DATA "M4K1: Return to PRE RUN menu"
6715 DATA "M4K2: Load Tunnel Conditions"
6720 DATA "M4K3: Save Tunnel Conditions"
6725 DATA "M4K4: Print Tunnel Conditions"
6730 DATA "M4K5: Enter Tunnel Condition Data"
6735 DATA "M4K6: Enter Tunnel Condition Names"
6740 DATA "M4K7: Enter Tunnel Condition Units"
6745 DATA "M4K8: Enter Tunnel Condition Images"
6750 DATA "M3K8: Traverse"
6755 DATA "M5K1: Return to PRE RUN menu"
6760 DATA "M5K2: View & Set TCS8 Positions"
6765 DATA "M5K3: View & Set TCS8 Units"
6770 DATA "M5K4: View & Set TCS8 Revolution"
6775 DATA "M5K5: View & Set TCS8 Velocity"
6780 DATA "M5K6: View & Set TCS8 Acceleration"
6785 DATA "M5K8: Recalc & Replot"
6790 DATA "M1K3: Post Run (Dump Graphics)"
6795 DATA "M1K4: Set Auto Move Positions"
6800 DATA "M1K5: Move traverse"
6805 DATA "M1K6: Take data"
6810 DATA "M1K7: Auto move and take"
6815 DATA "M1K8: Display Histograms"
6820 SUBEND
6825 Menu_disp: SUB Menu_disp(Menu,Menu$(*))
6830 ! Description:
6835 ! This subprogram displays the current menu at the top of the CRT.
6840 ! Variables:
6845 ! Menu Used as an index to the string array Menu$(*).
6850 ! Key Used as an index to the string array Menu$(*).
6855 ! Menu$(*) String array where each element describes its corresponding menu subroutine's function.
6860 COM /Color1/ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta
6865 COM /Color2/ White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
6870 PRINTER IS CRT
6875 PRINT PEN Blue ! Print the menu using blue text.
6880 PRINT CHR$(128);CHR$(129); ! Turn on inverse video.
6885 IF Menu=0 THEN Menu=1
6890 FOR Key=1 TO 8

```

```

6895         Menu$(Menu,Key)=Menu$(Menu,Key)&RPTS(" ",50-LEN(Menu$(Menu,Key)))
6900         PRINT TABXY(2,Key);" ";Menu$(Menu,Key)[3]
6905         NEXT Key
6910         PRINT CHR$(128);                ! Turn off inverse video.
6915         PRINT PEN Black                  ! Set printing color to black.
6920     SUBEND
6925 Menu_status: SUB Menu_status(Menu,Key,Pen,Menu$(*))
6930     ! Description:
6935     !     This subprogram displays the current menu selection in red or blue text. The red text
6940     !     style indicates that the subroutine for the current menu selection is busy. The blue text
6945     !     style indicates that the subroutine for the current menu selection is has completed.
6950     ! Variables:
6955     !     Menu      Indicates which of the menus has been selected as the current menu.
6960     !     Key       Indicates which one of eight menu subroutines in the menu is to be executed.
6965     !     Pen       Indicates Busy/Ready Status. Busy: Pen=Red. Ready: Pen=Blue.
6970     !     Menu$(*) String array. Each element describes its corresponding menu subroutine's function.
6975     COM /Color1/ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta
6980     COM /Color2/ White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
6985     PRINT PEN Pen
6990     PRINTER IS CRT
6995     PRINT CHR$(128);CHR$(129);                ! Turn on inverse video.
7000     PRINT TABXY(2,Key);" ";Menu$(Menu,Key)[3];CHR$(128) ! Print menu selection & turn off inverse video.
7005     PRINT PEN Black                            ! Set printing color to black.
7010     WAIT .1
7015     SUBEND
7020 Enter: !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
7025 Enter_value: SUB Enter_value(Name$,Value,Image$)
7030     ! Description:
7035     !     This subprogram displays the current value of a variable and then has the user enter its new
7040     !     value. The old value will be kept if the RETURN key is pressed and no data is entered.
7045     ! Variables:
7050     !     Name$      Name of the variable.
7055     !     Image$     Image format of the variable. Used for printing the variable with a format.
7060     !     Value      Contains the initial value and then the updated value for the variable.
7065     IF Name$="Date" OR Name$="Time" THEN SUBEXIT ! Turn on inverse video.
7070     DISP CHR$(129);                             ! Display name and old value for the variable.
7075     DISP USING 7080;Name$
7080     IMAGE #,"Old ",K,"="
7085     IF Image$<>" " THEN DISP USING #,"&Image$;Value
7090     IF Image$="" THEN DISP USING "#,K";Value
7095     DISP USING 7100;Name$
7100     IMAGE #,"      Enter new ",K
7105     INPUT " ? ",Value                          ! The user enters the new value here.
7110     DISP CHR$(128);                             ! Turn off inverse video.
7115     SUBEND
7120 Enter_string: SUB Enter_string(Name$,Value$,Image$)
7125     ! Description:
7130     !     This subprogram displays the current value of a string variable and then has the user enter its
7135     !     new value. The old value will be kept if the RETURN key is pressed and no data is entered.
7140     ! Variables:
7145     !     Name$      Name of the variable.
7150     !     Value$     Contains the initial value and then the updated value for the string variable.
7155     DISP CHR$(129);                             ! Turn on inverse video.
7160     DISP USING 7165;Name$                       ! Display name and old value for the string.
7165     IMAGE #,"Old ",K,"="
7170     DISP USING "#,&Image$;Value$"
7175     DISP USING 7180;Name$
7180     IMAGE #,"      Enter new ",K
7185     INPUT " ? ",Value$                         ! The user enters the new string value here.
7190     DISP CHR$(128);                             ! Turn off inverse video.
7195     SUBEND
7200 Array: !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
7205 Array_init: SUB Array_init(Name$(*),Array(*),Image$(*),Units$(*))
7210     ! Description:
7215     !     This subprogram reads in default data for each of the variable's names, values, image formats,
7220     !     and units. These variables include, but are not limited to, the tunnel conditions, laser
7225     !     parameters, graph scales, traverse positions, and coordinate system transformation matrices.
7230     ! Variables:
7235     !     Array(*)   Array of tunnel conditions, laser parameters, graph scales, etc.
7240     !     Name$(*)   Names for the variables in Array(*).
7245     !     Image$(*)  Image formats for the variables in Array(*).
7250     !     Units$(*)  Units for the variables in Array(*).
7255     !     X          Used as an index to the above arrays and string arrays.
7260     !     Y          Used as an index to the above arrays and string arrays.
7265     !     Before     Number of digits before the decimal point in the image format.
7270     !     After      Number of digits after the decimal point in the image format.
7275     ON ERROR GOTO 7365
7280     READ Y
7285     FOR X=1 TO SIZE(Name$,2)
7290         READ Name$(Y,X),Array(Y,X),Image$(Y,X),Units$(Y,X)

```

```

7295      SELECT Image$(Y,X)
7300      CASE "0"
7305          Image$(Y,X) = "9D"
7310      CASE "1" TO "7"
7315          After = VAL(Image$(Y,X))
7320          Before = 8 - After
7325          Image$(Y,X) = VAL$(Before) & "D." & VAL$(After) & "D"
7330      CASE "K"
7335      CASE "N"
7340      CASE ELSE
7345          Image$(Y,X) = "9D"
7350      END SELECT
7355  NEXT X
7360  GOTO 7280
7365  SUBEXIT
7370  ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
7375  DATA 1, Date , 0,0,"", Mach , 7.0,4,"", STemp , 0,0,"R", Tt Gain , 100,0,""
7380  DATA 2, Time , 0,0,"", Re/Ft , 0,0,"/Ft", TTemp , 0,0,"R", Alpha1 , 0,4,"°"
7385  DATA 3, Run , 0,2,"", Uedge , 1,4,m/s, Tt , 0,3,mv, Alpha2 , 0,4,"°"
7390  DATA 4, File , 0,0,"", Uinf , 1,4,m/s, Tt (raw), 0,3,v , Alpha3 , 0,4,"°"
7395  ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
7400  DATA 11, Xtun1 , 0,4,in , Xtun2 , 0,4,in , Xmod1 , 0,4,in , Xmod2 , 0,4,in
7405  DATA 12, Ytun1 , 0,4,in , Ytun2 , 0,4,in , Ymod1 , 0,4,in , Ymod2 , 0,4,in
7410  DATA 13, Ztun1 , 0,4,in , Ztun2 , 0,4,in , Zmod1 , 0,4,in , Zmod2 , 0,4,in
7415  DATA 14, Atun1 , 0,4,in , Atun2 , 0,4,in , Amod1 , 0,4,in , Amod2 , 0,4,in
7420  ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
7425  DATA 31, UBeamSpc, .3125,3,in , VBeamSpc, .3125,3,in , WBeamSpc, .3125,3,in , "" , 0,0,""
7430  DATA 32, UFocLen, 30.00,3,in , VFocLen, 30.00,3,in , WFocLen, 30.00,3,in , "" , 0,0,""
7435  DATA 33, UBeamSep, 0.000,3,"°", VBeamSep, 0.000,3,"°", WBeamSep, 0.000,3,"°", "" , 0,0,""
7440  DATA 34, UWaveLen, 514.5,3,nm , VWaveLen, 488.0,3,nm , WWaveLen, 476.5,3,nm , "" , 0,0,""
7445  DATA 35, UFrngSpc, 00.00,3,um , VFrngSpc, 00.00,3,um , WFrngSpc, 00.00,3,um , "" , 0,0,""
7450  DATA 36, Ubrag , 40.00,4,MHz, Vbrag , 40.00,4,MHz, Wbrag , 40.00,4,MHz, "" , 0,0,""
7455  DATA 37, Umix , 0.00,4,MHz, Vmix , 0.00,4,MHz, WMix , 0.00,4,MHz, "" , 0,0,""
7460  DATA 38, UmeaSgn , -1,0,"", VmeaSgn , +1,0,"", WmeaSgn , +1,0,"", "" , 0,0,""
7465  DATA 39, UbrgSgn , +1,0,"", VbrgSgn , -1,0,"", WbrgSgn , -1,0,"", "" , 0,0,""
7470  DATA 40, UmixSgn , -1,0,"", VmixSgn , +1,0,"", WmixSgn , +1,0,"", "" , 0,0,""
7475  DATA 41, U coin , 1,0,"", V coin , 1,0,"", W coin , 0,0,"", "" , 0,0,""
7480  DATA 42, UFreqMin, 8,4,MHz, VFreqMin, 25,4,MHz, WFreqMin, -99,4,MHz, "" , 0,0,""
7485  DATA 43, UFreqMax, 32,4,MHz, VFreqMax, 55,4,MHz, WFreqMax, 99,4,MHz, "" , 0,0,""
7490  ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
7495  DATA 51, Nreads , 1000,0,"", Atime , 5,6,s , AText , 12,0,"", Paxis , 2,0,""
7500  DATA 52, Nsam , 1000,0,"", Ctime , 1E-2,6,s , CText , 7,0,"", Clip , 1,0,""
7505  DATA 53, "" , 0,0,"", "" , 0,0,"", "" , 0,0,"", Nose , 139,1,cm
7510  ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
7515  DATA 61, Xmin1 , 0.00,0,"", Xmax1 , 60.00,0,"", Ymin1 , 0,0,"", Ymax1 , 100,0,""
7520  DATA 62, Xmin2 , 0.00,0,"", Xmax2 , 60.00,0,"", Ymin2 , 0,0,"", Ymax2 , 100,0,""
7525  DATA 63, Xmin3 , 0.00,0,"", Xmax3 , 60.00,0,"", Ymin3 , 0,0,"", Ymax3 , 100,0,""
7530  DATA 64, Xmin5 , -5.00,0,"", Xmax5 , 5.00,0,"", Ymin5 , 0,0,"", Ymax5 , 1000,0,""
7535  DATA 65, Xmin5 , -5.00,0,"", Xmax5 , 5.00,0,"", Ymin5 , 0,0,"", Ymax5 , 1000,0,""
7540  DATA 66, Xmin6 , -0.50,1,"", Xmax6 , 1.50,1,"", Ymin6 , 0.00,2,"", Ymax6 , 4.00,2,""
7545  DATA 67, Xmin7 , -5.00,0,"", Xmax7 , 5.00,0,"", Ymin7 , 0.00,2,"", Ymax7 , 4.00,2,""
7550  DATA 68, Xmin8 , 0,1,"", Xmax8 , 2000,1,"", Ymin8 , 0.00,2,"", Ymax8 , 4.00,2,""
7555  DATA 69, Xmin9 , -1.50,1,"", Xmax9 , 1.50,1,"", Ymin9 , 0.00,2,"", Ymax9 , 4.00,2,""
7560  ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
7565  DATA 71, Xmin1 , 835,0,pxl, Xmax1 , 1235,0,pxl, Ymin1 , 725,0,pxl, Ymax1 , 825,0,pxl
7570  DATA 72, Xmin2 , 835,0,pxl, Xmax2 , 1235,0,pxl, Ymin2 , 585,0,pxl, Ymax2 , 685,0,pxl
7575  DATA 73, Xmin3 , 835,0,pxl, Xmax3 , 1235,0,pxl, Ymin3 , 445,0,pxl, Ymax3 , 545,0,pxl
7580  DATA 74, Xmin4 , 835,0,pxl, Xmax4 , 1235,0,pxl, Ymin4 , 305,0,pxl, Ymax4 , 405,0,pxl
7585  DATA 75, Xmin5 , 835,0,pxl, Xmax5 , 1235,0,pxl, Ymin5 , 165,0,pxl, Ymax5 , 265,0,pxl
7590  DATA 76, Xmin6 , 75,0,pxl, Xmax6 , 325,0,pxl, Ymin6 , 525,0,pxl, Ymax6 , 825,0,pxl
7595  DATA 77, Xmin7 , 425,0,pxl, Xmax7 , 675,0,pxl, Ymin7 , 525,0,pxl, Ymax7 , 825,0,pxl
7600  DATA 78, Xmin8 , 75,0,pxl, Xmax8 , 325,0,pxl, Ymin8 , 165,0,pxl, Ymax8 , 465,0,pxl
7605  DATA 79, Xmin9 , 425,0,pxl, Xmax9 , 675,0,pxl, Ymin9 , 165,0,pxl, Ymax9 , 465,0,pxl
7610  ! Y *****X=1***** *****X=2***** *****X=3***** *****X=4*****
7615  DATA 81, Xdiv1 , 6,0,"", Ydiv1 , 5,0,"", Xdiv6 , 4,0,"", Ydiv6 , 8,0,""
7620  DATA 82, Xdiv2 , 6,0,"", Ydiv2 , 5,0,"", Xdiv7 , 10,0,"", Ydiv7 , 8,0,""
7625  DATA 83, Xdiv3 , 6,0,"", Ydiv3 , 5,0,"", Xdiv8 , 8,0,"", Ydiv8 , 8,0,""
7630  DATA 84, Xdiv4 , 10,0,"", Ydiv4 , 5,0,"", Xdiv9 , 6,0,"", Ydiv9 , 8,0,""
7635  DATA 85, Xdiv5 , 10,0,"", Ydiv5 , 5,0,"", "" , 0,0,"", "" , 0,0,""
7640  SUBEND
7645  Array_print: SUB Array_print(Array(*),Name$(*),Image$(*),Units$(*))
7650      ! Description:
7655      ! This subprogram prints the values of each of the variables with their names, image formats, and
7660      ! units. These variables include, but are not limited to, the tunnel conditions, laser
7665      ! parameters, and graph scales.
7670      ! Variables:
7675      ! Array(*) Array of tunnel conditions, laser parameters, graph scales, etc.
7680      ! Name$(*) Names for the variables in Array(*).
7685      ! Image$(*) Image formats for the variables in Array(*).
7690      ! Units$(*) Units for the variables in Array(*).

```

```

7695      !      X      Used as in index to the above arrays and string arrays.
7700      !      Y      Used as in index to the above arrays and string arrays.
7705  PRINT USING "#,5/"
7710  FOR Y=1 TO SIZE(Array,1)
7715      MAT SEARCH Array(Y,"),#LOC(<>0);L1
7720      MAT SEARCH Name$(Y,"),#LOC(<>"");L2
7725      IF L1+L2=0 AND L3=0 THEN 7845
7730      L3=L1+L2
7735      PRINT USING "#,28X"
7740      FOR X=1 TO SIZE(Array,2)
7745          SELECT Name$(Y,X)
7750          CASE "" ! If the variable has no name, then print just blanks.
7755              PRINT USING "#,28X"
7760          CASE "Date" ! Use a special printing format for printing the date.
7765              LS=DATE$(Array(Y,X))
7770              LS=LS[1,2]&LS[4,6]&LS[8,11]
7775              PRINT USING "#,8A,A,9A,X,3A,6X";TRIMS(Name$(Y,X)),"=",LS,Units$(Y,X)
7780          CASE "Time" ! Use a special printing format for printing the time.
7785              LS="&TIMES(Array(Y,X))
7790              PRINT USING "#,8A,A,9A,X,3A,6X";TRIMS(Name$(Y,X)),"=",LS,Units$(Y,X)
7795          CASE ELSE ! All others use a standard format.
7800              IF Image$(Y,X)="" THEN Image$(Y,X)="9D"
7805              ON ERROR GOTO 7820
7810              PRINT USING "#,8A,A,&Image$(Y,X)&","X,3A,6X";TRIMS(Name$(Y,X)),"=",Array(Y,X),Units$(Y,X)
7815              GOTO 7830
7820              OFF ERROR
7825              PRINT USING "#,8A,A,K,X,3A,6X";TRIMS(Name$(Y,X)),"=",Array(Y,X),Units$(Y,X)
7830          END SELECT
7835      NEXT X
7840      PRINT
7845  NEXT Y
7850  SUBEND
7855  Change: !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
7860  Change: SUB Change(Type$,Array(*),Name$(*),Image$(*),Units$(*))
7865      ! Description:
7870      ! This subprogram displays on the CRT the values of each of the variables with their names,
7875      ! image formats, and units. The user can select one of the variables and enter a new value,
7880      ! name, image format, or units. The user selects the particular variable by using the
7885      ! left, right, up, and down cursor keys. The selected variable will appear in inverse video.
7890      ! When it is not selected, it will appear in normal text. When the user has selected the
7895      ! appropriate variable he should then press the "Select" key on the keyboard. Then, depending on
7900      ! the value of Type$, he will be asked to enter a new value, name, image format, or units. To
7905      ! exit the change variables mode press the "Escape" key.
7910      ! There are three types of data that are passed to the subprogram. The first type of data
7915      ! include, but are not limited to, the tunnel conditions, laser parameters, and graph scales.
7920      ! With this first type the user is allowed to enter new variable values, names, image formats, and
7925      ! units. The second type of data are the "Auto Move and Take" data. These data are for the pre
7930      ! programmed traverse positions used in a profile scan. The third type of data are the "View and
7935      ! Set TCS8 parameters" data acquired from and then sent back to the TCS8.
7940      ! Variables:
7945      ! Array(*) Array whose values, names, image formats, or units are to be modified.
7950      ! Name$(*) Names for the variables in Array(*).
7955      ! Image$(*) Image formats for the variables in Array(*).
7960      ! Units$(*) Units for the variables in Array(*).
7965      ! Type$ Indicates which type of data is to be entered.
7970      ! Type$="VALUES" has the user enter a new value for the selected variable.
7975      ! Type$="NAMES" has the user enter a new name for the selected variable.
7980      ! Type$="IMAGES" has the user enter a new image format for the selected variable.
7985      ! Type$="UNITS" has the user enter a new units for the selected variable.
7990      ! X,X1,X2 Used as in index to the above arrays and string arrays.
7995      ! Y,Y1,Y2 Used as in index to the above arrays and string arrays.
8000  GRAPHICS OFF ! Turn off the graphics contents of the CRT.
8005  PRINTER IS CRT ! Direct printed output to the CRT.
8010  FOR Y=1 TO SIZE(Array,1)
8015      ! Search Array(*) for section containing variables.
8020      FOR Y1=Y TO SIZE(Array,1)
8025          FOR X=1 TO SIZE(Array,2)
8030              IF Name$(Y1,X)<>" THEN 8055
8035          NEXT X
8040      NEXT Y1
8045      CLEAR SCREEN ! If no more variables are found in Array(*), the Clear the CRT display and exit.
8050      SUBEXIT
8055      ! Search Array(*) for section empty of variables.
8060      FOR Y2=Y1 TO SIZE(Array,1)
8065          FOR X=1 TO SIZE(Array,2)
8070              IF Name$(Y2,X)<>" THEN 8085
8075          NEXT X
8080          GOTO 8090
8085      NEXT Y2
8090      ! Find the length of the following empty section.

```

```

8095         FOR Y2=Y2 TO SIZE(Array,1)
8100             FOR X=1 TO SIZE(Array,2)
8105                 IF Name$(Y2,X) <> "" THEN 8120
8110                     NEXT X
8115             NEXT Y2
8120             Y2=Y2-1
8125             ! Clear the CRT and then display the section contain variables and the following empty section.
8130             CLEAR SCREEN
8135             CALL Display(Type$,Y1,Y2,Array(*),Name$(*),Image$(*),Units$(*))
8140             Done=0
8145             X=1
8150             Y=Y1
8155             ON KBD ALL,15 GOSUB Kbd
8160 Wait:      IF NOT Done THEN Wait      ! The program will wait hear until a key is pressed on the keyboard.
8165             OFF KBD
8170             CLEAR SCREEN
8175             Y=Y2
8180         NEXT Y
8185         GRAPHICS ON ! Turn the graphic part of the CRT back on.
8190         SUBEXIT
8195 Kbd:        ! This subroutine will be called when one of the cursor, select, etc. keys is pressed.
8200             CALL Update(Type$,X,Y,Y1,Y2,Done,Array(*),Name$(*),Image$(*),Units$(*))
8205             RETURN
8210     SUBEND
8215 Display:    SUB Display(Type$,Y1,Y2,Array(*),Name$(*),Image$(*),Units$(*))
8220             ! Description:
8225             ! This subprogram displays on the CRT the values of each of variables with their names, image
8230             ! formats, and units.
8235             ! Variables:
8240             ! Array(*)      Array whose values, names, image formats, or units are to be modified.
8245             ! Name$(*)      Names for the variables in Array(*).
8250             ! Image$(*)     Image formats for the variables in Array(*).
8255             ! Units$(*)    Units for the variables in Array(*).
8260             ! Type$        Indicates which type of data is to be entered.
8265             !               Type$="VALUES" has the user enter a new value for the selected variable.
8270             !               Type$="NAMES"  has the user enter a new name for the selected variable.
8275             !               Type$="IMAGES" has the user enter a new image format for the selected variable.
8280             !               Type$="UNITS"  has the user enter a new units for the selected variable.
8285             ! X,X1,X2      Used as in index to the above arrays and string arrays.
8290             ! Y,Y1,Y2      Used as in index to the above arrays and string arrays.
8295             FOR Y=Y1 TO Y2
8300                 FOR X=1 TO SIZE(Array,2)
8305                     CALL Select(Type$,X,Y,Y1,Y2,0,Array(*),Name$(*),Image$(*),Units$(*))
8310                 NEXT X
8315             NEXT Y
8320             CALL Select(Type$,1,Y1,Y1,Y2,1,Array(*),Name$(*),Image$(*),Units$(*))
8325         SUBEND
8330 Select:     SUB Select(Type$,X,Y,Y1,Y2,C,Array(*),Name$(*),Image$(*),Units$(*))
8335             ! Description:
8340             ! This subprogram displays on the CRT the value of one variable along with its names, image
8345             ! format, and units.
8350             ! Variables:
8355             ! Array(*)      Array whose values, names, image formats, or units are to be modified.
8360             ! Name$(*)      Names for the variables in Array(*).
8365             ! Image$(*)     Image formats for the variables in Array(*).
8370             ! Units$(*)    Units for the variables in Array(*).
8375             ! Type$        Indicates which type of data are to be entered.
8380             !               Type$="VALUES" has the user enter a new value for the selected variable.
8385             !               Type$="NAMES"  has the user enter a new name for the selected variable.
8390             !               Type$="IMAGES" has the user enter a new image format for the selected variable.
8395             !               Type$="UNITS"  has the user enter a new units for the selected variable.
8400             ! X            Used as in index to the above arrays and string arrays.
8405             ! Y,Y1,Y2      Used as in index to the above arrays and string arrays.
8410             PRINT CHR$(128+C);TABXY(26*X-24,15+Y-Y1+1);      ! If C=0 then normal. If C=1 then inverse video.
8415             PRINT RPT$( " ",23);TABXY(26*X-24,15+Y-Y1+1);
8420             IF Name$(Y,X)="" AND Array(Y,X)=0 THEN 8520
8425             Img$=Image$(Y,X)
8430             Unt$=Units$(Y,X)
8435             IF Image$(Y,X)="" THEN Img$="K"
8440             IF Units$(Y,X)="" THEN Unt$=" "
8445             SELECT Type$
8450             CASE "VALUES"      ! If Type$="VALUES" then display the variable's value.
8455                 SELECT Name$(Y,X)
8460                 CASE "Date"
8465                 CASE "Time"
8470                 CASE ELSE
8475                     PRINT USING "#,10A,A,"&Img$&","X,3A";Name$(Y,X),":",Array(Y,X),Unt$
8480                 END SELECT
8485             CASE "NAMES"      ! If Type$="NAMES" then display the variable's name.
8490                 PRINT USING "#,10A,A,8A";Name$(Y,X),":",Name$(Y,X)

```

```

8495 CASE "UNITS" ! If Type$="UNITS" then display the variable's units.
8500 PRINT USING "#,10A,A,8A";Name$(Y,X),":",Units$(Y,X)
8505 CASE "IMAGES" ! If Type$="IMAGES" then display the variable's image format.
8510 PRINT USING "#,10A,A,8A";Name$(Y,X),":",Image$(Y,X)
8515 END SELECT
8520 PRINT CHR$(128); ! Turn off inverse video printing.
8525 SUBEND
8530 Update: SUB Update(Type$,X,Y,Y1,Y2,Done,Array(*),Name$(*),Image$(*),Units$(*))
8535 ! Description:
8540 ! This subprogram scrolls through the variables displayed on the CRT and has the user enter
8545 ! updated values. The user can select one of the variables and enter a new value, name, image
8550 ! format, or units. The user selects the particular variable by using the left, right, up, down
8555 ! cursor keys. This subprogram will only have been called after a keyboard key has been pressed.
8560 ! If a cursor key has been pressed then the previously selected variable will be redisplayed in
8565 ! normal text and the new selected variable will appear in inverse video text. When the user has
8570 ! selected the appropriate variable he will have pressed the "Select" key on the keyboard. Then,
8575 ! depending on the value of the Type$ he will be asked to enter a new value, name, image format,
8580 ! or units. To exit the change variables mode the user will have pressed the "Escape" key.
8585 ! Variables:
8590 ! Array(*) Array of tunnel conditions, laser parameters, graph scales, etc.
8595 ! Name$(*) Names for the variables in Array(*).
8600 ! Image$(*) Image formats for the variables in Array(*).
8605 ! Units$(*) Units for the variables in Array(*).
8610 ! Type$ Indicates which type of data is to be entered.
8615 ! Type$="VALUES" has the user enter a new value for the selected variable.
8620 ! Type$="NAMES" has the user enter a new name for the selected variable.
8625 ! Type$="IMAGES" has the user enter a new image format for the selected variable.
8630 ! Type$="UNITS" has the user enter a new units for the selected variable.
8635 ! X Used as in index to the above arrays and string arrays.
8640 ! Y,Y1,Y2 Used as in index to the above arrays and string arrays.
8645 DISABLE ! Disable the keyboard.
8650 KS=KBDS ! Get the key pressed from the keyboards buffer.
8655 IF KS="" THEN 8885
8660 SELECT NUM(KS[1,1])
8665 CASE 27 ! ESC key pressed.
8670 Done=1
8675 CASE 255
8680 CALL Select(Type$,X,Y,Y1,Y2,0,Array(*),Name$(*),Image$(*),Units$(*))
8685 SELECT NUM(KS[2,2])
8690 CASE 73,80 ! Break or Stop key pressed.
8695 PAUSE
8700 CASE 124 ! Menu
8705 Done=1
8710 CASE 38 ! Select key pressed.
8715 CALL Select(Type$,X,Y,Y1,Y2,1,Array(*),Name$(*),Image$(*),Units$(*))
8720 SELECT Type$
8725 CASE "VALUES"
8730 IF Name$(Y,X)="" THEN CALL Enter_string("Name for "&Name$(Y,X),Name$(Y,X),"K")
8735 IF Image$(Y,X)="" THEN CALL Enter_string("Image for "&Name$(Y,X),Image$(Y,X),"K")
8740 CALL Enter_value(Name$(Y,X),Array(Y,X),Image$(Y,X))
8745 CASE "NAMES"
8750 CALL Enter_string("Name for "&Name$(Y,X),Name$(Y,X),"K")
8755 CASE "UNITS"
8760 CALL Enter_string("Units for "&Name$(Y,X),Units$(Y,X),"K")
8765 CASE "IMAGES"
8770 CALL Enter_string("Image for "&Name$(Y,X),Image$(Y,X),"K")
8775 END SELECT
8780 CALL Select(Type$,X,Y,Y1,Y2,0,Array(*),Name$(*),Image$(*),Units$(*))
8785 IF X=SIZE(Array,2) THEN Y=Y+1
8790 X=X+1 ! Left key pressed.
8795 CASE 60
8800 X=X-1 ! Right key pressed.
8805 CASE 62
8810 X=X+1 ! Up key pressed.
8815 CASE 94
8820 Y=Y-1 ! Down key pressed.
8825 CASE 86
8830 Y=Y+1 ! First key pressed.
8835 CASE 92
8840 X=1
8845 Y=1
8850 END SELECT
8855 X=(X-1) MOD SIZE(Array,2)+1
8860 Y=(Y-1+1-1) MOD (Y2-Y1+1)+Y1
8865 IF X<1 THEN X=SIZE(Array,2)
8870 IF Y<Y1 THEN Y=Y2
8875 CALL Select(Type$,X,Y,Y1,Y2,1,Array(*),Name$(*),Image$(*),Units$(*))
8880 END SELECT
8885 ENABLE
8890 SUBEXIT

```



```

8895 SUBEND
8900 Misc: !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
8905 Convert2words: SUB Convert2words(Real,INTEGER High,Low)
8910 ! Description:
8915 ! This subprogram converts a single real precision variable into two 16 bit words. The initial
8920 ! real precision variables is converted in to a 32 bit integer and then separated into high and
8925 ! low 16 bit integers. The most significant 16 bits will be in the "High" variable while the
8930 ! least significant 16 bits will be placed the the "Low" variable. The main purpose of this
8935 ! subprogram is to provide a means to send a 32 bit integer to the LVDAS over the 16 bit high
8940 ! speed interface.
8945 ! Variables:
8950 ! Real Initial real precision value for the variable.
8955 ! Hex$ Hex value of "Real". String length will be 8 bytes for 32 bits.
8960 ! High Most significant 16 bits of integerized "Real".
8965 ! Low Least significant 16 bits of integerized "Real".
8970 Hex$=DVAL$(Real,16)
8975 High=IVAL(Hex$[1,4],16)
8980 Low=IVAL(Hex$[5,8],16)
8985 SUBEND
8990 Error: SUB Error
8995 ! Description:
9000 ! This subprogram will print an error message when ever a program error occurs. The error message
9005 ! will be displayed at the top of the CRT and also printed on the printers paper. Such errors
9010 ! might occur when data to be printed will not fit in the image formats. Other errors will also
9015 ! generate a displayed and printed error message.
9020 BEEP
9025 DISP ERRMS
9030 OUTPUT PRT;ERRMS
9035 PRT=VAL(SYSTEMS("PRINTER IS"))
9040 PRINTER IS CRT
9045 PRINT TABXY(95,1);ERRMS
9050 PRINTER IS PRT
9055 ERROR SUBEXIT
9060 SUBEND
9065 Scale: SUB Scale(G)
9070 ! Description:
9075 ! This subprogram selects one of nine histogram or profile plots. The plot's area of the CRT is
9080 ! selected and scaled to the appropriate scales.
9085 OPTION BASE 1
9090 COM /Graph1/ Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*)
9095 VIEWPORT Vwprt(G,1)/10.23,Vwprt(G,2)/10.23,Vwprt(G,3)/10.23,Vwprt(G,4)/10.23
9100 WINDOW Wndw(G,1),Wndw(G,2),Wndw(G,3),Wndw(G,4)
9105 SUBEND
9110 Table: !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
9115 Table: SUB Table(Table(*))
9120 ! Description:
9125 ! This subprogram is used to create a lookup table array. The lookup table array facilitates
9130 ! the rapid conversion of raw encoded Macrodyne data into a usable frequency. Once the table
9135 ! has been filled, then the raw Macrodyne data can be used as an index to the table array.
9140 ! Variables:
9145 ! Table(*) Lookup table of frequencies.
9150 ! Mantissa(*) The 10 bit mantissa part of the raw Macrodyne data (0..1023).
9155 ! Fringes The 1 bit Fringe Count part of the raw Macrodyne data (0:16, 1:8 fringes).
9160 ! Exponent The 4 bit Exponent part of the raw Macrodyne data.
9165 ! Time(*) An array of measurement times for a given number of Fringes and Exponent.
9170 ! Freq(*) An array of measured frequencies for a given number of Fringes and Exponent.
9175 ! Bin Used to index Mantissa(*).
9180 ! Min Used as a subrange index for Table(*).
9185 ! Max Used as a subrange index for Table(*).
9190 OPTION BASE 1
9195 REAL Mantissa(0:1023),Time(0:1023),Freq(0:1023)
9200 ! If the last entry in the table in not zero then the table has already been created.
9205 IF Table(32766) THEN SUBEXIT
9210 FOR Bin=0 TO 1023 ! Fill Mantissa array.
9215 Mantissa(Bin)=Bin
9220 NEXT Bin
9225 Mantissa(0)=1
9230 Min=0
9235 FOR Fringes=0 TO 1 ! 0 indicates 16 fringes while 1 indicates 8 fringes.
9240 FOR Exponent=0 TO 15
9245 Max=Min+1023
9250 IF Max=32767 THEN ! Maximum size of an array is 32766.
9255 Max=32766
9260 REDIM Mantissa(0:1022),Time(0:1022),Freq(0:1022)
9265 END IF
9270 DISP Fringes,Exponent
9275 MAT Time= Mantissa*(2^(Exponent-1))/500000000 ! Use this line with new macrodynes.
9280 !MAT Time= Mantissa*(2^(Exponent-3))/500000000 ! Use this line with old macrodynes.
9285 MAT Freq= (2^(4-Fringes))/Time
9290 MAT Freq= Freq/(1000000)

```

```

9295             MAT Table(Min:Max)= Freq
9300             Min=Min+1024
9305             NEXT Exponent
9310             NEXT Fringes
9315             SUBEND
9320 Lvdas:      !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
9325 Lvdas_init: SUB Lvdas_init(@Gpio)
9330             ! Description:
9335             ! This subprogram is used to initialize the HP98622-66501 Rev B 16-bit General Purpose
9340             ! Input Output (GPIO) interface. The subprogram also opens the LVDAS path on the HP computer
9345             ! for command and data transfer. The I/O path is given the name "@Lvdas". Data transferred
9350             ! from the HP to the LVDAS will use the "OUTPUT @Lvdas" statement. Data transferred to the HP
9355             ! from LVDAS will use the "ENTER @Lvdas" statement.
9360             ! The I/O path has a select code of 12 and is initialized to perform unformatted word
9365             ! transfers without any end of line designations. The DIP switches on the HP98622-66501 Rev B
9370             ! printed circuit board need to be set as shown below:
9375             ! DIP switches for INT LVL : Bit1=0 Bit0=0
9380             ! DIP switches for Select Code : Bit4=0 Bit3=1 Bit2=1 Bit1=0 Bit0=0
9385             ! DIP switches for DI15to08 clk: RDY =1 BSY =0 RD =1
9390             ! DIP switches for DI07to00 clk: RDY =1 BSY =0 RD =1
9395             ! DIP switches for Hndsk Levels: DOUT=0 DIN =0 HSHK=1 PSTS=0 PFLG=0 PCTL=1
9400             ASSIGN @Gpio TO 12;WORD,FORMAT OFF,EOL ""
9405             OUTPUT @Gpio USING "$,AA";"HP"
9410             SUBEND
9415 Lvdas_take: SUB Lvdas_take(@Lvdas,Atime,Ctime,INTEGER At_exp,Ct_exp,Cmask,Nsam)
9420             ! Description:
9425             ! This subprogram samples the two analog, three digital, and two external trigger channels
9430             ! from the LVDAS. The HP sends a "CS" to sample the LVDAS data with coincidence. Following the
9435             ! "CS" the HP sends the LVDAS an additional eight words to specify the acquisition and
9440             ! coincidence times, the inter-arrival and coincidence time exponents, the coincidence mask, and
9445             ! the number of desired samples. After the desired number of samples is acquired or the desired
9450             ! acquisition time expires then the LVDAS sends to the HP an updated number of samples (Nsam).
9455             ! The updated Nsam may be less than the original Nsam if the desired acquisition time expires
9460             ! before the desired Nsam samples are realized.
9465             ! Variables:
9470             ! Atime The maximum desired acquisition time (seconds).
9475             ! Ctime The maximum desired coincidence time (seconds).
9480             ! At1 The upper word of integer of 10000000*Atime.
9485             ! At2 The lower word of integer of 10000000*Atime.
9490             ! Ct1 The upper word of integer of 10000000*Ctime.
9495             ! Ct2 The lower word of integer of 10000000*Ctime.
9500             ! At_exp Exponent for inter-arrival times.
9505             ! Ct_exp Exponent for coincidence times.
9510             ! Nsam Number of desired samples.
9515             ! Cmask Coincidence Mask for U,V,W selection.
9520             ! Raw(*) Array of raw data acquired LVDAS data.
9525             OPTION BASE 1
9530             COM /Data1/ REAL Table(*),INTEGER Raw(*),Valid(*)
9535             INTEGER At1,At2,Ct1,Ct2
9540             DISP "Taking Data"
9545             CALL Convert2words(Atime*10000000,At1,At2)
9550             CALL Convert2words(Ctime*10000000,Ct1,Ct2)
9555             OUTPUT @Lvdas USING "AA,8(W)";"CS",At1,At2,Ct1,Ct2,At_exp,Ct_exp,Cmask,Nsam
9560             ENTER @Lvdas USING "$,W";Nsam
9565             IF Nsam=0 THEN SUBEXIT
9570             REDIM Raw(1:Nsam,1:10)
9575             ENTER @Lvdas USING "$,W";Raw(*)
9580             SUBEND
9585 Lvdas_sample: SUB Lvdas_sample(@Lvdas,Channel,Table(*),REAL Vave,Vsdv,Tave,Tsdv)
9590             ! Description:
9595             ! This subprogram samples one of the two analog, three digital, or two external trigger channels
9600             ! from the LVDAS. The HP sends the "DT","SC","RM", and "ET" commands to the LVDAS. The disable
9605             ! timer "DT" command tells the LVDAS to disable the LVDAS's internal timer interrupts. This
9610             ! prevents the LVDAS front panel displays from being updated but it also ensures that the data
9615             ! sampling will occur uninterrupted and at a maximum data rate. The sample channel "SC" tells the
9620             ! LVDAS to sample the specified channel and return 1000 data samples. Inter-arrival times are
9625             ! also returned. The read memory "RM" command reads back the data. The enable timer "ET"
9630             ! command enables the LVDAS's internal timer interrupts so that the front panel displays are
9635             ! updated.
9640             ! Variables:
9645             ! Channel Specifies one of the two analog, three digital, or two external trigger channels.
9650             ! Channel=0: Specifies the U digital channel.
9655             ! Channel=1: Specifies the V digital channel.
9660             ! Channel=2: Specifies the W digital channel.
9665             ! Channel=3: Specifies the A analog channel.
9670             ! Channel=4: Specifies the B analog channel.
9675             ! Channel=5: Specifies the External Trigger Timer channel.
9680             ! Channel=6: Specifies the Inter-arrival Timer channel.
9685             ! Data(*) Array of raw analog or digital data with inter-arrival time data.
9690             ! Data(*,1) Upper word of inter-arrival time data.

```

```

9695      Data(*,2)  Lower word of inter-arrival time data.
9700      Data(*,3)  Channel number for the data sampled.
9705      Data(*,4)  Data of the channel sampled.
9710      V(*)       Array of data for the channel sampled.
9715      Vv(*)      Squares of the V data array.
9720      T(*)       Array of inter-arrival times for the channel sampled.
9725      Tt(*)      Squares of the T inter-arrival time array.
9730      Vave       Average value of the channel's data.
9735      Vsdv       Standard deviation of the channel's data.
9740      Tave       Average value of the channel's inter-arrival time data.
9745      Tsdv       Standard deviation of the channel's inter-arrival time data.
9750      OPTION BASE 1
9755      INTEGER Data(1000,4),V(1000),Vv(1000),T(1000),Tt(1000)
9760      OUTPUT @Lvdas USING "#,AA";"DT"
9765      OUTPUT @Lvdas USING "#,AA,W";"SC",Channel+1      ! LVDAS expects to see 1 to 7, not 0 to 6.
9770      OUTPUT @Lvdas USING "AA";"RM"
9775      OUTPUT @Lvdas USING "W,W";IVAL("08F0",16),IVAL("0000",16)
9780      OUTPUT @Lvdas USING "W,W";IVAL("08F0",16),IVAL("1F3F",16)
9785      ENTER @Lvdas USING "#,W";Data(*)
9790      OUTPUT @Lvdas USING "#,AA";"ET"
9795      N=SIZE(Data,1)
9800      Channel=Data(1,3)
9805      SELECT Channel
9810      CASE 0,1,2      ! Convert raw digital data to frequencies.
9815          FOR I=1 TO N
9820              V(I)=Table(BINAND(32767,BINCMP(Data(I,4))))
9825          NEXT I
9830      CASE 3,4      ! Convert raw analog data to voltages.
9835          MAT V= Data(*,4)
9840          MAT V= V*(5/32768)
9845      CASE 5,6      ! The external trigger channels have no data.
9850          MAT V= (0)
9855      END SELECT
9860      MAT Vv= V . V
9865      MAT T= Data(*,2)
9870      MAT T= T/(100000000)
9875      MAT Tt= T . T
9880      Vave=SUM(V)/N
9885      Tave=SUM(T)/N
9890      Vsdv=SQR(ABS(SUM(Vv)/N-Vave*Vave))
9895      Tsdv=SQR(ABS(SUM(Tt)/N-Tave*Tave))
9900      MAT SEARCH Data(*,1),#LOC(<0);Bad1
9905      MAT SEARCH Data(*,2),#LOC(<0);Bad2
9910      SUBEXIT
9915      PRINT USING 9920;Channel,Vave,Vsdv,Tave,Tsdv,Bad1,Bad2
9920      IMAGE 4D,2(M8D.4D),2(M2D.6D),10X,2(5D)
9925      SUBEND
9930 Data:
9935 Data_reducel: SUB Data_reducel(INTEGER At_exp,Ct_exp,Nsam)
9940      ! Description:
9945      ! This subprogram separates the ten by Nsam Raw(*) data array into multiple one by Nsam
9950      ! arrays. The frequency arrays Ui,Vi,Wi are extracted from columns 6,7,8 of the Raw data array. The
9955      ! voltage arrays Ai & Bi are extracted from columns 9 & 10 of the Raw data array. The coincidence
9960      ! inter-arrival time array Ii is extracted from column 1 of the Raw data array. The validation word array
9965      ! time array Ci is extracted from column 2 of the Raw data array. The validation word array
9970      ! Valid is extracted from column 5 of the Raw data array. If i'th sample acquired contains
9975      ! valid data, then Valid(i) will be equal to one, and zero otherwise. All values for the Valid
9980      ! array are initially set to one by the LVDAS.
9985      ! The raw data from arrays Ui,Vi,Wi are converted into frequencies by using their initial
9990      ! values as indexes to the frequency look up table array Table(*). The raw data from arrays
9995      ! Ai & Bi are converted into voltages by multiplying their initial values by 5 volts over 2^15.
10000      ! The raw data from array Ii are converted into inter-arrival times by multiplying their initial
10005      ! values by 2^At_exp over 10 to get us. The raw data from array Ci are converted into
10010      ! coincidence times by multiplying their initial values by 2^Ct_exp over 10 to get us.
10015      ! Variables:
10020      ! Table(*)  Lookup table of frequencies.
10025      ! Raw(*)   Array of raw data acquired LVDAS data.
10030      ! Ui(*)   Array of extracted raw U frequency data.
10035      ! Vi(*)   Array of extracted raw V frequency data.
10040      ! Wi(*)   Array of extracted raw W frequency data.
10045      ! Ai(*)   Array of extracted raw A voltage data.
10050      ! Bi(*)   Array of extracted raw B voltage data.
10055      ! Ii(*)   Array of extracted raw inter-arrival time data.
10060      ! Ci(*)   Array of extracted raw coincidence time data.
10065      ! Valid(*) Array of extracted raw validation words.
10070      ! At_exp  Exponent of inter-arrival times.
10075      ! Ct_exp  Exponent of coincidence times.
10080      ! Nsam    Number of samples acquired.
10085      OPTION BASE 1
10090      COM /Data/ REAL Table(*),INTEGER Raw(*),Valid(*)

```

```

10095      COM /Data2/ REAL U1(*),V1(*),W1(*),A1(*),B1(*),I1(*),C1(*)
10100      REDIM U1(Nsam),V1(Nsam),W1(Nsam),A1(Nsam),B1(Nsam),I1(Nsam),C1(Nsam),Valid(Nsam)
10105      DISP "Reducing Data"
10110      MAT Valid= Raw(*,5)
10115      MAT I1= Raw(*,1)      ! Extract the inter-arrival times from the raw data array.
10120      MAT C1= Raw(*,2)      ! Extract the coincidence times from the raw data array.
10125      MAT U1= Raw(*,6)      ! Extract the instantaneous U velocities from the raw data array.
10130      MAT V1= Raw(*,7)      ! Extract the instantaneous V velocities from the raw data array.
10135      MAT W1= Raw(*,8)      ! Extract the instantaneous W velocities from the raw data array.
10140      MAT A1= Raw(*,9)      ! Extract the instantaneous A analog voltages from the raw data array.
10145      MAT B1= Raw(*,10)     ! Extract the instantaneous B analog voltages from the raw data array.
10150      FOR K=1 TO Nsam
10155          U1(K)=Table(U1(K)) ! The raw data of U1 is used to index the frequency lookup table.
10160          V1(K)=Table(V1(K)) ! The raw data of V1 is used to index the frequency lookup table.
10165          W1(K)=Table(W1(K)) ! The raw data of W1 is used to index the frequency lookup table.
10170      NEXT K
10175      MAT A1= A1*(5/32768)    ! The raw data for A1 is converted into a voltage (+/- 5 volts).
10180      MAT B1= B1*(5/32768)    ! The raw data for B1 is converted into a voltage (+/- 5 volts).
10185      MAT I1= I1*(2^At_exp/10) ! The raw data for I1 is converted into the inter-arrival time.
10190      MAT C1= C1*(2^Ct_exp/10) ! The raw data for C1 is converted into the coincidence time.
10195      SUBEND
10200 Data_reduce2: SUB Data_reduce2(Array(*))
10205      ! Description:
10210      ! This subprogram takes the frequency values from the arrays U1,V1,W1 and replaces them with
10215      ! velocities after doing the frequency to velocity conversion.
10220      ! Variables:
10225      ! Array(*)      An array containing relevant LDV laser and tunnel condition parameters
10230      ! Frng_spc(*)    Fringe Spacings extracted from Array(*).
10235      ! Brq_frq(*)     Bragg Frequencies extracted from Array(*).
10240      ! Mix_frq(*)     Mixing Freqs. extracted from Array(*).
10245      ! Mea_sgn(*)     Measured Freq's. Signs extracted from Array(*)
10250      ! Brq_sgn(*)     Bragg Freq's. Signs extracted from Array(*).
10255      ! Mix_sgn(*)     Mixing Freq's. Signs extracted from Array(*).
10260      ! U1(*)          Array of instantaneous U data.
10265      ! V1(*)          Array of instantaneous V data.
10270      ! W1(*)          Array of instantaneous W data.
10275      ! Equations:
10280      ! The following equations are used to convert the frequencies to velocities
10285      ! Velocity = Fs * Ftotal
10290      ! Ftotal = MeaSgn*Fmeas+BrqSgn*Fbrag+MixSgn*Fmix
10295      OPTION BASE 1
10300      COM /Data1/ REAL Table(0:32766),INTEGER Raw(*),Valid(*)
10305      COM /Data2/ REAL U1(*),V1(*),W1(*),A1(*),B1(*),I1(*),C1(*)
10310      DIM Frng_spc(3),Brq_frq(3),Mix_frq(3),Mea_sgn(3),Brq_sgn(3),Mix_sgn(3)
10315      DISP "Converting Data"
10320      MAT Frng_spc= Array(35,1:3)
10325      MAT Brq_frq= Array(36,1:3)
10330      MAT Mix_frq= Array(37,1:3)
10335      MAT Mea_sgn= Array(38,1:3)
10340      MAT Brq_sgn= Array(39,1:3)
10345      MAT Mix_sgn= Array(40,1:3)
10350      MAT U1= U1*(Mea_sgn(1))
10355      MAT V1= V1*(Mea_sgn(2))
10360      !MAT W1= W1*(Mea_sgn(3))
10365      MAT U1= U1+(Brq_sgn(1)*Brq_frq(1)+Mix_sgn(1)*Mix_frq(1))
10370      MAT V1= V1+(Brq_sgn(2)*Brq_frq(2)+Mix_sgn(2)*Mix_frq(2))
10375      !MAT W1= W1+(Brq_sgn(3)*Brq_frq(3)+Mix_sgn(3)*Mix_frq(3))
10380      MAT U1= U1*(Frng_spc(1))
10385      MAT V1= V1*(Frng_spc(2))
10390      !MAT W1= W1*(Frng_spc(3))
10395      MAT W1= (0)
10400      SUBEND
10405 Data_reduce3: SUB Data_reduce3(Gain)
10410      ! Description:
10415      ! This subprogram takes the voltage values from the array A1 and replaces them with the total
10420      ! temperature after doing the voltage to temperature conversion.
10425      ! Variables:
10430      ! A1(*)          Array of instantaneous A data.
10435      ! Nsam           Number of acquired samples.
10440      ! N              Exponent for the terms in the polynomial equations.
10445      ! An             Coefficients for the terms in the polynomial equations.
10450      ! Gain           Gain for the analog channels voltage.
10455      ! Mv(*)          Array of gained raw voltages converted to millivolts.
10460      ! Mvn(*)         Array of Mv(*) values raised to the power of N.
10465      ! Amvn(*)        Array of Mvn(*) values multiplied by the polynomial coefficients An.
10470      ! Sum(*)         Summation of the terms of polynomial equation.
10475      ! Equations:
10480      ! The following equations are used to convert the voltages to temperatures.
10485      ! Temp=A7*A1^7 + A6*A1^6 + ... + A0*A1^0 + 460
10490      DISP "Converting Data"

```

```

10495      OPTION BASE 1
10500      COM /Data2/ REAL U1(*),V1(*),W1(*),A1(*),B1(*),I1(*),C1(*)
10505      DIM Mv(1000),Mvn(1000),Amvn(1000),Sum(1000)
10510      Nsam=SIZE(A1,1)
10515      REDIM Mv(Nsam),Mvn(Nsam),Amvn(Nsam),Sum(Nsam)
10520      MAT Mv= A1*(1000/Gain)      ! Tt_mv=Tt_raw/Gain*1000
10525      MAT Sum= (0)
10530      MAT Mvn= (1)
10535      FOR N=0 TO 7
10540          READ An
10545          MAT Amvn= (An)*Mvn
10550          MAT Sum= Sum+Amvn
10555          MAT Mvn= Mvn . Mv
10560      NEXT N
10565      MAT Ai= Sum+(460)
10570      SUBEXIT
10575      !      A0,      A1,      A2,      A3,      A4,      A5,      A6,      A7
10580      DATA 150,257.10163,-28.16138,6.064559,-.792687,.05708673,-.002103462,.00003110036
10585      SUBEND
10590 Data_reduce4: SUB Data_reduce4(Mach,Mv,Ts,Tt)
10595      ! Description:
10600      ! This subprogram takes the analog voltage and converts it to a temperature.
10605      ! Variables:
10610      ! Mach      Mach number.
10615      ! Tt      Total Temperature in degrees Rankine.
10620      ! Ts      Stagnation Temperature in degrees Rankine.
10625      ! N      Exponent for the terms in the polynomial equations.
10630      ! An      Coefficients for the terms in the polynomial equations.
10635      ! Mv      The gained raw voltage converted to millivolts.
10640      ! Equations:
10645      ! The following equations are used to convert the voltages to temperatures.
10650      ! Temp=A7*A1^7 + A6*A1^6 + ... + A0*A1^0 + 460
10655      Tt=0
10660      FOR N=0 TO 7
10665          READ An
10670          Tt=Tt+An*Mv^N
10675      NEXT N
10680      Tt=Tt+460
10685      Ts=.09259*Tt
10690      IF Mach<>7 THEN BEEP
10695      IF Mach<>7 THEN PAUSE
10700      SUBEXIT
10705      !      A0,      A1,      A2,      A3,      A4,      A5,      A6,      A7
10710      DATA 150,257.10163,-28.16138,6.064559,-.792687,.05708673,-.002103462,.00003110036
10715      SUBEND
10720 Data_clip: SUB Data_clip(INTEGER Nsam,REAL Umin,Umax,Vmin,Vmax,Wmin,Wmax)
10725      ! Description:
10730      ! This subprogram compares each of the instantaneous U,V,W frequencies with user selectable
10735      ! minimum and maximum frequencies. If the instantaneous value is less than the desired
10740      ! minimum, then the validation word is set to zero. Also, if the instantaneous value is
10745      ! greater than the desired maximum, then the validation word is set to zero. The setting of the
10750      ! validation words to zero will have the net effect of discarding the data samples from the data
10755      ! set. In other words, the data are weighted as zero for the average, sdv, shear stress, and
10760      ! cross correlation calculations.
10765      ! Variables:
10770      ! Nsam      Number of samples acquired.
10775      ! U1(*)      Array of instantaneous U frequencies (MHz).
10780      ! V1(*)      Array of instantaneous V frequencies (MHz).
10785      ! W1(*)      Array of instantaneous W frequencies (MHz).
10790      ! Valid(*)   Array of sample validation words.
10795      ! Umin      The minimum acceptable U frequency (MHz).
10800      ! Umax      The maximum acceptable U frequency (MHz).
10805      ! Vmin      The minimum acceptable V frequency (MHz).
10810      ! Vmax      The maximum acceptable V frequency (MHz).
10815      ! Wmin      The minimum acceptable W frequency (MHz).
10820      ! Wmax      The maximum acceptable W frequency (MHz).
10825      OPTION BASE 1
10830      COM /Data1/ REAL Table(0:32766),INTEGER Raw(*),Valid(*)
10835      COM /Data2/ REAL U1(*),V1(*),W1(*),A1(*),B1(*),I1(*),C1(*)
10840      DISP "Clipping Histograms"
10845      FOR K=1 TO Nsam
10850          MAT SEARCH U1(*),LOC(<Umin);L,K
10855          IF L<Nsam THEN Valid(L)=0
10860          K=L
10865      NEXT K
10870      FOR K=1 TO Nsam
10875          MAT SEARCH U1(*),LOC(>Umax);L,K
10880          IF L<Nsam THEN Valid(L)=0
10885          K=L
10890      NEXT K

```

```

10895     FOR K=1 TO Nsam
10900         MAT SEARCH Vi(*),LOC(<Vmin);L,K
10905         IF L<Nsam THEN Valid(L)=0
10910         K=L
10915     NEXT K
10920     FOR K=1 TO Nsam
10925         MAT SEARCH Vi(*),LOC(>Vmax);L,K
10930         IF L<Nsam THEN Valid(L)=0
10935         K=L
10940     NEXT K
10945     !FOR K=1 TO Nsam
10950     !     MAT SEARCH Wi(*),LOC(<Wmin);L,K
10955     !     IF L<Nsam THEN Valid(L)=0
10960     !     K=L
10965     !NEXT K
10970     !FOR K=1 TO Nsam
10975     !     MAT SEARCH Wi(*),LOC(>Wmax);L,K
10980     !     IF L<Nsam THEN Valid(L)=0
10985     !     K=L
10990     !NEXT K
10995 SUBEND
11000 Data_sum: SUB Data_sum(INTEGER Nsam)
11005     ! Description:
11010     !     This subprogram performs the summations on the instantaneous LDV and analog data. Data
11015     !     will be weighted as zero in the summations if the value of the validation word is set to zero.
11020     !     Intermediate arrays will be made so that summations of the products of the LDV and analog data
11025     !     can be determined.
11030     ! Variables:
11035     !     Nsam      Number of samples acquired.
11040     !     Valid(*)  Array of sample validation words.
11045     !     Ui(*)     Array of instantaneous U frequency or velocity samples.
11050     !     Vi(*)     Array of instantaneous V frequency or velocity samples.
11055     !     Wi(*)     Array of instantaneous W frequency or velocity samples.
11060     !     Ai(*)     Array of instantaneous A voltage samples.
11065     !     Bi(*)     Array of instantaneous B voltage samples.
11070     !     Ii(*)     Array of inter-arrival times.
11075     !     Ci(*)     Array of coincidence times.
11080     !     Puu(*)    Instantaneous product of the instantaneous Ui & Ui.
11085     !     Pvv(*)    Instantaneous product of the instantaneous Vi & Vi.
11090     !     Pww(*)    Instantaneous product of the instantaneous Wi & Wi.
11095     !     Paa(*)    Instantaneous product of the instantaneous Ai & Ai.
11100     !     Pbb(*)    Instantaneous product of the instantaneous Bi & Bi.
11105     !     Pii(*)    Instantaneous product of the instantaneous Ii & Ii.
11110     !     Pcc(*)    Instantaneous product of the instantaneous Ci & Ci.
11115     !     Puv(*)    Instantaneous product of the instantaneous Ui & Vi.
11120     !     Pvw(*)    Instantaneous product of the instantaneous Vi & Wi.
11125     !     Pwu(*)    Instantaneous product of the instantaneous Wi & Ui.
11130     !     Pab(*)    Instantaneous product of the instantaneous Ai & Bi.
11135     !     Pua(*)    Instantaneous product of the instantaneous Ui & Ai.
11140     !     Pva(*)    Instantaneous product of the instantaneous Vi & Ai.
11145     !     Pwa(*)    Instantaneous product of the instantaneous Wi & Ai.
11150     !     Sumu     Summation of the array Ui.
11155     !     Sumv     Summation of the array Vi.
11160     !     Sumw     Summation of the array Wi.
11165     !     Suma     Summation of the array Ai.
11170     !     Sumb     Summation of the array Bi.
11175     !     Sumi     Summation of the array Ii.
11180     !     Sumc     Summation of the array Ci.
11185     !     Sumuu    Summation of the array Puu.
11190     !     Sumvv    Summation of the array Pvv.
11195     !     Sumww    Summation of the array Pww.
11200     !     Sumaa    Summation of the array Paa.
11205     !     Sumbb    Summation of the array Pbb.
11210     !     Sumii    Summation of the array Pii.
11215     !     Sumcc    Summation of the array Pcc.
11220     !     Sumuv    Summation of the array Puv.
11225     !     Sumvw    Summation of the array Pvw.
11230     !     Sumwu    Summation of the array Pwu.
11235     !     Sumab    Summation of the array Pab.
11240     !     Sumua    Summation of the array Pua.
11245     !     Sumva    Summation of the array Pva.
11250     !     Sumwa    Summation of the array Pwa.
11255     !     Suml     Number of valid samples acquired.
11260     OPTION BASE 1
11265     COM /Data1/ REAL Table(0:32766),INTEGER Raw(*),Valid(*)
11270     COM /Data2/ REAL Ui(*),Vi(*),Wi(*),Ai(*),Bi(*),Ii(*),Ci(*)
11275     COM /Data3/ REAL Puu(*),Pvv(*),Pww(*),Paa(*),Pbb(*),Pii(*),Pcc(*)
11280     COM /Data4/ REAL Puv(*),Pvw(*),Pwu(*),Pab(*),Pua(*),Pva(*),Pwa(*)
11285     COM /Sum1/ REAL Sumu,Sumv,Sumw,Suma,Sumb,Sumi,Sumc,Suml
11290     COM /Sum2/ REAL Sumuu,Sumvv,Sumww,Sumaa,Sumbb,Sumii,Sumcc

```

```

11295 COM /Sum3/ REAL Sumuv, Sumvw, Sumwu, Sumab, Sumua, Sumva, Sumwa
11300 REDIM Puu(Nsam), Pvv(Nsam), Pww(Nsam), Paa(Nsam), Pbb(Nsam), Pii(Nsam), Pcc(Nsam)
11305 REDIM Puv(Nsam), Pvw(Nsam), Pwu(Nsam), Pab(Nsam), Pua(Nsam), Pva(Nsam), Pwa(Nsam)
11310 DISP "Summing Data"
11315 MAT Ui= Ui . Valid ! Ui(I) is the instantaneous velocity Ui(I).
11320 MAT Vi= Vi . Valid ! Vi(I) is the instantaneous velocity Vi(I).
11325 MAT Wi= Wi . Valid ! Wi(I) is the instantaneous velocity Wi(I).
11330 MAT Ai= Ai . Valid ! Ai(I) is the instantaneous Channel #1 Analog Voltage Ai(I).
11335 MAT Bi= Bi . Valid ! Bi(I) is the instantaneous Channel #2 Analog Voltage Bi(I).
11340 MAT Ii= Ii . Valid ! Ii(I) is the inter-arrival time Ii(I).
11345 MAT Ci= Ci . Valid ! Ci(I) is the coincidence time Ci(I).
11350 MAT Puu= Ui . Ui ! Puu(I) is the square of the instantaneous velocity Ui(I).
11355 MAT Pvv= Vi . Vi ! Pvv(I) is the square of the instantaneous velocity Vi(I).
11360 MAT Pww= Wi . Wi ! Pww(I) is the square of the instantaneous velocity Wi(I).
11365 MAT Paa= Ai . Ai ! Paa(I) is the square if the instantaneous Analog Voltage Ai(I).
11370 MAT Pbb= Bi . Bi ! Pbb(I) is the square if the instantaneous Analog Voltage Bi(I).
11375 MAT Pii= Ii . Ii ! Pii(I) is the square if the inter-arrival time Ii(I).
11380 MAT Pcc= Ci . Ci ! Pcc(I) is the square if the coincidence time Ci(I).
11385 MAT Puv= Ui . Vi ! Puv(I) is the product of Ui(I) and Vi(I).
11390 MAT Pvw= Vi . Wi ! Pvw(I) is the product of Vi(I) and Wi(I).
11395 MAT Pwu= Wi . Ui ! Pwu(I) is the product of Wi(I) and Ui(I).
11400 MAT Pab= Ai . Bi ! Pab(I) is the product of Ai(I) and Bi(I).
11405 MAT Pua= Ui . Ai ! Pua(I) is the product of Ui(I) and Ai(I).
11410 MAT Pva= Vi . Ai ! Pva(I) is the product of Vi(I) and Ai(I).
11415 MAT Pwa= Wi . Ai ! Pwa(I) is the product of Wi(I) and Ai(I).
11420 Sumu=SUM(Ui) ! Sumu is the summation of Ui(I).
11425 Sumv=SUM(Vi) ! Sumv is the summation of Vi(I).
11430 Sumw=SUM(Wi) ! Sumw is the summation of Wi(I).
11435 Suma=SUM(Ai) ! Suma is the summation of Ai(I).
11440 Sumb=SUM(Bi) ! Sumb is the summation of Bi(I).
11445 Sumi=SUM(Ii) ! Sumi is the summation of Ii(I).
11450 Sumc=SUM(Ci) ! Sumc is the summation of Ci(I).
11455 Sumuu=SUM(Puu) ! Sumuu is the summation of Ui(I)*Ui(I).
11460 Sumvv=SUM(Pvv) ! Sumvv is the summation of Vi(I)*Vi(I).
11465 Sumww=SUM(Pww) ! Sumww is the summation of Wi(I)*Wi(I).
11470 Sumaa=SUM(Paa) ! Sumaa is the summation of Ai(I)*Ai(I).
11475 Sumbb=SUM(Pbb) ! Sumbb is the summation of Bi(I)*Bi(I).
11480 Sumii=SUM(Pii) ! Sumii is the summation of Ii(I)*Ii(I).
11485 Sumcc=SUM(Pcc) ! Sumcc is the summation of Ci(I)*Ci(I).
11490 Sumuv=SUM(Puv) ! Sumuv is the summation of Ui(I)*Vi(I).
11495 Sumvw=SUM(Pvw) ! Sumvw is the summation of Vi(I)*Wi(I).
11500 Sumwu=SUM(Pwu) ! Sumwu is the summation of Wi(I)*Ui(I).
11505 Sumab=SUM(Pab) ! Sumab is the summation of Ai(I)*Bi(I).
11510 Sumua=SUM(Pua) ! Sumua is the summation of Ui(I)*Ai(I).
11515 Sumva=SUM(Pva) ! Sumva is the summation of Vi(I)*Ai(I).
11520 Sumwa=SUM(Pwa) ! Sumwa is the summation of Wi(I)*Ai(I).
11525 Suml=SUM(Valid) ! Suml is the number of valid samples.
11530 SUBEND
11535 Data_calc: SUB Data_calc
11540 ! Description:
11545 ! This subprogram uses the summations on the instantaneous LDV and analog data as well as the
11550 ! summations of the products of the LDV and analog data. The subprogram takes these summations
11555 ! and calculates the averages, standard deviations, and shear stresses.
11560 ! Variables:
11565 ! Suml Number of valid samples acquired.
11570 ! Sumu Summation of the array Ui.
11575 ! Sumv Summation of the array Vi.
11580 ! Sumw Summation of the array Wi.
11585 ! Suma Summation of the array Ai.
11590 ! Sumb Summation of the array Bi.
11595 ! Sumi Summation of the array Ii.
11600 ! Sumc Summation of the array Ci.
11605 ! Sumuu Summation of the array Puu.
11610 ! Sumvv Summation of the array Pvv.
11615 ! Sumww Summation of the array Pww.
11620 ! Sumaa Summation of the array Paa.
11625 ! Sumbb Summation of the array Pbb.
11630 ! Sumii Summation of the array Pii.
11635 ! Sumcc Summation of the array Pcc.
11640 ! Sumuv Summation of the array Puv.
11645 ! Sumvw Summation of the array Pvw.
11650 ! Sumwu Summation of the array Pwu.
11655 ! Sumab Summation of the array Pab.
11660 ! Sumua Summation of the array Pua.
11665 ! Sumva Summation of the array Pva.
11670 ! Sumwa Summation of the array Pwa.
11675 ! N Number of valid samples acquired.
11680 ! U Average U frequency or velocity.
11685 ! V Average V frequency or velocity.
11690 ! W Average W frequency or velocity.

```

```

11695      !      A      Average A voltage.
11700      !      B      Average B voltage.
11705      !      I      Average inter-arrival time.
11710      !      C      Average coincidence time.
11715      !      U1     Standard deviation for U frequency or velocity.
11720      !      V1     Standard deviation for V frequency or velocity.
11725      !      W1     Standard deviation for W frequency or velocity.
11730      !      A1     Standard deviation for A voltage.
11735      !      B1     Standard deviation for B voltage.
11740      !      I1     Standard deviation for inter-arrival time.
11745      !      C1     Standard deviation for coincidence time.
11750      !      U1v1   Velocity:Velocity Shear Stress.
11755      !      V1w1   Velocity:Velocity Shear Stress.
11760      !      W1u1   Velocity:Velocity Shear Stress.
11765      !      Alb1   Voltage :Voltage Cross Correlation.
11770      !      U1a1   Velocity:Voltage Cross Correlation.
11775      !      V1a1   Velocity:Voltage Cross Correlation.
11780      !      W1a1   Velocity:Voltage Cross Correlation.
11785      COM /Sum1/ REAL Sumu, Sumv, Sumw, Suma, Sumb, Sumi, Sumc, Suml
11790      COM /Sum2/ REAL Sumuu, Sumvv, Sumww, Sumaa, Sumbb, Sumii, Sumcc
11795      COM /Sum3/ REAL Sumuv, Sumvw, Sumwu, Sumab, Sumua, Sumva, Sumwa
11800      COM /Reduced/ N,U,V,W,A,B,I,C,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,Alb1,U1a1,V1a1,W1a1
11805      DISP "Calculating Results"
11810      N=Sum1
11815      IF N>0 THEN
11820          U=Sumu/N
11825          V=Sumv/N
11830          W=Sumw/N
11835          A=Suma/N
11840          B=Sumb/N
11845          I=Sumi/N
11850          C=Sumc/N
11855          U1=SQR(ABS(Sumuu/N-U*U))
11860          V1=SQR(ABS(Sumvv/N-V*V))
11865          W1=SQR(ABS(Sumww/N-W*W))
11870          A1=SQR(ABS(Sumaa/N-A*A))
11875          B1=SQR(ABS(Sumbb/N-B*B))
11880          I1=SQR(ABS(Sumii/N-I*I))
11885          C1=SQR(ABS(Sumcc/N-C*C))
11890          U1v1=Sumuv/N-U*V
11895          V1w1=Sumvw/N-V*W
11900          W1u1=Sumwu/N-W*U
11905          Alb1=Sumab/N-A*B
11910          U1a1=Sumua/N-U*A
11915          V1a1=Sumva/N-V*A
11920          W1a1=Sumwa/N-W*A
11925      ELSE
11930          U=0
11935          V=0
11940          W=0
11945          A=0
11950          B=0
11955          I=0
11960          C=0
11965          U1=0
11970          V1=0
11975          W1=0
11980          A1=0
11985          B1=0
11990          I1=0
11995          C1=0
12000          U1v1=0
12005          V1w1=0
12010          W1u1=0
12015          Alb1=0
12020          U1a1=0
12025          V1a1=0
12030          W1a1=0
12035      END IF
12040      SUBEND
12045      Data_trnsfrm: SUB Data_trnsfrm(REAL K3x3(*),U,V,W,U1,V1,W1,U1v1,V1w1,W1u1,U1a1,V1a1,W1a1)
12050          !      Description:
12055          !          This subprogram performs a coordinate system transformation on the averages, standard
12060          !          deviations, and shear stresses. The coordinate system transformation to be applied is passed
12065          !          through the "K3X3" array. If a TUNNEL to MODEL coordinate system transformation is to be
12070          !          performed, then the array "Tun2mod" array will be passed to the "K3X3" array.
12075          !          NOTE: This sub-program performs a three dimensional coordinate system transformation on
12080          !          averages, standard deviations, shear stresses, and cross correlations. It performs this
12085          !          transformation for averages, standard deviations, shear stresses, and cross correlations that
12090          !          include one or more of the velocities U,V, or W. The delivered system is a two component

```



```

12095      ! system. Therefore, the third component W will have been set to be equal to zero. Other terms
12100      ! containing W are also set to zero by the main program. (W=W1=U1w1=W1v1=W1a1=0).
12105      ! Variables:
12110      ! U      Average U velocity.
12115      ! V      Average V velocity.
12120      ! W      Average W velocity.
12125      ! U1     Standard deviation for U velocity.
12130      ! V1     Standard deviation for V velocity.
12135      ! W1     Standard deviation for W velocity.
12140      ! U1v1   Velocity:Velocity Normal Stress.
12145      ! U1v1   Velocity:Velocity Shear Stress.
12150      ! U1w1   Velocity:Velocity Shear Stress.
12155      ! V1v1   Velocity:Velocity Shear Stress.
12160      ! V1v1   Velocity:Velocity Normal Stress.
12165      ! V1w1   Velocity:Velocity Shear Stress.
12170      ! W1v1   Velocity:Velocity Shear Stress.
12175      ! W1v1   Velocity:Velocity Shear Stress.
12180      ! W1w1   Velocity:Velocity Normal Stress.
12185      ! U1a1   Velocity:Voltage Cross Correlation.
12190      ! V1a1   Velocity:Voltage Cross Correlation.
12195      ! W1a1   Velocity:Voltage Cross Correlation.
12200      ! R(*)   Original U,V,W.
12205      ! F(*)   Original U1a1,V1a1,W1a1.
12210      ! P(*)   Original stress terms U1v1,U1w1,...,W1w1.
12215      ! K3X3   Coordinate system transformation matrix for average and Velocity:Voltage cross
12220      !          correlation conversions.
12225      ! K9X9   Coordinate system transformation matrix for Velocity:Velocity normal and shear
12230      !          stress conversions.
12235      ! S(*)   Transformed U,V,W.
12240      ! H(*)   Transformed U1a1,V1a1,W1a1.
12245      ! Q(*)   Transformed stress terms U1v1,U1w1,...,W1w1.
12250      ! OPTION BASE 1
12255      REAL R(3),S(3),F(3),H(3),P(9),Q(9),K9x9(9,9)
12260      DISP "Transforming Results"
12265      ! Calculate U1v1,V1v1,W1w1 using U1,V1,W1.
12270      U1v1=U1*U1
12275      V1v1=V1*V1
12280      W1w1=W1*W1
12285      ! Set U1w1,V1w1,W1v1 equal to W1v1,U1v1,V1w1.
12290      U1w1=W1v1
12295      V1w1=U1v1
12300      W1v1=V1w1
12305      ! Fill the matrix R with U,V,W.
12310      R(1)=U
12315      R(2)=V
12320      R(3)=W
12325      ! Fill the matrix F with U1a1,V1a1,W1a1.
12330      F(1)=U1a1
12335      F(2)=V1a1
12340      F(3)=W1a1
12345      ! Fill the matrix P with U1v1,U1w1,V1v1,V1w1,W1v1,W1w1.
12350      P(1)=U1v1
12355      P(2)=U1w1
12360      P(3)=U1w1
12365      P(4)=V1v1
12370      P(5)=V1w1
12375      P(6)=V1w1
12380      P(7)=W1v1
12385      P(8)=W1v1
12390      P(9)=W1w1
12395      ! Define the matrix K9x9 using products of the elements from then matrix K3x3.
12400      FOR X=1 TO 9
12405          FOR Y=1 TO 9
12410              Y1=((Y-1) DIV 3)+1
12415              X1=((X-1) DIV 3)+1
12420              Y2=((Y-1) MOD 3)+1
12425              X2=((X-1) MOD 3)+1
12430              K9x9(Y,X)=K3x3(Y1,X1)*K3x3(Y2,X2)
12435          NEXT Y
12440      NEXT X
12445      ! Transform matrix R to S using K3x3.
12450      MAT S= K3x3*R
12455      ! Transform matrix F to H using K3x3.
12460      MAT H= K3x3*F
12465      ! Transform matrix P to Q using K9x9.
12470      MAT Q= K9x9*P
12475      ! Extract the transformed U,V,W from the matrix S.
12480      U=S(1)
12485      V=S(2)
12490      W=S(3)

```

```

12495      ! Extract the transformed U1a1,V1a1,W1a1 from the matrix H.
12500      U1a1=H(1)
12505      V1a1=H(2)
12510      W1a1=H(3)
12515      ! Extract the transformed U1u1,U1v1,U1w1,V1u1,V1v1,V1w1,W1u1,W1v1,W1w1 from the matrix Q.
12520      U1u1=Q(1)
12525      U1v1=Q(2)
12530      U1w1=Q(3)
12535      V1u1=Q(4)
12540      V1v1=Q(5)
12545      V1w1=Q(6)
12550      W1u1=Q(7)
12555      W1v1=Q(8)
12560      W1w1=Q(9)
12565      ! Calculate U1,V1,W1 using U1u1,V1v1,W1w1.
12570      U1=SQR(ABS(U1u1))
12575      V1=SQR(ABS(V1v1))
12580      W1=SQR(ABS(W1w1))
12585      ! Return transformed U,V,W,U1,V1,W1,U1v1,V1w1,W1u1,U1a1,V1a1,W1a1 to main program.
12590      SUBEND
12595 Data_print1: SUB Data_print1(Run,File,Pos(*),CS)
12600      ! Description:
12605      ! This subprogram prints the averages, standard deviations, shear stresses, and cross
12610      ! correlations in tabular form. This subprogram prints the reduced velocity data when their
12615      ! units are in frequency (MHz).
12620      ! Variables:
12625      ! U      Average U frequency (MHz).
12630      ! V      Average V frequency (MHz).
12635      ! W      Average W frequency (MHz).
12640      ! A      Average A voltage.
12645      ! B      Average B voltage.
12650      ! I      Average inter-arrival time (us).
12655      ! C      Average coincidence time (us).
12660      ! U1     Standard deviation for U frequencies (MHz).
12665      ! V1     Standard deviation for V frequencies (MHz).
12670      ! W1     Standard deviation for W frequencies (MHz).
12675      ! A1     Standard deviation for A voltages.
12680      ! B1     Standard deviation for B voltages.
12685      ! I1     Standard deviation for inter-arrival times (us).
12690      ! C1     Standard deviation for coincidence times (us).
12695      ! U1v1   Velocity:Velocity Shear Stress.
12700      ! V1w1   Velocity:Velocity Shear Stress.
12705      ! W1u1   Velocity:Velocity Shear Stress.
12710      ! A1b1   Voltage :Voltage Cross Correlation.
12715      ! U1a1   Velocity:Voltage Cross Correlation.
12720      ! V1a1   Velocity:Voltage Cross Correlation.
12725      ! W1a1   Velocity:Voltage Cross Correlation.
12730      ! Axis   Indicates one of the three axes X,Y,Z being traversed.
12735      ! Pos(*)  Current Traverse Positions.
12740      ! N      Number of valid samples acquired.
12745      ! CS     Indicates units and/or coordinate system of data printed.
12750      OPTION BASE 1
12755      COM /Reduced/ N,U,V,W,A,B,I,C,U1,V1,W1,A1,B1,I1,C1,U1v1,V1w1,W1u1,A1b1,U1a1,V1a1,W1a1
12760      DISP "Printing Results"
12765      ON ERROR CALL Error
12770      PRINTER IS PRT;WIDTH 144
12775      LS=CHR$(NUM("X")+Axis-1)&"="
12780      PRINT USING 12820;"Xtun=",Pos(1),"in      U=",U,"MHz      U'V'=",U1v1,"
          U'A'=",U1a1,"      CT'=",C,"us"
12785      PRINT USING 12825;"Ytun=",Pos(2),"in      V=",V,"MHz      V'W'=",V1w1,"
          V'A'=",V1a1,"      IAT'=",I,"us"
12790      PRINT USING 12830;"Ztun=",Pos(3),"in      W=",W,"MHz      W'U'=",W1u1,"
          W'A'=",W1a1,"      CT'=",C1,"us"
12795      PRINT USING 12835;"Run =",Run,"      A=",A,"v      A'=",A1,"v      ",
          IAT'=",I1,"us"
12800      PRINT USING 12840;"File=",File,"      B=",B,"v      B'=",B1,"v      A'B'=",A1b1,"
          N      ="N,""
12805      PRINT
12810      PRINTER IS CRT
12815      OFF ERROR
12820      IMAGE      8X, K,3D.4D,      K,5D.3D,      K,5D.3D,      K,8D.2D,      K,6D.5D,      K,9D,K
12825      IMAGE      8X, K,3D.4D,      K,5D.3D,      K,5D.3D,      K,8D.2D,      K,6D.5D,      K,9D,K
12830      IMAGE      8X, K,3D.4D,      K,5D.3D,      K,5D.3D,      K,8D.2D,      K,6D.5D,      K,9D,K
12835      IMAGE      8X, K,5D.2D,      K,5D.3D,      K,5D.3D,      K,11X ,      K,12X ,      K,9D,K
12840      IMAGE      8X, K,5D.2D,      K,5D.3D,      K,5D.3D,      K,8D.2D,      K,12X ,      K,9D,K
12845      SUBEND
12850 Data_print2: SUB Data_print2(Run,File,Pos(*),CS)
12855      ! Description:
12860      ! This subprogram prints the averages, standard deviations, and shear stresses, and cross
12865      ! correlations in tabular form. This subprogram prints the reduced velocity data when their

```

```

12870      !           units are in m/s.
12875      ! Variables:
12880      !     U           Average U velocity.
12885      !     V           Average V velocity.
12890      !     W           Average W velocity.
12895      !     A           Average A voltage.
12900      !     B           Average B voltage.
12905      !     I           Average inter-arrival time.
12910      !     C           Average coincidence time.
12915      !     U1          Standard deviation for U velocities (m/s).
12920      !     V1          Standard deviation for V velocities (m/s)..
12925      !     W1          Standard deviation for W velocities (m/s)..
12930      !     A1          Standard deviation for A voltages.
12935      !     B1          Standard deviation for B voltages.
12940      !     I1          Standard deviation for inter-arrival times (us).
12945      !     C1          Standard deviation for coincidence times (us).
12950      !     Ulvl        Velocity:Velocity Shear Stress.
12955      !     V1wl        Velocity:Velocity Shear Stress.
12960      !     W1ul        Velocity:Velocity Shear Stress.
12965      !     Albl        Voltage :Voltage Cross Correlation.
12970      !     U1al        Velocity:Voltage Cross Correlation.
12975      !     V1al        Velocity:Voltage Cross Correlation.
12980      !     W1al        Velocity:Voltage Cross Correlation.
12985      !     Axis        Indicates one of the three axes X,Y,Z being traversed.
12990      !     Pos(*)       Current Traverse Positions.
12995      !     N            Number of valid samples acquired.
13000      !     CS          Indicates units and/or coordinate system of data printed.
13005      OPTION BASE 1
13010      COM /Reduced/ N,U,V,W,A,B,I,C,U1,V1,W1,A1,B1,I1,C1,Ulvl,V1wl,W1ul,Albl,U1al,V1al,W1al
13015      DISP "Printing Results"
13020      ON ERROR CALL Error
13025      PRINTER IS PRT;WIDTH 144
13030      LS=CHRS(NUM("X")+Axis-1)&"="
13035      PRINT USING 13070;"Xmod=",Pos(1),"in      U=",U,"m/s      U'=",U1,"m/s      U'V'=",Ulvl,"
13040      U'A'=",U1al,"      CT'=",C,"us"
13045      PRINT USING 13075;"Ymod=",Pos(2),"in      V=",V,"m/s      V'=",V1,"m/s      V'W'=",V1wl,"
13050      V'A'=",V1al,"      IAT'=",I,"us"
13055      PRINT USING 13080;"Zmod=",Pos(3),"in      W=",W,"m/s      W'=",W1,"m/s      W'U'=",W1ul,"
13060      W'A'=",W1al,"      CT'=",C1,"us"
13065      PRINT USING 13085;"Run =",Run,"      A=",A,"v      A'=",A1,"v      ", "      "
13070      IAT'=",I1,"us"
13075      PRINT USING 13090;"File=",File,"      B=",B,"v      B'=",B1,"v      A'B'=",Albl,"      ", "
13080      N      ="N,"
13085      PRINTER IS CRT
13090      OFF ERROR
13095      IMAGE      8X, K,3D.4D,      K,5D.3D,      K,5D.3D,      K,8D.2D,      K,6D.5D,      K,9D,K
13100      IMAGE      8X, K,3D.4D,      K,5D.3D,      K,5D.3D,      K,8D.2D,      K,6D.5D,      K,9D,K
13105      IMAGE      8X, K,3D.4D,      K,5D.3D,      K,5D.3D,      K,8D.2D,      K,6D.5D,      K,9D,K
13110      IMAGE      8X, K,5D.2D,      K,5D.3D,      K,5D.3D,      K,11X ,      K,12X ,      K,9D,K
13115      IMAGE      8X, K,5D.2D,      K,5D.3D,      K,5D.3D,      K,8D.2D,      K,12X ,      K,9D,K
13120      SUBEND
13125      SUB Data_plot (Array(*),Symbols(*),Plot,Sy,Y,X)
13130      ! Description:
13135      ! This subprogram plots the averages, standard deviations, and shear stresses in the 4 profile
13140      ! plots on the CRT. This subprogram will typically be called up to 4 times for each of the
13145      ! four profile plots. The first profile plot will contain the average velocities and their
13150      ! standard deviations normalized by Uedge. The second profile plot will contain the average
13155      ! voltages and their standard deviations for the two analog channels. The third profile plot
13160      ! will contain the average temperature and its standard deviation. The forth and last profile
13165      ! plot will contain the velocity shear stress terms. Data points outside the plot boundaries
13170      ! will be plotted at the plot boundary.
13175      ! Variables:
13180      ! Array(*)      Array containing the plot positions and scales.
13185      ! Symbols(*)   Array of Symbol arrays. Each symbol array contains a distinct geometric symbol.
13190      ! Plot         Indicates which plot that the data X will be plotted against Y in.
13195      ! Y           Vertical position of the normalized data points in the plot.
13200      ! X           Horizontal position of the data point.
13205      ! Wndw(*)      Array containing the plot's scales.
13210      ! Vwprt(*)     Array containing the plot's CRT position.
13215      ! Symbol(*)    Array containing a distinct geometric symbol.
13220      ! Sy          Specifies which distinct geometric symbol is to be used.
13225      ! Noc         Specifies the number of coordinates that make up the distinct geometric symbol.
13230      OPTION BASE 1
13235      COM /Color1/ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta
13240      COM /Color2/ White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
13245      DIM Wndw(4),Vwprt(4),Symbol(20,3)
13250      DISP "Plotting Results"
13255      MAT Wndw= Array(60+Plot,*)
13260      MAT Vwprt= Array(70+Plot,*)
13265      Noc=Symbols(Sy,0,1)

```

```

13245 REDIM Symbol(Noc,3)
13250 MAT Symbol= Symbols(Sy,1:Noc,*)
13255 SELECT Sy
13260 CASE 1 ! The symbol chosen is a square with black edges and filled with red.
13265 PEN 16*Black
13270 AREA PEN 16*Red
13275 CASE 2 ! The symbol chosen is a octagon with black edges and filled with yellow.
13280 PEN 16*Black
13285 AREA PEN 16*Yellow
13290 CASE 3 ! The symbol chosen is a diamond with black edges and filled with green.
13295 PEN 16*Black
13300 AREA PEN 16*Green
13305 CASE 4 ! The symbol chosen is a triangle with black edges and filled with blue.
13310 PEN 16*Black
13315 AREA PEN 16*Blue
13320 END SELECT
13325 Xm=MIN(MAX(X,Wndw(1)),Wndw(2)) ! If X is out of bounds then set X to the edge of the graph.
13330 Ym=MIN(MAX(Y,Wndw(3)),Wndw(4)) ! If Y is out of bounds then set Y to the edge of the graph.
13335 LOG 5
13340 MOVE Xm,Ym
13345 SYMBOL Symbol(*),FILL,EDGE ! This draws the selected symbol.
13350 SUBEND
13355 Tcs8:
13360 Tcs8init: SUB Tcs8init(@Tcs8)
13365 ! Description:
13370 ! This subprogram is used to initialize this computer's internal RS232 serial interface.
13375 ! The subprogram also opens the TCS8 path on the Hewlett Packard series 9000 model 3XX computer
13380 ! for command and data transfer. The I/O path is given the name "@Tcs8". Data transferred
13385 ! from the HP to the TCS8 will use the "OUTPUT @Tcs8" statement. Data transferred to the HP
13390 ! from TCS8 will use the "ENTER @Tcs8" statement.
13395 ! The I/O path has a select code of 9 and is initialized to perform unformatted byte
13400 ! transfers without any end of line designations.
13405 REAL I(1:8),C(1:8)
13410 ASSIGN @Tcs8 TO 9;BYTE,FORMAT OFF,EOL ""
13415 CONTROL 9,0;1 ! Reset interface.
13420 CONTROL 9,3;9600 ! Select a baud rate of 9600.
13425 CONTROL 9,4;31 ! Select even parity, enable parity, 2 stop bits, 8 bits per character.
13430 CONTROL 9,12;IVAL("EF",16) ! Enable Carrier Detect. Disable Data Set Ready. Disable Clear To Send.
13435 CONTROL 9,13;9600 ! Default baud rate of 9600.
13440 CONTROL 9,14;31 ! Default character format: Even parity enabled, 2 stop, 8 bits/ char.
13445 SUBEND
13450 Tcs8set: SUB Tcs8set(C$,@Tcs8)
13455 ! Description:
13460 ! This subprogram allows the user to view and then set the various initialization parameters
13465 ! of each channel of the TCS8. These parameters are the current position, counts per inch,
13470 ! counts per revolution, motor velocity, motor acceleration, plus and minus limit switches,
13475 ! home switch, and motor stall indication. All of these parameters can be viewed and set except
13480 ! the limit and home switches and the stall indication. They can only be viewed.
13485 ! Variables:
13490 ! Command$ A TCS8 command string which indicates which parameter we want to view & set.
13495 ! View(*) Array of old TCS8 parameters viewed (received from TCS8). One for each channel.
13500 ! Set(*) Array of new TCS8 parameters to be set (sent to TCS8). One for each channel.
13505 ! Name$(*) String array of TCS8 parameter names.
13510 ! Image$(*) String array of image formats.
13515 ! Units$(*) String array of units.
13520 ! Channel Indicates the TCS8 channel number. Used to index the above arrays.
13525 OPTION BASE 1
13530 DIM View(8,1),Set(8,2),Name$(8,1)[10],Image$(8,1)[10],Units$(8,1)[10]
13535 OUTPUT @Tcs8 USING "K,/";"V"&C$&"0" ! Tell the TCS8 we want to View a parameter.
13540 ENTER @Tcs8 USING "8(K)";View(*) ! Enter the parameter specified by Command$.
13545 ! Initialize the Name$,Image$,Units$ and Set arrays.
13550 DATA X1,X2,Y1,Y2,Z1,Z2,A1,A2
13555 READ Name$(*)
13560 MAT Image$= ("6D.4D")
13565 FOR Channel=1 TO 8
13570 Set(Channel,1)=Channel
13575 SELECT C$
13580 CASE "P" ! Command$="P" indicates we want to view the encoder Positions in inches.
13585 Name$(Channel,1)=Name$(Channel,1)&" (pos)"
13590 Units$(Channel,1)="in"
13595 CASE "U" ! Command$="U" indicates we want to view the Units in counts per inch.
13600 Name$(Channel,1)=Name$(Channel,1)&" (cpi)"
13605 Units$(Channel,1)="cnt"
13610 CASE "R" ! Command$="R" indicates we want to view the number counts per Revolution.
13615 Name$(Channel,1)=Name$(Channel,1)&" (cpr)"
13620 Units$(Channel,1)="cnt"
13625 CASE "V" ! Command$="V" indicates we want to view the Velocity in revolution per second.
13630 Name$(Channel,1)=Name$(Channel,1)&" (vel)"
13635 Units$(Channel,1)="rev"
13640 CASE "A" ! Command$="A" indicates we want to view the Acceleration in revolution per second^2.

```

```

13645         Name$(Channel,1)=Name$(Channel,1)&" (acc)"
13650         Units$(Channel,1)="rev"
13655     CASE "+" ! Command$="+ " indicates we want to view the current + direction limit switches.
13660         Name$(Channel,1)=Name$(Channel,1)&" (+LS)"
13665         Units$(Channel,1)=" "
13670     CASE "-" ! Command$="- " indicates we want to view the current - direction limit switches.
13675         Name$(Channel,1)=Name$(Channel,1)&" (-LS)"
13680         Units$(Channel,1)=" "
13685     CASE "S" ! Command$="S" indicates we want to view the current motor Stall indication status.
13690         Name$(Channel,1)=Name$(Channel,1)&" (STALL)"
13695         Units$(Channel,1)=" "
13700     CASE "H" ! Command$="H" indicates we want to view the current Home limit switches.
13705         Name$(Channel,1)=Name$(Channel,1)&" (HS)"
13710         Units$(Channel,1)=" "
13715     END SELECT
13720 NEXT Channel
13725 ! The "Change" subprogram allows the user to see and then change the values of the viewed parameters.
13730 CALL Change("VALUES",View(*),Name$(*),Image$(*),Units$(*))
13735 ! The "Set" parameters command is now sent to the TCS8.
13740 SELECT CS
13745 CASE "P","U","R","V","A"
13750     MAT Set(*,2)= View(*,1)
13755     OUTPUT @Tcs8 USING 13760;"S"&CS,Set(*)
13760     IMAGE K,8(D,";",M6D.4D,""),/
13765 END SELECT
13770 SUBEND
13775 Tcs8read: SUB Tcs8read(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*))
13780 ! Description:
13785 ! This subprogram reads the current TCS8 positions. The 8 positions are read in TUNNEL
13790 ! coordinates with the units being in inches. Four of the eight positions (X1,Y1,Z1,A1) which
13795 ! are the transmitting side traverse positions are entered into the Tun1 array. The other four
13800 ! positions (X2,Y2,Z2,A2) which are the receiving side traverse positions are entered into the
13805 ! Tun2 array. The Tun1 & Tun2 arrays are converted from TUNNEL to MODEL coordinates.
13810 ! The current updated positions in the two coordinate systems are printed on the top of the
13815 ! CRT. They are also returned to the main program. The auxiliary channels A1 & A2 are not used.
13820 ! They can be used in the future to position probes such as hot wires and pitot tubes.
13825 ! Variables:
13830 ! Tun1(*) TCS8 transmitting side traverse positions (X1,Y1,Z1,A1) in TUNNEL coordinates.
13835 ! Tun2(*) TCS8 receiving side traverse positions (X2,Y2,Z2,A2) in TUNNEL coordinates.
13840 ! Mod1(*) TCS8 transmitting side traverse positions in MODEL coordinates.
13845 ! Mod2(*) TCS8 receiving side traverse positions in MODEL coordinates.
13850 ! Tun2mod(*) Coordinate system transformation matrix for converting TUNNEL to MODEL.
13855 ! Mod2tun(*) Coordinate system transformation matrix for converting MODEL to TUNNEL.
13860 COM /Color1/ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta
13865 COM /Color2/ White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
13870 OUTPUT @Tcs8 USING "K,/";"VPO"
13875 ENTER @Tcs8 USING "8(K)";Tun1(1),Tun2(1),Tun1(2),Tun2(2),Tun1(3),Tun2(3),Tun1(4),Tun2(4)
13880 REDIM Tun1(1:3),Tun2(1:3),Mod1(1:3),Mod2(1:3)
13885 MAT Mod1= Tun2mod*Tun1
13890 MAT Mod2= Tun2mod*Tun2
13895 REDIM Tun1(1:4),Tun2(1:4),Mod1(1:4),Mod2(1:4)
13900 CALL Tcs8print(Tun1(*),Tun2(*),Mod1(*),Mod2(*))
13905 SUBEND
13910 Tcs8print: SUB Tcs8print(Tun1(*),Tun2(*),Mod1(*),Mod2(*))
13915 ! Description:
13920 ! This subprogram prints the current updated TCS8 positions at the top of the CRT. The
13925 ! positions are printed in TUNNEL and MODEL coordinates for each side (Tx & Rx).
13930 ! Variables:
13935 ! Tun1(*) TCS8 transmitting side traverse positions (X1,Y1,Z1,A1) in TUNNEL coordinates.
13940 ! Tun2(*) TCS8 receiving side traverse positions (X2,Y2,Z2,A2) in TUNNEL coordinates.
13945 ! Mod1(*) TCS8 transmitting side traverse positions in MODEL coordinates.
13950 ! Mod2(*) TCS8 receiving side traverse positions in MODEL coordinates.
13955 COM /Color1/ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta
13960 COM /Color2/ White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
13965 PRINT PEN Red ! Print the traverse positions with red text.
13970 PRINT CHR$(128);CHR$(129); ! Print using inverse video text.
13975 PRINT TABXY(52,1);"
13980 PRINT TABXY(52,2);" TUN1 TUN2 MOD1 MOD2 "
13985 PRINT TABXY(52,3);"
13990 PRINT TABXY(52,4);"
13995 PRINT USING "#,K,4(M3D.4D),X";" X:",Tun1(1),Tun2(1),Mod1(1),Mod2(1)
14000 PRINT TABXY(52,5);"
14005 PRINT USING "#,K,4(M3D.4D),X";" Y:",Tun1(2),Tun2(2),Mod1(2),Mod2(2)
14010 PRINT TABXY(52,6);"
14015 PRINT USING "#,K,4(M3D.4D),X";" Z:",Tun1(3),Tun2(3),Mod1(3),Mod2(3)
14020 PRINT TABXY(52,7);"
14025 PRINT USING "#,K,4(M3D.4D),X";" A:",Tun1(4),Tun2(4),Mod1(4),Mod2(4)
14030 PRINT TABXY(52,8);"
14035 PRINT CHR$(128); ! Turn off inverse video.
14040 PRINT PEN Black ! Set printing color to black.

```

```

14045      SUBEND
14050 Tcs8move: SUB Tcs8move(@Tcs8,Tun1(*),Tun2(*),Mod1(*),Mod2(*),Tun2mod(*),Mod2tun(*),Side$,Coor$,Mode$,K,Movement)
14055      ! Description:
14060      ! This subprogram allows for the movement of the probe volume and collecting optics in one of
14065      ! two coordinate systems. The two coordinate systems implemented are the TUNNEL and the MODEL
14070      ! coordinate systems. Two movements modes are available. The first movement mode makes moves
14075      ! relative to the current position. The second movement mode makes moves to an absolute fixed
14080      ! position. Both the transmitting side and receiving side traverses can be moved in tandem
14085      ! or separately.
14090      ! Variables:
14095      ! Tun1(*)      TCS8 transmitting side traverse positions (X1,Y1,Z1,A1) in TUNNEL coordinates.
14100      ! Tun2(*)      TCS8 receiving side traverse positions (X2,Y2,Z2,A2) in TUNNEL coordinates.
14105      ! Mod1(*)      TCS8 transmitting side traverse positions in MODEL coordinates.
14110      ! Mod2(*)      TCS8 receiving side traverse positions in MODEL coordinates.
14115      ! Tun2mod(*)   Coordinate system transformation matrix for converting TUNNEL to MODEL.
14120      ! Mod2tun(*)   Coordinate system transformation matrix for converting MODEL to TUNNEL.
14125      ! Side$       Indicates which sides are to be moved:
14130      !               Tx      : Transmitting side only.
14135      !               Rx      : Receiving side only.
14140      !               Tx & Rx : Both sides together.
14145      ! Coor$        Indicates which coordinate system the movement is to be made in:
14150      !               TUNNEL  : TUNNEL coordinates.
14155      !               MODEL   : MODEL coordinates.
14160      ! Mode$        Indicates which movement mode is to be completed:
14165      !               RELATIVE: Movements are relative to current positions.
14170      !               ABSOLUTE: Movements are to absolute positions.
14175      ! K            Indicates which axis of the four axes is to be moved.
14180      ! Movement     Indicates the desired movement for the selected axis.
14185      ! I(*)         Array of viewed TCS8 "Initialized" parameters.
14190      ! C(*)         Array of viewed TCS8 "Currents On" parameters.
14195      OPTION BASE 1
14200      DIM LS(100)
14205      ! If all of the channels have not yet been initialized, then do so now.
14210      REAL Move(8,2),I(8),C(8)
14215      OUTPUT @Tcs8 USING "K,/";"VIO"
14220      ENTER @Tcs8 USING "8(K)";I(*)
14225      IF SUM(I)>8 THEN OUTPUT @Tcs8 USING "K,/";"SIO"
14230      ! If all of the channels do not have their currents turned on, then do so now.
14235      OUTPUT @Tcs8 USING "K,/";"VCO"
14240      ENTER @Tcs8 USING "8(K)";C(*)
14245      IF SUM(C)>8 THEN OUTPUT @Tcs8 USING "K,/";"SCO:1,"
14250      ! If the movement mode is to be RELATIVE, then clear all of the previously read positions.
14255      IF Mode$="RELATIVE" THEN
14260          MAT Tun1= (0)
14265          MAT Tun2= (0)
14270          MAT Mod1= (0)
14275          MAT Mod2= (0)
14280      END IF
14285      ! Set the new Tun1(*) and Tun2(*) position arrays.
14290      SELECT Coor$
14295      CASE "MODEL"
14300          Mod1(K)=Movement
14305          Mod2(K)=Movement
14310          REDIM Tun1(1:3),Tun2(1:3),Mod1(1:3),Mod2(1:3)
14315          IF POS(Side$,"Tx") THEN MAT Tun1= Mod2tun*Mod1
14320          IF POS(Side$,"Rx") THEN MAT Tun2= Mod2tun*Mod2
14325          REDIM Tun1(1:4),Tun2(1:4),Mod1(1:4),Mod2(1:4)
14330      CASE "TUNNEL"
14335          IF POS(Side$,"Tx") THEN Tun1(K)=Movement
14340          IF POS(Side$,"Rx") THEN Tun2(K)=Movement
14345      END SELECT
14350      ! File the move array.
14355      FOR Channel=1 TO 8
14360          Move(Channel,1)=Channel
14365      NEXT Channel
14370      Move(1,2)=Tun1(1)
14375      Move(2,2)=Tun2(1)
14380      Move(3,2)=Tun1(2)
14385      Move(4,2)=Tun2(2)
14390      Move(5,2)=Tun1(3)
14395      Move(6,2)=Tun2(3)
14400      Move(7,2)=Tun1(4)
14405      Move(8,2)=Tun2(4)
14410      ! Initiate the start of the move.
14415      IF Mode$="ABSOLUTE" THEN OUTPUT @Tcs8 USING 14425;"MA",Move(*)
14420      IF Mode$="RELATIVE" THEN OUTPUT @Tcs8 USING 14425;"MR",Move(*)
14425      IMAGE K,8(D,";",S2D.5D,""),/
14430      ! The TCS8 will return the new updated positions only after the move is complete.
14435      ENTER @Tcs8 USING "8(K)";Tun1(1),Tun2(1),Tun1(2),Tun2(2),Tun1(3),Tun2(3),Tun1(4),Tun2(4)
14440      ! Turn off the motor drive currents.

```

```

14445 OUTPUT @Tcs8 USING *K,/;"SCO:0,"
14450 SUBEND
14455 !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
14460 Ctm: SUB Ctm(Alpha(*),Tun2mod(*),Mod2tun(*) )
14465 ! Description:
14470 ! This subprogram computes directly the MODEL to TUNNEL coordinate system transformation
14475 ! matrix "Mod2tun(*)". However, the desired coordinate system transformation matrix "Tun2mod" is
14480 ! required. It is the matrix inverse of "Mod2tun".
14485 ! Variables:
14490 ! Alpha(*) Angles of attack, yaw, and roll.
14495 ! T1(*) Partial coordinate system transformation matrix for converting from MODEL to
14500 ! TUNNEL coordinates. Takes into account a model at angle of attack.
14505 ! T2(*) Partial coordinate system transformation matrix for converting from MODEL to
14510 ! TUNNEL coordinates. Takes into account a model at angle of yaw.
14515 ! T3(*) Partial coordinate system transformation matrix for converting from MODEL to
14520 ! TUNNEL coordinates. Takes into account a model at angle of roll.
14525 ! Mod2tun(*) Coordinate system transformation matrix for converting from MODEL to TUNNEL.
14530 ! Tun2mod(*) Coordinate system transformation matrix for converting from TUNNEL to MODEL.
14535 OPTION BASE 1
14540 REAL T1(3,3),T2(3,3),T3(3,3),Temp(3,3)
14545 ! Define 1st coordinate transformation matrix for Mod2tun.
14550 ! Rotation in the x-y plane about the z-axis.
14555 ! Used when model is at an angle of attack.
14560 T1(1,1)=COS(Alpha(1))
14565 T1(1,2)=SIN(Alpha(1))
14570 T1(1,3)=0
14575 T1(2,1)=-SIN(Alpha(1))
14580 T1(2,2)=COS(Alpha(1))
14585 T1(2,3)=0
14590 T1(3,1)=0
14595 T1(3,2)=0
14600 T1(3,3)=1
14605 ! Define 2nd coordinate transformation matrix for Mod2tun.
14610 ! Rotation in the x-z plane about the y-axis.
14615 ! Used when model is at an angle of yaw.
14620 T2(1,1)=COS(Alpha(2))
14625 T2(1,2)=0
14630 T2(1,3)=-SIN(Alpha(2))
14635 T2(2,1)=0
14640 T2(2,2)=1
14645 T2(2,3)=0
14650 T2(3,1)=SIN(Alpha(2))
14655 T2(3,2)=0
14660 T2(3,3)=COS(Alpha(2))
14665 ! Define 3rd coordinate transformation matrix for Mod2tun.
14670 ! Rotation in the y-z plane about the x-axis.
14675 ! Used when model is at an angle of roll.
14680 T3(1,1)=1
14685 T3(1,2)=0
14690 T3(1,3)=0
14695 T3(2,1)=0
14700 T3(2,2)=COS(Alpha(3))
14705 T3(2,3)=SIN(Alpha(3))
14710 T3(3,1)=0
14715 T3(3,2)=-SIN(Alpha(3))
14720 T3(3,3)=COS(Alpha(3))
14725 ! Mod2tun converts MODEL coordinates to TUNNEL coordinates.
14730 MAT Temp= T2*T1
14735 MAT Mod2tun= T3*Temp
14740 ! Tun2mod converts TUNNEL coordinates to MODEL coordinates.
14745 MAT Tun2mod= INV(Mod2tun)
14750 SUBEND
14755 Color: !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!
14760 Crt_init: SUB Crt_init
14765 ! Description:
14770 ! This subprogram initializes the CRT as the plotting device and clears both the alpha
14775 ! numerics and graphics part of the CRT. The color map for both of the alpha numeric printing
14780 ! plains and the graphics drawing plains are defined here.
14785 COM /Color1/ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta
14790 COM /Color2/ White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
14795 CALL Color ! Define the color maps for the alpha numeric and the graphics plains.
14800 !CALL Map ! Draw the color map.
14805 !CALL Dump ! Dump the color map to the printer.
14810 PRINTER IS CRT ! Select the CRT as the printing device.
14815 PRINTALL IS CRT ! Send ERROR and DISP messages to CRT.
14820 KEY LABELS OFF ! Hide the special function key labels for f1..f8.
14825 CLEAR SCREEN ! Clear the alpha numeric printing plains of the CRT.
14830 GCLEAR ! Clear the graphics drawing plains of the CRT.
14835 SUBEND
14840 Color: SUB Color

```

```

14845      ! Description:
14850      !   This subprogram defines the color map for both alpha numeric printing and graphics drawing.
14855      !   Four of eight plains are dedicated to alpha numerics to provide for sixteen colors. The other
14860      !   four plains are dedicated to graphics to provide for sixteen colors.
14865      COM /Color1/ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta
14870      COM /Color2/ White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
14875      DIM Map(255,2)
14880      INTEGER Gmask(1)
14885      READ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta,White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
14890      DATA 0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15
14895      PLOTTER IS CRT,"INTERNAL";COLOR MAP      ! Select the CRT as the plotting device.
14900      CONTROL CRT,14;3
14905      SET PEN 0 INTENSITY 1,1,1      ! Set pen 0 equal to clear (white).
14910      SET PEN 1 INTENSITY 0,0,0      ! Set pen 1 equal to black.
14915      SET PEN 8 INTENSITY 1,1,1      ! Set pen 8 equal to white.
14920      SET PEN 14 INTENSITY 26/30,16/30,8/30      ! Set pen 14 equal to brown.
14925      SET PEN 15 INTENSITY .6,.6,.6      ! Set pen 15 equal to gray.
14930      GESCAPE CRT,2;Map(*)      ! Read RGB intensity for pens 0 to 255.
14935      Gmask(0)=IVAL("11110000",2)      ! Define graphics write enable mask.
14940      Gmask(1)=IVAL("11110000",2)      ! Define graphics display enable mask.
14945      GESCAPE CRT,7,Gmask(*)      ! Set graphics write & display enable masks.
14950      SET ALPHA MASK IVAL("00001111",2)      ! Set alpha write enable mask.
14955      SET DISPLAY MASK IVAL("00001111",2)      ! Set alpha display enable mask.
14960      GESCAPE CRT,4      ! Select normal dominant writing mode.
14965      GCLEAR      ! Clear the graphics screen.
14970      CLEAR SCREEN      ! Clear the alpha screen.
14975      GRAPHICS ON      ! Turn graphics on.
14980      ! Copy the alpha colors to the graph colors (use the same 16 alpha colors for graph colors.)
14985      FOR Alpha=0 TO 15
14990          FOR Graph=0 TO 15
14995              Pen=16*Graph+Alpha      ! Define pen number for Alpha:Graph combination.
15000              Color=Graph*(Alpha=0)+Alpha*(Alpha>0)      ! Choose the color for the
15005              IF Alpha=Graph THEN Color=Black*(Alpha>1)      ! Alpha:Graph combination.
15010              MAT Map(Pen,*)= Map(Color,*)      ! Get the RGB intensities for the color.
15015              SET PEN Pen INTENSITY Map(Pen,0),Map(Pen,1),Map(Pen,2)      ! Set the RGB for the pen.
15020          NEXT Graph
15025      NEXT Alpha
15030      ! AREA PEN White      ! Select white for area fills.
15035      ! PEN Black      ! Select black for line drawing and labeling.
15040      ALPHA PEN Black      ! Select black for printing.
15045      KEY LABELS PEN Blue      ! Select blue for special function key labels.
15050      PRINT PEN Black      ! Select black for printing.
15055      KEY LABELS OFF      ! Hide the special function key labels for f1..f8.
15060      SUBEND
15065      Map:
15070      SUB Map
15075      ! Description:
15080      !   This subprogram displays the color map on the CRT. The sixteen colors for the alpha plains are
15085      !   superimposed on top of the graphics plains to show the dominance interaction of alpha and
15090      !   graphics colors being printed and drawn on top of each other.
15095      COM /Color1/ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta
15100      COM /Color2/ White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
15105      VIEWPORT 25/10.23,(25+16*4*10)/10.23,210/10.23,850/10.23
15110      WINDOW 0,16,16,0
15115      Pen=0
15120      FOR Alpha=0 TO 15
15125          FOR Graph=0 TO 15
15130              AREA PEN 16*Graph+Alpha
15135              PEN 16*Black
15140              MOVE Alpha,Graph
15145              RECTANGLE 1,1,FILL,EDGE
15150              PRINT PEN Alpha
15155              PRINT TABXY(4+4*Alpha,10+2*Graph);
15160              PRINT USING "ZZZ";Pen
15165              Pen=Pen+1
15170          NEXT Graph
15175      NEXT Alpha
15180      ALPHA PEN Black
15185      KEY LABELS PEN Blue
15190      PRINT PEN Black
15195      SUBEND
15195      Dump:
15200      SUB Dump
15205      ! Description:
15210      !   This subprogram dumps the graphics contents of the CRT to the printer. This facilitates
15215      !   the printing of the histogram and profile plots. The CSUB binary subprogram is used to
15220      !   transfer the colored plots to the color paint jet printer.
15225      OUTPUT PRT USING "#,@"
15230      IF NOT (INMEM("Gdump_colored")) THEN LOADSUB ALL FROM "CDUMP6"
15235      IF NOT (INMEM("Bstore")) THEN LOADSUB ALL FROM "BPL0T6"
15240      IF NOT (INMEM("Bload")) THEN LOADSUB ALL FROM "BPL0T6"
      !OUTPUT PRT USING "#,5/"

```



```

15245      !CALL Gdump_colored(CRT,PRT,"NORMAL",180,"OFF","DITHER")
15250      !CALL Gdump_colored(CRT,PRT,"ROTATE",90,"ON","ERRDIF")
15255      CALL Gdump_colored(CRT,PRT,"NORMAL",180,"ON","DITHER")
15260      SUBEND
15265 Read_symbols: SUB Read_symbols(Symbols(*))
15270      ! Description:
15275      !       This subprogram defines 5 geometric symbols to be used with the SYMBOL statement. The
15280      !       symbols provided are as follows: Square,Octagon,Diamond, and Triangles (upwards & downwards
15285      !       pointing triangles). All of the symbols have a dot added to their center.
15290      ! Variables:
15295      !       Symbols(*) Array of Symbol arrays. Each symbol arrays contains a distinct geometric symbol.
15300      !       Symbol(*) Array of coordinates which when connected produce a distinct geometric symbol.
15305      !       Dot(*) Array of coordinates which produce a dot. The dot symbol is added to all symbols.
15310      !       Noc       The number of coordinates in a symbol.
15315      !       S         Used to index the Symbols array.
15320      OPTION BASE 1
15325      REAL Symbol(20,3),Dot(2,3)
15330      READ Dot(*)
15335      FOR S=1 TO 5
15340          READ Noc
15345          REDIM Symbol(Noc,3)
15350          READ Symbol(*)
15355          MAT Symbols(S,1:Noc,*)= Symbol
15360          MAT Symbols(S,Noc+1:Noc+2,*)= Dot
15365          Symbols(S,0,1)=Noc+2
15370      NEXT S
15375 Dot:
15380 Square:
15385 Octagon:
15390 Diamond:
15395 Utriangle:
15400 Dtriangle:
15405      SUBEND
15410 Graph:
15415 Setup_graph: SUB Setup_graph(Array(*),Image$(*),Paxis,Symbols(*))
15420      ! Description:
15425      !       This subprogram sets up nine empty plots on the CRT screen. Four plots are profile plots
15430      !       while the other five plots are histogram plots. The profile and histogram plots provided are
15435      !       as follows:
15440      !
15445      !
15450      !
15455      !
15460      !
15465      !
15470      !
15475      !
15480      !
15485      ! Variables:
15490      !       Array(*) Array containing the plot positions and scales.
15495      !       Image$(*) String array containing image formats for the axes labeling.
15500      !       Wndw(*) Array containing the plot's scales.
15505      !       Vwprt(*) Array containing the plot's CRT position.
15510      !       Xdiv(*) Array containing the number of X divisions for the plot's X axis.
15515      !       Ydiv(*) Array containing the number of Y divisions for the plot's Y axis.
15520      !       Xlabel$(*) String array containing labels for the X axis.
15525      !       Ylabel$(*) String array containing labels for the Y axis.
15530      !       Title$(*) String array containing labels for the Plots.
15535      !       Ximage$(*) String array containing image formats for the X axis labeling.
15540      !       Yimage$(*) String array containing image formats for the Y axis labeling.
15545      !       Legend$(*) String array containing labels for each symbol in a profile plot.
15550      !       Symbols(*) Array of Symbol arrays. Each symbol arrays contains a distinct geometric symbol.
15555      !       G         Used as an index to the above arrays. Specifies one of nine plots.
15560      !       I         Used as an index to the Legend$ array.
15565      OPTION BASE 1
15570      COM /Graph1/ Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*)
15575      COM /Graph2/ Titles(*),Ximage$(*),Yimage$(*),Legend$(*)
15580      COM /Color1/ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta
15585      COM /Color2/ White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
15590      MAT Wndw= Array(61:69,*)
15595      MAT Vwprt= Array(71:79,*)
15600      MAT Xdiv(1:5)= Array(81:85,1)
15605      MAT Xdiv(6:9)= Array(81:84,3)
15610      MAT Ydiv(1:5)= Array(81:85,2)
15615      MAT Ydiv(6:9)= Array(81:84,4)
15620      MAT Ximage$= Image$(61:69,1)
15625      MAT Yimage$= Image$(61:69,3)
15630      FOR G=1 TO 9
15635          READ G,Xlabel$(G)

```

```

15640         FOR I=1 TO SIZE(Legend$,2)
15645             READ Legend$(G,I)
15650         NEXT I
15655         SELECT G
15660         CASE 1 TO 5
15665             Ylabel$(G)=""
15670         CASE 6 TO 9
15675             Ylabel$(G)=CHR$(NUM("X")+Paxis-1)
15680         END SELECT
15685         CALL Set_up(G,Symbols(*))
15690     NEXT G
15695     SUBEXIT
15700     ! G, X axis Label
15705     ! Symbol #1...5 labels
15710     DATA 1, ""
15715     DATA 2, ""
15720     DATA 3, ""
15725     DATA 4, ""
15730     DATA 5, ""
15735     DATA 6, "U,V,U',V' /Uinf"
15740     DATA 7, "A,B,A',B' volts"
15745     DATA 8, "Tt:dR Uinf:m/s Uedge:m/s"
15750     DATA 9, "Shear Stress Terms / Uinf^2"
15755     SUBEND
15755     Set_up: SUB Set_up(G,Symbols(*))
15760     ! Description:
15765     ! This subprogram clears and then redraws one of nine empty plots on the CRT screen.
15770     ! Variables:
15775     ! Wndw(*) Array containing the plot's scales.
15780     ! Vwprt(*) Array containing the plot's CRT position.
15785     ! Xdiv(*) Array containing the number of X divisions for the plot's X axis.
15790     ! Ydiv(*) Array containing the number of Y divisions for the plot's Y axis.
15795     ! Xlabel$(*) String array containing labels for the X axis.
15800     ! Ylabel$(*) String array containing labels for the Y axis.
15805     ! Title$(*) String array containing labels for the Plots.
15810     ! Ximage$(*) String array containing image formats for the X axis labeling.
15815     ! Yimage$(*) String array containing image formats for the Y axis labeling.
15820     ! Legend$(*) String array containing labels for each symbol in a profile plot.
15825     ! Symbols(*) Array of Symbol arrays. Each symbol arrays contains a distinct geometric symbol.
15830     ! G Used as an index to the above arrays. Specifies one of nine plots.
15835     OPTION BASE 1
15840     COM /Graph1/ Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*)
15845     COM /Graph2/ Title$(*),Ximage$(*),Yimage$(*),Legend$(*)
15850     COM /Color1/ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta
15855     COM /Color2/ White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
15860     DIM L$(80)
15865     ON ERROR CALL Error
15870     CSIZE 100*15/1023 ! Select a character labeling size of 15 pixels high.
15875     ! Define the values for the left,right,bottom,top ends of the horizontal and vertical scales.
15880     Xmin=Wndw(G,1)
15885     Xmax=Wndw(G,2)
15890     Ymin=Wndw(G,3)
15895     Ymax=Wndw(G,4)
15900     ! Define the values for the left,right,bottom,top pixel locations for the plot.
15905     Xpix1=Vwprt(G,1)
15910     Xpix2=Vwprt(G,2)
15915     Ypix1=Vwprt(G,3)
15920     Ypix2=Vwprt(G,4)
15925     ! Define the step size between grid lines, axis tick marks, and axis labels.
15930     Xstep=(Xmax-Xmin)/Xdiv(G)
15935     Ystep=(Ymax-Ymin)/Ydiv(G)
15940     ! Define the amount of scale X and Y which equals the size of one pixel (picture element).
15945     Xpixel=(Xmax-Xmin)/(Xpix2-Xpix1)
15950     Ypixel=(Ymax-Ymin)/(Ypix2-Ypix1)
15955     ! Clear the plots back ground & plot area and also draw the plots borders, grids, and axes.
15960     AREA PEN 16*White
15965     !GOSUB Clear_screen
15970     AREA PEN 16*White
15975     GOSUB Back_ground
15980     AREA PEN 16*White
15985     GOSUB Plot_area
15990     GOSUB Scale
15995     PEN 16*Blue
16000     GOSUB Axes
16005     GOSUB Grid
16010     GOSUB Scale
16015     PEN 16*Black
16020     CLIP OFF
16025     ! Draw the X and Y axis labels.
16030     GOSUB Ylabel
16035     GOSUB Xlabel

```

```

16040      ! Create a legend to define which symbol is used with which data.
16045      CALL Legend(G,Symbols(*))
16050      OFF ERROR
16055      SUBEXIT
16060 Clear_screen:  ! This subroutine fills the entire CRT screen with the specified color.
16065      VIEWPORT 0/10.23,1279/10.23,0/10.23,1023/10.23
16070      WINDOW -1.E+9,1.E+9,-1.E+9,1.E+9
16075      MOVE 0,0
16080      WINDOW 0,1279,0,1023
16085      MOVE 0,0
16090      RECTANGLE 1279,1023,FILL
16095      RETURN
16100 Back_ground:  ! This subroutine clears the plot's background.
16105      VIEWPORT (Xpix1-80)/10.23,(Xpix2+15)/10.23,(Ypix1-33)/10.23,(Ypix2+10)/10.23
16110      WINDOW -1.E+9,1.E+9,-1.E+9,1.E+9
16115      MOVE 0,0
16120      WINDOW Xmin,Xmax,Ymin,Ymax
16125      MOVE Xmin,Ymin
16130      RECTANGLE (Xmax-Xmin),(Ymax-Ymin),FILL
16135      RETURN
16140 Plot_area:  ! This subroutine selects part of the CRT plot area and give it scales for the X and Y axes.
16145      VIEWPORT Xpix1/10.23,Xpix2/10.23,Ypix1/10.23,Ypix2/10.23
16150      WINDOW -1.E+9,1.E+9,-1.E+9,1.E+9
16155      MOVE 0,0
16160      WINDOW Xmin,Xmax,Ymin,Ymax
16165      MOVE Xmin,Ymin
16170      RECTANGLE (Xmax-Xmin),(Ymax-Ymin),FILL
16175      RETURN
16180 Axes:  ! This subroutine draws the plot's X and Y axes.
16185      VIEWPORT (Xpix1-1)/10.23,(Xpix2+1)/10.23,(Ypix1-6)/10.23,(Ypix1-1)/10.23
16190      WINDOW Xmin,Xmax,1,0
16195      AXES Xstep,2,Xmin,0,1,1,1
16200      VIEWPORT (Xpix1-1)/10.23,(Xpix2+1)/10.23,(Ypix2+1)/10.23,(Ypix2+6)/10.23
16205      WINDOW Xmin,Xmax,0,1
16210      AXES Xstep,2,Xmin,0,1,1,1
16215      VIEWPORT (Xpix1-6)/10.23,(Xpix1-1)/10.23,(Ypix1-1)/10.23,(Ypix2+1)/10.23
16220      WINDOW 1,0,Ymin,Ymax
16225      AXES 2,Ystep,0,Ymin,1,1,1
16230      VIEWPORT (Xpix2+1)/10.23,(Xpix2+6)/10.23,(Ypix1-1)/10.23,(Ypix2+1)/10.23
16235      WINDOW 0,1,Ymin,Ymax
16240      AXES 2,Ystep,0,Ymin,1,1,1
16245      RETURN
16250 Grid:  ! This subroutine draws the plot's X and Y grid lines.
16255      VIEWPORT (Xpix1-1)/10.23,(Xpix2+1)/10.23,(Ypix1-1)/10.23,(Ypix2+1)/10.23
16260      WINDOW Xmin,Xmax,Ymin,Ymax
16265      LINE TYPE 4
16270      GRID Xstep,Ystep,Xmin,Ymin
16275      LINE TYPE 1
16280      RETURN
16285 Scale:  ! This subroutine labels the X axis and also names the X axis.
16290      VIEWPORT Xpix1/10.23,Xpix2/10.23,Ypix1/10.23,Ypix2/10.23
16295      WINDOW Xmin,Xmax,Ymin,Ymax
16295      RETURN
16300 Xlabel:  ! This subroutine labels the X axis and also names the X axis.
16305      LOG 5
16310      FOR X=Xmin TO Xmax+Xstep/100 STEP Xstep
16315          MOVE X,Ymin-14*Ypixel
16320          OUTPUT L$ USING Ximage$(G);X
16325          LABEL TRIMS(L$)
16330      NEXT X
16335      MOVE (Xmin+Xmax)/2,Ymin-27*Ypixel
16340      !LABEL Xlabel$(G)
16345      RETURN
16350      IF G=8 THEN
16355          LOG 5
16360          DEG
16365          LDIR 45
16370          MOVE (Xmin+Xmax)/2,(Ymax+Ymin)/2
16375          CSIZE 100*100/1023
16380          LABEL "VOID"
16385          CSIZE 100*15/1023
16390          LDIR 0
16395      END IF
16400      RETURN
16405 Ylabel:  ! This subroutine labels the Y axis and also names the Y axis.
16410      LOG 8
16415      Len=0
16420      FOR Y=Ymin TO Ymax+Ystep/100 STEP Ystep
16425          MOVE Xmin-7*Ypixel,Y
16430          OUTPUT L$ USING Yimage$(G);Y
16435          LABEL TRIMS(L$)

```

```

16440      Len=MAX(Len,LEN(TRIMS(LS)))
16445      NEXT Y
16450      LOG 5
16455      MOVE Xmin-(18+7*Len)*Xpixel,(Ymin+Ymax)/2      !+20*Ypixel
16460      LABEL Ylabel$(G)
16465      RETURN
16470  SUBEND
16475  Legend: SUB Legend(G,Symbols(*))
16480      ! Description:
16485      ! This subprogram produces a legend within one of the nine plots on the CRT screen.
16490      ! Variables:
16495      ! Wndw(*)      Array containing the plot's scales.
16500      ! Vwprt(*)     Array containing the plot's CRT position.
16505      ! Xdiv(*)      Array containing the number of X divisions for the plot's X axis.
16510      ! Ydiv(*)      Array containing the number of Y divisions for the plot's Y axis.
16515      ! Xlabel$(*)   String array containing labels for the X axis.
16520      ! Ylabel$(*)   String array containing labels for the Y axis.
16525      ! Title$(*)    String array containing labels for the Plots.
16530      ! Ximage$(*)   String array containing image formats for the X axis labeling.
16535      ! Yimage$(*)   String array containing image formats for the Y axis labeling.
16540      ! Legend$(*)   String array containing labels for each symbol in a profile plot.
16545      ! Symbols(*)   Array of Symbol arrays. Each symbol arrays contains a distinct geometric symbol.
16550      ! Symbol(*)    Array of coordinates which when connected produce a distinct geometric symbol.
16555      ! G           Used as an index to the above arrays. Specifies one of nine plots.
16560      ! S           Used to index the Legend$ array.
16565      ! Noc         The number of coordinates in a symbol.
16570      ! Len         Total Length of all Legend$ array elements.
16575      OPTION BASE 1
16580      COM /Graph1/ Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*)
16585      COM /Graph2/ Title$(*),Ximage$(*),Yimage$(*),Legend$(*)
16590      COM /Color1/ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta
16595      COM /Color2/ White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
16600      DIM Symbol(20,3)
16605      VIEWPORT Vwprt(G,1)/10.23,Vwprt(G,2)/10.23,Vwprt(G,3)/10.23,Vwprt(G,4)/10.23
16610      WINDOW Vwprt(G,1),Vwprt(G,2),Vwprt(G,3),Vwprt(G,4)
16615      CLIP OFF
16620      CSIZE 100*15/1023      ! Select a character labeling size of 15 pixels high.
16625      LOG 2
16630      ! Calculate the total length of all of the symbol labels.
16635      Len=0
16640      FOR S=1 TO SIZE(Legend$,2)
16645      IF LEN(Legend$(G,S)) THEN
16650      Len=Len+LEN(TRIMS(Legend$(G,S)))+2.2
16655      END IF
16660      NEXT S
16665      X=(Vwprt(G,1)+Vwprt(G,2))/2
16670      Y=(Vwprt(G,3)+Vwprt(G,4))/2
16675      MOVE X,Y
16680      X=(Vwprt(G,1)+Vwprt(G,2))/2-5*Len+10
16685      Y=Vwprt(G,3)-28
16690      ! For each symbol put up a sample symbol and its label.
16695      FOR S=1 TO SIZE(Legend$,2)
16700      IF LEN(Legend$(G,S))=0 THEN 16825
16705      Noc=Symbols(S,0,1)
16710      REDIM Symbol(Noc,3)
16715      MAT Symbol= Symbols(S,1:Noc,*)
16720      ! Define the colors for symbol filling and edge drawing.
16725      SELECT S
16730      CASE 1
16735      AREA PEN 16*Red
16740      PEN 16*Black
16745      CASE 2
16750      AREA PEN 16*Yellow
16755      PEN 16*Black
16760      CASE 3
16765      AREA PEN 16*Green
16770      PEN 16*Black
16775      CASE 4
16780      AREA PEN 16*Blue
16785      PEN 16*Black
16790      END SELECT
16795      MOVE X,Y      ! Move to the place of next symbol.
16800      SYMBOL Symbol(*),FILL,EDGE      ! Draw the next symbol.
16805      X=X+12      ! Move the X placement to the right 12 pixels.
16810      MOVE X,Y-1      ! Move to the place of next label.
16815      LABEL Legend$(G,S)      ! Draw the next label.
16820      X=X+10*LEN(Legend$(G,S))+10      ! Move the X placement to the right 10+10*Len pixels
16825      NEXT S
16830  SUBEND
16835  Histo: !!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!!

```

```

16840 Rt_histo: SUB Rt_histo(@Lvdas,Symbols(*),Repeat,Kbd$)
16845 ! Description:
16850 ! This subprogram plots real time histograms within five of the nine plots on the CRT screen.
16855 ! The histogram data are acquired from the LVDAS over a specified acquisition time.
16860 ! Variables Defined in Main Program:
16865 ! Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),Legend$(*)
16870 ! Local Variables:
16875 ! Histo(*) Array of bin numbers, old histogram bin heights, and new histogram bin heights.
16880 ! Nbins Number of bins in the Histo(*).
16885 ! Bin 2^Bin is the bin width of individual histogram vertical bars.
16890 ! Min Minimum value for histogram. Left side of histogram scale.
16895 ! Max Maximum value for histogram. right side of histogram scale.
16900 ! F1 Upper 16bits of integerized Min.
16905 ! F2 Lower 16bits of integerized Min.
16910 ! A1 Upper 16bits of integerized histogram acquisition time.
16915 ! A2 Lower 16bits of integerized histogram acquisition time.
16920 ! Nnew Number of samples in the most up to date histogram.
16925 ! Nold Number of samples in the previous histogram.
16930 ! N(*) Number of samples for each histogram of the five separate channels.
16935 ! Channel Used to select the LVDAS channel that will be sampled for a histogram.
16940 ! Kw Converts Hz to MHz or raw data to volts.
16945 ! Ww Window width of each vertical histogram bar.
16950 ! Old Histogram height of previous histogram at a particular bin.
16955 ! New Histogram height of current histogram at a particular bin.
16960 ! X1 Horizontal position of histogram rectangle.
16965 ! Y1 Vertical position of histogram rectangle.
16970 ! X2 Horizontal width of histogram rectangle.
16975 ! Y2 Vertical width of histogram rectangle.
16980 ! I Used as an index to the Histo(*). Specifies one of Nbins bins.
16985 ! G Used as an index to the graphics arrays. Specifies one of nine plots.
16990 OPTION BASE 1
16995 COM /Graph1/ Wndw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*)
17000 COM /Graph2/ Title$(*),Ximage$(*),Yimage$(*),Legend$(*)
17005 COM /Color1/ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta
17010 COM /Color2/ White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
17015 INTEGER Histo(1000,3),Nplots,Nbins,F1,F2,A1,A2
17020 REAL Nnew,Nold,N(5)
17025 ! Clear all of the histogram data within the LVDAS.
17030 OUTPUT @Lvdas USING "AA";"CA"
17035 ! Draw new plots for the five histograms.
17040 FOR Channel=1 TO 5
17045 CALL Set_up(Channel,Symbols(*))
17050 NEXT Channel
17055 ! Calculate the acquisition time. 0.1*10000000 will give an acquisition of 0.1 seconds.
17060 CALL Convert2words(.1*10000000,A1,A2) ! Atime=.1 seconds
17065 ! Enable the keyboard to terminate histogram plotting.
17070 ON KBD GOSUB Hdone
17075 REPEAT
17080 FOR Channel=1 TO 5
17085 G=Channel
17090 SELECT Channel
17095 CASE 1,2 ! Channels 1,2,3 are for LDV frequency data.
17100 Kw=1000000 ! Converts Hz to MHz.
17105 Min=Kw*Wndw(G,1) ! Minimum frequency for left histogram scale.
17110 Max=Kw*Wndw(G,2) ! Maximum frequency for right histogram scale.
17115 Bin=INT(LGT((Max-Min)/100)/LGT(2))+1 ! 2^Bin is the window width of each vertical bar.
17120 Ww=2^Bin ! Window width of each vertical histogram bar.
17125 CALL Convert2words(Min,F1,F2)
17130 CASE 4 ! Channels 4,5 are for analog voltage data.
17135 Kw=32768/5 ! Converts raw data to volts.
17140 Min=Kw*Wndw(G,1) ! Minimum voltage for left histogram scale.
17145 Max=Kw*Wndw(G,2) ! Maximum voltage for right histogram scale.
17150 Bin=INT(LGT((Max-Min)/100)/LGT(2))+1 ! 2^Bin is the window width of each vertical bar.
17155 Ww=2^Bin ! Window width of each vertical histogram bar.
17160 CALL Convert2words(Min,F1,F2)
17165 CASE ELSE
17170 GOTO 17350
17175 END SELECT
17180 Hsend: ! Tell the LVDAS to Take a Histogram.
17185 OUTPUT @Lvdas USING "AA,6(W)";"TH",F1,F2,Bin,A1,A2,Channel
17190 Henter: ! Enter number of bins in the histogram.
17195 ENTER @Lvdas USING "#,W";Nbins
17200 ! Redimension the Histo(*) and the enter the histogram data.
17205 IF Nbins>0 THEN
17210 REDIM Histo(Nbins,3)
17215 ENTER @Lvdas USING "#,W";Histo(*)
17220 END IF
17225 ! Enter the number of samples for the previous and current histogram.
17230 ENTER @Lvdas USING "#,W";Nnew,Nold
17235 Hplot: ! Scale part of the CRT for the histogram plotting.

```

```

17240 VIEWPORT Vwprt (G,1)/10.23,Vwprt (G,2)/10.23,Vwprt (G,3)/10.23,Vwprt (G,4)/10.23
17245 WINDOW Kw*Wdw(G,1),Kw*Wdw(G,2),Wdw(G,3),Wdw(G,4)
17250 Xpixel=Kw*(Wdw(Channel,2)-Wdw(Channel,1))/(Vwprt (Channel,2)-Vwprt (Channel,1))
17255 N1=N(Channel)
17260 N2=N(Channel)-Nold+Nnew
17265 N(Channel)=N(Channel)-Nold+Nnew
17270 PEN 16*Aqua ! Select the pen for the histogram bars edge.
17275 AREA PEN 16*Aqua ! Select the pen for the histogram bars fill.
17280 FOR I=1 TO Nbins
17285 Old=MIN(Histo(I,3),Wdw(Channel,4))
17290 New=MIN(Histo(I,2),Wdw(Channel,4))
17295 X1=Histo(I,1)*Ww*Min ! Calculate histogram bar horizontal position.
17300 X2=Ww ! Calculate histogram bar horizontal width.
17305 Y1=Old ! Calculate histogram bar vertical position.
17310 Y2=New-Old ! Calculate histogram bar vertical width.
17315 IF X1<Kw*Wdw(G,1) THEN X1=Kw*Wdw(G,1) ! If X1<Xmin then set X1=Xmin
17320 IF X1>Kw*Wdw(G,2)-X2 THEN X1=Kw*Wdw(G,2)-X2 ! If X1>Xmax then set X1=Xmax
17325 MOVE X1,Y1
17330 CONTROL CRT,14;6 ! Change to complimentary drawing mode.
17335 RECTANGLE X2-Xpixel,Y2,FILL,EDGE ! Draw the rectangle representing one bar of the bargraph.
17340 CONTROL CRT,14;3 ! Switch back to dominant drawing mode.
17345 NEXT I
17350 NEXT Channel
17355 Kbd$=KBDS$
17360 UNTIL Kbd$<>"** OR NOT Repeat ! Quit if any key on the keyboard has been pressed.
17365 SUBEXIT
17370 Hdone: Done=1
17375 RETURN
17380 SUBEND
17385 Pt_histo: SUB Pt_histo(Symbols(*),Run,File,Pos,INTEGER Nsam)
17390 ! Description:
17395 ! This subprogram plots post time histograms within five of the nine plots on the CRT screen.
17400 ! The histogram data are acquired from the LVDAS over a specified acquisition time.
17405 ! Variables Defined in Main Program:
17410 ! Wdw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*),Title$(*),Ximage$(*),Yimage$(*),Legend$(*)
17415 ! U1(*),V1(*),Wi(*),Ai(*),Bi(*)
17420 ! Local Variables:
17425 ! Histo(*) Array of histogram bin heights indexed by bin number.
17430 ! Data(*) Array of instantaneous U,V,W velocity or A,B voltage data.
17435 ! Nsam Number of samples acquired.
17440 ! Xmin Minimum value for histogram. Left side of histogram scale.
17445 ! Xmax Maximum value for histogram. right side of histogram scale.
17450 ! Xwin Window width of each vertical histogram bar.
17455 ! K Used as an index to the above arrays.
17460 ! L Used as an index to the Histo(*). Specifies one of 100 bins.
17465 ! Xpixel Horizontal length of one picture on the CRT in scale units.
17470 ! Channel Selects one of the 5 channels of U1(*),V1(*),Wi(*),Ai(*),Bi(*) data.
17475 ! G Used as an index to the graphics arrays. Specifies one of nine plots.
17480 OPTION BASE 1
17485 COM /Data2/ REAL U1(1000),V1(1000),Wi(1000),Ai(1000),Bi(1000),I1(1000),C1(1000)
17490 COM /Graph1/ Wdw(*),Vwprt(*),Xdiv(*),Ydiv(*),Xlabel$(*),Ylabel$(*)
17495 COM /Graph2/ Title$(*),Ximage$(*),Yimage$(*),Legend$(*)
17500 COM /Color1/ Clear,Black,Red,Yellow,Green,Cyan,Blue,Magenta
17505 COM /Color2/ White,Olive,Aqua,Royal,Maroon,Brick,Brown,Gray
17510 INTEGER Histo(0:100)
17515 REAL Data(1000)
17520 REDIM Data(Nsam)
17525 FOR Channel=1 TO 5
17530 ! Fill the data array with U1(*),V1(*),Wi(*),Ai(*), or Bi(*) depending on Channel.
17535 G=Channel
17540 IF Channel=1 THEN MAT Data= U1
17545 IF Channel=2 THEN MAT Data= V1
17550 IF Channel=3 THEN MAT Data= Wi
17555 IF Channel=4 THEN MAT Data= Ai
17560 IF Channel=5 THEN MAT Data= Bi
17565 ! Draw a new empty histogram plot.
17570 CALL Set_up(Channel,Symbols(*))
17575 Hsort: Xmin=Wdw(Channel,1)
17580 Xmax=Wdw(Channel,2)
17585 Xwin=(Xmax-Xmin)/100
17590 ! Sort the data into a histogram.
17595 MAT Data= Data-(Xmin)
17600 MAT Data= Data/((Xmax-Xmin)/100)
17605 MAT Histo= (0)
17610 FOR K=1 TO Nsam
17615 L=MAX(MIN(Data(K),100),0)
17620 Histo(L)=Histo(L)+1
17625 NEXT K
17630 Hplot: ! Scale part of the CRT for histogram plotting.
17635 VIEWPORT Vwprt (G,1)/10.23,Vwprt (G,2)/10.23,Vwprt (G,3)/10.23,Vwprt (G,4)/10.23

```

```

17640      WINDOW 0,100,Wndw(G,3),Wndw(G,4)
17645      Xpixel=(100-0)/(Vwprt(Channel,2)-Vwprt(Channel,1))
17650      ! Draw the histogram.
17655      FOR K=0 TO 100
17660          IF Histo(K) THEN
17665              MOVE K-.5,0
17670              AREA PEN 16*Green
17675              RECTANGLE 1-Xpixel,Histo(K),FILL
17680          END IF
17685      NEXT K
17690      NEXT Channel
17695      SUBEXIT
17700  SUBEND
17705  CSUB Gdump_colored(From_ds,To_ds,OPTIONAL Rotates$,INTEGER Resolution,Background$,Algorithm$)
17710  CSUB Bload(INTEGER A(*),Xpixels,Ypixels,OPTIONAL INTEGER Rule,REAL Xstart,Ystart)
17715  CSUB Bstore(INTEGER A(*),Xpixels,Ypixels,OPTIONAL INTEGER Rule,REAL Xstart,Ystart)
DONE

```

APPENDIX C

DATA REDUCTION AND COORDINATE SYSTEM TRANSFORMATION EQUATIONS.

APPENDIX C

Data Reduction and Coordinate System Transformation Equations.

| SECTION | <u>TABLE OF CONTENTS</u> | Page |
|----------------|--|-------------|
| 1. | Introduction | 2 |
| 2. | List of Variables | 4 |
| 3. | Laser to Tunnel Coordinate System Transformation | 6 |
| 4. | Tunnel to Model Coordinate System Transformation | 10 |
| 5. | Data Reduction Equations | 13 |
| 6. | Coordinate System Transformation Equations | 18 |
| 7. | Proofs for Coordinate System Transformation Equations | 20 |
| 8. | Matrix Notation for Coordinate System Transformation Equations | 24 |

| FIGURE | <u>LIST OF FIGURES</u> | Page |
|---------------|---|-------------|
| 1a. | U Velocity Laser Beam Pair Orientation (Violet) | 6 |
| 1b. | V Velocity Laser Beam Pair Orientation (Green) | 7 |
| 1c. | W Velocity Laser Beam Pair Orientation (Blue) | 7 |

1. Introduction

The purpose of this write-up is to describe the data reduction performed on the raw data acquired from the Laser Velocimeter Data Acquisition System. The digital Macrodyne data are converted from raw 16bit integer words into frequencies. The frequency results are in turn converted into particle velocities. The analog data are converted from raw two's complement 16bit integer words into voltages. Example types for the analog data might originate from such sources as temperature probes, laser fluorescence anemometers, hot wire anemometers, etc.

Section 2 contains a list of variables that are used throughout this write-up. A brief description of each variable is provided along with the corresponding variable name that is used in the software program. (NOTE: This chapter has been written for three component LDV systems. The delivered system is a two component system. Therefore, the third component W of the three components U,V,W is not measured.)

Velocities are measured in "Laser Coordinates" directly. That is, the measured velocity of each component is parallel to a vector which is orthogonal to the fringe planes in the probe volume. These vectors may or may not be parallel to the tunnel coordinate system. If they are not, then it is desirable to convert the velocities from "Laser Coordinates" to "Tunnel Coordinates." In other words, a coordinate system transformation needs to be applied to the measured velocities to obtain velocities in tunnel coordinates. Section 3 describes how this laser to tunnel coordinate system transformation is performed. (NOTE: The delivered system is a two component system whose laser beam pairs have been orientated orthogonally to the wind tunnel's X,Y,Z axes. Therefore, velocities measured in "Laser Coordinates" will be equal to velocities transformed to "Tunnel Coordinates". For this reason, the the "Laser" to "Tunnel" coordinate system transformations have not been included in the "3.5'HWT" data acquisition programs listed in Appendixes A & B.)

In some cases it is preferred to perform an additional coordinate system transformation

to obtain velocities in “Model Coordinates.” For example, if the model is at an angle of attack, then the model’s coordinate system would be at rotation with respect to the tunnel’s coordinate system. Other model attitude angles in addition to the angle of attack, such as roll and yaw, can be used determine the transformation required to convert from tunnel to model coordinates. Section 4 describes how this tunnel to model coordinate system transformation is performed.

Section 5 contains the equations that are used to calculate the average, standard deviation, as well as normal and shear stress terms for the velocity and voltage data. Equations are included for both the original and transformed to coordinate systems. Normal text is used in the equations for variables that represent the original coordinate system while italicized text is used for variables that represent the transformed to coordinate system.

Section 6 contains the equations that are used to convert average, standard deviation, as well as normal and shear stress terms from the original to the transformed to coordinate system.

Section 7 contains proofs demonstrating that we can perform the coordinate system transformations on the reduced averaged data without having to perform the transformation on the instantaneous values. This saves costly run time because there are typically thousands of instantaneous values that contribute one averaged value.

Section 8 shows how the equations of section 6 can be represented in matrix notation. The matrix notation for the coordinate system transformation is an elegant way to show the multitude of complex equations in compact and concise format.

2. List of Variables

The following is a list of variables that are used throughout this write-up. A brief description of each variable is provided along with the corresponding variable name that is used in the software program. Normal text style (not italicized) indicates the original coordinate system while italicized text style indicates a transformed to coordinate system.


| <u>Original</u> | <u>Transformed</u> | <u>Description</u> | <u>Variable</u> |
|-------------------|--------------------|--------------------------------|-----------------|
| U_i | U_i | Instantaneous U velocity. | U_i (I) |
| V_i | V_i | Instantaneous V velocity. | V_i (I) |
| W_i | W_i | Instantaneous W velocity. | W_i (I) |
| A_i | A_i | Instantaneous A voltage. | A_i (I) |
| B_i | B_i | Instantaneous B voltage. | B_i (I) |
| \bar{U} | \bar{U} | Average U velocity. | U |
| \bar{V} | \bar{V} | Average V velocity. | V |
| \bar{W} | \bar{W} | Average W velocity. | W |
| \bar{A} | \bar{A} | Average A voltage. | A |
| \bar{B} | \bar{B} | Average B voltage. | B |
| U' | U' | U velocity standard deviation. | U_1 |
| V' | V' | V velocity standard deviation. | V_1 |
| W' | W' | W velocity standard deviation. | W_1 |
| A' | A' | A voltage standard deviation. | A_1 |
| B' | B' | B voltage standard deviation. | B_1 |
| $\overline{A'A'}$ | $\overline{A'A'}$ | A-A normal stress term. | A_1a_1 |
| $\overline{B'B'}$ | $\overline{B'B'}$ | B-B normal stress term. | B_1b_1 |
| $\overline{A'B'}$ | $\overline{A'B'}$ | A-B shear stress term. | A_1b_1 |
| $\overline{U'A'}$ | $\overline{U'A'}$ | U-A shear stress term. | U_1a_1 |
| $\overline{V'A'}$ | $\overline{V'A'}$ | V-A shear stress term. | V_1a_1 |
| $\overline{W'A'}$ | $\overline{W'A'}$ | W-A shear stress term. | W_1a_1 |

| <u>Original</u> | <u>Transformed</u> | <u>Description</u> | <u>Variable</u> |
|-------------------------|---------------------|---|-----------------|
| $\overline{U'U'}$ | $\overline{U''U''}$ | U-U normal stress term. | U1u1 |
| $\overline{U'V'}$ | $\overline{U''V''}$ | U-V shear stress term. | U1v1 |
| $\overline{U'W'}$ | $\overline{U''W''}$ | U-W shear stress term. | U1w1 |
| $\overline{V'U'}$ | $\overline{V''U''}$ | V-U shear stress term. | V1u1 |
| $\overline{V'V'}$ | $\overline{V''V''}$ | V-V normal stress term. | V1v1 |
| $\overline{V'W'}$ | $\overline{V''W''}$ | V-W shear stress term. | V1w1 |
| $\overline{W'U'}$ | $\overline{W''U''}$ | W-U shear stress term. | W1u1 |
| $\overline{W'V'}$ | $\overline{W''V''}$ | W-V shear stress term. | W1v1 |
| $\overline{W'W'}$ | $\overline{W''W''}$ | W-W normal stress term. | W1w1 |
| X | | Tunnel or Model X axis | X |
| Y | | Tunnel or Model Y axis | Y |
| Z | | Tunnel or Model Z axis | Z |
| \emptyset_{UX} | | Angle between Laser U and Tunnel X | ThetaAU |
| \emptyset_{UY} | | Angle between Laser U and Tunnel Y | ThetaAV |
| \emptyset_{UZ} | | Angle between Laser U and Tunnel Z | ThetaAW |
| \emptyset_{VX} | | Angle between Laser V and Tunnel X | ThetaBU |
| \emptyset_{VY} | | Angle between Laser V and Tunnel Y | ThetaBV |
| \emptyset_{VZ} | | Angle between Laser V and Tunnel Z | ThetaBW |
| \emptyset_{WX} | | Angle between Laser W and Tunnel X | ThetaCU |
| \emptyset_{WY} | | Angle between Laser W and Tunnel Y | ThetaCV |
| \emptyset_{WZ} | | Angle between Laser W and Tunnel Z | ThetaCW |
| α_1 | | Model angle of attack | Alpha (1) |
| α_2 | | Model angle of roll | Alpha (2) |
| α_3 | | Model angle of yaw | Alpha (3) |
| J _{3x3} | | Coordinate system transformation matrix | --- |
| K _{3x3} | | Coordinate system transformation matrix | --- |
| K _{9x9} | | Coordinate system transformation matrix | --- |


3. Laser to Tunnel Coordinate System Transformation

Velocities are measured in "Laser Coordinates" directly. That is, the measured velocity of each component is parallel to a vector which is orthogonal to the fringe planes in the probe volume. These vectors may or may not be parallel to the tunnel coordinate system. If they are not, then it is desirable to convert the velocities from "Laser Coordinates" to "Tunnel Coordinates." In other words, a coordinate system transformation needs to be applied to the measured velocities to obtain velocities in tunnel coordinates. This section describes how this laser to tunnel coordinate system transformation is performed.


Tunnel Coordinate have the axes label as X,Y,Z while velocities measured in laser coordinates typically named U,V,W. The direction of each of the measured velocities in laser coordinates can be defined in terms of the angle it is off of the tunnel coordinate axes. The three angles ϕ_{UX} , ϕ_{UY} , ϕ_{UZ} define the angular relationship between measured U velocities in laser coordinates and the axes X,Y,Z of the tunnel coordinate system (Fig. 1a).



The three angles θ_{VX} , θ_{VY} , θ_{VZ} define the angular relationship between measured V velocities in laser coordinates and the axes X,Y,Z of the tunnel coordinate system (Fig. 1b)



The three angles θ_{WX} , θ_{WY} , θ_{WZ} define the angular relationship between measured W velocities in laser coordinates and the axes X,Y,Z of the tunnel coordinate system (Fig. 1c)



When a particle travels through the probe volume, its velocity is measured as U,V,W in

laser coordinates. However, it is desired to have these velocities (U,V,W) transformed to tunnel coordinate velocities (U,V,W). Each of tunnel coordinate velocities U,V,W would be parallel to its X,Y,Z tunnel axis. The laser coordinate velocities U,V,W can be defined in terms of the tunnel coordinate velocities using the follow equations:

$$\begin{aligned} U &= U \cos(\theta_{UX}) + V \cos(\theta_{UY}) + W \cos(\theta_{UZ}) \\ V &= U \cos(\theta_{VX}) + V \cos(\theta_{VY}) + W \cos(\theta_{VZ}) \\ W &= U \cos(\theta_{WX}) + V \cos(\theta_{WY}) + W \cos(\theta_{WZ}) \end{aligned}$$

The coefficients of the of the three equations above can be used to define the coordinate transformation matrix $\mathbf{J}_{3 \times 3}$ as shown below:

$$\begin{aligned} J_{11} &= \cos(\theta_{UX}) & J_{12} &= \cos(\theta_{UY}) & J_{13} &= \cos(\theta_{UZ}) \\ J_{21} &= \cos(\theta_{VX}) & J_{22} &= \cos(\theta_{VY}) & J_{23} &= \cos(\theta_{VZ}) \\ J_{31} &= \cos(\theta_{WX}) & J_{32} &= \cos(\theta_{WY}) & J_{33} &= \cos(\theta_{WZ}) \end{aligned}$$

$$\mathbf{J}_{3 \times 3} = \begin{bmatrix} J_{11} & J_{12} & J_{13} \\ J_{21} & J_{22} & J_{23} \\ J_{31} & J_{32} & J_{33} \end{bmatrix} = \begin{bmatrix} \cos(\theta_{UX}) & \cos(\theta_{UY}) & \cos(\theta_{UZ}) \\ \cos(\theta_{VX}) & \cos(\theta_{VY}) & \cos(\theta_{VZ}) \\ \cos(\theta_{WX}) & \cos(\theta_{WY}) & \cos(\theta_{WZ}) \end{bmatrix}$$

The coordinate transformation matrix $\mathbf{J}_{3 \times 3}$ can be used to convert tunnel coordinate velocities to laser coordinate velocities.

$$\begin{bmatrix} U \\ V \\ W \end{bmatrix} = \begin{bmatrix} J_{11} & J_{12} & J_{13} \\ J_{21} & J_{22} & J_{23} \\ J_{31} & J_{32} & J_{33} \end{bmatrix} \times \begin{bmatrix} U \\ V \\ W \end{bmatrix}$$

$$\begin{aligned} U &= J_{11}U + J_{12}V + J_{13}W \\ V &= J_{21}U + J_{22}V + J_{23}W \\ W &= J_{31}U + J_{32}V + J_{33}W \end{aligned}$$

However, we need to perform just the opposite coordinate transformation. The coordinate transformation matrix $\mathbf{K}_{3 \times 3}$ is defined as inverse of transformation matrix $\mathbf{J}_{3 \times 3}$.

$$\mathbf{K}_{3 \times 3} = \mathbf{J}_{3 \times 3}^{-1} = \begin{bmatrix} J_{11} & J_{12} & J_{13} \\ J_{21} & J_{22} & J_{23} \\ J_{31} & J_{32} & J_{33} \end{bmatrix}^{-1} = \begin{bmatrix} K_{11} & K_{12} & K_{13} \\ K_{21} & K_{22} & K_{23} \\ K_{31} & K_{32} & K_{33} \end{bmatrix}$$

The coordinate transformation matrix $\mathbf{K}_{3 \times 3}$ can be used to convert laser coordinate velocities to tunnel coordinate velocities.

$$\begin{bmatrix} U \\ V \\ W \end{bmatrix} = \begin{bmatrix} K_{11} & K_{12} & K_{13} \\ K_{21} & K_{22} & K_{23} \\ K_{31} & K_{32} & K_{33} \end{bmatrix} \times \begin{bmatrix} U \\ V \\ W \end{bmatrix}$$

$$U = K_{11}U + K_{12}V + K_{13}W$$

$$V = K_{21}U + K_{22}V + K_{23}W$$

$$W = K_{31}U + K_{32}V + K_{33}W$$

4. Tunnel to Model Coordinate System Transformation

In some cases it is preferred to perform an additional coordinate system transformation to obtain velocities in "Model Coordinates." For example, if the model is at an angle of attack, then the model's coordinate system would be at rotation with respect to the tunnel's coordinate system. Other model attitude angles in addition to the angle of attack, such as roll and yaw, can be used determine the transformation required to convert from tunnel to model coordinates. This section describes how this tunnel to model coordinate system transformation is performed.

The angle of attack, roll, and yaw angles are defined as follows:

α_1 angle of attack

α_2 angle of roll

α_3 angle of yaw

Velocities that are calculated in tunnel coordinates U, V, W can be transformed to model coordinate velocities (U, V, W). Each of the tunnel coordinate velocities U, V, W are parallel to the tunnel's X, Y, Z axes. Each of them can be transformed to model coordinates where each of the model coordinate velocities U, V, W would be parallel to the model's X, Y, Z axes.

If the model were at angle of attack ($\alpha_1 \neq 0$), then the tunnel coordinate velocities can be defined in terms of the model coordinate velocities using the follow equations:

$$\begin{aligned} U &= + \cos(\alpha_1) U + 0 V + \sin(\alpha_1) W \\ V &= + 0 U + 1 V + 0 W \\ W &= - \sin(\alpha_1) U + 0 V + \cos(\alpha_1) W \end{aligned}$$

If the model were at angle of roll ($\alpha_2 \neq 0$), then the tunnel coordinate velocities can be defined in terms of the model coordinate velocities using the follow equations:

$$\begin{aligned} U &= + 1 U + 0 V + 0 W \\ V &= + 0 U + \cos(\alpha_2) V - \sin(\alpha_2) W \\ W &= + 0 U + \sin(\alpha_2) V + \cos(\alpha_2) W \end{aligned}$$

If the model were at angle of yaw ($\alpha_3 \neq 0$), then the tunnel coordinate velocities can be defined

in terms of the model coordinate velocities using the follow equations:

$$\begin{aligned} U &= + \cos(\alpha_3) U - \sin(\alpha_3) V + W \\ V &= + \sin(\alpha_3) U + \cos(\alpha_3) V + 0 W \\ W &= + U + 0 V + W \end{aligned}$$

The model to tunnel coordinate system transformation matrices for each of the three sets of equations are defined as follows:

$$\mathbf{J}_{\alpha_1} = \begin{bmatrix} +\cos(\alpha_1) & 0 & +\sin(\alpha_1) \\ 0 & 1 & 0 \\ -\sin(\alpha_1) & 0 & +\cos(\alpha_1) \end{bmatrix}$$

$$\mathbf{J}_{\alpha_2} = \begin{bmatrix} 1 & 0 & 0 \\ 0 & +\cos(\alpha_2) & -\sin(\alpha_2) \\ 0 & +\sin(\alpha_2) & +\cos(\alpha_2) \end{bmatrix}$$

$$\mathbf{J}_{\alpha_3} = \begin{bmatrix} +\cos(\alpha_3) & -\sin(\alpha_3) & 0 \\ +\sin(\alpha_3) & +\cos(\alpha_3) & 0 \\ 0 & 0 & 1 \end{bmatrix}$$

If the angles of attack, roll, and yaw are used in combination, then an equivalent model to tunnel coordinate system transformation matrix can be obtained by computing the cross products of the three individual transformations.

$$\mathbf{J}_{3 \times 3} = \mathbf{J}_{\alpha_3} \times \mathbf{J}_{\alpha_2} \times \mathbf{J}_{\alpha_1}$$

$$\mathbf{J}_{3 \times 3} = \begin{bmatrix} J_{11} & J_{12} & J_{13} \\ J_{21} & J_{22} & J_{23} \\ J_{31} & J_{32} & J_{33} \end{bmatrix}$$

$$\mathbf{J}_{3 \times 3} = \begin{bmatrix} +\cos(\alpha_3) & -\sin(\alpha_3) & 0 \\ +\sin(\alpha_3) & +\cos(\alpha_3) & 0 \\ 0 & 0 & 1 \end{bmatrix} \times \begin{bmatrix} 1 & 0 & 0 \\ 0 & +\cos(\alpha_2) & -\sin(\alpha_2) \\ 0 & +\sin(\alpha_2) & +\cos(\alpha_2) \end{bmatrix} \times \begin{bmatrix} +\cos(\alpha_1) & 0 & +\sin(\alpha_1) \\ 0 & 1 & 0 \\ -\sin(\alpha_1) & 0 & +\cos(\alpha_1) \end{bmatrix}$$

The coordinate transformation matrix $\mathbf{J}_{3 \times 3}$ can be used to convert model coordinate velocities to tunnel coordinate velocities.

$$\begin{bmatrix} U \\ V \\ W \end{bmatrix} = \begin{bmatrix} J_{11} & J_{12} & J_{13} \\ J_{21} & J_{22} & J_{23} \\ J_{31} & J_{32} & J_{33} \end{bmatrix} \times \begin{bmatrix} U \\ V \\ W \end{bmatrix}$$

$$U = J_{11}U + J_{12}V + J_{13}W$$

$$V = J_{21}U + J_{22}V + J_{23}W$$

$$W = J_{31}U + J_{32}V + J_{33}W$$

However, we need to perform just the opposite coordinate transformation. The coordinate transformation matrix $\mathbf{K}_{3 \times 3}$ is defined as inverse of transformation matrix $\mathbf{J}_{3 \times 3}$.

$$\mathbf{K}_{3 \times 3} = \mathbf{J}_{3 \times 3}^{-1} = \begin{bmatrix} J_{11} & J_{12} & J_{13} \\ J_{21} & J_{22} & J_{23} \\ J_{31} & J_{32} & J_{33} \end{bmatrix}^{-1} = \begin{bmatrix} K_{11} & K_{12} & K_{13} \\ K_{21} & K_{22} & K_{23} \\ K_{31} & K_{32} & K_{33} \end{bmatrix}$$

The coordinate transformation matrix $\mathbf{K}_{3 \times 3}$ can be used to convert tunnel coordinate velocities to model coordinate velocities.

$$\begin{bmatrix} U \\ V \\ W \end{bmatrix} = \begin{bmatrix} K_{11} & K_{12} & K_{13} \\ K_{21} & K_{22} & K_{23} \\ K_{31} & K_{32} & K_{33} \end{bmatrix} \times \begin{bmatrix} U \\ V \\ W \end{bmatrix}$$

$$U = K_{11}U + K_{12}V + K_{13}W$$

$$V = K_{21}U + K_{22}V + K_{23}W$$

$$W = K_{31}U + K_{32}V + K_{33}W$$

5. Data Reduction Equations

This section contains the equations that are used to calculate the average, standard deviation as well as normal and shear stress terms for the velocity and voltage data. Equations are included for both the original and transformed to coordinate systems. Normal text is used in the equations for variables that represent the original coordinate system while italicized text is used for variables that represent the transformed to coordinate system.

The following equations are used to calculate the velocity and voltage averages:

$$\begin{aligned}\bar{U} &= \frac{\sum_{i=1}^N U_i}{N} & \bar{V} &= \frac{\sum_{i=1}^N V_i}{N} & \bar{W} &= \frac{\sum_{i=1}^N W_i}{N} & \bar{A} &= \frac{\sum_{i=1}^N A_i}{N} & \bar{B} &= \frac{\sum_{i=1}^N B_i}{N} \\ \bar{U} &= \frac{\sum_{i=1}^N U_i}{N} & \bar{V} &= \frac{\sum_{i=1}^N V_i}{N} & \bar{W} &= \frac{\sum_{i=1}^N W_i}{N} & \bar{A} &= \frac{\sum_{i=1}^N A_i}{N} & \bar{B} &= \frac{\sum_{i=1}^N B_i}{N}\end{aligned}$$

The following equations are used to calculate their standard deviations:

$$\begin{aligned}U' &= \sqrt{\frac{\sum_{i=1}^N (U_i - \bar{U})^2}{N}} = \sqrt{\frac{\sum_{i=1}^N U_i^2}{N} - \bar{U}^2} & U' &= \sqrt{\frac{\sum_{i=1}^N (U_i - \bar{U})^2}{N}} = \sqrt{\frac{\sum_{i=1}^N U_i^2}{N} - \bar{U}^2} \\ V' &= \sqrt{\frac{\sum_{i=1}^N (V_i - \bar{V})^2}{N}} = \sqrt{\frac{\sum_{i=1}^N V_i^2}{N} - \bar{V}^2} & V' &= \sqrt{\frac{\sum_{i=1}^N (V_i - \bar{V})^2}{N}} = \sqrt{\frac{\sum_{i=1}^N V_i^2}{N} - \bar{V}^2} \\ W' &= \sqrt{\frac{\sum_{i=1}^N (W_i - \bar{W})^2}{N}} = \sqrt{\frac{\sum_{i=1}^N W_i^2}{N} - \bar{W}^2} & W' &= \sqrt{\frac{\sum_{i=1}^N (W_i - \bar{W})^2}{N}} = \sqrt{\frac{\sum_{i=1}^N W_i^2}{N} - \bar{W}^2} \\ A' &= \sqrt{\frac{\sum_{i=1}^N (A_i - \bar{A})^2}{N}} = \sqrt{\frac{\sum_{i=1}^N A_i^2}{N} - \bar{A}^2} & A' &= \sqrt{\frac{\sum_{i=1}^N (A_i - \bar{A})^2}{N}} = \sqrt{\frac{\sum_{i=1}^N A_i^2}{N} - \bar{A}^2} \\ B' &= \sqrt{\frac{\sum_{i=1}^N (B_i - \bar{B})^2}{N}} = \sqrt{\frac{\sum_{i=1}^N B_i^2}{N} - \bar{B}^2} & B' &= \sqrt{\frac{\sum_{i=1}^N (B_i - \bar{B})^2}{N}} = \sqrt{\frac{\sum_{i=1}^N B_i^2}{N} - \bar{B}^2}\end{aligned}$$

The following equations are used to calculate the normal and shear stress terms for all of the relevant velocity:velocity, velocity:voltage, and voltage:voltage combinations:

$$\overline{U'U'} = \frac{\sum_{i=1}^N (U_i - \bar{U})(U_i - \bar{U})}{N} = \frac{\sum_{i=1}^N U_i U_i}{N} - \bar{U} \bar{U} = \frac{\sum_{i=1}^N U_i^2}{N} - \bar{U}^2$$

$$\overline{U'V'} = \frac{\sum_{i=1}^N (U_i - \bar{U})(V_i - \bar{V})}{N} = \frac{\sum_{i=1}^N U_i V_i}{N} - \bar{U} \bar{V}$$

$$\overline{U'W'} = \frac{\sum_{i=1}^N (U_i - \bar{U})(W_i - \bar{W})}{N} = \frac{\sum_{i=1}^N U_i W_i}{N} - \bar{U} \bar{W}$$

$$\overline{V'U'} = \frac{\sum_{i=1}^N (V_i - \bar{V})(U_i - \bar{U})}{N} = \frac{\sum_{i=1}^N V_i U_i}{N} - \bar{V} \bar{U}$$

$$\overline{V'V'} = \frac{\sum_{i=1}^N (V_i - \bar{V})(V_i - \bar{V})}{N} = \frac{\sum_{i=1}^N V_i V_i}{N} - \bar{V} \bar{V} = \frac{\sum_{i=1}^N V_i^2}{N} - \bar{V}^2$$

$$\overline{V'W'} = \frac{\sum_{i=1}^N (V_i - \bar{V})(W_i - \bar{W})}{N} = \frac{\sum_{i=1}^N V_i W_i}{N} - \bar{V} \bar{W}$$

$$\overline{W'U'} = \frac{\sum_{i=1}^N (W_i - \bar{W})(U_i - \bar{U})}{N} = \frac{\sum_{i=1}^N W_i U_i}{N} - \bar{W} \bar{U}$$

$$\overline{W'V'} = \frac{\sum_{i=1}^N (W_i - \bar{W})(V_i - \bar{V})}{N} = \frac{\sum_{i=1}^N W_i V_i}{N} - \bar{W} \bar{V}$$

$$\overline{W'W'} = \frac{\sum_{i=1}^N (W_i - \bar{W})(W_i - \bar{W})}{N} = \frac{\sum_{i=1}^N W_i W_i}{N} - \bar{W} \bar{W} = \frac{\sum_{i=1}^N W_i^2}{N} - \bar{W}^2$$

$$\overline{U'U'} = \frac{\sum_{i=1}^N (U_i - \bar{U})(U_i - \bar{U})}{N} = \frac{\sum_{i=1}^N U_i U_i}{N} - \bar{U} \bar{U} = \frac{\sum_{i=1}^N U_i^2}{N} - \bar{U}^2$$

$$\overline{U'V'} = \frac{\sum_{i=1}^N (U_i - \bar{U})(V_i - \bar{V})}{N} = \frac{\sum_{i=1}^N U_i V_i}{N} - \bar{U} \bar{V}$$

$$\overline{U'W'} = \frac{\sum_{i=1}^N (U_i - \bar{U})(W_i - \bar{W})}{N} = \frac{\sum_{i=1}^N U_i W_i}{N} - \bar{U} \bar{W}$$

$$\overline{V'U'} = \frac{\sum_{i=1}^N (V_i - \bar{V})(U_i - \bar{U})}{N} = \frac{\sum_{i=1}^N V_i U_i}{N} - \bar{V} \bar{U}$$

$$\overline{V'V'} = \frac{\sum_{i=1}^N (V_i - \bar{V})(V_i - \bar{V})}{N} = \frac{\sum_{i=1}^N V_i V_i}{N} - \bar{V} \bar{V} = \frac{\sum_{i=1}^N V_i^2}{N} - \bar{V}^2$$

$$\overline{V'W'} = \frac{\sum_{i=1}^N (V_i - \bar{V})(W_i - \bar{W})}{N} = \frac{\sum_{i=1}^N V_i W_i}{N} - \bar{V} \bar{W}$$

$$\overline{W'U'} = \frac{\sum_{i=1}^N (W_i - \bar{W})(U_i - \bar{U})}{N} = \frac{\sum_{i=1}^N W_i U_i}{N} - \bar{W} \bar{U}$$

$$\overline{W'V'} = \frac{\sum_{i=1}^N (W_i - \bar{W})(V_i - \bar{V})}{N} = \frac{\sum_{i=1}^N W_i V_i}{N} - \bar{W} \bar{V}$$

$$\overline{W'W'} = \frac{\sum_{i=1}^N (W_i - \bar{W})(W_i - \bar{W})}{N} = \frac{\sum_{i=1}^N W_i W_i}{N} - \bar{W} \bar{W} = \frac{\sum_{i=1}^N W_i^2}{N} - \bar{W}^2$$

$$\overline{A'A'} = \frac{\sum_{i=1}^N (A_i - \bar{A})(A_i - \bar{A})}{N} = \frac{\sum_{i=1}^N A_i A_i}{N} - \bar{A} \bar{A} = \frac{\sum_{i=1}^N A_i^2}{N} - \bar{A}^2$$

$$\overline{B'B'} = \frac{\sum_{i=1}^N (B_i - \bar{B})(B_i - \bar{B})}{N} = \frac{\sum_{i=1}^N B_i B_i}{N} - \bar{B} \bar{B} = \frac{\sum_{i=1}^N B_i^2}{N} - \bar{B}^2$$

$$\overline{A'B'} = \frac{\sum_{i=1}^N (A_i - \bar{A})(B_i - \bar{B})}{N} = \frac{\sum_{i=1}^N A_i B_i}{N} - \bar{A} \bar{B}$$

$$\overline{U'A'} = \frac{\sum_{i=1}^N (U_i - \bar{U})(A_i - \bar{A})}{N} = \frac{\sum_{i=1}^N U_i A_i}{N} - \bar{U} \bar{A}$$

$$\overline{V'A'} = \frac{\sum_{i=1}^N (V_i - \bar{V})(A_i - \bar{A})}{N} = \frac{\sum_{i=1}^N V_i A_i}{N} - \bar{V} \bar{A}$$

$$\overline{W'A'} = \frac{\sum_{i=1}^N (W_i - \bar{W})(A_i - \bar{A})}{N} = \frac{\sum_{i=1}^N W_i A_i}{N} - \bar{W} \bar{A}$$

$$\overline{A'A'} = \frac{\sum_{i=1}^N (A_i - \bar{A})(A_i - \bar{A})}{N} = \frac{\sum_{i=1}^N A_i A_i}{N} - \bar{A} \bar{A} = \frac{\sum_{i=1}^N A_i^2}{N} - \bar{A}^2$$

$$\overline{B'B'} = \frac{\sum_{i=1}^N (B_i - \bar{B})(B_i - \bar{B})}{N} = \frac{\sum_{i=1}^N B_i B_i}{N} - \bar{B} \bar{B} = \frac{\sum_{i=1}^N B_i^2}{N} - \bar{B}^2$$

$$\overline{A'B'} = \frac{\sum_{i=1}^N (A_i - \bar{A})(B_i - \bar{B})}{N} = \frac{\sum_{i=1}^N A_i B_i}{N} - \bar{A} \bar{B}$$

$$\overline{U'A'} = \frac{\sum_{i=1}^N (U_i - \bar{U})(A_i - \bar{A})}{N} = \frac{\sum_{i=1}^N U_i A_i}{N} - \bar{U} \bar{A}$$

$$\overline{V'A'} = \frac{\sum_{i=1}^N (V_i - \bar{V})(A_i - \bar{A})}{N} = \frac{\sum_{i=1}^N V_i A_i}{N} - \bar{V} \bar{A}$$

$$\overline{W'A'} = \frac{\sum_{i=1}^N (W_i - \bar{W})(A_i - \bar{A})}{N} = \frac{\sum_{i=1}^N W_i A_i}{N} - \bar{W} \bar{A}$$

6. Coordinate System Transformation Equations

The following equations are used to convert the instantaneous as well as the average velocities from the original to the transformed to coordinate system. The instantaneous and average values for the voltage data are the same in either coordinate systems:

$$U_i = K_{11}U_i + K_{12}V_i + K_{13}W_i \quad \bar{U} = K_{11}\bar{U} + K_{12}\bar{V} + K_{13}\bar{W}$$

$$V_i = K_{21}U_i + K_{22}V_i + K_{23}W_i \quad \bar{V} = K_{21}\bar{U} + K_{22}\bar{V} + K_{23}\bar{W}$$

$$W_i = K_{31}U_i + K_{32}V_i + K_{33}W_i \quad \bar{W} = K_{31}\bar{U} + K_{32}\bar{V} + K_{33}\bar{W}$$

$$A_i = A_i \quad \bar{A} = \bar{A}$$

$$B_i = B_i \quad \bar{B} = \bar{B}$$

The following equations are used to convert the velocity:voltage and voltage:voltage normal and shear stress terms from the original to the transformed to coordinate system.

The voltage:voltage normal and shear stress terms are the same in either coordinate systems:

$$\overline{U'A'} = K_{11}\overline{U'A'} + K_{12}\overline{V'A'} + K_{13}\overline{W'A'}$$

$$\overline{V'A'} = K_{21}\overline{U'A'} + K_{22}\overline{V'A'} + K_{23}\overline{W'A'}$$

$$\overline{W'A'} = K_{31}\overline{U'A'} + K_{32}\overline{V'A'} + K_{33}\overline{W'A'}$$

$$\overline{A'A'} = \overline{A'A'}$$

$$\overline{B'B'} = \overline{B'B'}$$

$$\overline{A'B'} = \overline{A'B'}$$

The following equations are used to convert the velocity:velocity normal and shear stress terms from the original to the transformed to coordinate system:

$$\begin{aligned}\overline{U'U'} &= K_{11}K_{11}\overline{U'U'} + K_{11}K_{12}\overline{U'V'} + K_{11}K_{13}\overline{U'W'} \\ &+ K_{12}K_{11}\overline{V'U'} + K_{12}K_{12}\overline{V'V'} + K_{12}K_{13}\overline{V'W'} \\ &+ K_{13}K_{11}\overline{W'U'} + K_{13}K_{12}\overline{W'V'} + K_{13}K_{13}\overline{W'W'}\end{aligned}$$

$$\begin{aligned}\overline{U'V'} &= K_{11}K_{21}\overline{U'U'} + K_{11}K_{22}\overline{U'V'} + K_{11}K_{23}\overline{U'W'} \\ &+ K_{12}K_{21}\overline{V'U'} + K_{12}K_{22}\overline{V'V'} + K_{12}K_{23}\overline{V'W'} \\ &+ K_{13}K_{21}\overline{W'U'} + K_{13}K_{22}\overline{W'V'} + K_{13}K_{23}\overline{W'W'}\end{aligned}$$

$$\begin{aligned}\overline{U'W'} &= K_{11}K_{31}\overline{U'U'} + K_{11}K_{32}\overline{U'V'} + K_{11}K_{33}\overline{U'W'} \\ &+ K_{12}K_{31}\overline{V'U'} + K_{12}K_{32}\overline{V'V'} + K_{12}K_{33}\overline{V'W'} \\ &+ K_{13}K_{31}\overline{W'U'} + K_{13}K_{32}\overline{W'V'} + K_{13}K_{33}\overline{W'W'}\end{aligned}$$

$$\begin{aligned}\overline{V'U'} &= K_{21}K_{11}\overline{U'U'} + K_{21}K_{12}\overline{U'V'} + K_{21}K_{13}\overline{U'W'} \\ &+ K_{22}K_{11}\overline{V'U'} + K_{22}K_{12}\overline{V'V'} + K_{22}K_{13}\overline{V'W'} \\ &+ K_{23}K_{11}\overline{W'U'} + K_{23}K_{12}\overline{W'V'} + K_{23}K_{13}\overline{W'W'}\end{aligned}$$

$$\begin{aligned}\overline{V'V'} &= K_{21}K_{21}\overline{U'U'} + K_{21}K_{22}\overline{U'V'} + K_{21}K_{23}\overline{U'W'} \\ &+ K_{22}K_{21}\overline{V'U'} + K_{22}K_{22}\overline{V'V'} + K_{22}K_{23}\overline{V'W'} \\ &+ K_{23}K_{21}\overline{W'U'} + K_{23}K_{22}\overline{W'V'} + K_{23}K_{23}\overline{W'W'}\end{aligned}$$

$$\begin{aligned}\overline{V'W'} &= K_{21}K_{31}\overline{U'U'} + K_{21}K_{32}\overline{U'V'} + K_{21}K_{33}\overline{U'W'} \\ &+ K_{22}K_{31}\overline{V'U'} + K_{22}K_{32}\overline{V'V'} + K_{22}K_{33}\overline{V'W'} \\ &+ K_{23}K_{31}\overline{W'U'} + K_{23}K_{32}\overline{W'V'} + K_{23}K_{33}\overline{W'W'}\end{aligned}$$

$$\begin{aligned}\overline{W'U'} &= K_{31}K_{11}\overline{U'U'} + K_{31}K_{12}\overline{U'V'} + K_{31}K_{13}\overline{U'W'} \\ &+ K_{32}K_{11}\overline{V'U'} + K_{32}K_{12}\overline{V'V'} + K_{32}K_{13}\overline{V'W'} \\ &+ K_{33}K_{11}\overline{W'U'} + K_{33}K_{12}\overline{W'V'} + K_{33}K_{13}\overline{W'W'}\end{aligned}$$

$$\begin{aligned}\overline{W'V'} &= K_{31}K_{21}\overline{U'U'} + K_{31}K_{22}\overline{U'V'} + K_{31}K_{23}\overline{U'W'} \\ &+ K_{32}K_{21}\overline{V'U'} + K_{32}K_{22}\overline{V'V'} + K_{32}K_{23}\overline{V'W'} \\ &+ K_{33}K_{21}\overline{W'U'} + K_{33}K_{22}\overline{W'V'} + K_{33}K_{23}\overline{W'W'}\end{aligned}$$

$$\begin{aligned}\overline{W'W'} &= K_{31}K_{31}\overline{U'U'} + K_{31}K_{32}\overline{U'V'} + K_{31}K_{33}\overline{U'W'} \\ &+ K_{32}K_{31}\overline{V'U'} + K_{32}K_{32}\overline{V'V'} + K_{32}K_{33}\overline{V'W'} \\ &+ K_{33}K_{31}\overline{W'U'} + K_{33}K_{32}\overline{W'V'} + K_{33}K_{33}\overline{W'W'}\end{aligned}$$

7. Proofs for Coordinate System Transformation Equations

This section contains proofs demonstrating that we can perform the coordinate system transformations on the reduced averaged data without having to perform the transformation on the instantaneous values. This saves costly run time because there are typically thousands of instantaneous values that contribute one averaged value.

The following equations show how the average velocities from the original coordinate system can be used along with the coordinate transformation matrix to provide velocities in the new transformed to coordinate system:

$$\bar{U} = \frac{\sum_{i=1}^N U_i}{N}$$

$$U_i = K_{11}U_i + K_{12}V_i + K_{13}W_i$$

$$\bar{U} = \frac{\sum_{i=1}^N (K_{11}U_i + K_{12}V_i + K_{13}W_i)}{N}$$

$$\bar{U} = \frac{\sum_{i=1}^N K_{11}U_i}{N} + \frac{\sum_{i=1}^N K_{12}V_i}{N} + \frac{\sum_{i=1}^N K_{13}W_i}{N}$$

$$\bar{U} = K_{11} \frac{\sum_{i=1}^N U_i}{N} + K_{12} \frac{\sum_{i=1}^N V_i}{N} + K_{13} \frac{\sum_{i=1}^N W_i}{N}$$

$$\bar{U} = K_{11}\bar{U} + K_{12}\bar{V} + K_{13}\bar{W}$$

With similar proofs we can show that the following equations apply:

$$\bar{V} = K_{21}\bar{U} + K_{22}\bar{V} + K_{23}\bar{W}$$

$$\bar{W} = K_{31}\bar{U} + K_{32}\bar{V} + K_{33}\bar{W}$$

The following equations show how the velocity:velocity normal and shear stress terms from the original coordinate system can be used along with the coordinate transformation matrix to provide velocity:velocity normal stress terms in the new transformed to coordinate system:

$$\overline{U'U'} = \frac{\sum_{i=1}^N (U_i - \bar{U})(U_i - \bar{U})}{N} = \frac{\sum_{i=1}^N U_i U_i}{N} - \bar{U} \bar{U}$$

$$U_i = K_{11}U_i + K_{12}V_i + K_{13}W_i$$

$$\bar{U} = K_{11}\bar{U} + K_{12}\bar{V} + K_{13}\bar{W}$$

$$\overline{U'U'} = \frac{\sum_{i=1}^N (K_{11}U_i + K_{12}V_i + K_{13}W_i)(K_{11}U_i + K_{12}V_i + K_{13}W_i)}{N} - (K_{11}\bar{U} + K_{12}\bar{V} + K_{13}\bar{W})(K_{11}\bar{U} + K_{12}\bar{V} + K_{13}\bar{W})$$

$$\overline{U'U'} = \frac{\sum_{i=1}^N \begin{pmatrix} K_{11}K_{11}U_iU_i + K_{11}K_{12}U_iV_i + K_{11}K_{13}U_iW_i \\ + K_{12}K_{11}V_iU_i + K_{12}K_{12}V_iV_i + K_{12}K_{13}V_iW_i \\ + K_{13}K_{11}W_iU_i + K_{13}K_{12}W_iV_i + K_{13}K_{13}W_iW_i \end{pmatrix}}{N} - \begin{pmatrix} K_{11}K_{11}\bar{U}\bar{U} + K_{11}K_{12}\bar{U}\bar{V} + K_{11}K_{13}\bar{U}\bar{W} \\ + K_{12}K_{11}\bar{V}\bar{U} + K_{12}K_{12}\bar{V}\bar{V} + K_{12}K_{13}\bar{V}\bar{W} \\ + K_{13}K_{11}\bar{W}\bar{U} + K_{13}K_{12}\bar{W}\bar{V} + K_{13}K_{13}\bar{W}\bar{W} \end{pmatrix}$$

$$\begin{aligned} \overline{U'U'} &= K_{11}K_{11} \left(\frac{\sum_{i=1}^N U_i U_i}{N} - \bar{U} \bar{U} \right) + K_{11}K_{12} \left(\frac{\sum_{i=1}^N U_i V_i}{N} - \bar{U} \bar{V} \right) + K_{11}K_{13} \left(\frac{\sum_{i=1}^N U_i W_i}{N} - \bar{U} \bar{W} \right) \\ &+ K_{12}K_{11} \left(\frac{\sum_{i=1}^N V_i U_i}{N} - \bar{V} \bar{U} \right) + K_{12}K_{12} \left(\frac{\sum_{i=1}^N V_i V_i}{N} - \bar{V} \bar{V} \right) + K_{12}K_{13} \left(\frac{\sum_{i=1}^N V_i W_i}{N} - \bar{V} \bar{W} \right) \\ &+ K_{13}K_{11} \left(\frac{\sum_{i=1}^N W_i U_i}{N} - \bar{W} \bar{U} \right) + K_{13}K_{12} \left(\frac{\sum_{i=1}^N W_i V_i}{N} - \bar{W} \bar{V} \right) + K_{13}K_{13} \left(\frac{\sum_{i=1}^N W_i W_i}{N} - \bar{W} \bar{W} \right) \end{aligned}$$

$$\begin{aligned} \overline{U'U'} &= K_{11}K_{11}\overline{U'U'} + K_{11}K_{12}\overline{U'V'} + K_{11}K_{13}\overline{U'W'} \\ &+ K_{12}K_{11}\overline{V'U'} + K_{12}K_{12}\overline{V'V'} + K_{12}K_{13}\overline{V'W'} \\ &+ K_{13}K_{11}\overline{W'U'} + K_{13}K_{12}\overline{W'V'} + K_{13}K_{13}\overline{W'W'} \end{aligned}$$

The following equations show how the velocity:velocity normal and shear stress terms from the original coordinate system can be used along with the coordinate transformation matrix to provide velocity:velocity shear stress terms in the new transformed to coordinate system:

$$\overline{U'V'} = \frac{\sum_{i=1}^N (U_i - \bar{U})(V_i - \bar{V})}{N} = \frac{\sum_{i=1}^N U_i V_i}{N} - \bar{U} \bar{V}$$

$$U_i = K_{11}U_i + K_{12}V_i + K_{13}W_i$$

$$\bar{U} = K_{11}\bar{U} + K_{12}\bar{V} + K_{13}\bar{W}$$

$$V_i = K_{21}U_i + K_{22}V_i + K_{23}W_i$$

$$\bar{V} = K_{21}\bar{U} + K_{22}\bar{V} + K_{23}\bar{W}$$

$$\overline{U'V'} = \frac{\sum_{i=1}^N (K_{11}U_i + K_{12}V_i + K_{13}W_i)(K_{21}U_i + K_{22}V_i + K_{23}W_i)}{N} - (K_{11}\bar{U} + K_{12}\bar{V} + K_{13}\bar{W})(K_{21}\bar{U} + K_{22}\bar{V} + K_{23}\bar{W})$$

$$\overline{U'V'} = \frac{\sum_{i=1}^N \left(\begin{array}{l} K_{11}K_{21}U_iU_i + K_{11}K_{22}U_iV_i + K_{11}K_{23}U_iW_i \\ + K_{12}K_{21}V_iU_i + K_{12}K_{22}V_iV_i + K_{12}K_{23}V_iW_i \\ + K_{13}K_{21}W_iU_i + K_{13}K_{22}W_iV_i + K_{13}K_{23}W_iW_i \end{array} \right)}{N} - \left(\begin{array}{l} K_{11}K_{21}\bar{U}\bar{U} + K_{11}K_{22}\bar{U}\bar{V} + K_{11}K_{23}\bar{U}\bar{W} \\ + K_{12}K_{21}\bar{V}\bar{U} + K_{12}K_{22}\bar{V}\bar{V} + K_{12}K_{23}\bar{V}\bar{W} \\ + K_{13}K_{21}\bar{W}\bar{U} + K_{13}K_{22}\bar{W}\bar{V} + K_{13}K_{23}\bar{W}\bar{W} \end{array} \right)$$

$$\begin{aligned} \overline{U'V'} &= K_{11}K_{21} \left(\frac{\sum_{i=1}^N U_iU_i}{N} - \bar{U}\bar{U} \right) + K_{11}K_{22} \left(\frac{\sum_{i=1}^N U_iV_i}{N} - \bar{U}\bar{V} \right) + K_{11}K_{23} \left(\frac{\sum_{i=1}^N U_iW_i}{N} - \bar{U}\bar{W} \right) \\ &+ K_{12}K_{21} \left(\frac{\sum_{i=1}^N V_iU_i}{N} - \bar{V}\bar{U} \right) + K_{12}K_{22} \left(\frac{\sum_{i=1}^N V_iV_i}{N} - \bar{V}\bar{V} \right) + K_{12}K_{23} \left(\frac{\sum_{i=1}^N V_iW_i}{N} - \bar{V}\bar{W} \right) \\ &+ K_{13}K_{21} \left(\frac{\sum_{i=1}^N W_iU_i}{N} - \bar{W}\bar{U} \right) + K_{13}K_{22} \left(\frac{\sum_{i=1}^N W_iV_i}{N} - \bar{W}\bar{V} \right) + K_{13}K_{23} \left(\frac{\sum_{i=1}^N W_iW_i}{N} - \bar{W}\bar{W} \right) \end{aligned}$$

$$\begin{aligned} \overline{U'V'} &= K_{11}K_{21}\overline{U'U'} + K_{11}K_{22}\overline{U'V'} + K_{11}K_{23}\overline{U'W'} \\ &+ K_{12}K_{21}\overline{V'U'} + K_{12}K_{22}\overline{V'V'} + K_{12}K_{23}\overline{V'W'} \\ &+ K_{13}K_{21}\overline{W'U'} + K_{13}K_{22}\overline{W'V'} + K_{13}K_{23}\overline{W'W'} \end{aligned}$$

The following equations show how the velocity:voltage shear stress terms from the original coordinate system can be used along with the coordinate transformation matrix to provide velocity:voltage shear stress terms in the new transformed to coordinate system:

$$\overline{U'A'} = \frac{\sum_{i=1}^N (U_i - \bar{U})(A_i - \bar{A})}{N} = \frac{\sum_{i=1}^N U_i A_i}{N} - \bar{U} \bar{A}$$

$$U_i = K_{11}U_i + K_{12}V_i + K_{13}W_i$$

$$\bar{U} = K_{11}\bar{U} + K_{12}\bar{V} + K_{13}\bar{W}$$

$$A_i = A_i$$

$$\bar{A} = \bar{A}$$

$$\overline{U'A'} = \frac{\sum_{i=1}^N (K_{11}U_i + K_{12}V_i + K_{13}W_i)(A_i)}{N} - (K_{11}\bar{U} + K_{12}\bar{V} + K_{13}\bar{W})(\bar{A})$$

$$\overline{U'A'} = \frac{\sum_{i=1}^N (K_{11}U_i A_i + K_{12}V_i A_i + K_{13}W_i A_i)}{N} - (K_{11}\bar{U} \bar{A} + K_{12}\bar{V} \bar{A} + K_{13}\bar{W} \bar{A})$$

$$\overline{U'A'} = K_{11} \left(\frac{\sum_{i=1}^N U_i A_i}{N} - \bar{U} \bar{A} \right) + K_{12} \left(\frac{\sum_{i=1}^N V_i A_i}{N} - \bar{V} \bar{A} \right) + K_{13} \left(\frac{\sum_{i=1}^N W_i A_i}{N} - \bar{W} \bar{A} \right)$$

$$\overline{U'A'} = K_{11}\overline{U'A'} + K_{12}\overline{V'A'} + K_{13}\overline{W'A'}$$

8. Matrix Notation for Coordinate System Transformation Equations

This section shows how the equations of section 1.6 can be represented in matrix notation. The matrix notation for the coordinate system transformation is an elegant way to show the multitude of complex equations in compact and concise format. The rest of this page contains various matrix definitions:

$$\mathbf{J}_{3 \times 3} = \begin{bmatrix} J_{11} & J_{12} & J_{13} \\ J_{21} & J_{22} & J_{23} \\ J_{31} & J_{32} & J_{33} \end{bmatrix} \quad \mathbf{K}_{3 \times 3} = \begin{bmatrix} K_{11} & K_{12} & K_{13} \\ K_{21} & K_{22} & K_{23} \\ K_{31} & K_{32} & K_{33} \end{bmatrix}$$

$$\mathbf{S} = \begin{bmatrix} \bar{U} \\ \bar{V} \\ \bar{W} \end{bmatrix} \quad \mathbf{S}_i = \begin{bmatrix} U_i \\ V_i \\ W_i \end{bmatrix} \quad \mathbf{F} = \begin{bmatrix} \frac{\bar{U}'A'}{\bar{V}'A'} \\ \frac{\bar{V}'A'}{\bar{W}'A'} \end{bmatrix} \quad \mathbf{G} = \begin{bmatrix} \frac{\bar{U}'B'}{\bar{V}'B'} \\ \frac{\bar{V}'B'}{\bar{W}'B'} \end{bmatrix}$$

$$\mathbf{R} = \begin{bmatrix} \bar{U} \\ \bar{V} \\ \bar{W} \end{bmatrix} \quad \mathbf{R}_i = \begin{bmatrix} U_i \\ V_i \\ W_i \end{bmatrix} \quad \mathbf{H} = \begin{bmatrix} \frac{\bar{U}'A'}{\bar{V}'A'} \\ \frac{\bar{V}'A'}{\bar{W}'A'} \end{bmatrix} \quad \mathbf{I} = \begin{bmatrix} \frac{\bar{U}'B'}{\bar{V}'B'} \\ \frac{\bar{V}'B'}{\bar{W}'B'} \end{bmatrix}$$

$$\mathbf{P} = \begin{bmatrix} \frac{\bar{U}'U'}{\bar{U}'V'} \\ \frac{\bar{U}'V'}{\bar{U}'W'} \\ \frac{\bar{U}'W'}{\bar{V}'U'} \\ \frac{\bar{V}'U'}{\bar{V}'V'} \\ \frac{\bar{V}'V'}{\bar{V}'W'} \\ \frac{\bar{V}'W'}{\bar{W}'U'} \\ \frac{\bar{W}'U'}{\bar{W}'V'} \\ \frac{\bar{W}'V'}{\bar{W}'W'} \end{bmatrix} \quad \mathbf{Q} = \begin{bmatrix} \frac{\bar{U}'U'}{\bar{U}'V'} \\ \frac{\bar{U}'V'}{\bar{U}'W'} \\ \frac{\bar{U}'W'}{\bar{V}'U'} \\ \frac{\bar{V}'U'}{\bar{V}'V'} \\ \frac{\bar{V}'V'}{\bar{V}'W'} \\ \frac{\bar{V}'W'}{\bar{W}'U'} \\ \frac{\bar{W}'U'}{\bar{W}'V'} \\ \frac{\bar{W}'V'}{\bar{W}'W'} \end{bmatrix}$$

$$\mathbf{K}_{9 \times 9} = \begin{bmatrix} K_{11}K_{11} & K_{11}K_{12} & K_{11}K_{13} & K_{12}K_{11} & K_{12}K_{12} & K_{12}K_{13} & K_{13}K_{11} & K_{13}K_{12} & K_{13}K_{13} \\ K_{11}K_{21} & K_{11}K_{22} & K_{11}K_{23} & K_{12}K_{21} & K_{12}K_{22} & K_{12}K_{23} & K_{13}K_{21} & K_{13}K_{22} & K_{13}K_{23} \\ K_{11}K_{31} & K_{11}K_{32} & K_{11}K_{33} & K_{12}K_{31} & K_{12}K_{32} & K_{12}K_{33} & K_{13}K_{31} & K_{13}K_{32} & K_{13}K_{33} \\ K_{21}K_{11} & K_{21}K_{12} & K_{21}K_{13} & K_{22}K_{11} & K_{22}K_{12} & K_{22}K_{13} & K_{23}K_{11} & K_{23}K_{12} & K_{23}K_{13} \\ K_{21}K_{21} & K_{21}K_{22} & K_{21}K_{23} & K_{22}K_{21} & K_{22}K_{22} & K_{22}K_{23} & K_{23}K_{21} & K_{23}K_{22} & K_{23}K_{23} \\ K_{21}K_{31} & K_{21}K_{32} & K_{21}K_{33} & K_{22}K_{31} & K_{22}K_{32} & K_{22}K_{33} & K_{23}K_{31} & K_{23}K_{32} & K_{23}K_{33} \\ K_{31}K_{11} & K_{31}K_{12} & K_{31}K_{13} & K_{32}K_{11} & K_{32}K_{12} & K_{32}K_{13} & K_{33}K_{11} & K_{33}K_{12} & K_{33}K_{13} \\ K_{31}K_{21} & K_{31}K_{22} & K_{31}K_{23} & K_{32}K_{21} & K_{32}K_{22} & K_{32}K_{23} & K_{33}K_{21} & K_{33}K_{22} & K_{33}K_{23} \\ K_{31}K_{31} & K_{31}K_{32} & K_{31}K_{33} & K_{32}K_{31} & K_{32}K_{32} & K_{32}K_{33} & K_{33}K_{31} & K_{33}K_{32} & K_{33}K_{33} \end{bmatrix}$$

This page consolidates all of the coordinate transformation equations in matrix notation.

$$\mathbf{S} = \mathbf{K}_{3 \times 3} \times \mathbf{R}$$

$$\begin{bmatrix} \bar{U} \\ \bar{V} \\ \bar{W} \end{bmatrix} = \begin{bmatrix} K_{11} & K_{12} & K_{13} \\ K_{21} & K_{22} & K_{23} \\ K_{31} & K_{32} & K_{33} \end{bmatrix} \times \begin{bmatrix} \bar{U} \\ \bar{V} \\ \bar{W} \end{bmatrix}$$

$$\mathbf{S}_i = \mathbf{K}_{3 \times 3} \times \mathbf{R}_i$$

$$\begin{bmatrix} U_i \\ V_i \\ W_i \end{bmatrix} = \begin{bmatrix} K_{11} & K_{12} & K_{13} \\ K_{21} & K_{22} & K_{23} \\ K_{31} & K_{32} & K_{33} \end{bmatrix} \times \begin{bmatrix} U_i \\ V_i \\ W_i \end{bmatrix}$$

$$\mathbf{H} = \mathbf{K}_{3 \times 3} \times \mathbf{F}$$

$$\begin{bmatrix} \bar{U}'A' \\ \bar{V}'A' \\ \bar{W}'A' \end{bmatrix} = \begin{bmatrix} K_{11} & K_{12} & K_{13} \\ K_{21} & K_{22} & K_{23} \\ K_{31} & K_{32} & K_{33} \end{bmatrix} \times \begin{bmatrix} \bar{U}'A' \\ \bar{V}'A' \\ \bar{W}'A' \end{bmatrix}$$

$$\mathbf{I} = \mathbf{K}_{3 \times 3} \times \mathbf{G}$$

$$\begin{bmatrix} \bar{U}'B' \\ \bar{V}'B' \\ \bar{W}'B' \end{bmatrix} = \begin{bmatrix} K_{11} & K_{12} & K_{13} \\ K_{21} & K_{22} & K_{23} \\ K_{31} & K_{32} & K_{33} \end{bmatrix} \times \begin{bmatrix} \bar{U}'B' \\ \bar{V}'B' \\ \bar{W}'B' \end{bmatrix}$$

$$\mathbf{Q} = \mathbf{K}_{9 \times 9} \times \mathbf{P}$$

$$\begin{bmatrix} \bar{U}'U' \\ \bar{U}'V' \\ \bar{U}'W' \\ \bar{V}'U' \\ \bar{V}'V' \\ \bar{V}'W' \\ \bar{W}'U' \\ \bar{W}'V' \\ \bar{W}'W' \end{bmatrix} = \begin{bmatrix} K_{11}K_{11} & K_{11}K_{12} & K_{11}K_{13} & K_{12}K_{11} & K_{12}K_{12} & K_{12}K_{13} & K_{13}K_{11} & K_{13}K_{12} & K_{13}K_{13} \\ K_{11}K_{21} & K_{11}K_{22} & K_{11}K_{23} & K_{12}K_{21} & K_{12}K_{22} & K_{12}K_{23} & K_{13}K_{21} & K_{13}K_{22} & K_{13}K_{23} \\ K_{11}K_{31} & K_{11}K_{32} & K_{11}K_{33} & K_{12}K_{31} & K_{12}K_{32} & K_{12}K_{33} & K_{13}K_{31} & K_{13}K_{32} & K_{13}K_{33} \\ K_{21}K_{11} & K_{21}K_{12} & K_{21}K_{13} & K_{22}K_{11} & K_{22}K_{12} & K_{22}K_{13} & K_{23}K_{11} & K_{23}K_{12} & K_{23}K_{13} \\ K_{21}K_{21} & K_{21}K_{22} & K_{21}K_{23} & K_{22}K_{21} & K_{22}K_{22} & K_{22}K_{23} & K_{23}K_{21} & K_{23}K_{22} & K_{23}K_{23} \\ K_{21}K_{31} & K_{21}K_{32} & K_{21}K_{33} & K_{22}K_{31} & K_{22}K_{32} & K_{22}K_{33} & K_{23}K_{31} & K_{23}K_{32} & K_{23}K_{33} \\ K_{31}K_{11} & K_{31}K_{12} & K_{31}K_{13} & K_{32}K_{11} & K_{32}K_{12} & K_{32}K_{13} & K_{33}K_{11} & K_{33}K_{12} & K_{33}K_{13} \\ K_{31}K_{21} & K_{31}K_{22} & K_{31}K_{23} & K_{32}K_{21} & K_{32}K_{22} & K_{32}K_{23} & K_{33}K_{21} & K_{33}K_{22} & K_{33}K_{23} \\ K_{31}K_{31} & K_{31}K_{32} & K_{31}K_{33} & K_{32}K_{31} & K_{32}K_{32} & K_{32}K_{33} & K_{33}K_{31} & K_{33}K_{32} & K_{33}K_{33} \end{bmatrix} \times \begin{bmatrix} \bar{U}'U' \\ \bar{U}'V' \\ \bar{U}'W' \\ \bar{V}'U' \\ \bar{V}'V' \\ \bar{V}'W' \\ \bar{W}'U' \\ \bar{W}'V' \\ \bar{W}'W' \end{bmatrix}$$

